


Instructions for Use

CHEETAH[®] XPLORE







FSX009x







		3
EN	INSTRUCTIONS FOR USE	5
DE	GEBRAUCHSANWEISUNG	10
FR	NOTICE D'UTILISATION	16
ES	INSTRUCCIONES PARA EL USO	21
IT	ISTRUZIONI PER L'USO	26
NO	BRUKSANVISNING	31
DA	BRUGSANVISNING	36
SV	BRUKSANVISNING	41
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	46
FI	KÄYTTÖOHJEET	52
NL	GEBRUIKSAANWIJZING	57
PT	INSTRUÇÕES DE UTILIZAÇÃO	62
PL	INSTRUKCJA UŻYTKOWANIA	67
CS	NÁVOD K POUŽITÍ	72
TR	KULLANIM TALIMATLARI	77
RU	ИНСТРУКЦИЯ ПО ИСПОЛЬЗОВАНИЮ	82
JA	取扱説明書	88
ZH	中文说明书	93
KO	사용 설명서	98

DIRECT LAMINATION

ENGLISH

MEDICAL CONDITIONS AND CLINICAL CLAIMS

This product is intended to compensate for loss/limitation of mobility due to a lost lower limb. It is designed and tested for everyday use while also allowing the user to partake in various sport activities.

It provides increased balance confidence throughout an active lifestyle while also providing stability in standing.

The device should only be fitted and aligned by a qualified health care professional.

INDICATION FOR USE

The foot is a non-invasive, single patient, reusable prosthetic device. The Cheetah Xplore is tested and validated for exoskeletal prosthetic fittings of lower extremities amputations.

TARGET POPULATION

Suitable for moderate to high active users
Maximum patient weight: 147 kg / 324 lbs
Minimum patient weight: 45 kg / 99 lbs

This product is designed and tested based on single patient use, and is not recommended for multiple patient use.

PRODUCT PARAMETERS (Figure 1)

- Weight: 646 g.
- Maximum height of the posterior brim: 521 mm
- Minimum height of the posterior brim: 348 mm
- Minimum clearance from the socket to the foot: 32 mm
- Foot size: 24-30

CONTRAINDICATIONS:

N/A.

ALIGNMENT (Figure 1)

BENCH ALIGNMENT

Sagittal plane

- Determine vertical socket position in alignment jig.
- Adjust to appropriate heel height.
- Introduce appropriate socket angles flexion.
- The load line should fall at the junction of the posterior and the middle third portions as shown in **(Figure 1)**.

Coronal plane

- Introduce appropriate socket abduction /adduction using alignment jig.
- Determine appropriate internal or external socket rotation.

NOTE: The recommendations for bench alignment are used as starting point.

CATEGORY SELECTION CHART

Please refer to the selection charts below to determine the appropriate stiffness required according to Össur recommendations.

WEIGHT KG	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130	131-147
WEIGHT LBS	99-115	116-130	131-150	151-170	171-194	195-221	222-256	257-287	288-324
Moderate Impact Level	1	2	3	4	5	6	7	8	9
High Impact Level	2	3	4	5	6	7	8	9	N/A

- Grind foot to length, with help from the grinding sticker placed on the top of the foot plate (**Figure 2**).

NOTE: When using a foot cover the size of the foot needs to be 2 cm shorter when grinding the toe to length. This requirement has to be taken into account when grinding the sole plate.

CAUTION: Wear dust mask and gloves when grinding and gluing the device.

LAMINATION

The foot may be connected with either lamination connector or direct lamination.

LAMINATION CONNECTOR (Figure 3)

- Lamination connector is laminated directly into the posterior socket wall.
- Placement of the lamination connector is recommended in the middle 1/3 of the socket for easier height adjustment.
- Make sure the socket surface is flat over the lamination connector to ensure good load distribution on the foot module when connected. Bolts go through blade into lamination connector in socket. Apply Loctite 410 to all threads and torque to following specifications:

Category	Bolt size	Torque
1-9	9,5 mm (3/8")	45 Nm (33 ft-lbs)

DIRECT LAMINATION

Before you start you need a well fitting conventional prosthesis. Either an existing socket or alternatively a prosthesis with a check socket on conventional prosthesis can be used.

1. Laminate the socket (Figure 4)

- a. Laminate the socket using carbon braid. Add extra layers of carbon in the upper posterior surface where the foot will be attached.
- b. Roughen the upper posterior surface to improve the adhesion to the socket when the foot module is attached (shown by grey area).

2. Trace around the original walking prosthesis (Figure 5)

- a. Ensure that the prosthetic foot is resting firmly on the shoes insole and fasten the laces tightly.
- b. Lay the original walking prosthesis over the grid on the A1 alignment sheet (included)
- c. Mark the lowest points of the anterior and posterior brims as a reference.
- d. With the shoe on, trace a line around the entire prosthesis. Then remove the shoe and trace the outline of the foot cover, ensuring the prosthesis remains aligned with the original tracing.

3. Fill in the patient information table on the alignment sheet (Figure 6)

- a. The socket adduction can be observed by looking at the original walking prosthesis in the frontal plane.
- b. Use the angle guide on the A1 alignment sheet to help measure the foot rotation in relation to the socket.

4. Attach Sole and Heel counter (Figure 7)

- a. Glue the heel counter onto the upper surface of the heel, aligning the center line with the split in the carbon heel.
- b. Glue the sole onto the underside of the foot, again aligning the center line with the splits in the carbon fiber.
- c. Grind the sole and heel into a shape that gives a good fit in the patient's shoe. Use the sole template on the back of the product box.

5. Cut the pylon to length (Figure 8)

- a. Align the new socket and the foot with the tracing of the original prosthesis.
- b. Mark the pylon at the height of the posterior brim of the socket using the drawing of the original prosthesis.
- c. Cut the pylon at this mark and round off the end, removing any sharp edges.

6. Mark the length of the attachment surface on the foot pylon (Figure 9)

- a. Mark the pylon at 90mm from the top (cut end).
- b. Roughen all the surfaces of the pylon above this mark (shown by grey area).

7. Align the foot and the socket with the tracing of the original walking prosthesis (Figure 10)

- a. The sole should be lined up with the outline of the foot cover sole.
- b. The new socket should be positioned over the tracing of the original socket with additional 5-7° of knee extension in the socket.
- c. Use a foam spacer/s (or other spacer) placed between the pylon and the socket to maintain the correct extension. Spacers should be placed at the end of the roughened attachment surface. Ensure the minimum clearance of 5 mm is maintained between the distal posterior area of the socket and the foot pylon.
- d. Any spacer used at the top of the pylon can have maximum thickness of 5 mm.
- e. Fasten the foam to the foot using tape.
- f. Mark a clear height on the socket at height of a foam spacer on the pylon.

8. Align the socket according to the patient information (Figure 11)

- a. Put the shoe on the foot and place it upright in the work surface.
- b. Hold the socket in position on the pylon according to:
 - The height
 - The recorded socket adduction
 - The recorded foot rotation in relation to socket (using angle guides on A1 sheet)
- c. The foam spacers will hold the socket in the correct extension.

9. Mark the attachment surface onto the socket (Figure 12)

- a. Holding the socket in the desired position, trace around the end of the pylon onto the socket above the height mark.

10. Fix the foot to the socket (Figure 13)

- a. Use Composite 1 adhesive that is provided with the foot kit, read all information on the adhesive packaging.
NOTE: Adhesive working time is 1 minute.
- b. Cut the tip of the adhesive mixer at 16 mm from the end.
- c. Apply the adhesive generously onto the marked attachment surface on the socket.
- d. Attach the foot onto the marked attachment area and make any small adjustment to ensure the alignment is correct.
- e. Smooth the surface of any adhesive that has been pushed out around the outside of the pylon.
- f. Hold the socket in position for 1 minute, then leave for a further 5 minutes for the adhesive to set.

11. Apply adhesive around the outside of the pylon (Figure 14)

- a. Apply adhesive generously onto the outside roughed surface of the pylon and spread outwards onto the roughened surface of the socket.
- b. Smooth the finish of the adhesive and allow to set for 5 minutes.

12. Test alignment before final lamination (Figure 15)

- a. Remove any foam spacers and tape before testing the alignment on the patient.
- b. If needed, minor adjustments can be made by shaping the underside of the sole.
- c. If the alignment needs further adjustment, grind/cut the foot off the socket and make the desired changes before repeating steps 10 and 11.
- d. If the alignment is good, proceed to laminate over the attachment area onto the socket using carbon braid and a hand lay up.

NOTE: The mixing sprout for the adhesive is a single use. When reusing the adhesive you need to replace it with a new mixing sprout.

WATERPROOF

A waterproof device allows use in a wet and/or humid environment and temporary (<30 mins.) submersion of up to 1m in fresh water. Exposure to fresh water under these circumstances shall have no harmful effect. Dry thoroughly after contact with fresh water. Fresh water: Includes tap water. Excludes salt and chlorinated water.

The components used with the foot need to be waterproof as well.

CAUTION: If the user notices any change of function or alignment of the device, he should contact his healthcare practitioner.

USE ENVIRONMENT


- Environmental operating temperature: -10°C (+50°F) to 40° C (104°F).
- Operating humidity: 10%-95% Relative humidity.
- Shipping and storage temperature: -20°C (-4°F) to 70°C (158° F).

LIABILITY

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

COMPLIANCE

This component has been tested according to ISO 10328 standard to two million load cycles. Depending on the amputee's activity this corresponds to a duration of use of two to three years. We recommend carrying out regular yearly safety checks

In the standard mentioned, test levels (P) are assigned to a certain maximal body masses (m in kg). In some cases, which are marked with, no test level is assigned to the product related maximal body mass. In these cases, the test loads

have been adapted adequately on the basis of the specified load level.

Category Össur high activity						
Category	Weight (Kg)	Lable text				
1	52	P3	ISO 10328 -	P3	52	kg
2	59	P3	ISO 10328 -	P3	59	kg
3	68	P4	ISO 10328 -	P4	68	kg
4	77	P4	ISO 10328 -	P4	77	kg
5	88	P4	ISO 10328 -	P5	88	kg
6	100	P5	ISO 10328 -	P5	100	kg
7	116	P5	ISO 10328 -	P6	116	kg
8	130	P6	ISO 10328 -	P6	130	kg
9	147	P6	ISO 10328 -	P6	147	kg

- This manual is intended for use by a certified prosthetist.

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EN – Caution: Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standard is achieved only when Össur products are used with other recommended Össur components. If un-usual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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