





Instructions for Use





---

ICEFORM<sup>®</sup> LINER

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	9
FR	Notice d'utilisation	14
ES	Instrucciones para el uso	18
PT	Instruções de Utilização	22

## INTRODUCTION

Thank you for choosing a high quality liner from Össur. We can assure you that you have made the right choice. Please carefully review and follow these instructions before fitting a patient with the liner or accessory products. Remember that informing the patient concerning proper handling and care of the liner, as set forth in these instructions, is important to achieve patient comfort and maximum product life.

## INDICATION FOR USE

The ICEFORM liner is a suspension device for low to moderate active below knee amputees, available as locking or cushion liner, to suspend the residual limb to the prosthetic leg.

## CONTRAINDICATIONS FOR USE

Össur does not recommend prescribing an ICEFORM liner for amputees that are not capable to use the liner correctly according to instructions, for example due to cognitive disability.

**Note** : Össur does not recommend prescribing an ICEFORM liner for amputees that do not have sufficient hand strength to don the liner properly and do not have caregivers to assist with donning.

## PRODUCT SELECTION AND FITTING

### *Liner sizing*

To determine the correct ICEFORM size:

1. Measure the circumference of the residual limb at 4 and 30 cm from the distal end of the stump (**Figure 1**).
2. Choose the ICEFORM size, which corresponds to the measurement, or the one which is immediately below (**Figure 2**).
3. Never choose a liner size, which exceeds the measurement, or round up.

**Example** : If the residual limb measures 10 inches (26cm) at the indicated distal location and more than 14 inches (36cm) at the indicated proximal location, than choose ICEFORM size M+.

Selecting the correct size is very important. If the liner is too tight, it may cause pistoning and/or numbness and blistering. If the liner is too loose, increased perspiration and/or movement of the limb inside the liner may occur. Perspiration and movement can lead to blisters and rashes. If any of the above symptoms or any other indication of improper fit are noted, the amputee should be instructed to contact his/her physician or prosthetist immediately.

## DISTAL ATTACHMENTS

Your ICEFORM locking liner will have fixed distal attachment, note that cushion liners have no distal attachments.

ICEFORM cushion liners should be used with the ICEFORM sleeve.

## DISTAL ATTACHMENT LOCKING LINERS (**Figure 3**)

When using an ICEFORM liner with a fixed distal attachment and an attachment pin, caution must be taken to properly tighten the attachment pin to the fixed distal attachment. Always apply Loctite

410/411 or equivalent compound to the threads of the Össur attachment pin (used with Icelock) and torque to 4 Nm (3ft.lbs). Össur does not recommend the use of pins without a flange with ICEFORM locking liners (**Figure 3**).

If excessive movement or product wear is detected between the attachment pin and the distal attachment at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her physician or prosthetist.

#### **TRIMMING AN ICEFORM LINER (Figure 4)**

ICEFORM liners may be trimmed for amputee preference and to increase range of motion. The ICEROSS CurveMaster is specifically designed for this function and to produce a smooth, rounded edge, which will minimize the risk of skin irritation and/or tears of the liner.

**Warning :** Do not cut the liner below the proximal trim-lines of the socket. Excessive trimming may reduce the natural suction between the liner and limb and/or compromise suspension.

#### **PROSTHETIC SOCKS**

If the amputee experiences insecurity or decreased contact with the socket due to volume fluctuations one or more ICEROSS Socks can be added to restore the socket fit.

#### **CUSTOM SHAPE CAPTURE**

##### ***Shape Capture***

Due to the material properties of the liner the unique shape of the limb will be automatically captured during normal use. To accelerate shape capture “*Cold capture*” or “*Hot capture*” techniques can be utilized.

Both techniques require a cast/scan of the residual limb with no liner and a subsequent positive model.

##### ***Cold Capture***

1. Add a protective layer to the clean positive model.
2. Apply ICEFORM liner to positive model.
3. Leave for a minimum time of one hour.

##### ***Hot Capture***

1. Add a protective layer to the clean dry positive model.
2. Place positive model into preheated oven at 65°C (150°F) for 60 minutes.
3. Apply ICEFORM liner to positive model and place into the oven.
4. After maximum of 10 minutes remove from oven.
5. Allow to cool.

**Note :** Make sure the positive model and protective layer is made out of material that can withstand 65°C (150°F).

**Note :** The ICEFORM liner must not be put into oven unless it is on a positive model.

**Note :** Direct physical contact with heating elements will permanently damage the ICEFORM liner and void the warranty.

**Note** : Use heat protective gloves at all times during this process.

## **DONNING THE ICEFORM LINER**

**Note** : Do NOT apply lotion to the residual limb immediately prior to donning the liner. The residual limb should be clean and dry, without any residue of lotions since that can damage the liner. Cover any open wounds or damaged skin with a bandage or some other appropriate covering to prevent direct contact between the wound and the liner.

1. The liner should be oriented with the logo on the front of the residuum (**Figure 5**).
2. Grip the top of the liner from the inside and slide it over the hand until the inner surface is fully exposed, taking care not to damage it with fingernails (**Figure 6**).
3. Be very careful that the inside of the liner is clean, dry and free from any foreign objects that can cause skin irritation.
4. After exposing as much of the distal end of the liner as possible, position it against the residual limb and with light compression roll upward onto the limb (**Figure 7**).
5. Roll the liner all the way up the limb, taking care not to damage it with fingernails (**Figure 8**). Do not tug or pull. Check that no air pockets are present.
6. For locking liners it is important to align the pin on the distal end properly such that it does not cause discomfort when the socket has been donned (**Figure 9**).

**Warning** : Failure to adhere to the donning instructions may compromise comfort and/or suspension.

## **DAILY CARE**

Proper product care is important. The liner needs to be washed every day, following use. Remove the liner, turn the liner inside out (**Figure 6**) and wash it with pH balanced, 100% fragrance- and dye-free soap. After cleaning, rinse the liner thoroughly with luke warm water (**Figure 10**) and pat dry on both sides with a lint-free cloth.

**Note** : Excessive heat may damage the liner. Take care not to scrub the surface of the liner, since that can roughen the surface which can cause skin irritation.

An ICEFORM liner can be used directly after cleaning. Always return the liner to its neutral state, with the textile cover facing out, as soon as possible after cleaning and do not store it inside-out since that can damage the liner.

ICEFORM liners should always be checked for damage or wear. Any damage may weaken the effectiveness of the liner and should be reported to the amputee's prosthetist immediately. Misuse may result in potential loss of suspension.

Daily cleaning of the residual limb is also essential. We recommend use of a mild liquid, pH balanced, 100% fragrance- and dye-free soap. If dry skin is noted, apply a pH balanced, 100% fragrance- and dye-free lotion to nourish and soften the skin.

**Note** : Care must be taken by both the prosthetist and the amputee not to expose a liner to glass or carbon fibers or other foreign particles. Such substances can become embedded in the liner material causing aggravation of the skin. Washing alone may not be sufficient to eliminate the problem. The liner should be returned to the prosthetist if it is inadvertently exposed to foreign substances or chemicals.

**Note** : Many common household or bath products, including soaps, deodorants, perfumes, aerosol or alcohol sprays or abrasive cleaners may cause or contribute to skin irritation and the patient should be so advised.

#### **WARRANTY**

Össur offers a 6 months limited warranty on the ICEFORM Liner for the indicated period from the invoice date.

#### **PERSPIRATION**

If a use of a liner results in excessive perspiration such that stability or suspension is compromised, it is recommended to remove the liner and dry off both residual limb and liner.

#### **STORAGE**

Store the liner with the fabric cover facing out in a cool and dry place. For extended periods of storage put liner in a clean plastic bag.

For further information, please contact Össur or your local Össur distributor. You can also visit the Össur web page at [www.ossur.com](http://www.ossur.com).

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**EN – Caution:** Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standard is achieved only when Össur products are used with other recommended Össur components. If un-usual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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