
Instructions for Use

K2 SENSATION®

Product number:

K2P0xyz(BR)

K2PUxyz(BR)

		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	6
FR	Notice d'utilisation	9
ES	Instrucciones para el uso	11
IT	Istruzioni per l'uso	14
DA	Brugsanvisning	17
SV	Bruksanvisning	19
EL	Οδηγίες χρήσης	21
NL	Gebruiksaanwijzing	24
PT	Instruções de Utilização	27
JA	取扱説明書	30
ZH	中文说明书	32
KO	사용 설명서	34

1


21 - 22 cm - 56 mm (23/16")
23 - 28 cm - 62 mm (27/16")

10 mm (3/8")


90°

BENCH ALIGNMENT


2


3

4

5

ENGLISH

K2 SENSATION® FOOT SYSTEM (Figure 1)

Build height:

- Pyramid adapter 56-62 mm (2 3/16" - 2 7/16").
- Maximum user weight: 136 kg (299 lbs).
- Available in 10 mm (3/8") heel height.
- Available in sizes 21-28 cm.

NOTE: Select foot category according to impact level.

BENCH ALIGNMENT (Figure 2 & 3)

- Fit foot with selected cover.
- Adjust to appropriate heel height (using the shoe).
- Introduce appropriate socket angles flexion/ extension and abduction/adduction.
- Divide the foot cover length into 3 equal portions.
- The load line should fall at the junction of the posterior and the middle third portions as shown in the diagram.

DYNAMIC ALIGNMENT

The heel stores energy after initial contact, slowly releasing it at midstance. The heel to toe action can be influenced by:

- A-P positioning of foot
- Dorsi-Plantarflexion

CATEGORY SELECTION CHART

Please refer to the selection charts below to determine the appropriate stiffness required according to Össur recommendations.

WEIGHT KG	36-44	45-61	62-80	81-97	98-114	115-136
WEIGHT LBS	79-97	98-134	135-176	177-214	215-251	252-300
21	0	1	N/A	N/A	N/A	N/A
22	0	1	N/A	N/A	N/A	N/A
23	0	1	2	N/A	N/A	N/A
24	0	1	2	N/A	N/A	N/A
25	N/A	1	2	3	4	5
26	N/A	1	2	3	4	5
27	N/A	N/A	2	3	4	5
28	N/A	N/A	2	3	4	5

HEEL STIFFNESS

Heel too soft

Symptoms

- Foot will come to flat position too early (amputee feels he/she is sinking into a hole).
- Extra energy is required to climb up over the toe.
- Toe will feel too stiff.
- Knee may hyperextend.

Actions

- Shift socket anteriorly.
- Dorsi-flex the foot.

Heel too hard

Symptoms

- Rapid heel to toe movement.
- Amputee has poor control of his prosthesis.
- Toe will feel too stiff - minimal energy return.
- Knee may become unstable.

Action

- Shift socket posteriorly.
- Plantar-flex the foot.

PROTECTIVE FILM ON PYRAMID

The protective film protects the pyramid from scratching. Remove the film when the fitting is complete.

SHOE HORN

When removing and installing the cover use shoehorn to avoid damage to foot or cover (Figure 4 & 5).

LIABILITY

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

COMPLIANCE

This component has been tested according to ISO 10328 standard to two million load cycles. Depending on the amputee’s activity this corresponds to a duration of use of two to three years. We recommend carrying out regular yearly safety checks

ISO 10328 - “P” - “m”kg *)

*) Body mass limit not to be exceeded!
For specific conditions and limitations of use see
manufacturer’s written instructions on intended use!

In the standard mentioned, test levels (P) are assigned to a certain maximal body masses (m in kg). In some cases, which are marked with, no test level is assigned to the product related maximal body mass.

In these cases, the test loads have been adapted adequately on the basis of the specified load level.

ISO 10328 - Lable		
Category	Weight (Kg)	Label text
0	44	ISO 10328–P3–44–kg
1	61	ISO 10328–P3–61–kg
2	80	ISO 10328–P4–80–kg
3	97	ISO 10328–P5–97–kg
4	114	ISO 10328–P5–114–kg
5	136	ISO 10328–P6–136–kg

- This manual is intended for use by a certified prosthetist.

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EN – Caution: Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standard is achieved only when Össur products are used with other recommended Össur components. If un-usual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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