
DK	Transfer Belt	2
GB/US	Transfer Belt	5
DE	Transfer Brustgurt	9
SE	Förflyttnings bälte	12
FR	Transfer Belt	15
IT	Cintura di trasferimento	18
NO	Transfer Belt	21
BR	Bolsa de transferência	24

Manual – Version 3.00

## DK... TRANSFER BELT

Vers. 3.00

Varenr.:  
284141 S  
284151 M  
284161 L

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Vers. 3.00

Item nos:  
284141 S  
284151 M  
284161 L

<b>1.00</b>	<b>Purpose and use</b> .....	<b>5</b>
1.01	Manufacturer .....	5
1.02	Intended purpose .....	5
1.03	Area of use .....	5
1.04	Conditions of use .....	5
1.05	Important/Precautions .....	5
1.06	UHF RFID tag .....	5
1.07	Use .....	5
<b>2.00</b>	<b>Maintenance</b> .....	<b>6</b>
2.01	Cleaning .....	6
2.02	The owner's daily maintenance duty .....	6
2.03	Disposal of belts .....	6
<b>3.00</b>	<b>Service and lifetime</b> .....	<b>6</b>
3.01	Safety/service inspections .....	6
3.02	Lifetime .....	7
<b>4.00</b>	<b>Technical specifications</b> .....	<b>7</b>
<b>5.00</b>	<b>EC-Declaration of conformity</b> .....	<b>7</b>
<b>6.00</b>	<b>Environmental policy statement – V. Guldmann A/S</b> .	<b>7</b>
<b>7.00</b>	<b>Warranty and service conditions</b> .....	<b>7</b>
A.	Warranty .....	7
B.	Service or Repair .....	8
	<b>Placing the belt</b> .....	<b>27</b>

### 1.00 Purpose and use

#### 1.01 Manufacturer

V. Guldmann A/S  
Graham Bells Vej 21-23A  
DK-8200 Aarhus N  
Tel. + 45 8741 3100  
www.guldmann.com

#### 1.02 Intended purpose

The belt is intended for lifting or supporting a person or body parts of a person.

#### 1.03 Area of use

The belt is used with the GTP1 Transfer Platform to activate users from sitting to standing position at hospitals, nursing homes, institutions and individual housing.

The belt is intended to support a standing person on the GTP1. The user must actively participate and also be able to communicate well.

### 1.04 Conditions of use

The use of the belt is subject to the following:

- The user can balance in a sitting position and has strength in the lower and upper body and can communicate
- The belt is used by trained staff or persons who have been instructed in the use of the belt in question.
- The correct size of belt is used.
- The maximum nominal load of 205 kg must not be exceeded.
- The helper pays attention to the well-being of the user when using the belt.
- The belt is used in connection with the GTP1 Transfer Platform from Guldmann.

### 1.05 Important/Precautions

- Read the instructions carefully before using the belt.
- The belt's maximum load must never be exceeded.
- The belt may only be used to lift a person.
- Before a belt is used, it must be examined according to point 2.02.
- Never use a belt that is too big for the user.
- Possible repairs must only be made by the manufacturer.
- Any serious incident that occurred in relation to this device should be reported to the manufacturer and the local competent authority.

### 1.06 UHF RFID tag

This product features a passive UHF RFID tag integrated in the product label. The RFID tag can be used for asset management & tracking purposes.

The RFID tag is readable with equipment compliant to EPC global UHF Class 1 Gen ISO 18000-63.

### 1.07 Use

*If there is any doubt about the selection or use of a belt, please contact your supplier.*

#### Important!

**Plan the move. Never leave the user standing in the belt unattended. Before the user stands up be sure to check that the user cannot get trapped and that the belt will not get tangled on the bed, wheelchair or other things.**

Be careful with any tubes or wires attached to the user and / or equipment.

Guldmann shall not be liable for faults or accidents due to incorrect use of the belt, or for reasons of inadequate attention on the part of the carer or user. If the belt is used in combination with products that are not manufactured by Guldmann, a risk assessment must be made by qualified staff.

#### Caution!

**Be careful when attaching the straps of the belt on the Transfer Platform. Check that the straps have been correctly placed in the hook.**

*Applying the belt, see page 27.*

## 2.00 Maintenance

### 2.01 Cleaning

Normal washing at the indicated temperature

Do not use bleaching agent

Tumble-drying at low temperature

Do not iron

### 2.02 The owner's daily maintenance duty

Check the belt for wear and damage before use according to the following checklist which is not intended to represent all potential inspection steps. Potential damage may vary. Judgment of inspector/site prevails.

#### Belt inspection checklist

Before using a Guldmann belt / sling / accessory check the following:

#### Is the belt clean?

Follow facility specific infection control procedure.

#### Is the belt's label present, legible and complete?

Missing, illegible or incomplete belt label(s) could make identification of appropriate size of the belt, function of belt, and or weight limit capacity of the belt impossible.

#### Are the lifting straps and stitches intact?

- Look for broken or worn stitches
- Look for knots in straps
- Look for tears or fraying of straps
- Look for snags or punctures or holes
- Look for any particles in fabric or straps

#### Is the fabric intact?

- Look for abnormal wear patterns, excessive wear, abrasive evidence
- Look for cuts or frayed fabric
- Look for unusual or significant discoloration
- Look for snags, punctures, tears, holes
- Look for frayed or insecure seams
- Look for any acid / caustic / thermal burns
- Look for changes in material consistency, e.g. increased stiffness
- Look for any imbedded particles

#### Has the shape of the belt been altered, made shorter or longer in relation to the original size using knots, needles, tape or other methods?

#### Conclusion

If the belt suffers from one or more of the above mentioned conditions then it must be taken out of service regardless of the weight of the person to be lifted/supported.

### 2.03 Disposal of belts

Belts are disposed of by incineration. By proper incineration polyester will be degraded to carbon dioxide and water.

## 3.00 Service and lifetime

### 3.01 Safety/service inspections

In accordance with international standard EN/ISO 10535 "Hoist for the transfer of disabled persons – Requirements and test methods" an inspection *must* be performed every 6-month ac-

ording to the following instructions, which is not intended to represent all potential inspection steps. Potential damage may vary. Judgment of inspector/site prevails.

#### Safe Operating Practices with Belts

Considerations for damaged or defective belts and taking them out of service:

Withdraw the belt from service if one or more of the following conditions are present:

- chemical or caustic burns
- melting or charring of any part of the belt
- snags, punctures, tears or cuts
- broken or worn stitches
- missing, illegible or incomplete tag
- knots in any part of the belt
- abrasion
- other visible damage that causes doubt as to the strength of the belt

Belt inspection is done for the protection of the user, the caregiver, and the overall hospital site safety. A belt inspection system has additional benefit. Systematic belt inspection will assist in the identification of damage trends, potentially leading to cost effective suggestions and results. The inspection process can also help to identify inventory duplicity in certain belt types and sizes.

#### Belt inspection system

Development of a specific procedure and program for the inspection of belts at your facility is your best safeguard. Consider employing a three part system of inspection. Belts that are removed from service and are not capable of repair should be disposed of so they are unfit for any future use and can not find a way back into active inventory.

#### 1) Initial

This level of inspection is done at the time that the belt is received into your facility. The inspector should ensure that no damage has occurred during transit, and also verify that the belt work load limits match those contained in the manufacturer's catalogue. If your facility documents the belt inspection process through written inspection records, the paper trail should begin at this stage.

#### 2) Frequent

The frequent level of inspection should be done by the belt user before each use. The belt should be examined and removed from service if damage is detected. The belt user should also determine that the belt is proper for the user conditions, care task required and the required weight capacity.

#### 3) Periodic

Your facility might want to consider implementing a program for a periodic level of inspection at regular intervals. The interval should be based upon the frequency of use, severity of the service cycle and information derived through the inspection process. Recommendations to prevent damage and enhance service life could be made by staff that perform the periodic inspections. If written inspection records are maintained, they should always reference the unique belt identification number, and be updated to record the condition of the belt. Not intended to represent all potential inspection steps or all potential aspects of product management program. Judgment of inspector/site prevails.

#### Belt inspection technique

The belt inspection procedure should be thorough, systematic and consistent; both visual and "hands on" inspection techniques are recommended. Certain forms of damage are far more discernable through hands-on inspection, than by visual inspection. For example, fabric stiffness, crushed webbing, as well as, thinning fabric can be identified through tactile inspection. Visual inspection alone may not reveal all forms of belt damage. Once signs of damage have been identified, do not downgrade the work load limit of the belt, with the intent of continuing to use it, but at limited capacity or frequency. This is sometimes done to get more service life out of a

damaged belt. The operating rule and standard should be: intact = use; damage = do not use.

Consider the practice of documenting belt inspections through written inspection records. The documentation should include information such as: the name of manufacturer, the belt stock number, width and length, the unique belt identification number (important in differentiating similar belts), as well as the condition of the belt.

Other important information might also include the date it was received or put into use at your facility and any special features (if applicable). A beneficial outcome of an inspection program would be the realization of repetitive forms of damage and the analysis that would lead to specific recommendations.

**Sample visual examples of synthetic belt damage \*)**


*Chemical/caustic burns*

*Broken stitching*


*Crushed / Frayed webbing*

*Knots*

*Melting / Charring*

\*) sample visual images not intended to represent all types of potential damage

**3.02 Lifetime**

The life of the belt is individual and depends on how it is used, washed etc. Before use the belt must be examined according to description in section 2.02 and if it does not meet the inspection requirements, it must be discarded if necessary.

**4.00 Technical specifications**

Lifting capacity, SWL . . . . . 175 kg / 385 lbs  
 Material . . . . . Polyester

**5.00 EC-Declaration of conformity**

The product is manufactured in compliance with regulation (EU) 2017/745 of the European parliament and of the Council of 5 April 2017, as medical device class 1.

**6.00 Environmental policy statement – V. Guldmann A/S**

Guldmann is continuously working towards ensuring that the company's impact on the environment, locally and globally, is reduced to a minimum.

It is Guldmann's goal to:

- Comply with the current environmental legislation (e.g. WEEE and REACH directives)
- Ensure that we, at the widest possible range, use RoHS compliant materials and components
- Ensure that our products do not have an unnecessary negative impact on the environment regarding use, recirculation or disposal
- Ensure that our products contribute to a positive working environment in the places they are utilised

Inspections are made annually by the Department for Nature and Environment from the Municipality of Aarhus using the Danish Environmental Protection Act, section 42 as a reference.

**7.00 Warranty and service conditions**

**A. Warranty**

Guldmann warrants its equipment is free from material defects under normal use, and will perform substantially in accordance with the specifications set forth in documentation provided with the equipment.

This express warranty shall be in effect for one year from the date of original purchase and installation (the "Warranty Period"). If a valid claim is made during the Warranty Period for malfunction or equipment defect, Guldmann will repair or replace the equipment at no additional cost to you. Guldmann retains sole discretion as to whether the equipment will be repaired or replaced.

The warranty does not cover any part of the equipment that has been subject to damage or abuse by the user or others. The warranty does not cover any part of the equipment that has been altered or changed in any way by the user or others. Guldmann does not warrant that the lifting device functions will meet your requirements, be uninterrupted or error free.

The warranty set forth is in lieu of all other express and implied warranties, whether oral, written or implied, and the remedies set forth above are your sole and exclusive remedies. Only an authorized officer of Guldmann may make modifications to this warranty, or additional warranties binding on Guldmann. Accordingly, additional statements such as advertising or presentations, whether oral or written, do not constitute warranties by Guldmann.

This warranty shall be null and void if the equipment is operated and maintained in any manner inconsistent with its intended use or the instructions provided with the product. Further, in order for the warranty to remain in effect for the full Warranty Period, all service to the equipment must be provided by a Guldmann certified technician. Any parts or components repaired or replaced by a Guldmann certified technician will be guaranteed for the remainder of the Warranty Period.

**Only for USA**

This warranty shall be null and void if the equipment is operated and maintained in any manner inconsistent with its intended use or the instructions provided with the product. Further, in order for the warranty to remain in effect for the full Warranty Period, all service to the equipment must be provided by a Guldmann Certified Technician. A Guldmann Certified Technician is a technician who has successfully completed Guldmann Service Training, and who holds a valid Service Training Certificate from Guldmann, and is in possession of a valid password to access Guldmann's Service and Information Console (SIC). A Guldmann Service Training Certificate and SIC password are valid for three years (only USA) from the date the technician is first certified. Thereafter, the techni-

cian must undergo re-certification training to obtain a new valid certificate and password. Any parts or components repaired or replaced by a Guldmann Certified Technician will be guaranteed for the remainder of the Warranty Period. In the event the warranty is rendered null and void, the purchaser shall indemnify and hold Guldmann harmless of and from any and all claims or liability arising as a result of equipment malfunction or misuse.

---

**B. Service or Repair**

Contact Guldmann Repair for an authorization to return any defective item during the Warranty Period. You will be provided with a return authorization number and address for returning the item for warranty service or replacement. Do not return items to Guldmann under warranty without receiving a Return Authorization Number.

If mailing the item, pack it carefully in a sturdy carton to prevent damage. Include your Return Authorization Number, a brief description of the problem and your return address and phone number. Guldmann does not assume the risk of loss or damage while in transit, so it is recommended you insure the package.

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**Guldmann GmbH**
**Guldmann Sverige AB**
**Guldmann Srl**
**Guldmann Inc.**