
Instruction for Use

---

OP4 KNEE

		5
EN	Instructions for Use	7
DE	Gebrauchsanweisung	11
FR	Notice d'utilisation	15
ES	Instrucciones para el uso	19
IT	Istruzioni per l'uso	23
DA	Brugsanvisning	27
NL	Gebruiksaanwijzing	31
PT	Instruções de Utilização	35
RU	Инструкция по использованию	39

## **ALIGNMENT**

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




## PRODUCT DESCRIPTION

The OP4 Knee is a monocentric knee joint with a stable aluminium frame construction and a brake. It is used by moderately active users, mobility 2 + (3), and is approved for a maximum body weight of 100 kg. The swing phase is reliably controlled by the powerful pneumatic system. This enables users with reduced muscle activity to benefit from this control.

## INDICATIONS FOR USE

The OP4 Knee is designed for patients with transfemoral and knee disarticulation amputation.

## CONTRAINDICATIONS FOR USE

N/A.

## TECHNICAL DATA

Material	AL/CA
Amputation level	Above knee
Maximum body weight	100 kg
Total fitted height	approx. 198 mm
Effective fitted height	approx. 146 mm
Fitted height	approx. 26 mm
Weight	approx. 680 g
Prox. connection	Pyramid
Dist. connection	Integrated 30 mm tube receiver
Swing phase	Pneumatic
Stance phase	Brake (load-dependent)
Axes	Monocentric
Flexion angle without shaft	approx. 145°

## ADJUSTMENT

The knee joint OP4 Knee is delivered in full working order in the factory setting.

Carefully start to test the function of the joint with the user at the dynamic trial fitting and check whether this setting suits the user's needs. If not, please go through the following steps.

### *Stance phase stability*

#### **Adjusting the load regulating screw L (13)**

The load regulating screw L (13) is for adjusting the brake under load i.e. as weight is exerted. Use a 5 mm hexagonal key for adjustment.

When setting up the joint individually with fine adjustments turn the screw in small steps of not more than about ¼ of a turn each time. Check the result.

#### **Brake effect inadequate**

- Turn the load regulating screw L (13) to the left
- The brake responds at a lower load

**Brake effect too strong/Brake sticks at the end of the stance phase**

- Turn the load regulating screw L (13) to the right
- The brake responds at a higher load

**Note:** The load regulating screw L (13) must not be screwed beyond the brake body otherwise the function of the knee joint can not be guaranteed.

**Readjusting the brake unit after prolonged use**

The service screw S (29) need not be adjusted under standard conditions. After a prolonged period of use, play may develop in the joint which can be readadjusted using the service screw S (29). To adjust the screw remove the protective cap (30) and use a 3 mm hexagonal key.

**When setting up the joint with fine adjustments turn the screw in small steps of a maximum of 1/8 of a turn each time. Check the result.**

**Joint develops play after prolonged period of use**

- Turn the service screw S (29) to the right
- Play in the joint is reduced




**Swing phase control**

You can set the extension and flexion separately to suit the needs of the user.

The pneumatic system is delivered in the basic setting. In the basic setting both valves are two turns open. You can make the adjustments from this basic setting to suit the individual gait.

In order to achieve a comfortable gait first adjust the flexion. Then harmonise the gait by adjusting the extension.

**When setting up the joint with fine adjustments turn the screw in small steps of about 1/4 of a turn each time. Check the result.**

	<p><b>Too much flexion</b> Close flexion valve F Turn the flexion valve to the right with a 3 mm hexagonal key.</p>
	<p><b>Too little flexion</b> Open flexion valve F Turn the flexion valve to the left with a 3 mm hexagonal key.</p>
	<p><b>Abrupt extension stop</b> Close extension valve E Turn the extension valve to the right with a 3 mm hexagonal key.</p>
	<p><b>Extension too slow</b> Open extension valve E Turn the extension valve to the left with a 3 mm hexagonal key.</p>

**Note:** Flexion and extension must be possible at all settings.

**Note:** It may be necessary to readjust the brake again after setting the swing phase.

### **Resetting to the basic setting**

In order to restore the joint to the basic setting please go through the following steps:

#### **Basic setting of the brake**

Close the load regulating screw L(13) clockwise until you feel a slight resistance. Open the load regulating screw 3 full turns anticlockwise to restore the basic setting.

#### **Basic setting of the pneumatic system**

Close both valves (extension valve/flexion valve) clockwise until you feel a slight resistance. Open both valves by 2 full turns anticlockwise to restore the basic setting.

### **SAFETY WARNINGS**

- Adapters used on the upper connection must always be tightened onto the respective modular part to the specified torque. Use a torque wrench.
- Adapters (e.g. tube adapters) used on the lower connection must be inserted down to the end stop in the tube receiver of the knee joint.
- Tighten the screw (7) on the tube receiver to a torque of 16 Nm and secure with Loctite. Use a torque wrench (5 mm Allen).
- Do not loosen any screws apart from those described above.
- Do not use any lubricants (e.g. talc, silicone spray, oil) otherwise the brake function can not be guaranteed.
- The load regulating screw L (13) must not be screwed beyond the brake body.

**Observe all the above points, otherwise the guarantee becomes null and void.**

### **MAINTENANCE**

The joint must be serviced at least every 6 months!

Inspect:

- the alignment
- the screw connections
- the suitability of the patient (e.g. max. body weight, degree of mobility)
- loss of lubricant
- damage to the joint and the adapter
- soiling of the bushes
- joint play AP + ML
- extension stop (12)
- soiling with talc

### **CARE**

Clean the joint with a soft cloth moistened with a little mild benzene. Do not use any more aggressive cleaning agents because these can damage the seals and bushes.

Do not use compressed air for cleaning. Compressed air can force dirt into the seals and bushes. This can lead to premature damage and wear.

### **LIABILITY**

In case of damage: Össur can only consider complaints accompanied by a copy of the delivery note or the Össur invoice together with a detailed description of the reasons for returning the product.

A manufacturer can only be held liable for the failure of his own fittings. The manufacturer can only be held liable beyond this when it can be proved that his fittings were causally responsible for the damage to or loss of function of fittings from other manufacturers.

### **CE MARKING/COMPATIBILITY**

All the modular components have been tested in compliance with EC Directive 93/42/EEC and bear the CE marking.

### **General conditions**

This medical aid is recommended for use by one patient only. Product liability pursuant to the medical Devices Act becomes null and void if it is used for treating more than one patient. Össur products are guaranteed compatible with modular components from others if the following points are observed:

The modular components are used only in compliance with their intended purpose (see Instructions for Use).

- If modular components with different max. weight limits are used, the weight limit of the weakest component applies.
- The use of tested individual components with the CE marking does not release the prosthetist from his obligation to check as well as possible the combination of fittings for their suitability, correct assembly and safety.
- Should there be any indications that a given combination of fittings may not assure the required degree of safety, the fittings must not be combined.
- The prostheses must be manufactured in compliance with the generally acknowledged professional rules of the orthopaedic appliances trade.
- If the patient's suitability for wearing the prosthesis changes (e.g. weight, activity, etc.) the prosthesis must be checked without delay.
- If the prosthesis has been exposed to unusually high stress, e.g. a fall, it must be inspected in a prosthetic workshop immediately for possible damage.
- Safety-relevant precautions for individual fittings (e.g. maintenance intervals, see Instructions for Use) must be observed.

### **Combination of different fittings**

- The above general conditions apply in addition to:
- Only use fittings that comply with the requirements of ISO 10328
- Only use fittings that comply with the requirements of the medical Devices Act
- The fittings must be suitable for each other and fit the standard pyramid or 30 mm system.

### **USEFUL LIFE**

The maximum useful life of the Össur OP4 is 36 months.

### **REPLACEMENT PARTS/ACCESSORIES**

No.	Name	Article No.
7	Screw for tube receiver	1.970.952
12	Extension stop OP4 Knee (new version starting with serial number 0712xxxx)	1.970.954
13	Load regulating screw OP4 Knee	1.970.953
14	Spring	1.970.951

Only the listed replacement parts and accessories are available.  
To have a joint inspected please send the joint to Össur for service.  
We will provide you with a service joint for the duration of the repair.

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**EN – Caution:** Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standards is achieved only when Össur products are used with other recommended Össur components. If an unusual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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