SEGUFX[®]-Hand Restraint



Item 2204 / Item 2304 / Item 2404

Consisting of: 2 SEGUFIX[®]-Hand Straps 2 SEGUFIX[®]-Locks (key not included)

Component part: item 2104 (for one hand, locks and keys not included)



Item 2224 / Item 2324 / Item 2424

Consisting of:

2 SEGUFIX[®]-Hand Straps extended

4 SEGUFIX[®]-Locks

1 SEGUFIX[®]-Key

Sizes: S/M/L/XL

Component part: item 2124 (for one hand, locks and keys not included)

Sizes: S/M/L/XL

Versions: basic version / r





Consisting of: 2 SEGUFIX[®]-Hand Straps

extended Velcro® (with buckle lock) Component part: item 2824 (for one hand) Size: M



The SEGUFIX[®]-Hand Restraint (yellow label) is part of the SEGUFIX®-Bed Restraint System and made from a cotton/ravon mix.

3 Versions> SEGUFIX [®] -Locks:				
		Je Jo		
1206 / 1208	1306 / 1308	1406 / 1408		
Item 2204	ltem 2304	Item 2404		
Item 2224	Item 2324	Item 2424		

Table of Contents	Page
SEGUFIX®-Hand Restraint with lock – Instructions	2 - 4
SEGUFIX®-Hand Restraint extended Velcro® – Instructions	4 - 5
Overview table: all versions, sizes, measurements and weights	5
SEGUFIX [®] -Locks and Keys	6 - 7
SEGUFIX®-Safety Information	8 - 11
Instructions for Care, Imprint	12

Versions – Item no. / identification			
Item no. / Version Product features			
	Basic version		
/r Eyelet strap reinforced for higher demands			

3 Locking Systems: Item 2204 / 2304 / 2404 Item 2224 / 2324 / 2424

Intended use:

- Adjustable restraint of one or both hands in bed
- Physically restraining patients (e.g. pathology or disability that requires restraint in a defined position in bed to restrict patients from harming themselves or care givers, e.g. anxiety disorders)

With Velcro[®] fastener / Buckle lock **Item 2724**

Intended use:

- Positioning of one or both hands in bed
- Also compatible with ambulance stretchers
- Not a restraint (note: some jurisdictions consider Velcro[®] products as a restraint, should the patient not be able to remove them without help)

Item 2204 / 2304 / 2404 SEGUEIX[®]-Hand Restraint

SEGUFIX^{*}-Hand Restraint



Item 2224 / 2324 / 2424 SEGUFIX®-Hand Restraint extended



Item 2724 SEGUFIX®-Hand Restraint extended Velcro®



Instructions for use item 2204 / 2304 / 2404 and 2104

Locks and keys see pages 6 and 7

Important: SEGUFIX[®]-Standard required. Install the SEGUFIX-Standard on the bed and fasten it to the patient. For detailed instructions please refer to the respective instructions.

How to attach the Hand Strap to the wide bed belt of the SEGUFIX[®]-Standard:

- Insert the pin of the lock (1) through the eyelet of the pin flap on the padded cuff of the hand strap (A).
- On the wide bed belt of the SEGUFIX[®]-Standard is a loop (B). Insert the padded cuff of the hand strap (A) through that loop, with the padding facing up and pointing towards the patient. Pull the cuff through, until the loop (B) is positioned between the metal buckle (2) and the pin of the lock (1) of the hand strap.



How to fasten the Hand Strap to the wrist:

- Wrap the padded cuff of the hand strap (A) snugly around the patient's wrist. The Velcro® fastener on the cuff prevents slipping after closing.
- Insert the eyelet strap (C) through the metal buckle (2) and pull it tightly back over the metal buckle.
- Lay the evelet strap (C) tightly with an evelet over the pin of the lock (1) and close the lock.
- Secure the evelet (4) to the lock (3) of the attachment strap of the Standard.
- Fasten the second hand strap in the same manner if needed.
- The cuff of the hand strap should be snug, but should not restrict blood circulation.
- Pull on each lock to test for a secure hold. Then lock bed side rails in an up position.

Alternatives for attaching the Hand Strap:

- Wrap the padded cuff of the hand strap (A) snugly around the patient's wrist, without using the loop (B) on the wide bed belt of the SEGUFIX®-Standard.
- Using the evelet (4), the hand strap can now be attached to any lock. This allows for the hand or arm to be positioned according to indication.

How to release the restraint from the patient:

See page 7 on how to open the respective locking systems.

- Open the locks (1, 3 and 5) of the hand strap and the waist belt, and pull the straps of the pins.
- Remove the hand strap from the patient and waist belt.

Instructions for use item 2224 / 2324 / 2424 and 2124

Locks and keys see pages 6 and 7

SEGUFIX[®]-Standard is **not required**.

How to attach the SEGUFIX[®]-Hand Restraint extended:

- Insert the pin of the lock (1 and 3, see page 4) through the respective eyelet of the pin flap.
- Fasten the padded cuff of the hand strap as shown on page 2 (fig. 2-5, see above), without using the loop of the SEGUFIX[®]-Standard.











- Direct the extended attachment strap (A) from inside out around the bed frame or the frame of the folding part of the bed. Attach the extended attachment strap in such manner that it cannot slide along the bed frame.
- Lay one of the eyelets at the end of the attachment strap over the pin of the lock (3). The many eyelets allow for gradually adjusting the restraint.
- Close the lock (3).
- Fasten the second hand strap in the same manner if needed.
- The cuff of the hand strap should be snug but should not restrict blood circulation.
- Pull on each lock to test for a secure hold. Then lock bed side rails in an up position.

How to release the restraint from the patient and bed:

See page 7 on how to open the respective locking systems.

- Open the locks (1 and 3) of the hand strap and pull the strap from the pin.
- Remove hand strap from patient and bed.

Follow SEGUFIX[®]-Safety Information and Warnings on pages 8 - 11

Instructions for use item 2724


SEGUFIX[®]-Standard is **not required**.

How to attach the SEGUFIX[®]-Hand Restraint extended Velcro[®]:

Note: not a restraint. It is possible for the patient to open the Velcro $^{\circ}$ fastener and buckle lock.

- Wrap the padded cuff with the padded part facing up snugly around the patient's wrist. The Velcro[®] fastener on the cuff prevents slipping after closing.
- Insert the strap (A) of the hand strap through the both metal loops of the buckle lock (1), back through the first loop, and pull tight, so that the buckle lock is fastened as in fig. 2.
- Direct the extended attachment strap (B) from inside out around the bed frame or the frame of the folding part of the bed. Attach the extended attachment strap in such manner that it cannot slide along the bed frame.
- For closing insert the end of the extended attachment strap (B) through the buckle lock (2) and pull the strap tight.
- Fasten the second hand strap in the same manner if needed.
- The cuff of the hand strap should be snug, but should not restrict blood circulation.
- Pull on each buckle lock to test for a secure hold. Then lock bed side rails in an up position.







How to release the hand strap from the patient and the bed:

- Undo the buckle lock (1) and remove the hand strap from the patient.
- Undo the buckle lock (2) of the extended attachment strap and remove the strap from the patient.

Follow SEGUFIX[®]-Safety Information and Warnings on pages 8 - 11

WARNING: Read the Safety Information prior to using the SEGUFIX[®]-Hand Restraint.

- The SEGUFIX[®]-Hand Restraint may only be used in connection with the SEGUFIX[®]-Standard. The extended versions of the Hand Restraint have longer straps for direct attachment to the bed frame, and may be used without the SEGUFIX[®]-Standard.
- The SEGUFIX[®]-Hand Restraint extended with Velcro[®] fastener is not a restraint. It is possible for the patient to open the Velcro[®] fastener and buckle lock. Note: some jurisdictions consider Velcro[®] products as a restraint, should the patient not be able to remove them without help.
- Note that attachment straps and moving mechanics of the bed or mattress may interfere with each other, causing damage, loss of function and injury.
- For more product info visit us online at www.segufix.com or www.segufix-shop.com

SEGUFIX [®] -Hand Restraint					Width of strap		Length of strap					
Size	Item	kg	Item	kg	Item	kg	Compo- nent part	kg	mm	inches	mm	inches
S	2204	0,33	2304	0,37	2404	0,37	2104	0,15	45	1.8	620	24.4
М	2204	0,38	2304	0,42	2404	0,42	2104	0,17	45	1.8	770	30.1
L	2204	0,41	2304	0,45	2404	0,45	2104	0,18	45	1.8	810	31.9
XL	2204	0,44	2304	0,48	2404	0,48	2104	0,21	45	1.8	1100	43.3
S	2204/r	0,36	2304/r	0,40	2404/r	0,40	2104/r	0,16	45	1.8	620	24.4
М	2204/r	0,40	2304/r	0,44	2404/r	0,44	2104/r	0,18	45	1.8	770	30.1
L	2204/r	0,42	2304/r	0,46	2404/r	0,46	2104/r	0,19	45	1.8	810	31.9
XL	2204/r	0,45	2304/r	0,49	2404/r	0,49	2104/r	0,21	45	1.8	1100	43.3
S	2224	0,62	2324	0,70	2424	0,70	2124	0,27	45	1.8	1360	53.5
М	2224	0,67	2324	0,75	2424	0,75	2124	0,29	45	1.8	1360	53.5
L	2224	0,72	2324	0,80	2424	0,80	2124	0,30	45	1.8	1490	58.7
S	2224/r	0,66	2324/r	0,75	2424/r	0,75	2124/r	0,28	45	1.8	1360	53.5
М	2224/r	0,71	2324/r	0,79	2424/r	0,79	2124/r	0,30	45	1.8	1360	53.5
L	2224/r	0,77	2324/r	0,85	2424/r	0,85	2124/r	0,33	45	1.8	1490	58.7
XL	2224/r	0,83	2324/r	0,91	2424/r	0,91	2124/r	0,36	45	1.8	1640	64.6
М	2724	0,48	-	-	-	-	2824	0,24	45	1.8	1420	55.9

Wrist girth				
Size	mm	inches		
S	120 - 190	4.7 - 7.5		
М	160 - 240	6.3 - 9.5		
L	220 - 290	8.7 - 11.4		
XL	260 - 340	10.2 - 13.4		

Tables: approximate measurements – about 8 % shrinkage possible (products are not pre-washed)

SEGUF X[®]-Locks and Keys

∎ �	SEGUFIX [®] -Magnetic locking system black		
Item	Consisting of:		
1206	 Magnetic Button black / segufix.com Metal Pin 		
1208	1 Magnetic Key red		
1204	1 Magnetic Button black / segufix.com		
1205	1 Metal Pin		
1209	1 Magnetic Key red with key ring		

Push locking system

Consisting of:

1 Push Button

1 Push Lock Kev

1 Push Button

1 Push Pin

1 Push Pin

Item

1406

1408

1404

1405

	SEGUFIX [®] -Patent Twist locking system green		
Item	Consisting of:		
1306	 Patent Twist Button green Square Pin 		
1308	1 Patent Twist Key yellow		
1304	1 Patent Twist Button green		
1305	1 Square Pin		

Intended use of the Lock: Closing of the SEGUFIX[®] systems.

Intended use of the Key: Quick release of the locks.

Caution: Each locking system only works with its own component parts and is not interchangeable with other locking systems.

CE

WARNING: Read the Safety Information prior to using the SEGUFIX[®]-Locks and Keys.

- Do not submerge locks in any type of liquid. This can cause malfunction.
- Locks and keys must not be laundered.
- Test locks for proper function before each use.
- Examine the head of the pin for damage (wear and tear, rounding off or deformation).
- When closing the lock make sure that the button locks securely onto the pin. You should not be able to pull the button off of the pin.
- When opening locks, ensure movement between the upper and lower points of contact (item 1206). Therefore each lock can hold up to a maximum of 4 belt layers (exception: products with "r" or "f" in the item number have thicker belts. These allow for a maximum of 3 belt layers!). When using the Patent Twist Lock (item no. 1306) use no more then 4 belt layers. When using the Push Lock (item no. 1406) use no more then 3 belt layers.
- If damaged locks fail to open, straps can be cut with scissors or belt cutters.
- For more product info visit us online at www.segufix.com or www.segufix-shop.com

Instructions for use: SEGUFIX[®]-Locks and Keys

Item 1206 / 1208 (SEGUFIX[®]-Magnetic locking system black) Item 1306 / 1308 (SEGUFIX[®]-Patent Twist locking system green) Item 1406 / 1408 (Push locking system)

How to open the Locks:

Item 1206 / 1208

- Set the red Magnetic Key with a light push on top of the black Magnetic Lock.
- Lift both the Magnetic Key and the Magnetic Button off of the Metal Pin (fig. 1).

Item 1406 / 1408

- Push the two pins of the Push Lock Key into the holes on top of the Push Button to unlock
- Lift both the Push Lock Key and the Push Button off of the Push Pin (fig. 2).

How to close the Locks:

- Set the respective button with the opening onto the pin so that it locks in place.
- Now test the lock for a secure hold: You should not be able to pull the button off of the pin.

Caution: Each locking system only works with its own component parts and is not interchangeable with other locking systems.

Follow SEGUFIX[®]-Safety Information and Warnings on pages 8 - 11

Item 1306 / 1308 Set the vellow Magnetic Twist Key on top of the green Magnetic Twist Lock.

 Unlock the Patent Twist Lock by turning the Patent

Twist Key, and then lift both the Patent Twist Key and the Patent Twist Button off of the Square Pin (fig. 3).





SEGUF X[®]-Safety Information

The following instructions for use are only applicable to original SEGUFIX[®] products. Do not use the SEGUFIX[®]-Restraint and Positioning System until you have read and understood the instructions for use.

Only use up-to-date instructions and keep them in a safe and readily available location. Product and instruction must match. Failure to use this product properly can cause serious injury or death. Your patient's safety depends on your using the SEGUFIX[®] system correctly.

1. Requirements for use

It is the user's responsibility to comply with all national laws. Always consider restraining protocols of your work place and current legislation. The instructions for use must be adhered to.

SEGUFIX[®] products may only be used by individuals with proper training and knowledge of the product. Legal provisions and doctor's orders must be strictly observed! The attending physician must ensure, that the patient's medical condition allows the application of the SEGUFIX[®] system.

2. Intended use

SEGUFIX[®] products may only be used for the medical purposes, which they were designed for.

3. Remove dangerous objects

Remove all objects from the reach of the patient, which could lead to injury or may damage the belts (e.g. glasses, jewelry and sharp objects).

4. Inspect products before use



Inspect before each use for broken stitches, torn, worn out or frayed parts and damaged eyelets. Damaged products could malfunction and cause injury. Do not use damaged or worn out products.

5. When opening, locks may not be under tension

When opening locks, ensure movement between the upper and lower points of contact (item 1206). Therefore each lock can hold up to a maximum of 4 belt layers (exception: products with "r" or "f" in the item number have thicker belts. These allow for a maximum of 3 belt layers!). When using the Patent Twist Lock (item 1306) use no more then 4 belt layers. When using the Push Lock (item 1406) use no more then 3 belt layers.

6. Always attach Side Positioning Straps

The SEGUFIX[®]-Standard must always be used with the Side Positioning Straps attached. Side Positioning Straps prevent the patient from helicoptering in bed or from getting over the edge of the bed causing compression of the chest, possibly resulting in death. In some cases it is possible that patients get over the edge of the bed even with the Side Positioning Straps attached and the side rails in the up position. In these cases additional restraint measures must be taken.

Important: Insert the Side Positioning Straps through the respective loops (according to instructions) before securing them to the bed belt.

7. Beds, chairs and stretchers

- Only use SEGUFIX[®] products on beds, stable chairs and stretchers that allow the product to be attached to, in accordance with the instructions.
- Products must be attached in such manner, that they cannot get out of place.
- Note that straps and moving mechanics of beds, stretchers or chairs (or alternating pressure pads and the like), may interfere with each other causing damage, loss of function and injury.
- The attachment points must be free of sharp edges.
- SEGUFIX[®] products do not prevent beds, chairs or stretchers from tipping over or from unintended movement.

Helpful hints:

European standards for beds: DIN EN 60601-2-52

European standards for wheelchairs: DIN EN 12182 and DIN EN 12183

8. All side rails must be in the up position



To prevent accidents always lock side rails in the up position (exception: 5-point-restraint).

<u>Caution</u>: Do not use with parted side rails except with a gap protector that effectively prevents the patient from slipping through the gap. An open gap could allow the patient to slip through even with the restraint fastened possibly causing compression of the chest.

9. Monitor patients in restraints



Patients in restraints require frequent monitoring. Aggressive, agitated or restless patients and patients in danger of aspiration require constant monitoring.

10. Test closure-systems prior to each use



Check SEGUFIX[®]-Magnetic Locks black/segufix.com, SEGUFIX[®]-Patent Twist Locks green, Push Locks, Snap Locks, Velcro[®] Fasteners and Buckle Locks for a secure hold.

Do not submerge locks in any type of liquid. This may cause malfunction. If damaged locks fail to open, straps can be cut with scissors or belt cutters. Check the head of the pin for signs of wear and tear (rounding off

or deformation).

11. Pacemaker Warning



A distance of 10 cm (4 inches) must be maintained between a patient's pacemaker and magnetic locks and keys to prevent the pacemaker from possible malfunction (alternatively use Push Locks). For possible interference of strong magnets with other patient implants please refer to manufacturer's instructions of these products.

12. Test products for correct attachment

Frequently check products for correct attachment. The straps must be secured to the frame of the bed (or the folding part of the frame of the bed) in such a way, that they cannot slide out of place.

13. Only use original SEGUFIX® products

For safety reasons $\mathsf{SEGUFIX}^{\texttt{B}}$ products must not be altered or be combined with third party products.

14. Use of SEGUFIX® products in vehicles



The SEGUFIX[®]-Transport Belt System does not replace the vehicle's legally required seat belt system. SEGUFIX[®]-Transport Belts may only be used as additional support or for restraint of a patient. Make sure that the patient can be rescued quickly in case of emergency (e.g. clearly visible and handy belt cutters).

15. Secure but not too tight

The restraints must be fastened tightly, but not interfere with breathing or blood circulation. Loosely fastened products may jeopardize the patient's safety. Frequently check products for correct position.

16. Storage

To prolong product life store belts and locking systems in dry and dark conditions at about 20 $^{\circ}\mathrm{C}.$

17. Flammability of SEGUFIX[®]-Bed Restraints



According to DIN EN 71-2 section 4.3 non-flammable. According to DIN 75200 the burning rate is zero. However, products may not come in contact with cigarettes or be exposed to open flame, as this could permanently damage the products, and this could lead to injury.

18. Correct sizing

Choosing the right size is of essential importance. Make sure that the size matches body and waist sizes (refer to sizing tables). Products that are too small or too large compromise the patient's comfort and safety.

19. Disposal

Carton and plastic wrapping can be disposed of in recyclable disposal, defective products in household garbage.

20. Declaration of conformity

For the declaration of conformity please visit: https://www.segufix.com/conformity_en.php

21. Expert recommandation to prevent accidents

(Prof. Dr.-Ing. U. Boenick, TU-Berlin)

- 1. Always use Side Positioning Straps
- 2. Lock non-parted side rails in the up position (exception: five-point-restraint)
- 3. The use of additional restraint elements may be necessary to eliminate bed egress.

22. Ensure prior to each use

- Test locks for a secure hold
- Visual inspection of belts, straps and eyelets
- Ensure sufficient medical indication
- Check legal requirements

Care Instructions

SEGUFIX[®]-Restraints (white) are made of a cotton/rayon mix and are washable in hot wash cycles up to 95 °C. Lower temperature and laundry nets will prolong product life. Avoid methods that might bend metal parts of the system (e.g. laundry line, dry pressing, rotary irons), eyelets may be damaged.

Not pre-washed, up to 8 % shrinkage. Store in a dry place.

Detergents: non-aggressive, no fabric softener, no bleach.

Dryer: delicate cycle at low temperature and preferably in the laundry net.

<u>Disinfection</u>: may be disinfected using chemothermal disinfection washing procedure. Overdosing can damage the product.

• SEGUFIX[®]-Transport Belts (coloured/black) are made of polyester and are washable in warm wash cycles up to 30 °C.

Avoid methods that might bend metal parts of the system (e.g. laundry line, dry pressing, rotary irons), eyelets may be damaged. Store in a dry place.

Detergents: non-aggressive, no fabric softener, no bleach.

Dryer: delicate cycle at low temperature and preferably in the laundry net.

- SEGUFIX[®]-Locks and metal pins: remove all locks and metal pins before washing or drying. Laundering will damage locks and possibly cause them to fail. Locks may be cleaned with a damp cloth. Never submerge locks in any type of liquid.
- Velcro[®] Fastener: fasten before washing and drying to prevent lint buildup. Occasionally use a stiff brush to remove lint from Velcro[®] Fastener. Even with normal use, Velcro[®] material has the tendency to collect lint. This can reduce the grip strength. Before each use, test grip strength of the Velcro[®] Fastener.
- **Durability:** depends on the intensity of use and is prolonged by gentle laundry procedures (detergent/laundry net).

Edition 15 · 2019-06

© by SEGUFIX®-Bandagen, 2019

The general business and delivery terms apply:

https://www.segufix.com/PDFs/AGBs/ Terms_and_Conditions.pdf

Product specifications subject to change. Deviations from pictures in colour and shape, errors and amendments excepted.

All rights reserved. Reprints or distribution using electronic systems/media, in its entirety or in parts, with prior explicit authorization in writing only. All requests concerning this must be directed at:

SEGUFIX[®]-Bandagen

Das Humane System GmbH & Co. KG

Allerbeeksring 33 · 21266 Jesteburg Germany

Phone: +49 (0) 41 83 500-0 Fax: +49 (0) 41 83 500-200

E-mail: export@segufix.de www.segufix.com www.segufix-shop.com