

Turn All 3CT CR

Place the Turn All 3CT CR system in the bed under the existing mattress, connect the pump and turn it on using the power button . Then select a tilt cycle of either 30, 60 or 90 minutes.

The system automatically places the patient in a 30-degree lateral position, alternating between the right and left sides in a cycle of 30, 60 or 90 minutes on each side.

There is a 5-degree counter-lift on the opposite side of the primary lift side, which provides the patient with security and prevents shear.

The system runs continuously until it is switched off.



Product name	Used for mattress size cm	Maximum user weight	Product number
Turn All 3CT CR	85 - 105 x 200	225	70226



1. The button  is used to start the system and is used again to turn off the system.
2. Then select the kip cycle of either ,  or  minutes.