

# EXTREME™ LINER

AKDT | AKDTHD | AKFR | AKFRHD

For Transfemoral and Active Transtibial  
Amputees – With or Without Distal Attachment



## DESCRIPTION

The ALPS® Extreme™ Liner is specifically designed for transfemoral and active transtibial amputees. It is available in GripGel or HD Gel to provide maximum comfort. The Extreme™ Liner has 80% less vertical stretch than other ALPS® gel liners, which provides more control and stability.

| PART #      | SUSPENSION | L-CODE |
|-------------|------------|--------|
| AKDT/AKDTHD | LOCKING    | L5673  |
| AKFR/AKFRHD | CUSHION    | L5679  |

Recommended for all activity levels.




## TOP FEATURES & BENEFITS




- ✓ Brown fabric offers limited vertical stretch to reduce movement of redundant tissue
- ✓ Wider distal end proportionate to the size of liner for ease of donning
- ✓ Slightly tacky feel reduces bunching in the popliteal region
- ✓ Available in 3mm and 6mm uniform

# CROSS SECTION DIAGRAM



## HOW TO SIZE

Measurements taken 6cm from the distal end



## SIZING CHART

| MEASURED SIZE (CM) | LOCKING<br>3.0MM or 6MM UNIFORM | CUSHION<br>3.0MM or 6MM UNIFORM |
|--------------------|---------------------------------|---------------------------------|
| 16-19 cm           | AKDT(HD) <b>16</b> - ( 3 or 6)  | AKFR(HD) <b>16</b> - ( 3 or 6)  |
| 20-23 cm           | AKDT(HD) <b>20</b> - ( 3 or 6)  | AKFR(HD) <b>20</b> - ( 3 or 6)  |
| 24-25 cm           | AKDT(HD) <b>24</b> - ( 3 or 6)  | AKFR(HD) <b>24</b> - ( 3 or 6)  |
| 26-27 cm           | AKDT(HD) <b>26</b> - ( 3 or 6)  | AKFR(HD) <b>26</b> - ( 3 or 6)  |
| 28-31 cm           | AKDT(HD) <b>28</b> - ( 3 or 6)  | AKFR(HD) <b>28</b> - ( 3 or 6)  |
| 32-34 cm           | AKDT(HD) <b>32</b> - ( 3 or 6)  | AKFR(HD) <b>32</b> - ( 3 or 6)  |
| 35-37 cm           | AKDT(HD) <b>35</b> - ( 3 or 6)  | AKFR(HD) <b>35</b> - ( 3 or 6)  |
| 38-43 cm           | AKDT(HD) <b>38</b> - ( 3 or 6)  | AKFR(HD) <b>38</b> - ( 3 or 6)  |
| 44-53 cm           | AKDT(HD) <b>44</b> - ( 3 or 6)  | AKFR(HD) <b>44</b> - ( 3 or 6)  |

\*Expired patent owned by ALPS and/or one of its subsidiaries



T: 727.528.8566 | 800.574.5426 | F: 727.528.8862  
[www.easyliner.com](http://www.easyliner.com) | [info@easyliner.com](mailto:info@easyliner.com)  
 2895 42nd Avenue North St. Petersburg, Florida 33714



CAN #10346 Form #100412 Rev K