


EASYSLEEVE™

FABRIC REINFORCED SFS


GripGel Suction Suspension Sleeve

ALPS® SFS has superior comfort provided by GripGel. The elasticity of the fabric allows the user excellent freedom of knee flexion. The SFS also achieves secure suspension by sealing with the skin without restricting circulation, while the GripGel sticks to the patients skin without causing shear force.

i Available in both 3mm and 6mm uniform thickness.




CROSS SECTION DIAGRAM



HOW TO SIZE

Measurements taken 10cm
above the patella



SIZING CHART

MEASURED SIZE (CM)	3MM OR 6MM UNIFORM
20-24 cm	SFS 16 - (3 or 6)
25-30 cm	SFS 20 - (3 or 6)
31-34 cm	SFS 24 - (3 or 6)
35-37 cm	SFS 26 - (3 or 6)
38-40 cm	SFS 28 - (3 or 6)
41-45 cm	SFS 32 - (3 or 6)
46-60 cm	SFS 38 - (3 or 6)
61-70 cm	SFS 44 - (3 or 6)



T: 727.528.8566 | 800.574.5426 | F: 727.528.8862
www.easyliner.com | info@easyliner.com
2895 42nd Avenue North St. Petersburg, Florida 33714

