

Adjusted by:

E-mail:

Adjusted for:

Tel:

Date:

Save the guide to the desktop and open it from there.

Total weight (g): 0

The weight changes based on the choices made below.

U3 has a ajustable rear axle in carbon fibre. To find a balance that suits you the rear axle can be adjusted into different positions. Note desired measure at "Settings". S3 has a seat height of 47 cm in the front and 43 at the rear, without cushion. To adjust the seat height you can chose other rearwheel sizes or a raising Kit with suitable castor fork from the List of options.

Max userweight: 100 kg for seat width 33-42 and 150kg for seat width 45 cm

Save the guide to the "desktop" and open it with Adobe Reader, for best function Design your wheelchair by clicking in the circles in this guide. Blue, underlined text is clickable to show info or image. BB = seatwidth, HH= backrest height

Original design

Alternative options

Other options

[info](#)

1) Chassis-colour

Black - Anodic Black

Optional colour - A higher price and

longer delivery time for optional colour, Fill in art.no. in the square

+ 310 + 13%

ΠΛΤ με ΦΠΑ
από 3400 + barcode

2) Seat width (BB)

art.no.

33 cm 36 cm 39 cm 42 cm 45 cm

Clear

Art.no include: chassis, seat, calf strap, backrest lock

3) Back rest

Standard
523BBHH VER2

Clear

Height (HH)

20 cm 35 cm
 25 cm 40 cm
 30 cm 45 cm

Clear

Tapered

525BBHH VER2

Curved back

524BBHH VER2

Other back rest

List of options art. _____

4) Rear axle

Carbon standard
33023BB

Clear

Carbon +15 mm

33023BC*

* seat width BB +1

Use this rear axle when chosing arm support at Accessories and/or increased seat height

Other rear axle or seat height

Increase the seat height +0 mm

List of options art. _____

5) Castor and fork

90 mm
art.no.

Clear

120 mm

art.no.

Ίδια τιμή
με 90mm

6) Footrest

Footrest Titanium
36915BB

Clear

Extended, reinforced

36006BB (w. 33-36) 36006BBB (w. 39-45)

Plastic plate for footrest

36010BB

+ 35 +13%

Other footrest

List of options art. _____

7) Brake

Brake S3
4630000

Clear

One hand brake

Art.no.

+ 100 + 13%

Other brake

List of options art. _____

8) Rear wheel - without pushrim

24" standard
art.no.

Clear

24" Spinergy X black spokes

art.no.

25" Spinergy X black spokes

art.no.

26" Spinergy X black spokes

art.no.

Rear wheel sizes

+ 450
+ 13%

Other rearwheels

List of options art. _____

Tire - choose tire

Schwalbe Right Run 24",25",26"

Clear

Rugged 24"

Schwalbe Marathon Plus 24",25",26"

+ 60 +13%

PU 24" (massive)

reTyre-tyre skins (does not fit PU tyres)

9) Push rim

Titanium
190DD11 (x2 / wheelchair)

Clear

DD = Wheeldiameter

Aluminum

+ 35 + 13%

1904404 x2 (04=24", 05=25", 06=26")

Tetragrip

19052DDBTx2

Paragrip

19051DDBT x2

Ίδια τιμή
με Titanium

Maxgrip 22", 24" (contains latex)

1903060 x2 (60=24", 80=22")

Accessories



Side guards

- [Side guards S3](#) [Video](#)
art.no.


Clear



Arm support

- [Arm support S3](#) [Video](#)
art.no.


+ 150 +13%



Side guard Carbon


- [Side guard Carbon](#)
art.no.

+ 280 + 13%



Push handles + 250 + 13%

- [Height adjustable S3](#) [Video](#)
4252060 VER2
- [Height adjustable S3 +10](#)
4252063 VER2
Max. userweight 100kg
- Clear



Push handle Foldable

- [Foldable S3](#) [Video](#)
4245000

+ 110 +13%



Push handle Pushbar

- [Pushbar S3](#)
42400BB VER2

Max. userweight 100kg



Anti-tip

- [Anti-tip std seat height](#)
4353000 [Video](#)

Av όχι -150 + 13%

[Other seat height](#)



Cushion Wedge cushion

- [2,5 cm](#)
572BB40
- [5 cm](#)
575BB40
- Clear
- [Wedge cushion](#)
57710BB

+ 110 + 13%



Spoke protectors

- [24" logo](#) 4500000
- [22" logo](#) 4500100
- [24" no logo](#) 4500011
- [22" no logo](#) 4500111
- Clear + 120 + 13%



Other

- [Tube protection front](#)
9900251 (2 pcs) > + 12 +13%
ΤΟ ΤΕΜΆΧΙΟ
- [Tube protection rear](#)
9900251 (1 pcs)
- [Crutch holder S3](#) _____ pcs
- 9900301 (x1/crutch)
- [Backrest wedge](#) _____ pcs
- 4710012

- [Table S3](#)
6810020



Chose art.no. for Freewheel according to the **footrest height** of the wheelchair

Read more about Freewheel on the [homepage](#)

- [Freewheel](#) Height 6-14 cm art.no. 2600102

Set-up adjustments



Back rest angle

Note measure (A) for desired back rest angle.
14-28 mm. Standard=20 mm

A= _____ mm



Rear axle adjustment

Note measure (A) for desired balance.
Standard=60 mm, at raised rear axle= 51 mm

A= _____ mm [Info - 26" rearwheel together with the S3 brake](#)



Footrest height

Note measure (A) for the footrest.

A= _____ mm



Transportation marking

- YES do mark the wheelchair with attachment points for transport in disabled service vehic

Art.no.9000152

NOTE! See the [manual](#) for more info.

Extra