

ΟΙ ΤΙΜΕΣ ΕΙΝΑΙ ΛΙΑΝΙΚΗΣ ΧΩΡΙΣ ΦΠΑ 13% ΚΑΙ ΔΕΝ ΣΥΜΠΕΡΙΛΑΜΒΑΝΟΥΝ ΕΞΟΔΑ ΑΠΟΣΤΟΛΗΣ

invoice / delivery address

Kom.:

☐ **Order**
☐ **Offer**

Rehab Car Seat **KIDSFLEX R129** NEW! Edition

Exportprices

Release: 2025

BASIC

	size 1	size 2	size 2XL
total height	74 cm	80 cm	90 cm
crown height	67 cm	73 cm	82 cm
back height	37 - 51 cm	41 - 57 cm	47 - 65 cm
seat width	31 cm	35 cm	40 cm
seat depth (adjustable)	28 - 33 cm	34 - 39 cm	40 - 46 cm
lateral support height (small)	26 - 36 cm	26 - 36 cm	34 - 44 cm
lateral support height (large)	30 - 40 cm	30 - 40 cm	38 - 48 cm
distance of lateral supports	16 - 28 cm	19 - 32 cm	24 - 37 cm
weight	8,2 kg	8,9 kg	12,3 kg
capacity	50 kg	50 kg	100 kg

common information (standard):



- seat depth adjustable
- headrest included (longer backrest)
- upholstery & cover
- backrest angle adjustable
- 5-point-harness with upholstery
- frame to fix a footrest, clamp bow or ISOFIX
- Crash tested according to **R129/03**
- **ISOFIX adapter** included
- **headrest KICK-BACK** included
- **trunk supports** small, fix, included
- **shoulder support** included


☐ size 1
black
(300000-E3)

3.250,00 €


☐ size 2
black
(310000-E3)

3.330,00 €


☐ size 2XL
black
(310000-E3-XL)

3.675,00 €

belt guide*:

- ☐ LEFT
- ☐ RIGHT

*please choose ONE of them

- ☐ wider hip
cushions (SW -4 cm)
(300850-E)
[€ 90,--]

fittings for different needs

head support „Kick-Back“
with aluminum frame
(free of charge)



☐ head fixation
(300630-E), [€ 375,-]



leg support

☐ left
(300350-E)
[€ 155,-]

☐ right
(300355-E)
[€ 155,-]




abduction block, removable

- ☐ size 1
(300400-E)
[€ 275,-]
- ☐ size 2
(310400-E)
[€ 285,-]
- ☐ size 2XL
(310400-E-XL)
[€ 300,-]




foot rest, height and angle adjustable

- ☐ size 1
(300450-E)
[€ 325,-]
- ☐ size 2, 2XL
(310450-E)
[€ 335,-]
- ☐ swing up
(300452-E)
[€ 190,-]



foot roll, height and angle adjustable

- ☐ size 1
(300455-E)
[€ 295,-]
- ☐ size 2, 2XL
(310455-E)
[€ 305,-]



trunk supports, small
(height: 11 cm, depth: 14 cm)

- ☒ fix
(free of charge)
- ☐ swing away
(300505-E3)
[€ 240,-]



trunk supports, large
(height: 21 cm, depth: 16 cm)

- ☐ fix
(300510-E3)
[€ 80,-]
- ☐ swing away
(300515-E3)
[€ 290,-]




shoulder supports
(free of charge)


 fittings frame

swivel base

- ☐ size 1 and 2
(300200-E)
[€ 425,-]
- ☐ RIGHT
- ☐ LEFT




fastener

(independent of the direction of rotation)

- ☐ inside
(standard)
- ☐ outside
(300220-E)

angular adjustment (0° - 10°)

- ☐ size 1
(300100-E)
[€ 220,-]
- ☐ size 2
(310100-E)
[€ 240,-]
- ☐ size 2XL
(310100-E-XL)
[€ 260,-]



extra angular adjustment (10° - 20°)

- ☐ size 1
(300105-E)
[€ 295,-]
- ☐ size 2
(310105-E)
[€ 275,-]
- ☐ size 2XL
(310105-E-XL)
[€ 310,-]



☐ Sandwich-Base (mounted under the base board) (300790-E), [€ 365,00]

Useful in the following situations:


- with very sporty shaped car seats, in order to compensate for the recess there. This results in a more stable position and more freedom of movement when using the swivel adapter.
- with car Isofix brackets that are very far inside to allow nevertheless contact with the car's backrest and the KIDSFLEX (the ISOFIX adapter sits in the sandwich base)
- approx. 1.5 kg additional weight

fittings for car installation

Fixation (one option included in delivery)


ISOFIX-Adapter

☐ size 1
(300575-E3)

☐ size 2 & 2XL
(310575-E3)



clamp bow

☐ size 1
(300550-E3)

☐ size 2 & 2XL
(310550-E3)


outrigger*

☐ (310453-E), [€ 190,-]

*300450-E
or 310450-E
necessary!

additional belt guide
(on the other side)

☐ (300075-E)
[€ 135,-]

 fittings security

security tray

☐ size 1
(300600-E)
[€ 455,-]

☐ size 2
(310600-E)
[€ 475,-]

☐ size 2XL
(310600-E-XL)
[€ 495,-]

Security pad for the buckle
Cover with snap fastener

☐ (300650-D)
[€ 165,-]

Security pad for the buckle
with HERNIK Security Key*


☐ (300655-D)
[€ 210,-]



*incl. belt cutter


protection mat
(to protect the original seats of the car)

☐ (300780-E)
[€ 85,-]

belt to connect the chest straps
(to improve the position of the belt in front of the breast)

☐ (300640-E)
[€ 125,-]

 individuality
☐ Individual Dimensions (310000-E3-SB), [€ 3.675,00 + Expense allowance]

The expense allowance is always calculated individually depending on the desired measurements
Therefore ALWAYS request a quotation.
(max. load: 75 kg; cover colour: grey-black)

☐ Patient dimensions:

a (seat width): _____

b (seat depth): _____

c (back height): _____

d (crown height): _____

e (lower leg length): _____

f (thorax width): _____

☐ Frame dimensions:

a (seat width): _____

b (seat depth): _____

c (back height): _____

d (crown height): _____

☐ fabric, per running metre
(300000-X), [€ 60,-]

☐ red (-B) _____ m

☐ blue (-D) _____ m

☐ black (-C) _____ m

☐ grey (-E) _____ m

☐ without cover
(300000-OP), [- € PR]

☐ integration of a seat shell
(300940-E)

☐ cover for a seat shell
(300810-E)