

ΟΙ ΤΙΜΕΣ ΕΙΝΑΙ ΛΙΑΝΙΚΗΣ ΧΩΡΙΣ ΦΠΑ 13% ΚΑΙ ΔΕΝ ΣΥΜΠΕΡΙΛΑΜΒΑΝΟΥΝ ΕΞΟΔΑ ΑΠΟΣΤΟΛΗΣ

invoice / delivery address

Reference:

☐ **Order**

☐ **Offer**

Rehab Car Seat *IPAI* **MINI**

Exportprices

Release: 2025


BASIC

Details:

IPAI-MINI	
SEAT WIDTH (pelvis)	25 cm
SEAT WIDTH (knee)	31 cm
SEAT DEPTH	26 cm
BACK HEIGHT	30 - 43 cm
CROWN HEIGHT	52 - 65 cm
HEIGHT LATERALS	17 - 31 cm
DISTANCE LATERALS	14 - 28 cm
angle	80° - 120°
total height	66 - 89 cm
total width	44 cm
total depth	44 cm
weight	ab 8 kg

common information (standard):

- crash tested according to ECE-R129/03
- 5-point-harness with upholstery, height adjustable
- belt mechanics installed into the front of the seat
- solid aluminium board, including the mechanism to recline the backrest stepless
- head rest height adjustable, including belt guides
- additional upholstery (20 mm) for the back & headrest
- ISOFIX adapter or clamp bow integrated in base board



☐ **IPAI-MINI red-black**
(480000-E1)



☐ **IPAI-MINI grey-black**
(480000-E2)



accessory that makes the base board **NECESSARY**


☐ **BASE BOARD**

Necessary for the options below!



ISOFIX

☐ (480000-BEI)



☐ (480000-BEK)

swivel base

☐ (480200-E)

☐ right

☐ left



☐ angular adjustment
(5° - 20°)
(480100-E)



Conversion set

☐ from ISOFIX adapter
to clampbow
(480020-E)

☐ from clampbow
to ISOFIX adapter
(480030-E)

foot rest

☐ (480450-E)




☐ swing up*
(480452-E)



*480450-E necessary!


☐ outrigger*
(4 0453-E)



*480450-E necessary!

foot roll

☐ (480455-E)



further options

protection mat

(to protect the original seats of the car)

☐ (480780-E)



Security pad for the buckle

☐ cover with snap fastener
(480650-E)



☐ with HERNIK Security Key*
(480655-E)



*incl. belt cutter

NEW!

fittings for stabilization

- ☐ abduction block, removable
(480400-E)




- ☐ head fixation
(480630-E)



- ☐ security tray
(480600-E)



- ☐ side supports, pair, adjust.
in height, deep, width
(480800-E)



- ☐ seatdepth adapter
(480850-E)



- ☐ lateral support
(480810-E)



-  belt to connect the chest straps
(480640-E)



Notes