

ΟΙ ΤΙΜΕΣ ΕΙΝΑΙ ΛΙΑΝΙΚΗΣ ΧΩΡΙΣ ΦΠΑ 13% ΚΑΙ ΔΕΝ ΣΥΜΠΕΡΙΛΑΜΒΑΝΟΥΝ ΕΞΟΔΑ ΑΠΟΣΤΟΛΗΣ

invoice / delivery address

Reference:

☐ **Order**

☐ **Offer**

Rehab Car Seat *IPAI* **LGT**

Exportprices

Release: 2025

BASIC

common information (standard):

- crash tested according to ECE-R129/03
- 5-point-harness with upholstery, height adjustable, can be tightened independently
- large pad between the buckle and the pelvis
- ISOFIX adapter included with one-hand-adjustment
- head rest height adjustable, including belt guides
- Side impact protection SICT, separately adjustable per side
- side supports, individually height-adjustable
- removable, breathable cover (washable at 30 °C)
- every accessory can be mounted in seconds, even later!



Details:

	IPAI-LGT
SEAT WIDTH (pelvis)	31 cm
SEAT WIDTH (knee)	34 cm
SEAT DEPTH	30 cm
BACK HEIGHT	41 - 54 cm
CROWN HEIGHT	59 - 74 cm
HEIGHT LATERALS	20 - 33 cm
DISTANCE LATERALS	14 - 28 cm
angle	85° - 160°
total height	68 - 88 cm
total width	51 cm
total depth	42 cm
weight	ab 7,0 kg

☐ **IPAI-LGT
black**
(460000-E)

accessories, that make the base plate necessary

☐ **basis plate**

Necessary for the following accessories
(460200-E)


☐ **abduction block**
(460400-E)

☐ **seatdepth adapter**
(460850-E)

foot rest
☐ (460450-E)

☐ **swing up***
(460452-E)


*460450-E necessary!

☐ **outrigger***
(460453-E)


*460450-E necessary!

further options

☐ **protection mat**
(to protect the seat of the car)
(460780-E)

☐ **head fixation**
(460630-E)

☐ **belt to connect the chest straps**
(460640-E)

☐ **security try**
(460600-E)

☐ **angular adjustment (10°)**
(460100-E),

Security pad for the buckle
☐ **cover with snap fastener** (460650-E),

☐ **with HERNIK Security Key*** (460655-E)


*incl. belt cutter