

REINFORCED SLEEVE RFS


Gel Suction Suspension Sleeve

ALPS® Reinforced Sleeve is constructed with grey knitted fabric, OptiGel, and an extra layer of internal fabric to deliver superior durability. This combination works together to provide the maximum comfort, control, and high-performance capabilities required for higher activity-level users.

i Available in 3mm uniform thickness.




CROSS SECTION DIAGRAM



HOW TO SIZE

Measurements taken 10cm
above the patella



SIZING CHART

MEASURED SIZE	3.0 MM UNIFORM
25-30 cm	RFS 20-3G
31-34 cm	RFS 24-3G
35-37 cm	RFS 26-3G
38-40 cm	RFS 28-3G
41-45 cm	RFS 32-3G
46-60 cm	RFS 38-3G
61-70 cm	RFS 44-3G