

Date .....

User name .....

Company stamp

Responsible .....

Folding frame wheelchair

Wheelchair measures

Seat width (LS)	<input type="checkbox"/> 33	<input type="checkbox"/> 36	<input type="checkbox"/> 39	<input type="checkbox"/> 42	<input type="checkbox"/> 45	<input type="checkbox"/> 48
The seat width is measured from the outer sides of the frame, while the distance between the two side guards is 1 cm wider.						
Seat depth (PS)	<input type="checkbox"/> 35	<input type="checkbox"/> 37,5	<input type="checkbox"/> 40	<input type="checkbox"/> 42,5	<input type="checkbox"/> 45	<input type="checkbox"/> 47,5
The real seat depth, from the front part of the seat cloth to the backrest tube will be 1,5 cm more than the declared one.						
Footplate distance (DP) .....	(Real distance between seat and footplate)					
<b>IMPORTANT NOTICE! From the footplate to the floor minimum 4cm.!! (Minimum Front height = DP + 4cm.)</b>						
Setting	☆ Standard <input type="checkbox"/> Prudential <input type="checkbox"/> Active					
Backrest Height (TS)	Std.Value 33,5 .....(min 24.5-max 47)					
camber	☆ 0° <input type="checkbox"/> 2°					

Front height STD 49	.....
Rear height STD 42	.....


Front frame

☆ Aluminium standard frame


Colour aluminium frame

- ☆ Bright silver
- ☐ Polished black
  - ☐ Opaque black
  - ☐ Bright blue
  - ☐ Bright anthracite
  - ☐ Bright orange
  - ☐ Bright fuchsia
  - ☐ Bright red
  - ☐ Bright lime
  - ☐ Bright blue china
  - ☐ Pearl white


Aluminium frame



☆



☐  
(max 5" caster)



☐  
(max 4" caster)

Backrest

☆ Adjustable velcro straps and padded transpiring upholstery

Handle tube

☆ Tubes with short handle

- ☐ SPORT tubes Without push handles (Lightest configuration)
- ☐ Folding handles



02958      200.00

Casters

☆ 4" Soft

- ☐ 3" Sport (Lightest configuration)
- ☐ 5" Soft
- ☐ 6" Soft (Available only with 100° angle)
- ☐ 6" Pneumatic (Available only with 100° angle)

- |  |                       |       |        |
|--|-----------------------|-------|--------|
| <input type="radio"/> 3" Active Soft (wide)    | (Aluminium alloy rim) | 06824 | 110.00 |
| <input type="radio"/> 4" Active Sport (narrow) | (Aluminium alloy rim) | 01398 | 110.00 |
| <input type="radio"/> 4" Active Soft (wide)    | (Aluminium alloy rim) | 03943 | 110.00 |
| <input type="radio"/> 5" Active Soft (wide)    | (Aluminium alloy rim) | 07518 | 110.00 |
| <input type="radio"/> 5" Pneumatic             |                       | 10844 | 50.00  |

Support front fork

☆ Standard support front fork



- ☐ Front fork support V Design 2.0  
Exclude the installation of : Quick release front forks



<b>Date</b> .....	<b>Order</b> .....	<b>User name</b> .....	<b>Address</b> .....
-------------------	--------------------	------------------------	----------------------

**Special Options**

<input type="checkbox"/> Rear wheels plate extension	01950	160.00
<input type="checkbox"/> lockingband/strip to close/fold the wheelchair	06833	20.00
<input type="checkbox"/> Drawstring to unlock the back with a singleaction	06834	110.00

**Colour set: Rear wheels and forks - To be matched only with: Rear wheels Ultra-light Progeo - Standard! (22 - 24)**

<input checked="" type="checkbox"/> Black	22" 24" 25" 26"		
<input type="checkbox"/> Bright silver	22" 24"	00833	110.00
<input type="checkbox"/> Bright anthracite	24"	00833	110.00
<input type="checkbox"/> Bright orange	22" 24"	00833	110.00
<input type="checkbox"/> Bright fuchsia	22" 24"	00833	110.00
<input type="checkbox"/> Bright blue	22" 24"	00833	110.00
<input type="checkbox"/> Bright red	22" 24"	00833	110.00
<input type="checkbox"/> Bright lime	22" 24"	00833	110.00
<input type="checkbox"/> Pearl white	22" 24"	00833	110.00

**Rear wheels**

<input checked="" type="checkbox"/> 24" Ultra-light Progeo - Standard!			
<input type="checkbox"/> 22" Ultra-light Progeo - Standard!			
<input type="checkbox"/> 25" Ultra-light Progeo - Standard!	(colour: black)		
<input type="checkbox"/> 26" Progeo - Standard!	(colour: black)		
<input type="checkbox"/> 24" ULTRA-LIGHT PROGEO KEY SPORT	(Lightest configuration)	04383	750.00
<input type="checkbox"/> 24" Ultra-light "Spinergy LX"		01685	750.00
<input type="checkbox"/> 24" Ultra-light "Spinergy LXL"	<input type="checkbox"/> Black hub and white spokes <input type="checkbox"/> Hub and spokes black	07846	750.00
<input type="checkbox"/> 24" Ultra-light "Spinergy Flex rim"	Integrated aluminium/rubber handrim	02270	1,200.00
<input type="checkbox"/> 25" Ultra-light "Spinergy LX"		02101	750.00
<input type="checkbox"/> 25" Ultra-light "Spinergy LXL"	<input type="checkbox"/> Black hub and white spokes	07847	750.00
<input type="checkbox"/> 26" Ultra-light "Spinergy LX"		02102	750.00
<input type="checkbox"/> 24" Off-road (aluminium handrim black)		01392	200.00

**Rear wheel tyre**

<input checked="" type="checkbox"/> High pressure tyre	(22"-24"-25"-26")	(Lightest configuration)	
<input type="checkbox"/> High pressure profile tyre	(24")		
<input type="checkbox"/> Marathon Plus	(24"-25")	with anti-puncture protection	01669      55.00
<input type="checkbox"/> Anti-puncture - Colour Black	(24")		11083      100.00
<input type="checkbox"/> Anti-puncture - Colour Grey	(24")		11083      100.00
<input type="checkbox"/> Anti-puncture - Colour Grey	(22")		04181      100.00


  

Handrim	Wheel distance <input checked="" type="checkbox"/> 3 cm	<input type="checkbox"/> 4 cm		
<input checked="" type="checkbox"/> Aluminium	(22"-24"-25"-26")	(Lightest configuration)		
<input type="checkbox"/> Aluminium black SOFT-TOUCH	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	09170	100.00
<input type="checkbox"/> Chromed Steel	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	01684	150.00
<input type="checkbox"/> Titanium	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00016	250.00
<input type="checkbox"/> Ultra Grepp	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00017	280.00
<input type="checkbox"/> Max Grepp	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00834	310.00
<input type="checkbox"/> Ergo Para	(24")		00632	310.00
<input type="checkbox"/> TOP FRICTION	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11023	260.00
<b>Oval countour with a rubber insert</b>				
<input type="checkbox"/> NOVA H	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11025	240.00
<b>Ergonomic black profile</b>				
<input type="checkbox"/> NOVA D	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11026	260.00
<b>Ergonomic profile with rubber coating (ideal for quadriplegics)</b>				
<input type="checkbox"/> SURGE LT	(25")	Oval countour with a rubber insert	06427	290.00
<input type="checkbox"/> Quad	(24")		00693	180.00
<input type="checkbox"/> Silicon cover for standard aluminium handrim		Available for measures (22" - 24" - 26")	04050	150.00

**Wheel locks**

<input checked="" type="checkbox"/> Bent lever push to lock		
<input type="checkbox"/> Straight lever push to lock		
<input type="checkbox"/> Straight lever pull to lock		
<input type="checkbox"/> Sport brake	(Lightest configuration)	06428      100.00



Date	Order	User name	Address
------	-------	-----------	---------

Automatic Closure footplate


☆ Footplate COMPACT, automatic Closure

(Lightest configuration)

Minimum footplate distance 35cm.

☐ Internal

☐ External



RAKE AND DEPTH ADJUSTABLE

Footplate - See reference tables N°3

Warning!! All flip-up footrests, need a stretch calf band and standard fork-support

Weight increase approx. 500gr.

☐ Plastic footrest. rake adjustable


Minimum footplate distance 35cm.

☆2/3

Interna

☐ 2/3

External




☐ One piece footplate - Flip-up!. rake and depth adjustable

Minimum footplate distance 35cm.

☐ Internal

☆2/3 Internal

☐ External



☐ One piece aluminium footplate - Flip-up!. rake and depth adjustable


Minimum footplate distance 35cm.

☐ Internal

☆2/3 Internal

☐ 2/3

External




☐ One piece Carbon fibre footplate, rake adjustable, swing away, flip-up.

Minimum footplate distance 35cm.

☆ Internal

☐ External



03291240.00


Mudguard

Mudguard - FIXED

☆ Carbon fibre mudguard

☐ Carbon fibre Straight sideguard

☐ Extra charge for : REMOVABLE VERSION




08921100.00

Armrests

☐ Tubular armrests

(Swing away, removable and height adjustable armrest)




00694220.00

Removable armrest with ABS mudguard


☐ Desk armrests PROGEO

(Removable and height adjustable)



09956200.00

☐ Removable armrest with carbon-fiber straight sideguard



09079180.00

Accessories

☐ Pad for front frame protectors

0705840.00

☐ Backrest-seat connection cloth

1024115.00

☐ Removable storage bag (with quick coupling)





09811110.00

Pag. 3 / 7



Rev. 2022.0

Date ..... Order ..... User name ..... Address .....

**Optional 1**

<input type="radio"/> Transit wheels			00836	130.00
<input type="radio"/> Anti-tipper with rapid hookup	(inside rotation to reduce the encumbrance)		00023	110.00
	<input type="radio"/> Right		00024	110.00
	<input type="radio"/> Left			
<input type="radio"/> Tipping aid	<input type="radio"/> Right		00025	70.00
	<input type="radio"/> Left		01046	70.00
<input type="radio"/> Crutches support	<input type="radio"/> Right		00838	70.00
	<input type="radio"/> Left		00837	70.00

**Optional 2**

<input type="radio"/> Tetraclick	(Quick release axle for quadriplegic)		03947	150.00
<input type="radio"/> Safety pelvic belt			03972	60.00
<input type="radio"/> Quick release front forks	Available with: Standard support front fork		00030	150.00
<input type="radio"/> Push handles for companion	(Removable and height adjustable)		01048	170.00
<input type="radio"/> Harness			04103	120.00
<input type="radio"/> Lumbar support	(to insert in the backrest upholstery)		00013	40.00
<input type="radio"/> Swing away lateral support	Not available with Posture backrest PHYSIO		02097	220.00
	<input type="radio"/> Right		02098	220.00
	<input type="radio"/> Left			
<input type="radio"/> Spoke guards Black	(22"-24")		01043	140.00
<input type="radio"/> Spoke guards Transparent	(22"-24")		00835	140.00
<input type="radio"/> Spoke guards Progeo design	(22"-24")		00019	160.00
<input type="radio"/> Spoke guards Dragon	(22"-24")		02106	160.00

**Additional Kit**

<input type="radio"/> 24" Off-road rear wheels	Off-Road tyre - aluminium handrim (black colour)	03378	450.00
<input type="radio"/> 24" FAT rear wheels	FAT tyre 3" - aluminium handrim (black colour)	09357	750.00

Date

Order








User name

Address

PHYSIO - Posture backrest - Std. Backrest upholstery not provided





Removable with quick release  
Suitable for tubes with diameter 18÷26mm  
Depth adjustment (4cm)  
Angle adjustment (±18°)

☐  Carbon shell  
not available for seat width 48

BACKREST HEIGHT		28		34		40	
Shell's shape		Lateral retaining wing A=high - R=regular					
7,5 cm		A	R	A	R	A	R
							
		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

840.00

☐  Aluminum shell

BACKREST HEIGHT		34	40	52
Shell's shape		Lateral retaining wing B=regular		
7,5 cm		R	R	R
				

810.00

Configuration DEEP with side support PAD (height adjustable)

☐ Aluminum lateral with padded, breathable upholstery.



Small (h 13cm)



☐ Right  
10272

25.00

☐ Left  
10272

25.00

Large (h 20cm)



☐ Right  
10273


30.00

☐ Left  
10273

30.00

SWING AWAY - Lateral support

Fixing and adjustment hardware complete with protection sleeve  
Alluminum lateral with padded, breathable upholstery.



☐ Large (18x12)

05647 150.00

☐ Small (14x9)

05643 150.00

☐ Lateral support extension

Quantity (number pcs)..... 10271 55.00

☐ Lateral support extension

Quantity (number pcs).....

☐ Small (14x9)

05642 150.00

☐ Large (18x12)

05646 150.00

☐ Headrest adjustable - Double ball mounting system

Available only with backrest Physio  
Not available for backrest Physio AIR carbon shell



☐ Large (34x14) 05648 200.00



☐ Small (28x14)

Date .....	Order .....	User name .....	Address .....
------------	-------------	-----------------	---------------

Note

Wheelchair : Progeo EGO - Rev. 2022.0	Prices in EURO - VAT excluded to be added to the price list	Send this page only if any optional have been chosen	
		<input type="checkbox"/> standard	(supplied if no requests are indicated)
		<input type="checkbox"/> on request	(without extra charge)
		<input type="checkbox"/> optional	(cost to be added at the basic price)

## REFERENCE TABLES

1	Front height			
With closed frame (95°) max 5" caster With closed frame (90°) max 5" caster With closed frame (85°) max 4" caster				
FORK	Casters			
	3"	4"	5"	6"
SMALL 	max 46	max 47	max 48	
MEDIUM 	max 50	max 51	max 52	max 54
The size of the fork will be defined only by Progeo in accordance with wheelchair final configuration.				

2	Rear height			
Rear wheels				
22"	24"	25"	26"	
min 35	min 38	min 39	min 40	
max 43	max 45	max 46	max 48	

3	Footplate option	Footplate distance
	Footplate COMPACT, automatic Closure	min 36
	Plastic footrest	min 36
	One piece footplate - Flip-up!	min 36
	One piece aluminium footplate - Flip-up!	min 36
	One piece Carbon fibre footplate	min 36
<b>IMPORTANT NOTICE! From the footplate to the floor minimum 4cm.!! (Minimum Front height = DP + 4cm.)</b>		
<b>NOTE! Lower footplate distances are possible only after our valuation</b>		
<b>Calculated with the best configuration and with the footplates positioned on inside.</b>		

4	Rear wheels	
	Model	Features
	22" Ultra-light Progeo - Standard!	(not available in anthracite color)
	24" Ultra-light Progeo - Standard!	(available for all colors)
	24" ULTRA-LIGHT PROGEO KEY SPORT	(silver hub, black rim and 16 spokes)
	24" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	24" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white spokes) - (black hub, rim and 18 spokes)
	24" Ultra-light "Spinergy Flex rim"	(silver hub, black rim and spokes) - (integrated aluminium/rubber handrim - encumbrance reduction up to 3cm. - all tyre)
	24" Off-road	(colour: black) - (aluminium handrim - tyre Off-road)
	25" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	25" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white spokes)
	26" Progeo - Standard!	(colour: black)
	26" Ultra-light "Spinergy LX"	(silver hub, black rim and 18 spokes)