

Date .....

User name .....

Responsible .....

Company stamp

Rigid frame wheelchair

### Wheelchair measures

Seat width (LS) 33 36 39 42 45 48  
 The seat width is measured from the outer sides of the frame, while the distance between the two side guards is 1.5 cm wider.

Seat depth (PS) 35 37,5 40 42,5 45 47,5 50  
 The real seat depth, from the front part of the seat cloth to the backrest tube will be 1,5 cm more than the declared one.


Footplate distance (DP) ..... **(Real distance between seat and footplate)**  
**IMPORTANT NOTICE! From the footplate to the floor minimum 4cm.!! (Minimum Front height = DP + 4cm.)**

Setting  Standard  Prudential  Active  Extreme

Backrest Height (TS) Std.Value 33,5 .....(min 24.5-max 47)

Camber  0°  2°  4°

### ☆ Standard rear adjustment




- Carbon fibre mudguard
- Carbon fibre Straight sideguard
- Removable version



08921 100.00

FULL CARBON rear adjustment plate with integrated mudguard  
 (Lightest configuration)



07646 800.00



- Version mudguard (AVAILABLE ONLY WITH REAR WHEELS 24")
- Version straight sideguard

### Front height - Tick the box with the required height

The possible difference in height between HIGH frame and LOW frame is 3cm (with same fork, caster and assembly hole).

	 SPORT fork (H95)	 Small size fork (H115)
<b>LOW FRAME 90°</b> Front height (min 44,5 - max 50) Footplate distance (min 33)  Footplate distance (min 21) with U-TUBE Aluminium footplate	Casters 3" <input type="checkbox"/> ~45 <input type="checkbox"/> ~46 4" <input type="checkbox"/> ~46 <input type="checkbox"/> ~47 5" <input type="checkbox"/> ~48,5	Casters 3" <input type="checkbox"/> ~45,5 <input type="checkbox"/> ~47 4" <input type="checkbox"/> ~46,5 <input type="checkbox"/> ~48 <input type="checkbox"/> ~49 5" <input type="checkbox"/> ~49 <input type="checkbox"/> ~50
<b>HIGH FRAME 90°</b> Front height (min 47,5 - max 53) Footplate distance (min 36)  Footplate distance (min 24) with U-TUBE Aluminium footplate	Casters 3" <input type="checkbox"/> ~48 <input type="checkbox"/> ~49 4" <input type="checkbox"/> ~49 <input type="checkbox"/> ~50 5" <input type="checkbox"/> ~51,5	Casters 3" <input type="checkbox"/> ~48,5 <input type="checkbox"/> ~50 4" <input type="checkbox"/> ~49,5 <input type="checkbox"/> ~51 <input type="checkbox"/> ~52 5" <input type="checkbox"/> ~52 <input type="checkbox"/> ~53

The real front-seat height could be slightly different from the chosen one considering it is affected by seat depth and rear seat height.

### Rear height .....

- Rear wheels 22" (min 36 - max 44)
- Rear wheels 24" (min 38 - max 46)
- Rear wheels 25" (min 39 - max 47)
- Rear wheels 26" (min 41 - max 48)

The difference between Front Height and Rear Height must be in the 3 - 12 cm range.

Date ..... Order ..... User name ..... Address .....

**Seat**

☆ Seat cloth with pocket (under seat)

Seat with tension straps 08883 **200.00**

Carbon rigid seat 08336 **400.00**

Pocket (under seat) [25x11x4,5] 08281 **60.00**

**CUSTOMIZE COLORS**

**1 Frame colour**

Pearl white  Opaque orange

Opaque black 08341 **60.00**

Bright silver

Bright red

Bright blue

**2 Fork colour**

☆ Opaque black

Bright silver  Bright blue

Bright anthracite  Bright orange

Pearl white  Bright fuchsia

Bright red  Bright lime

08342 **55.00**

**3 Colour rear wheels**

To be matched only with: Rear wheels Ultra-light Progeo - Standard! (22 - 24)

☆ Black


Bright silver  Bright orange

Bright anthracite (24")  Bright blue

Pearl white  Bright fuchsia

Bright red  Bright lime

04692 **55.00**



**Rear wheels axle**

☆ Aluminum rear wheels axle

Carbon fiber rear wheels axle (Lightest configuration) **07647 140.00**

**Backrest adjustment and fixing system**

☆ Backrest with double locking (upright and folded) with rake adjustable from -4° to +16°

Backrest with upright locking system ONLY (no locking when folded) with rake adjustable from -4° to +16°

Fixed back(non-folding back) with rake adjustable from -4° to +16°


**Backrest**

☆ Adjustable velcro straps and padded transpiring upholstery

☆ Aluminium backrest

Titanium backrest 06811 **160.00**

Carbon fibre backrest 04389 **600.00**









**AIR** Integrated posture backrest (NOT REMOVABLE) **840.00**

Carbon shell - Lateral retaining 7,5cm

Depth adjustment ±1cm. - Angle adjustment ±10°

not available for seat width 48

H 28		H 34		H 40	
Backrest Height (TS) min 34		Backrest Height (TS) min 38		Backrest Height (TS) min 44	
( ..... )		( ..... )		( ..... )	
Lateral retaining wing A=high - R=regular					
A	R	A	R	A	R
					
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>


**Installation on:**

☆ Backrest tubes


Aluminium backrest

Titanium backrest 06811 **160.00**

Carbon fibre backrest 04389 **600.00**



Date .....	Order .....	User name .....	Address .....
------------	-------------	-----------------	---------------

<b>Handle tube</b>			
☆ Tubes with standard handle			
<input type="checkbox"/> Tubes with short handle			
<input type="checkbox"/> SPORT tubes Without push handles			
<input type="radio"/> Carbon fibre tubes with short handle	Weight reduction approx. 110gr.	04390	120.00
<b>Not available with: Posture backrest PHYSIO - Posture backrest A I R</b>			
<input type="radio"/> Folding handles			02958 200.00












<b>Casters</b>			
☆ 4" Soft			
<input type="checkbox"/> 3" Sport			
<input type="checkbox"/> 5" Soft			
<input type="radio"/> 3" Active Soft (wide)	(Aluminium alloy rim)	06824	110.00
<input type="radio"/> 4" Active Sport (narrow)	(Aluminium alloy rim)	01398	110.00
<input type="radio"/> 4" Active Soft (wide)	(Aluminium alloy rim)	03943	110.00
<input type="radio"/> 5" Active Soft (wide)	(Aluminium alloy rim)	07518	110.00
<input type="radio"/> 5" Pneumatic		10844	50.00

<b>Rear wheels</b>			
☆ 24" Ultra-light Progeo - Standard!			
<input type="checkbox"/> 22" Ultra-light Progeo - Standard!			
<input type="checkbox"/> 25" Ultra-light Progeo - Standard!	(colour: black)		
<input type="checkbox"/> 26" Progeo - Standard!	(colour: black)		
<input type="radio"/> 24" ULTRA-LIGHT PROGEO KEY SPORT	Weight reduction approx. 450gr.	04383	750.00
<input type="radio"/> 24" Ultra-light "Spinergy LX"		01685	750.00
<input type="radio"/> 24" Ultra-light "Spinergy LXL"	<input type="checkbox"/> Black hub and white spokes <input type="checkbox"/> Hub and spokes black	07846	750.00
<input type="radio"/> 24" Ultra-light "Spinergy Flex rim"	<b>Integrated aluminium/rubber handrim</b>	02270	1,200.00
<input type="radio"/> 25" Ultra-light "Spinergy LX"		02101	750.00
<input type="radio"/> 25" Ultra-light "Spinergy LXL"	<input type="checkbox"/> Black hub and white spokes	07847	750.00
<input type="radio"/> 26" Ultra-light "Spinergy LX"		02102	750.00
<input type="radio"/> 24" Off-road (aluminium handrim black)		01392	200.00

<b>Rear wheel tyre</b>			
☆ High pressure tyre	(22"-24"-25"-26")		
<input type="checkbox"/> High pressure profile tyre	(24")		
<input type="radio"/> Marathon Plus	(24"-25")	with anti-puncture protection	01669 55.00
<input type="radio"/> Anti-puncture - Colour Grey	(22")		04181 100.00
<input type="radio"/> Anti-puncture - Colour Black	(24")		11083 100.00
<input type="radio"/> Anti-puncture - Colour Grey	(24")		11083 100.00

<b>Handrim</b>		<b>Wheel distance</b> ☆ 3 cm <input type="checkbox"/> 4 cm	
☆ Aluminium	(22"-24"-25"-26")		
<input type="radio"/> Aluminium black SOFT-TOUCH	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	09170 100.00
<input type="radio"/> Chromed Steel	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	01684 150.00
<input type="radio"/> Titanium	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00016 250.00
<input type="radio"/> Ultra Grepp	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00017 280.00
<input type="radio"/> Max Grepp	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00834 310.00
<input type="radio"/> Ergo Para	(24")		00632 310.00
<input type="radio"/> TOP FRICTION	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11023 260.00
Oval countour with a rubber insert			
<input type="radio"/> NOVA H	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11025 240.00
Ergonomic black profile			
<input type="radio"/> NOVA D	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11026 260.00
Ergonomic profile with rubber coating (ideal for quadriplegics)			
<input type="radio"/> SURGE LT	(25")	Oval countour with a rubber insert	06427 290.00
<input type="radio"/> Quad (24")	(24")		00693 180.00
<input type="radio"/> Silicon cover for standard aluminium handrim		Available for measures (22"- 24" - 26")	04050 150.00

<b>Wheel locks</b>			
☆ Bent lever push to lock			
<input type="checkbox"/> Straight lever push to lock			
<input type="checkbox"/> Straight lever pull to lock			
<input type="radio"/> Sport brake			
<input type="radio"/> Brake with extension lever push to lock			
<input type="radio"/> Brake with extension lever pull to lock			06428 100.00
			01949 60.00
			03372 60.00








Date .....	Order .....	User name .....	Address .....
<b>Footplate</b>			
☆ Aluminium footplate. Rake not adjustable			
☆ Aluminium footplate	<input type="checkbox"/> with plastic plate		
<input type="radio"/> One piece aluminium footplate - Flip-up! FLIP-BACK! - rake and depth adjustable			09420 200.00
<input type="radio"/> Aluminium footplate - rake and depth adjustable			04456 220.00
☆ Internal	<input type="checkbox"/> 2/3 Internal	<input type="checkbox"/> External	<input type="checkbox"/> 2/3 External
<input type="radio"/> U-TUBE - Aluminium footplate, rake and depth adjustable	<b>Low Frame minimum footplate distance 21 - High Frame minimum footplate distance 24</b>		
<input type="radio"/> One piece Carbon fibre footplate. Rake adjustable			04110 300.00
☆ Internal	<input type="checkbox"/> External		
<b>Armrests</b>			
<input type="radio"/> Tubular armrests	(Swing away, removable and height adjustable armrest)		00694 220.00
<input type="radio"/> Removable armrest with carbon-fiber straight sideguard	Available only with: Standard rear adjustment		09079 180.00
<b>Accessories</b>			
<input type="radio"/> Pad for front frame protectors			07058 40.00
<input type="radio"/> Backrest-seat connection cloth			10241 15.00
<input type="radio"/> Removable storage bag (with quick coupling)			09811 110.00
<b>Optional 1</b>			
<u>Only one optional devices for side can be chosen</u>			
<input type="radio"/> Transit wheels	Exclude the installation of Anti-tippers		03290 130.00
<input type="radio"/> Anti-tipper with rapid hookup	Exclude the installation of Transit wheels (inside rotation to reduce the encumbrance)		06821 110.00
	<input type="radio"/> Right		06822 110.00
	<input type="radio"/> Left		
<b>Optional 2</b>			
<input type="radio"/> Push handles for companion	(Removable and height adjustable)		03645 170.00
<input type="radio"/> Tetraclip	(Quick release axle for quadriplegic)		03947 150.00
<input type="radio"/> Safety pelvic belt			03972 60.00
<input type="radio"/> Lumbar support	(to insert in the backrest upholstery)		00013 40.00
<input type="radio"/> Swing away lateral support	<input type="radio"/> Right		02097 220.00
Not available with: Posture backrest PHYSIO - Posture backrest A I R	<input type="radio"/> Left		02098 220.00
<input type="radio"/> Spoke guards Black	(22"-24")		01043 140.00
<input type="radio"/> Spoke guards Transparent	(22"-24")		00835 140.00
<input type="radio"/> Spoke guards Progeo design	(22"-24")		00019 160.00
<input type="radio"/> Spoke guards Dragon	(22"-24")		02106 160.00
<b>Additional Kit</b>			
<input type="radio"/> 24" Off-road rear wheels	Off-Road tyre - aluminium handrim (black colour)		03378 450.00
<input type="radio"/> 24" FAT rear wheels	FAT tyre 3" - aluminium handrim (black colour)		09357 750.00

Date ..... Order ..... User name ..... Address .....

### PHYSIO - Posture backrest - Std. Backrest upholstery not provided





Removable with quick release  
 Suitable for tubes with diameter 18÷26mm  
 Depth adjustment (4cm)  
 Angle adjustment (±18°)

- AIR** Carbon shell  
not available for seat width 48

BACKREST HEIGHT	28	34	40			
Shell's shape 	Lateral retaining wing A=high - R=regular					
	A	R	A	R	A	R
						
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

840.00

- POSTURE** Aluminum shell

BACKREST HEIGHT	34	40	52			
Shell's shape 	Lateral retaining wing B=regular					
	R	R	R			
						
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			

810.00

### Configuration DEEP with side support PAD (height adjustable)

- Alluminum lateral with padded, breathable upholstery.



Small (h 13cm)



Right  
10272

25.00

Left  
10272

25.00

Large (h 20cm)



Right  
10273


30.00

Left  
10273

30.00

### SWING AWAY - Lateral support

Fixing and adjustment hardware complete with protection sleeve  
 Alluminum lateral with padded, breathable upholstery.



- Large (18x12)

05647 150.00

- Small (14x9)

05643 150.00

- Lateral support extension

Quantity (number pcs).....  
10271 55.00

- Lateral support extension

Quantity (number pcs).....

- Small (14x9)

05642 150.00

- Large (18x12)

05646 150.00

- Headrest adjustable - Double ball mounting system



Available only with backrest Physio  
 Not available for backrest Physio AIR carbon shell

- Large (34x14) 05648 200.00
- Small (28x14)

Date .....	Order .....	User name .....	Address .....
------------	-------------	-----------------	---------------

--

--

--

--

<b>Note</b>

Wheelchair : Progeo JOKER R2 - Rev. 2022.0	Prices in EURO - VAT excluded to be added to the price list	<b>Send this page only if any optional have been chosen</b>	
		<input type="checkbox"/> <b>standard</b> (supplied if no requests are indicated)	<input type="checkbox"/> <b>on request</b> (without extra charge)
		<input type="checkbox"/> <b>optional</b> (cost to be added at the basic price)	