

Please contact DHCare EU Export Customer Services on :-
Tel : +33 247 626980 ; Email : serviceclient_export@dhcare.com
 Or visit our website : www.dhcare.eu.com



Rea® Azalea® Max

Manual steel wheelchair with folding frame



Description:

- Seat tilt with maintenance of the center of gravity and the correct position of the leg rests: 0°/20° and of the backrest: 0°/30° by pneumatic cylinders.
- Adjustable legrests from 0° to 80° and width
- adjustable armrests.
- Frame colour: Pearl gray (RAL 0515).
- Weight: from 54 kg.
- Maximum user weight: 180 kg.
- 2-year warranty period excluding wear parts



Date :

Purchase Order :

Customer Code :

Distributor Stamp:


DHCare Unique Customization Service

Please specify the quote or special production number :




☐ Please make your selection with a cross in the box

☒ Standard feature

☐ Free of charge option

CODE		DESCRIPTION	NC : Non Compatible NC : 9999 = not compatible with the reference 9999	Euros	
BSA0010	☐	Rea® Azalea® Max			
PROPULSION TYPE					
BSA0030	☐	SELF PROPELLED: 24"/ 22" WHEELS WITH USER BRAKES	NC : 0524, 0545	◆	
BSA0040	☐	SELF- AND ASSISTANT PROPELLED: 24" / 22" WHEELS WITH HUB BRAKES	NC : 0522		
SEAT DIMENSIONS					
SEAT WIDTH					
BSA0050	☐	550 mm - Adjustable up to 600 mm		○	 0130
BSA0060	☐	610 mm - Adjustable up to 660 mm		○	
Seat width pre-adjustment					
BSA0070	☐	Supplied with adjustment + 50 mm (at 600 or 660 mm)		○	
Seat depth					
BSA0090	☐	Adjustable from 530 to 600 mm		◆	
Floor/seat height (measured without cushion)					
BSA0100	☐	400 mm	NC : 0290, 0495, 0516	○	
BSA0110	☐	450 mm	NC : 0480, 0514	◆	
SEAT TYPE					
BSA0120	☐	DEPTH ADJUSTABLE		◆	
SEAT CUSHION					
BSA0130	☐	MAX SEAT CUSHION		◆	
SEAT UPHOLSTERY					
BSA0155	☐	BLACK PU TR26 - MUTLI STRETCH PU COATED FABRIC		◆	
BACKREST					
BACKREST TYPE					
BSA0160	☐	BACKREST, TENSION ADJUSTABLE - PRE-ADJUSTED SPACE UPWARDS		◆	 0160
BACKREST COVER/ CUSHION					
BSA0170	☐	COVER, LATERAL - BACKREST UPHOLSTERY MANDATORY		◆	 0180  0170
BSA0180	☐	CUSHION, LAGUNA (LATERAL SUPPORT) - BACKREST UPHOLSTERY MANDATORY	NC : 0415	○	
BACKREST UPHOLSTERY					
BSA0215	☐	BLACK PU TR26 - MUTLI STRETCH PU COATED FABRIC		◆	
PUSH HANDLES					
BSA0220	☐	PUSH HANDLES BRACED		◆	 0220
BSA0230	☐	HIGH PUSH BAR			
SEAT TILT AND BACKREST REGULATION TYPE					
BSA0240	☐	ASSISTANT OPERATED GAS PISTON SEAT AND BACKREST TILT. - TILT LOCK MANDATORY		◆	
BSA0250	☐	ELECTRIC BY REMOTE CONTROL - COMPLIES WITH EN 12184 : 2014	NC : 1263, 1266		
BSA0605	☐	TILT LOCK (LOCKING OF PRE-SET ANGLE)		◆	

LEGRESTS / FOOTRESTS				
LEGRESTS / HANGERS TYPE				
BSA0260	☐	HANGER ANGLE ADJUSTABLE WITH KNEE PAD AND REINFORCED, SEPARATE CALF SUPPORT, PAIR		◆
FOOTRESTS TYPE				
BSA0280	☐	FOOTREST LOWER PART, DEPTH AND ANGLE ADJUSTABLE WITH SHORT TUBE AND FOOTBOARD CONVERTER, PAIR		○
BSA0290	☐	FOOTREST LOWER PART, DEPTH AND ANGLE ADJUSTABLE WITH LONG TUBE AND FOOTBOARD CONVERTER, PAIR		◆
ARMRESTS (Height adjustable, Autolock attachment included)				
ARMRESTS LH TYPE				
BSA0300	☐	REINFORCED, WIDE PAD, LH		◆
ARMRESTS RH TYPE				
BSA0320	☐	REINFORCED, WIDE PAD, RH		◆
SEATING OPTIONS				
ARMREST OPTIONS				
BSA0340	☐	SIDE SUPPORT COVER		
BSA0350	☐	CUSHION TO SIDE SUPPORT COVER		
SEAT OPTIONS				
BSA0360	☐	INCONTINENCE COVER		
BSA0365	☐	TILT SCALE		
BACKREST OPTIONS				
BSA0370	☐	HEADREST, SIDEWISE ADJUSTABLE	NC : 0400	○
BSA0390	☐	HEADREST WITH CHEEK SUPPORT, SIDEWISE ADJUSTABLE	NC : 0380, 0400	
BSA0380	☐	COVER FOR HEADREST		
BSA0400	☐	NECKREST		◆
OTHERS BACKREST OPTIONS				
BSA0415	☐	LATERAL POSITIONING PADS (WITH LOOP BAND BETWEEN COVER AND BACKREST UPHOLSTERY), PAIR		
BSA0420	☐	HANDLES CLINIC KIT FOR HEAD-/NECKREST		
BSA0430	☐	LUMBAR SUPPORT CUSHION		
BSA0440	☐	WEDGES (ONE PAIR)		
CHASSIS				
CHASSIS COLOUR				
BSA0441	☐	PEARL GREY		◆
BSA0442	☐	AZURITE		○
BSA0443	☐	SAND		○
BSA0444	☐	ELECTRIC GREEN		○
BSA0445	☐	HAPPY RED		○
CHASSIS SAFETY				
BSA0450	☐	CHASSIS EXTENDER, CASTORS 50 MM FORWARD.		
CASTOR ATTACHMENT TYPE				
BSA0460	☐	Fixed		◆
FRONT CASTORS				
BSA0480	☐	5,5"x1.3/4" - 140X45 BLACK SOFT		
BSA0495	☐	8"x2" - 200X45 BLACK SOLID TYRE		◆
REAR WHEELS SIZE				
REAR WHEELS SIZE				
BSA0514	☐	22" / 560 MM - QR REAR WHEEL ATTACHMENT		○
BSA0516	☐	24" / 600 MM - QR REAR WHEEL ATTACHMENT		◆
REAR WHEELS TYPE				
BSA0522	☐	REINFORCED (COLORLESS HUB, CLEAR ANODIZED RIM & STAINLESS STEEL SPOKES)		◆
BSA0524	☐	DRUM BRAKE (22/24" : COLORLESS HUB, CLEAR ANODIZED RIM & STAINLESS STEEL SPOKES / 16" : BLACK PLASTIC HUB & RIM)		○
TYRE TYPE				
BSA0535	☐	TYRES SOLID BLACK 1" 3/8		◆
HANDRIM TYPE				
BSA0540	☐	ALUMINIUM		◆
BSA0545	☐	PLASTIC COATED		
SPOKEGUARD				
BSA0550	☐	FROST		



SAFETY				
BRAKE TYPE				 0561  0590
BSA0560	☐	USER MANOUEVERED TYRE BRAKES		
BSA0561	☐	USER MANOUEVERED TYRE BRAKES, EXTENDED ARM		
SAFETY OPTIONS				 0570
BSA0570	☐	PELVIC BELT - PLASTIC BUCKLE		
BSA0590	☐	REFLECTORS KIT (WHEELS AND CHASSIS)		
BSA0595	☐	RETROFIT TRANSPORTATION BRACKET		
OPTIONS				
Think of Bodypoint positioning accessories: belts, harnesses, ankle supports.				
OTHER OPTIONS				
BSA0600	☐	TABLE WITH ELBOW SUPPORTS		 0600
RAMPS				
1654128	☐	Basic ramp 2 telescopic parts in aluminum and non-slip coating. Folded dimensions: 1210 x 194 mm. 8.8 kg per pair - maximum capacity: 270 kg		 1654128
1654129	☐	Quick2go Compact 3-part aluminum telescopic ramp with non-slip coating Folded dimensions: 870 x 237 mm 13.6 kg per pair - maximum capacity: 260 kg		 1654129
1654130	☐	Aluminium folding mono ramp Folded dimensions: 1824 x 714 mm 13.9 kg - maximum capacity: 300 kg		 1654130