

Rea® Azalea® Tall

Please contact DHCare EU Export Customer Services on :-
Tel : +33 247 626980 ; Email : serviceclient_export@dhcare.com
 Or visit our website : www.dhcare.eu.com



Rea® Azalea® Tall

Manual folding frame steel chair - for tall users



CE



Description:

- Seat tilt with maintenance of the center of gravity and the correct position of the legrests: +1°/+24° and of the backrest: 0°/30° by pneumatic jacks.
- Width-adjustable legrests.
- Heel straps.
- Frame colour: Pearl gray (RAL 0515).
- Weight: from 34 kg.
- Maximum user weight: 135 kg.
- 2-year warranty period excluding wear parts

Date :

Purchase Order :

Customer Code :




Distributor Stamp:

DHCare Unique Customization Service
 Please specify the quote or special production number :

☐ Please make your selection with a cross in the box

◆ Standard feature

○ Free of charge option

CODE		DESCRIPTION	NC : Non Compatible NC : 9999 = not compatible with the reference 9999	Euros		
BNH0019	☐	Rea [®] Azalea [®] Tall				
PROPULSION TYPE						
BNH0025	☐	SELF PROPELLED: 24"/ 22" WHEELS WITH USER BRAKES	NC : 1076	◆		
BNH0030	☐	SELF- AND ASSISTANT PROPELLED: 24" / 22" WHEELS WITH HUB BRAKES	NC : 1072, 1074			
SEAT						
SEAT DIMENSIONS						
BNH0040	☐	390 mm - Adjustable from 390 to 440 mm	NC : 0244, 0245, 0246, 0247	○		
BNH0050	☐	440 mm - Adjustable from 440 to 490 mm	NC : 0241, 0242, 0246, 0247	○		
BNH0060	☐	490 mm - Adjustable from 490 to 540 mm	NC : 0206, 0241, 0242, 0243, 0244	○		
Seat width pre-adjustable						
BNH0065	☐	Supplied with adjustment + 50 mm (at 440, 490 or 540 mm)		○		
Seat depth (SD)						
BNH0105	☐	Adjustable from 480 to 550 mm		◆		
Floor/seat height (measured without cushion)						
BNH0125	☐	500 mm		◆		
SEAT ANGLE						
BNH0127	☐	NEGATIVE SEAT ANGLE -1 TO +23º (STD +1 TO +24º)		○		
SEAT TYPE						
BNH0130	☐	DEPTH ADJUSTABLE		◆		
SEAT CUSHION						
BNH0165	☐	CUSHION FLOSHAPE		◆		
BNH0167	☐	CUSHION VICAIR ADJUSTER O2 (10 cm)	NC : 0191, 0201, 0202, 1225			
BNH0175	☐	CUSHION MISTRAL 2-LAYER	NC : 0191, 0201, 0202	○		
BNH0179	☐	CUSHION VICAIR VECTOR O2 100 mm	NC : 0191, 0201, 0202, 1225			
SEAT UPHOLSTERY						
BNH0191	☐	BLACK PLUSH TR35				
BNH0201	☐	BLACK PU TR26 - MUTLI STRETCH PU COATED FABRIC		◆		
BNH0202	☐	BLACK SPACER TR36		○		
BACKREST						
BACKREST WIDTH EXTENSION (PRE ADJUSTMENT)						
BNH0205	☐	+ 50 MM		○		
BNH0206	☐	+ 100 MM	NC : 0220	○		
BNH0207	☐	WITHOUT BACKREST WIDTH EXTENSION		◆		

BACKREST TYPE				
BNH0210	☐	FLEX 3 BACKREST, FIXED PLATE (WIDTH ADJUSTABLE)	NC : 0235, 0240, 0241 to 0247, 0341, 0343, 0905	☐
BNH0215	☐	BACKREST, PLATE "LAGUNA 2"	NC : 0235, 0240, 0241 to 0247, 0260, 0270, 0335, 0343, 0905	☒
BNH0220	☐	FLEX 3 BACKREST, TENSION ADJUSTABLE (WIDTH ADJUSTABLE)	NC : 0235, 0241 to 0247, 0341, 0343	☐
BNH0221	☐	MATRIX ELITE E2, (INCL. HARDWARE MOUNTED TO FRAME, BACKREST SUPPLIED IN BOX)	NC : 0240, 0250, 0260, 0270, 0311, 0321,0322, 0335, 0341, 0800, 0810, 0880, 0890, 0905, 0920	
BNH0223	☐	MATRIX ELITE DEEP E2, (INCL. HARDWARE MOUNTED TO FRAME, BACKREST SUPPLIED IN BOX)	NC : 0240, 0250, 0260, 0270, 0311, 0321, 0322, 0335, 0341, 0800, 0810, 0880, 0890, 0905, 0920	
BNH0225	☐	MATRIX POSTURE BACK, (INCL. HARDWARE MOUNTED TO FRAME, BACKREST SUPPLIED IN BOX)	NC : 0235, 0240, 0241 à 0247, 0250, 0260, 0270, 0311, 0321, 0322, 0335, 0341, 0800, 0810, 0880, 0890, 0905, 0920	
BACKREST HEIGHT MATRIX BACKRESTS				
BNH0233	☐	MATRIX HEIGHT 41 CM		☐
BNH0235	☐	MATRIX HEIGHT 46 CM		
BNH0234	☐	MATRIX HEIGHT 51 CM		☐
Matrix backrest width				
BNH0241	☐	360 mm - indicative bust 265 mm		☐
BNH0242	☐	380 mm - indicative bust 290 mm		☐
BNH0243	☐	410 mm - indicative bust 315 mm		☐
BNH0244	☐	430 mm - indicative bust 340 mm		☐
BNH0245	☐	460 mm - indicative bust 370 mm		☐
BNH0246	☐	480 mm - indicative bust 390 mm		☐
BNH0247	☐	510 mm - indicative bust 420 mm		☐
BACKREST COVER/ CUSHION				
BNH0240	☐	COVER, LATERAL		☐
BNH0250	☐	CUSHION, LAGUNA (LATERAL SUPPORT)	NC : 0905	☒
BNH0260	☐	CUSHION, MISTRAL 2 (WAIST SUPPORT)	NC : 0322, 0905	
BNH0270	☐	CUSHION, PASSAD 2 (SHOULDER SUPPORT)	NC : 0905	
BNH0276	☐	VICAIR MULTIFONCTION BACK	NC : 0305, 0311, 0321, 0322, 1225	
BNH0277	☐	VICAIR FLEX BACK	NC : 0305, 0311, 0321, 0322, 1225	
BACKREST UPHOLSTERY				
BNH0311	☐	BLACK PLUSH TR35		
BNH0321	☐	BLACK PU TR26		☒
BNH0322	☐	BLACK SPACER TR36		☐
PUSH HANDLES				
BNH0335	☐	HEIGHT ADJUSTABLE AND REMOVABLE PUSHHANDLES		☒
BNH0341				
BNH0343				
BNH0345	☐	ANGLE ADJUSTABLE PUSHBAR	NC : 1266	
SEAT TILT AND BACKREST REGULATION TYPE				
BNH0350	☐	ASSISTANT OPERATED GAS PISTON SEAT AND BACKREST TILT.		☒
BNH0360	☐	ELECTRIC BY REMOTE CONTROL - COMPLIES WITH EN 12184 : 2014	NC : 1263, 1266	
BNH1263	☐	TILT LOCK (LOCKING OF PRE-SET ANGLE)		☒
LEGRESTS / FOOTRESTS				
CENTRAL LEGRESTS TYPE				
BNH0370	☐	CENTRAL LEGREST, INCL CALFPAD AND FOOTREST.	NC : 0376, 0378, 0380, 0390, 0410, 0430, 0440, 0460, 0476, 0520, 0580, 0720, 0941, 0942, 0952, 0954, 0956	
LEGRESTS / HANGERS ADJUSTMENT				
BNH0375	☐	LOW POSITION, STANDARD		☒
BNH0376	☐	HIGH POSITION FOR LONGER LEGS		
BNH0378	☐	90° HANGER ATTACHMENT	NC : 0390, 0410, 0440, 0460	



0210



0215



0220



Matrix Elite E2 /
Elite Deep E2



0240



0250



0260



0270






0345





0370

Matrix Elite E2 /
Elite Deep E2

VICAIR®

LEGRESTS / HANGERS LH TYPE					 0380 / 0430  0390 / 0440  0410 / 0460	
BNH0380	☐	HANGER ANGLE ADJUSTABLE LH. WITH KNEE PAD AND CALF SUPPORT	NC : 0475, 0490, 0942	◆		
BNH0390	☐	HANGER FIXED 80 DEG, SWING IN SWING OUT - ALU 2003, LH		○		
BNH0410	☐	LEGREST AMPUTEE - ALU 2008 LH	NC : 0475, 0476, 0490, 0520, 0941, 0942			
LEGRESTS / HANGERS RH TYPE						
BNH0430	☐	HANGER ANGLE ADJUSTABLE RH. WITH KNEE PAD AND CALF SUPPORT	NC : 0475, 0550, 0942	◆	 0390 / 0440  0410 / 0460	
BNH0440	☐	HANGER FIXED 80 DEG, SWING IN SWING OUT - ALU 2003, RH		○		
BNH0460	☐	LEGREST AMPUTEE - ALU 2008 RH	NC : 0475, 0476, 0550, 0580, 0941, 0942			
FOOTRESTS						
FOOTRESTS ONE PIECE						 0476  0520 / 0580  0520 / 0580
BNH0475	☐	FOOTREST ONE PIECE DEPTH AND ANGLE ADJUSTABLE, NO WIDTH ADJUSTMENT - SHORT TUBE	NC : 0490, 0520, 0550, 0580, 0930, 0940, 0952, 0954, 0956			
BNH0476	☐	FOOTREST ONE PIECE DEPTH AND ANGLE ADJUSTABLE, NO WIDTH ADJUSTMENT - LONG TUBE	NC : 0520, 0580, 0952, 0954, 0956			
FOOTRESTS LH TYPE						
BNH0490	☐	FOOTREST LOWER PART, DEPTH AND ANGLE ADJUSTABLE WITH SHORT TUBE AND HEEL STRAP, LH		○	 0520 / 0580  0520 / 0580	
BNH0520	☐	FOOTREST LOWER PART, DEPTH AND ANGLE ADJUSTABLE WITH LONG TUBE AND HEEL STRAP, LH		◆		
FOOTRESTS RH TYPE						
BNH0550	☐	FOOTREST LOWER PART, DEPTH AND ANGLE ADJUSTABLE WITH SHORT TUBE AND HEEL STRAP, RH		○	 0520 / 0580  0520 / 0580	
BNH0580	☐	FOOTREST LOWER PART, DEPTH AND ANGLE ADJUSTABLE WITH LONG TUBE AND HEEL STRAP, RH		◆		
ARMRESTS (Height adjustable, Autolock attachment included)						
ARMRESTS LH TYPE						 0600 / 0630  0610 / 0640
BNH0600	☐	STANDARD PAD, LH	NC : 1232	○		
BNH0606	☐	CONCAVE GUTTER ADJUSTABLE IN HEIGHT AND DEPTH 70 MM WITHOUT TOOLS (COMPATIBLE WITH SHELF)	NC : 1232	◆		
BNH0610	☐	PADDED, LH.	NC : 0680, 0690, 1230, 1231, 1233, 1250	○		
BNH0615	☐	HEMIPLEGIC ARMRESTS, LH.	NC : 1230, 1231, 1233, 1240, 1245, 1250	○	 0610 / 0640  0615 / 0645	
BNH0620	☐	WITHOUT ARMREST, LH	NC : 0690, 1230, 1231, 1240, 1245, 1250	○		
ARMRESTS RH TYPE						
BNH0630	☐	STANDARD PAD, RH	NC : 1233	○	 0610 / 0640  0615 / 0645	
BNG0636	☐	CONCAVE GUTTER ADJUSTABLE IN HEIGHT AND DEPTH 70 MM WITHOUT TOOLS (COMPATIBLE WITH SHELF)	NC : 1232	◆		
BNH0640	☐	PADDED, RH.	NC : 0680, 0690, 1230, 1231, 1232, 1250	○		
BNH0645	☐	HEMIPLEGIC ARMRESTS, RH.	NC : 1230, 1231, 1232, 1240, 1245, 1250	○		
BNH0650	☐	WITHOUT ARMREST, RH	NC : 0690, 1230, 1231, 1240, 1245, 1250	○		
SEATING OPTIONS						
ARMREST OPTIONS						 0720
BNH0680	☐	CUSHION TO SIDE REST COVER				
BNH0690	☐	PAD, SIDE REST				
SEAT OPTIONS						
BNH0710	☐	INCONTINENCE COVER				
BNH0720	☐	ABDUCTION CUSHION				

BACKREST OPTIONS					
BNH0730	☐	HEADREST, SIDEWISE ADJUSTABLE - FIXED HARDWARE INCLUDED	NC : 0742, 0744, 0746, 0750, 0770	☐	 0730
BNH0750	☐	HEADREST WITH CHEEK SUPPORT, SIDEWISE ADJUSTABLE - FIXED HARDWARE INCLUDED	NC : 0740, 0742, 0744, 0746, 0770, 0780		 0750
BNH0731	☐	MATRIX ELAN STANDARD PAD (254 mm/10")	NC : 0740, 0770		
BNH0736	☐	MATRIX ELAN 4-POINT PAD (229W x 177H mm / 9"W x 7"H)	NC : 0740, 0770		 0736
BNH0742	☐	MATRIX ELAN HARDWARE			0731/ 0742
BNH0740	☐	COVER FOR HEADREST			
BNH0770	☐	NECKREST	NC : 0771	☐	 0770
BNH0782	☐	UNIVERSAL HEAD/NECKREST HOLDER			
BNH0790	☐	HANDLES CLINIC KIT FOR HEAD-/NECKREST			
TRUNK SUPPORT LH					
BNH0800	☐	TRUNK SUPPORT MULTI ADJUSTABLE WITH STANDARD ARM (LH)			0782
BNH0880	☐	TRUNK SUPPORT SWING-AWAY TR26 LEFT CPLT (V2025)& BRACKET, FIX BACK, PASSIVE		NEW	 0800 / 0810
TRUNK SUPPORT RH					
BNH0810	☐	TRUNK SUPPORT MULTI ADJUSTABLE WITH STANDARD ARM (RH)			 0880 / 0890
BNH0890	☐	TRUNK SUPPORT SWING-AWAY TR26 RIGHT CPLT (V2025)& BRACKET, FIX BACK, PASSIVE		NEW	
BACKREST OPTIONS (Following)					
BNH0905	☐	LATERAL POSITIONING PADS (WITH LOOP BAND BETWEEN COVER AND BACKREST UPHOLSTERY), PAIR			 0905
BNH0915	☐	TILT SCALE, BACKREST RECLINE			 0952 / 0954 / 0956
BNH0920	☐	WEDGES (ONE PAIR)			
LEGREST OPTIONS					
BNH0941	☐	COVER FOOTPLATE			
BNH0942	☐	CALF STRAP PADDED			
BNH0952	☐	BODYPOINT ANKLE HUGGERS LARGE	MOUNTING KIT INCLUDED		
BNH0954	☐	BODYPOINT ANKLE HUGGERS MEDIUM	MOUNTING KIT INCLUDED		
BNH0956	☐	BODYPOINT ANKLE HUGGERS SMALL	MOUNTING KIT INCLUDED		
CHASSIS					
CHASSIS COLOUR					
BNH0944	☐	PEARL GREY		☐	 0944
BNH0946	☐	AZURITE		☐	 0946
BNH0947	☐	SAND		☐	 0947
BNH0948	☐	ELECTRIC GREEN		☐	 0948
BNH0949	☐	HAPPY RED		☐	 0949
CASTOR ATTACHMENT TYPE					
BNH0960	☐	FIXED		☐	
FRONT CASTORS					
BNH1015	☐	8"x1.1/4" - 200X30 BLACK SOLID TYRE		☐	
BNH1022	☐	8"x1.1/4" - 200X32 PNEUMATIC BLACK		☐	
BNH1025	☐	8"x2" - 200X45 BLACK SOLID TYRE		☐	
REAR WHEELS					
REAR WHEELS SIZE					
BNH1062	☐	24" / 600 MM - QR REAR WHEEL ATTACHMENT		☐	
REAR WHEELS TYPE					
BNH1072	☐	STANDARD (BLACK HUB, COLORLESS RIM & ZINC PLATTED SPOKES)	NC : 1120	☐	
BHN1074	☐	REINFORCED (COLORLESS HUB, CLEAR ANODIZED RIM & STAINLESS STEEL SPOKES)		☐	
BHN1076	☐	DRUM BRAKE (22/24" : COLORLESS HUB, CLEAR ANODIZED RIM & STAINLESS STEEL SPOKES / 16" : BLACK PLASTIC HUB & RIM)	NC : 1086	☐	
TYRE TYPE					
BNH1082	☐	PNEUMATIC PROFILED BLACK 1" 3/8		☐	
BNH1086	☐	PNEUMATIC LOW PROFILE BLACK 1"		☐	
BNH1092	☐	TYRES SOLID BLACK 1" 3/8		☐	
HANDRIM TYPE					
BNH1100	☐	ALUMINIUM		☐	
BNH1120	☐	PLASTIC COATED		☐	
WHEEL OPTIONS					
SPOKEGUARD					
BNH1160	☐	FROST		☐	

SAFETY					
BNH1170	☐	USER MANOUEVERED TYRE BRAKES		◆	 1200  1220
BNH1180	☐	ONE ARM BRAKE LH			
BNH1185	☐	ONE ARM BRAKE RH			
BNH1200	☐	PELVIC BELT - PLASTIC BUCKLE			
BNH1201	☐	BODYPOINT BELT - 2 POINT CENTER-PULL PUSH-BUTTON MEDIUM			
BNH1202	☐	BODYPOINT BELT - 2 POINT CENTER-PULL PUSH-BUTTON LARGE			
BNH1203	☐	BODYPOINT BELT - 4 POINT CENTER-PULL PUSH-BUTTON MEDIUM			
BNH1204	☐	BODYPOINT BELT - 4 POINT CENTER-PULL PUSH-BUTTON LARGE			
BNH1206	☐	BODYPOINT PIVOTFIT HARNESS MEDIUM	ONLY WITH ELITE E2 /ELITE DEEP E2 BACKS OR WITH HEADREST OR NECKREST		 1201 / 1202  1203 / 1204  1206 / 1207
BNH1207	☐	BODYPOINT PIVOTFIT HARNESS LARGE	ONLY WITH ELITE E2 /ELITE DEEP E2 BACKS OR WITH HEADREST OR NECKREST		
OPTIONS					
Think of Bodypoint positioning accessories: belts, harnesses, ankle supports.					
OTHER OPTIONS					
TRAYS					
BNH1230	☐	TABLE TRAY INCL MOUNTING KIT	NC : 1237		 1230  1231
BNH1231	☐	TABLE TRAY INCL MOUNTING KIT, WITH ELBOW PADS	NC : 1237		
BNH1232	☐	HALF TABLE TRAY, CLEAR, FLIP UP, LH (ARMREST WITH STRAIGHT NARROW PAD INCL.)	NC : 1235		 1232 / 1233
BNH1233	☐	HALF TABLE TRAY, CLEAR, FLIP UP, RH (ARMREST WITH STRAIGHT NARROW PAD INCL.)	NC : 1235		
BNH1240	☐	SWING AWAY TABLE	NC : 1235, 1237, 1250		 1245
BNH1245	☐	SWING AWAY TABLE WITH ELBOW PADS	NC : 1235, 1237, 1250		
TABLE TRAY CUSHION					
BNH1235	☐	TABLE TRAY CUSHION			 1237
BNH1237	☐	PADDING			
OTHER OPTIONS					
BNH1250	☐	LOCKING DEVICE FOR TABLE TRAY			 1245
BNH1266	☐	MOUNTING BRACKET VIAMOBIL			
BNH1300	☐	HARNESS ATTACHMENT			 1237
BNH1310	☐	ATTACHMENT WHEN USING THE PELVIC BELT FOR POSITIONING			
BNH1320	☐	DRIP STAND			
RAMPS					
1654128	☐	Basic ramp 2 telescopic parts in aluminum and non-slip coating. Folded dimensions: 1210 x 194 mm. 8.8 kg per pair - maximum capacity: 270 kg			 1654128
1654129	☐	Quick2go Compact 3-part aluminum telescopic ramp with non-slip coating Folded dimensions: 870 x 237 mm 13.6 kg per pair - maximum capacity: 260 kg			 1654129
1654130	☐	Aluminium folding mono ramp Folded dimensions: 1824 x 714 mm 13.9 kg - maximum capacity: 300 kg			 1654130