

Cirrus G5



ORDER FORM

Errors and imissions excepted.

Billing Address Or Customer-ID			Shipping Address	Identical To Billing Address
C.-ID:			C.-ID:	
Name:			Name:	
Street:			Street:	
Town/PC:			Town/PC:	
E-Mail:			E-Mail:	
Phone:			Phone:	
Order Information			User Weight	
Order	Calculation	Built-4-Me [®]	User Weight Capacity: 135 kg	
Purchase Order No.:			User Weight: kg	
Art. No.: 35001000		Reimbursement Code: xxx	<h1>Base Price</h1> <h1>€,-</h1>	
Order Date:	Day	Month		

Please follow this measuring method:

The following method is recommended on the basis of extensive practical experience. Deviations are possible in individual instances.

Seat Width (SW):

- a) For the user: Allow space between sideguard and hip on each side of the chair (approx. 10 - 20 mm altogether).
- b) On the wheelchair: From outside seat frame tube to outside seat frame tube.

Seat Depth (SD):

- a) For the user: When sitting upright: leave space for 2 - 3 fingers between front edge of the seat upholstery and back of knee (approx. 30 - 60 mm). Care must be taken with adj. upholstery.
- b) On the wheelchair: From back tube to front seat upholstery.

Front Seat Height:

- a) For the user: When sitting with 90° knee angle: from back of knee to bottom of shoe sole.
- b) On the wheelchair: From front of seat tube to the ground. Min. 30 mm longer than user measurement. Remember to include the cushion height.

Rear Seat Height:

On the Wheelchair: From back of seat tube to the ground.

Note: Max. seat height difference rear to front is 110 mm

Back Height:

- a) For the user: From the seat, including cushion, to be no more than 10 - 20 mm below the lower shoulder blade.
- b) On the wheelchair: From back of seat tube to upper edge of backrest upholstery.

Note: The user's ability to sit upright has a major influence on the back height and it can therefore differ considerably

Back angle:

Important: The back height depends on the degree of disability. 0° is the angle 90° to the ground. Negative angles are <90° and positive angles are >90°.

Centre of Gravity:

On the wheelchair: From the back tube to the middle of the axle plate.

Technical Data

Max. Overall Height:	1100 mm
Max. Overall Length:	1150 mm
Overall Width:	SW + 240 mm
Weight:	48 kg

Product is crash tested according to ISO 7176-19 and therefore approved for transportation in a vehicle in combination with crash test option. For transport without a user, fix the chair safely in the vehicle (e.g. car or plane).

Seat Width

	SW39	SW42	SW45	SW48	SW51	SW54
mm	390	420	450	480	510	540

Seat Depth

	SD44	SD47	SD50			
mm	440 (390-490)	470 (420-520)	500 (450-550)			

Front Seat Height (1)

	FSH375	FSH40	FSH425	FSH45	FSH475	FSH50	FSH525
mm	375	400	425	450	475	500	525
5" / 16" (9)			-	-	-	-	-
6" / 16" (9)	-	-				-	-
8" / 16" (9)	-	-	-			-	-
5" / 20"			-	-	-	-	-
6" / 20"	-	-					-
8" / 20"	-	-	-				
5" / 22"	-		-	-	-	-	-
6" / 22"	-	-					-
8" / 22"	-	-	-				
5" / 24"	-	-	-	-	-	-	-
6" / 24"	-	-					-
8" / 24"	-	-	-				

Seat (1)

€

CIR020086	Seat, tilting (tilt range 0° - 33°)	-
	CIR020088 Locking mechanism tilt	-
CIR020089	Seat plate, fixed, +/-5 cm seat depth	-
	CIR020091 Seat plate, without pelvic hole (2)	-
CIR020090	Seat, sliding, +/-5 cm seat depth (3)	-
	CIR020092 Seat plate, with pelvic hole	-

Seat Cushions

€

CIR020093	Seat cushion, Trevira	-
CIR020094	Seat cushion, Nylon, PU coating	-
CIR020110	Omit seat cushion	-

SEDEO PRO Seat Cushions

CIR020073	SEDEO PRO, seat cushion, profiled, height 70 mm (10)(12)	-
CIR020075	SEDEO PRO, seat cushion, support, height 90 mm (10)	-
CIR020078	SEDEO PRO, seat cushion, stabilizer, height 90 mm (10)	-
CIR020076	SEDEO PRO, seat cushion, support memory, height 90 mm (11)	-

SEDEO PRO Seat Cushion Cover

CIR020518	Seat cushion cover, incontinence protection	-
CIR020080	Seat cushion cover, comfort	-

JAY Cushions

€

JAY020006	JAY Easy Visco F (4)	-
JAY020015	JAY Balance, Fluid (7)	-
JAY020020	JAY J2 Deep Contour (8)	-

(1) Seat height will be 25 mm higher with CIR020090 Seat, sliding, +/- 5 cm seat depth
(2) Not in combination with CIR020090 Seat, sliding, +/- 5 cm seat depth
(3) Not in combination with CIR020089 Seat plate, fixed, +/- 5 cm seat depth
(4) Medium skin protection

(7) High skin protection
(8) Very high skin protection, only for SW 390/450/510 mm

(9) Only in combination with CIR070008 Rear wheel, drum brakes, 36 spokes (silver), cross-wise spokes
(10) Not with SD44
(11) SD50 only
(12) Only in combination with Seat, sliding CIR020090

Seat cushions & Seat plate



// CIR020093



// CIR020073



// CIR020075



// CIR020076




// CIR020078



// JAY020006



// JAY020020




// JAY020009



// JAY020024



// CIR20092



// CIR20091

Seat Accessories €

CIR020070	Incontinence cover (1)	-
CIR090202	Abduction pommel, small, 70 x 60 x 100 mm	-
CIR090302	Abduction pommel, large, 90 x 90 x 100 mm	-

Positioning Belt with lock

CIR090135	Positioning belt, carseat buckle	-
CIR090156	Upholstery positioning belt, small	-
CIR090136	Upholstery positioning belt, large	-

JAY Pelvic Belts

JAY090071	2 Point Anterior Pelvic Belt 40 mm Push Button - Single Pull 190 mm Pad	,-
JAY090072	2 Point Anterior Pelvic Belt 40 mm Push Button - Single Pull 240 mm Pad	-
JAY090077	JAY Medial 4-Point Pelvic Belt, Push Button Buckle, Pad 140 x 40 mm	,-
JAY090078	JAY Medial 4-Point Pelvic Belt, Push Button Buckle, Pad 165 x 40 mm	,-
JAY090102	JAY V-Contour Pelvic Belt, Side Release Buckle, Pad, Medium	-
JAY090103	JAY V-Contour Pelvic Belt, Side Release Buckle, Pad, Large	-

Backrest €

Backrests

CIR030123	Backrest recline (0° - 35°), gas strut	-
	BH55 Back height 550 mm	-
	BH65 Back height 650 mm	-
CIR030134	Backrest, reclining, 0° - 35°, gas strut, for JAY, SEDEO and solid backs	-
CIR030124	Locking mechanism recline	-

Push Handles

CIR030126	Push bar	-
CIR030070	Push bar, angle-adjustable (2)	-
CIR030204	Push handles, height-adjustable	,-

Back Cushions

CIR030127	Back cushion, Trevira	,-
CIR030128	Back cushion, Nylon, PU coating	-
CIR030209	Omit back cushion	-

SEDEO LITE Back Cushions (4)

CIR030130	SEDEO LITE, back cushion, profiled, height 500 mm, comfort cover	-
CIR030131	SEDEO LITE, back cushion, profiled, adjustable straps, height 570 mm, comfort cover	-
CIR030132	SEDEO LITE, back cushion, deep contour, passive, adjustable straps, height 570 mm, comfort cover	-
CIR030133	SEDEO LITE, back cushion, deep contour, active, adjustable straps, height 570 mm, comfort cover	-

Headrest €

CIR030064	Headrest Standard, small, 139 x 210 x 40 mm	-
	CIR030137 Receiver, vertical and horizontal angle adjustable	-
	CIR030138 Cover, padded	-
CIR030066	Headrest, lateral guidance, small, 140 x 230 x 70	-
	CIR030137 Receiver, vertical and horizontal angle adjustable	-
	CIR030138 Cover, padded	-
CIR030135	Head-/neckrest, comfort, tension adjustable, Trevira	,-
	CIR030137 Receiver, vertical and horizontal angle adjustable	-
CIR030136	Head-/neckrest, comfort, tension adjustable, PU coated	,-
	CIR030137 Receiver, vertical and horizontal angle adjustable	-

(1) Not with SEDEO & JAY seat cushions

(2) Not in combination with CIR030134 Backrest, reclining, 0° - 35°, gas strut, for JAY and solid backs

(4) Only with CIR030134 Backrest, reclining, 0° - 35°, gas strut, for JAY, SEDEO LITE and solid backs

Headrest

Whitmyer Headrests

WHT030406	Whitmyer, Headrest Pad, Adjust-A-Plush, Large, (LxH) 460 x 127 mm (1)	
WHT030412	Whitmyer, Headrest Pad, Contoured Cradle, Medium (LxH) 250 x 158 mm (1)	
WHT030413	Whitmyer, Headrest Pad, Contoured Cradle, Large (LxH) 330 x 190 mm (1)	
WHT030414	Whitmyer, Headrest Mount, Cirrus G5, vertical and angle adjustable	

Seat accessories



// CIR090137/144/150/151/152



// JAY090071



// CIR090302/202



// JAY090077



// JAY090102

Backrest



// CIR30127



// CIR030126



// CIR030070



// CIR030204

(1) Including WHT030414 Whitmyer, Headrest Mount, Cirrus G5, vertical and angle adjustable

JAY J3 Backs (For complete product range see JAY order forms)

€

JAY J3 Back, Shallow Contour (SC)	Short (S)		Medium (M)		Tall (T)		
Upper Thoracic (UT)	JAY030012	420 mm	JAY030019	460 mm	JAY030025	500 mm	-
Shoulder Height (SH)	JAY030015	530 mm	JAY030022	570 mm	JAY030028	610 mm	-
JAY J3 Back, Mid Deep Contour (MDC)	Short (S)		Medium (M)		Tall (T)		
Upper Thoracic (UT)			JAY030042	460 mm			-
Shoulder Height (SH)	JAY030043	530 mm					-
JAY J3 Back, Deep Contour (DC)	Short (S)		Medium (M)		Tall (T)		
Upper Thoracic (UT)	JAY030014	420 mm	JAY030021	460 mm	JAY030027	500 mm	-
Shoulder Height (SH)	JAY030017	530 mm	JAY030024	570 mm	JAY030030	610 mm	-

JAY J3 Back Headrests

€

JAY030004	JAY Flip Back Headrest	-
JAY030006	JAY Std. Headrest	-
WHT030415	Whitmyer Headrest, COBRA XTRA, flip back (only with JAY J3 back)	-
WHT030416	Whitmyer Headrest, AXYS, Standard (only with JAY J3 back)	-

JAY J3 Lateral Thoracic Pads

€

JAY090024	JAY J3 Lateral Thoracic Pad, (LxH) 152 mm x 127 mm, Right	-
JAY090015	JAY J3 Lateral Thoracic Pad, (LxH) 152 mm x 127 mm, Left	-
JAY090025	JAY J3 Lateral Thoracic Pad, (LxH) 152 mm x 152 mm, Right	-
JAY090016	JAY J3 Lateral Thoracic Pad, (LxH) 152 mm x 152 mm, Left	-

Back Accessories

€

CIR030139	Positioning wedge kit	-
CIR030143	Lateral, thoracic support, height- and depth-adjustable, swing away, right	-
CIR030142	Lateral, thoracic support, height- and depth-adjustable, swing away, left	-
CIR030144	Cover, lateral thoracic support, padded (2)	-
JAY090087	Classic Style Anterior Trunk Support - Dynamic, Medium	-
JAY090088	Classic Style Anterior Trunk Support - Dynamic, Large	-
JAY090104	Contour Style Anterior Trunk Support - Dynamic, Medium	-
JAY090105	Contour Style Anterior Trunk Support - Dynamic, Large	-

Headrests



// CIR030064



// CIR030066



// CIR30135



// WHT030412



// WHT030406



// WHT030414

JAY J3 Backs



// Shallow Contour (SC)



// Deep Contour (DC)



// Mid Deep Contour (MDC)

Back accessories



// CIR030143



// JAY090087



// JAY090104

Electric back angle and seat tilt adjustment €

CIR020071	Seat tilt adjustment, powered	-
CIR030145	Back angle adjustment, powered	-
CIR020072	Lock for powered adjustment	-

Armrests €

CIR040042	Single-post sideguard, armpad (400 mm), tool height-adjustable, removable	-
	CIR040031 Armpad hemiplegia, integral foam, right	-
	CIR040032 Armpad hemiplegia, integral foam, left	-
	CIR020095 Seat width extension (+ 30 mm)	-
CIR040043	Single-post sideguard, armpad (440 mm), height-adjustable, reclining, flip back, removable (1)	-
	CIR040031 Arm pad hemiplegia, integral foam, right	-
	CIR040032 Arm pad hemiplegia, integral foam, left	-
CIR040000	Without sideguard	-

Armrest accessories

CIR040044	Cover arm pad, padded, PU (pair)	-
CIR040045	Cover, side panel, half-padded, PU (pair) (2)	-
CIR040046	Cover, side panel, complete padded, PU (pair) (2)	-

Tray tables

CIR090131	Tray table, transparent, slide on, width-adjustable 390 - 540 mm	-
	CIR090125 Partial padding for tray table (6)	-
CIR090139	Tray table, Half-Side, slide on, right (7)	-
CIR090140	Tray table, Half-Side, slide on, left (8)	-

Footrest & Footplates €

Hanger aluminium, angle adjustable (70°/80°/90°)

CIR050117	Hanger, aluminium, tool angle-adjustable 70°/80°/90°, swing-away, right	-
CIR050118	Hanger, aluminium, tool angle-adjustable 70°/80°/90°, swing-away, left	-

Hanger aluminium, flip-up (0° - 90°)

CIR050010	Hanger, aluminium, flip-up (0° - 90°), swing away, right	-
CIR050009	Hanger, aluminium, flip-up (0° - 90°), swing away, left	-

Amputee support

CIR050013	Amputee support, padded, right	-
CIR050014	Amputee support, padded, left	-

Cast support

CIR050094	Cast support, right	-
CIR050093	Cast support, left	-

Omit hanger and footplate

CIR050000	Omit hanger and footplate	-
-----------	---------------------------	---

Footplates

CIR050119	Footboard, divided, composite, angle-adjustable, flip-up (sideways), black	-
CIR050120	Footboard, divided, composite, depth- and angle-adjustable, flip-up (sideways), black	-
CIR050121	Footboard, platform, aluminium, angle-adjustable, black (3)	-

Footbox (3)

CIR050122	Footbox, padded, PU coated	-
-----------	----------------------------	---

Footrest accessories

CIR050144	Clinic kit for lower leg length adjustment, without tools	-
CIR050145	Padding for footboard, PU coated (3)	-
CIR050123	Calfpad, padded, not split (3)	-
CIR050124	Knee pad, PU coated (3)	-
CIR050134	Universal padding, PU coated (4)	-
CIR050133	Footplate extension, aluminium, black (5)	-

(1) Not with CIR030134 Backrest, reclining, 0° - 35°, gas strut, for JAY, SEDEO LITE & solid backs
 (2) Only with Single-post sideguard, armpad (440 mm), height-adj., reclining, flip back, removable
 (3) Only for CIR050009-10 Hanger, aluminium, flip-up (0° - 90°), swing away left/right

(5) Only CIR050119 Footboard, divided, composite, angle-adjustable, flip-up (sideways), black & CIR050120
 Footboard, divided, composite, depth- and angle-adjustable, flip-up (sideways), black
 (6) Not with Arm pad hemiplegia, integral foam

(4) E.g. for footplate, amputation support, leg support

Armrests & Accessories



// CIR090131



// CIR090140



// CIR040044



// CIR040045



// CIR040046



// CIR040042



// CIR040031

Footrests & Footplates



// CIR050117



// CIR050010



// CIR050013



// CIR050094



// CIR050119



// CIR050122



// CIR050145



// CIR050124



// CIR050134



// CIR050123

Center Of Gravity (COG)		€
COG040	-40 mm (Safe)	-
COG00	0 mm (Standard)	-
COG40	40 mm (Active)	-

Frame accessories

Tip assist (step tube)		€
CIR090003	Tip assist, pair	-

Anti-tip tubes swing away		€
CIR090006	Anti-tip tubes, swing away, pair	-
CIR090060	Crutch holder, right	-
CIR090039	Crutch holder, left	-

Castor Wheels

CIR080301	Castor, solid (hard), rim: composite, black	-
CIR080305	Castor, pneumatic, rim: composite, black (1)	-

Rear Wheels

CIR070073	Rear wheel, 36 spokes (silver), cross-wise spoked (2)	-
CIR070008	Rear wheel, drum brakes, 36 spokes (silver), cross-wise spokes (7)	-

Tyres

CIR070076	Tyre, 1", pneumatic, tread (3)	-
CIR070077	Tyre, 1 3/8", pneumatic, tread (4)	-
CIR070074	Tyre, 1", puncture proof, tread (3)	-
CIR070075	Tyre, 1 3/8", puncture proof, tread	-

Handrims

CIR070202	Handrims, aluminium, anodised, silver (4)	-
CIR070204	Handrims, stainless steel, silver (4)	-
CIR070207	Handrims, Supergrip® (4)	-
CIR070208	Handrims, Maxgrepp® (4) (5)	-
CIR070251	Handrim, cover, silicon (4) (6)	-

(1) Only for 6" & 8" castors

(2) Only for 20", 22" and 24" rear wheels

(3) Not with drumbrake, 16" wheel

(4) Not with 16" wheel

(5) Not with 20" wheel

(6) Not with CIR070207 Handrims, Supergrip® and CIR070208 Handrims, Maxgrepp®

(7) 16" rear wheel with composite rim

Frame accessories



// CIR090003



// CIR090006



// CIR090060

Rear wheels & Castors



// CIR070008



// CIR070073



// CIR080301



// CIR070208



// CIR070207

Spoke Guards (1) €

CIR070831	Spoke guards, transparent	-
-----------	---------------------------	---

Wheel locks (Choose only one option) €

CIR060009	Brake, pull to lock, brake lever long	-
-----------	---------------------------------------	---

One arm wheel lock

CIR060021	Brake, pull to lock, one-arm, right actuation	-
CIR060022	Brake, pull to lock, one-arm, left actuation	-

Accessories €

CIR090000	Tool kit	-
CIR090149	Air pump, powered	-
CIR090025	Air pump	-
CIR090030	Back pack	-
CIR090083	IV support with multiposition clamp	-
CIR090107	Oxygen bottle holder	-
CIR090148	Universal headrest fixation	-

Colours €

CIR100134	Matt-Anthracite	-
-----------	-----------------	---

Wheel locks



// CIR060009



// CIR060021

Accessories



// CIR090148



// CIR090083



// CIR090107

Comments for personalization


 built4me

Cirrus G5

€

Base price	<input type="text"/>
Total Of Additional Equipment	<input type="text"/>
Packaging And Transport	<input type="text"/>
Exclusive VAT	<input type="text"/>
VAT	<input type="text"/>
Total price	<input type="text"/>

Date

Signature

Day Month Year

**Sunrise Medical HCM**

Phone: +31 (0) 492 593 888

Fax: +31 (0) 492 546 497

E-Mail: customerservice@sunrisemedical.nlwww.SunriseMedical.eu