

Address

Fax +49 7261-7351-10

Shipping address

Berollka-aktiv Rollstuhltechnik GmbH
 Jahnstraße 16
 74889 Sinsheim - Germany
Phone: +49 7261-7351-0
 info@berollka.de
 www.berollka.de



Com.:

- Order
 Estimate

Basic Price 1860,- €

Price valid from February 1st 2022
 (until revoked)

Cruiser

mobility base for seating system
 for indoor and outdoor



Όλες οι αναγραφόμενες τιμές δεν περιλαμβάνουν μεταφορικά έξοδα και Φ.Π.Α

Standard equipment:

- rigid aluminium frame with abduction
- back with hand bow, angle adjustable, height adjustable and tiltable
- 140mm / 7" castors
- 12" rear wheels with drum brake

capacity up to 90kg

user weight + seating system: _____ kg

basic measure

frame		frame tube	width of seating frame A	width of frame front B	total width C
Cruiser mobility base for seating shells	frame size 1 <input type="checkbox"/> Code U2401	25mm	38 cm	51 cm	53 cm
	frame size 2 <input type="checkbox"/> Code U2402	25mm	44 cm	57 cm	59 cm

frame colour

technical data

- 25 RAL 9010 white
- 24 RAL 1028 melon yellow
- 27 mandarin met.
- 06 dormant red
- 43 RAL 4010 telemagenta met.
- 44 RAL 5015 sky blue
- 16 RAL 5002 light blue
- 17 dormant blue
- 45 RAL 6018 yellow green
- 19 silver metallic
- 20 grey metallic
- 22 RAL 9011 black mat
- 23 RAL 9005 black



met. = metallized

Seat tilt **Back incline**
-5 to 30° **90 to 120°**

Weight **min. seat height**
from 12 kg **44 cm**

concerning seat frame (trapezium adapter +2cm)





- H** total height: 95 - 114cm
- L1** total length wheels: max. 91,5cm
- L2** total length with trapezium bow: 101cm
- L2** total length with U-bow: 95cm
- T** Seat depth: 33 - 40cm




Because of the material different colours are possible

frame accessories



armrest

<p>Tip shackle</p> <p><input type="checkbox"/> right <input type="checkbox"/> left</p>  <p>Code U2411</p> <p><input type="checkbox"/> Piece 78,- €</p>	<p>Anti-tip wheel with adapter</p> <p>max. 2 part s possible</p> <p><input type="checkbox"/> right <input type="checkbox"/> left</p>  <p>Code U2412</p> <p><input type="checkbox"/> Piece 130,- €</p>	<p>Armrest holder with arm pads fixed at the back tube</p>  <p>(arm pads swing up)</p> <p>Code U2446</p> <p><input type="checkbox"/> Pair 268,- €</p>	<p>Armrest adjustable in height with raster, for mounting on the seating shell</p> <p><input type="checkbox"/> armrest long <input type="checkbox"/> armrest short</p>  <p>Code U2448 (self installation on the seating shell)</p> <p><input type="checkbox"/> Pair 198,- €</p>	<p>Armrest short adjustable in height and angle adjustable with raster, for mounting on the seating shell</p>  <p>Code U2447 (self installation on the seating shell)</p> <p><input type="checkbox"/> Pair 248,- €</p>
--	---	---	--	--





footrest

<p>Single footrest with knee angle adjustable</p>  <p>Code U2420</p> <p><input type="checkbox"/> Piece 258,- €</p>	<p>Single footrest with knee angle and gear segment adjustable</p>  <p>Code U2421</p> <p><input type="checkbox"/> Piece 318,- €</p>	<p>Footrest angle adjustable additional joint fo foot angle adjustable (for single footrest with knee angle CODE U2420 or U2421)</p>  <p>Code U2428</p> <p><input type="checkbox"/> Piece OE 95,- €</p>	<p>Footrest folding upwards (for single footrest with knee angle CODE U2420 and U2421)</p>  <p>Code U2429 not in combination with CODE U2428</p> <p><input type="checkbox"/> Piece OE 158,- €</p>	<p>Spasticity lock for footrest U2429 height retractable (only in combination with CODE U2429 possible)</p> <p>Code U2429S</p> <p><input type="checkbox"/> Set OE 65,- €</p>
--	---	---	---	---

footrest


<p>Single footrest alu. flat-profile, black with gear segment knee angle</p>  <p>Code U2422</p> <p><input type="checkbox"/> Piece 326,- €</p>	<p>Heel edge for singlefootrest (black)</p>  <p>Code U2422F</p> <p><input type="checkbox"/> Piece 128,- €</p>	<p>Double footrest alu. flat-profile black with gear segment knee angle</p>  <p>Code U2423</p> <p><input type="checkbox"/> Pair 372,- €</p>	<p>Footrest coloured: for footrest alu flat profile for CODE U2422, U2423, U2425 or CODE U2426 (for seating shell)</p> <p>colour: _____</p> <p>Code U2424 (see frame colour page 1)</p> <p><input type="checkbox"/> OE 71,- €</p>	<p>Double footrest with footrest holder detachable, alu. flat-profile, black from size 1</p>  <p>Code U2425</p> <p><input type="checkbox"/> Pair 380,- €</p>	<p>Double footrest with footrest holder detachable and height adjustable, alu. flat-profile, black from size 1</p>  <p>Code U2426</p> <p><input type="checkbox"/> Pair 598,- €</p>
---	---	---	--	---	--

back

<p>Back with handbow height adjustable, angle adjustable and tiltable (with 2 single brake levers)</p> <p><input type="checkbox"/> hand bow U-type at size 2 = standard <input type="checkbox"/> hand bow trapezium type at size 1 = standard</p>  <p>back angle: Back incline stepless fixed adjustable 90° - 120° (without gas spring)</p> <p><input checked="" type="checkbox"/> Code U2440</p> <p>Code U2441</p> <p><input type="checkbox"/> Pc. OE 113,- €</p>	<p>Back angle with gas spring stepless adjustable from 90° - 120°</p> <p>Code U2441</p> <p><input type="checkbox"/> Piece 246,- €</p>	<p>Back guide for seating shells with slide bush on back tube</p>  <p>Code U2442</p> <p><input type="checkbox"/> Stück 268,- €</p>	<p>Back guide for seating shells with slide bush on back tube and angle adjustable bracket</p>  <p>Code U2443</p> <p><input type="checkbox"/> Piece 312,- €</p>	<p>Back guide PD central with guide for seating shells with a physiological centre</p>  <p>Code U2444</p> <p><input type="checkbox"/> Pc. 68,- €</p>	<p>Headrest admission for back guide</p> <p>Code U2444KH</p>
--	--	--	--	--	---

seat


adapter


<p>Seating plate seat angle adjustable with gas spring</p> <p>Seat depth adjustable from 33 to 40cm</p>  <p>Code U2430</p> <p><input type="checkbox"/> Piece 180,- €</p>	<p>Seat plate longe (+10cm) seat angle adjustable with gas spring</p> <p>Seat depth +10cm from 43 to 50cm adjustable</p>  <p>Code U2435</p> <p><input type="checkbox"/> Set 286,- €</p>	<p>Trapezium adapter 2-pieces</p>  <p>Code U2431</p> <p><input type="checkbox"/> Piece 118,- €</p>	<p>Trapezium adapter upper part (for seating shell)</p>  <p>Code U2432</p> <p><input type="checkbox"/> Piece 169,- €</p>	<p>Trapezium adapter lower part alu. (for mobility base)</p>  <p>Code U2433</p> <p><input type="checkbox"/> Set 374,- €</p>	<p>HC-Trapezium adapter for seating shell, 2-Pieces alu. upper and lower part made of aluminium</p>  <p>Code U2439</p>
--	---	--	--	--	--

Sum

rear wheel			
Rear wheel 12" with drum brake, screw axle and pneumatic tyre (for companion) 7-spokes-rim black (grey tyred) wheel camber 0° <input type="checkbox"/> 12 inch Code U2460 without hand rim / with screw axle	Rear wheel 12" with drum brake, screw axle and pneumatic tyre with black tyre (for companion) 7 spokes rim black <input type="checkbox"/> 12 inch black cover Code U2461 without handrim / with screw axle	Rear wheel 16" with drum brake, with quick release axle and pneumatic tyre with black cover 16"x3/4 and black hub and rim (Black Edition) (for companion) <input type="checkbox"/> 16 inch black cover Code U2466	
Standard	<input type="checkbox"/> Pair OE 52,- €	<input type="checkbox"/> Pair OE 252,- €	

castor			
Castor with Design aluminium fork and screw axle <input type="checkbox"/> 140mm black soft (solid) Code U2481S <input type="checkbox"/> 140mm soft (solid) Code U2481 <input type="checkbox"/> 6" solid Code U2484 <input type="checkbox"/> 6" pneumatic Code U2483 <input type="checkbox"/> 7" solid small Code U2485	Castor with Design aluminium fork and screw axle 7" solid rubber small black tyre black plastic rim <input type="checkbox"/> Code U2486	Soft castor with Design aluminium fork and screw axle 6" soft castor black soft tyre black plastic rim <input type="checkbox"/> Code U2488	Design-Alu-fork coloured powder coated colour: _____ Code U2498 (see colours from page 1)
Standard	<input type="checkbox"/> Pair OE 56,- €	<input type="checkbox"/> Pair OE 65,- €	<input type="checkbox"/> Pair OE 88,- €

accessories		Restraint system for tranport in car		
Accessory box (box made of aluminium, black powder coated)  Code U2490	Fixation bow for accessory box  Code U2491	Restraint system for hook for transport with a mobility base for seating system in a car (according to ISO 7176-19) (without pelvic restraint and without buckle tongue) 4 piece transport brackets (2 in front and 2 back) Code U4499  additional you need an extended pelvic strap and matching retractors from specialist shop.	Pelvic restraint for the restraint system (for CODE U4499) Pelvic restraint for connecting the person with the restraint system Code U4498	Anchor plate system anchor plates for transport in cars (ISO 7176-19) Assembly kit for mobility base including mounting and integrated pelvic restraint  Code U4497
<input type="checkbox"/> Piece 587,- €	<input type="checkbox"/> Piece 88,- €	<input type="checkbox"/> Set net price 245,- €	<input type="checkbox"/> Piece net 78,- €	<input type="checkbox"/> Set net 440,- €

brake	accessories	
Fixing brake for rear wheel 12" or 16"  Code U2493 Note: Combination with Code U2495 not possible II <input type="checkbox"/> Pair OE 150,- €	Inlay made of metal, coloured powder coated (plugged in frame) colour: _____ (see colours on page 1) Code U2495 Combination with fixing brake Code U2493 not possible <input type="checkbox"/> Set 298,- €	

Cruiser mobility base for seating systems		basic price	€ 1860,-
For all orders apply the general terms and conditions.		Sum optinal extras	€
		Shipping and handling	€
		Sum	€
		Tax	€
		Total price	€
Purchaser _____	Telephone _____		

Comment:

Technical changes and print mistakes are reserved.
 VD2419E-14

Cruiser

accessories

Trike

traking extension
as steering and pushing aid



Cruiser

traking extension **TRIKE**
as steering and pushing aid

special holder lateral at frame
with hook closure
castor 10x2" black alu fork
arrangement with gear segment
with ramp as mounting aid

hook closure



capacity up to **60kg**
(wheelchair incl. payload)

capacity up to **90kg**
(wheelchair incl. payload)

CODE U2496

CODE U2497

Piece

689,- €

Piece

772,- €

TRIKE colour frame Edition "...-Color"

<p>Note: the TRIKE frame parts are made of aluminium</p>	<p>TRIKE- frame colour powder coated</p> <p>Colour : frame colour (see page 1)</p> <ul style="list-style-type: none"> - frame type right / left - cross tube - castor fork <p>Code 9908</p>	<p>Edition "Duo-Color"</p> <p>colour 1: _____ (castor fork)</p> <p>colour 2: _____ (frame type right / left, cross tube)</p> <p>Code 9909</p>	<p>Castor fork with brake for directional stability</p> <p>Code 9909F</p>
	<p>without extra charge <input type="checkbox"/></p> <p>OE 115,- €</p>	<p><input type="checkbox"/></p> <p>OE 125,- €</p>	<p><input type="checkbox"/> Piece</p> <p>OE 144,- €</p>

accessories

Note:

- Cruiser: - combination TRIKE with anchor plate system is possible !!
- combination TRIKE with Inlay (Code 2495) is possible !!

General Terms and Conditions of Business (GTC)

Validity

All sales are subject to our delivery and payment conditions, without any exception. Orders and other agreements will not be valid unless confirmed by us in writing. We expressly object to the customer's purchase conditions.

Quotations and contract conclusion

Quotations are always subject to change, unless they are expressly limited in terms of time. Orders require a written confirmation. After confirmation, prices quoted in orders are fixed but we reserve the right to make adequate price adjustments in the event of price increases for which we are not responsible (e.g. tax increases, wage increases, price changes on the part of our suppliers). Cost estimates for customized products and repairs are always non-binding. We reserve the right to make technical modifications. The contract will not be concluded before the order is confirmed.

Shipment

To the extent the customer has not made any specification, shipment will be made in usual packaging using an adequate method. Insurance will only be taken out upon request and at the expense of the customer.

Delivery will be made carriage forward. We are free to select a method of shipment. If the customer has special requests, e.g. express delivery, express freight, etc., the additional cost will also be borne by the recipient. The goods will always be shipped at the buyer's risk. In case of damage during transport, the recipient must apply for an ascertainment of the facts with the respective forwarding agent.

Prices

All prices quoted by BEROLLKA will be quoted ex works, excluding any applicable statutory value added tax, public charges and levies such as potential customs tariffs, and insurance and packaging costs that will be charged separately.

Minimum quantity surcharge

If a small order is placed with a net value below 20 euros, an additional lump-sum minimum quantity surcharge of 15.00 euros will additionally be charged to cover administrative costs.

Delivery periods

The delivery dates and periods stated are subject to change. Any delay in delivery will not entitle the customer to withdraw from the contract. The assertion of damage claims is ruled out. BEROLLKA can make partial deliveries for justified reasons.

Wrong orders / return of merchandise

Any goods erroneously ordered will only be taken back provided we are responsible for such an error and the merchandise is in impeccable condition. In the event that an item erroneously ordered or merchandise that was returned is taken back in an individual case as a gesture of goodwill, attention must be paid to the following:

1. It is necessary to state the article number, the designation of the article, the date and number of the delivery note and the invoice number.
2. The return shipment must be carriage paid.
3. The merchandise must correspond to the article number(s) listed in the delivery note. In no event will any individual parts removed from an article be taken back.
4. The returned merchandise must be in impeccable condition.
5. Processing fees amounting to 10% of the net merchandise value not, however, less than €35.- will be charged for any instance where goods are taken back.

Deliveries effected more than 3 months earlier as well as sewing pieces and customized products cannot be exchanged. A copy of the delivery note or of the invoice will serve as proof in this regard.

Warranty and notices of defects

Excluding any further claims, the seller is liable for defects of the merchandise in that the seller will mend or replace all those parts free of charge that have become defective or unusable within 6 months due to provable defects of material or workmanship, and of which the seller has been notified in writing immediately after knowledge was gained of the defect. Should longer periods of time be required for statutory reasons, the warranty period is agreed to be 12 months.

The customer must grant the seller the required time and opportunity for effecting delivery of the replacement parts. Replaced parts will become the property of the seller. The buyer must check the merchandise without undue delay after it has been delivered by the seller and, if a defect is found, notify the seller of the defect without undue delay. Should the buyer fail to notify the seller, the merchandise will be deemed accepted unless a defect occurs that could not have been found when the merchandise was checked.

In any event in which a complaint is brought to the notice of BEROLLKA, BEROLLKA invariably has the right to inspect and check the item delivered that is the subject matter of said complaint, for which the customer will grant BEROLLKA the required time and opportunity. BEROLLKA can also demand that the customer sends the item delivered that is the subject matter of a complaint back to BEROLLKA at the expense of BEROLLKA. Should a notice of defects given by the customer turn out to be wilfully or gross negligently unfounded, the customer will be obliged to compensate BEROLLKA for any and all expenses incurred in this connection. BEROLLKA will be entitled to remedy warranty defects at its own discretion either by rectification of the defect or by replacement of the defective part or of the entire item delivered at no cost for the customer. Should the customer itself repair an item that gave rise to a complaint, or should the customer endeavour to rectify the defects itself without prior written approval of the seller, the seller will not grant any claims for reimbursement of expenses (such as transport, road, labour or material costs) or warranty claims.

The following are excluded from the liability for damages:

damage arising from natural wear and tear, damage due to negligent and/or improper treatment, or violent damage. A revocation of the purchase contract, a reduction of the purchase price or damage claims of any kind are ruled out. The warranty obligation will cease if repairs and/or modifications or replacements of individual parts are performed by other parties than the seller or the delivering company.

Reservation of title

With regard to resellers: The merchandise delivered will remain our property until any and all invoices – also for different deliveries – have been fully paid. Our merchandise that is hence still our property may only be resold within the ordinary course of business. Amounts receivable from any resale are assigned to us. We must be notified without undue delay and on the reseller's own initiative of any disposal of said merchandise outside of the ordinary course of business, e.g. through pledging our property or amounts receivable assigned to us.

Terms of payment

Any invoices are payable within 30 days from the date of the invoice strictly net cash. A 2% discount will be granted for payments received within 8 days provided any and all payment claims under earlier deliveries have been settled. We are entitled to charge default interest if the 30-day payment deadline is exceeded. The interest rate will be the customary market rate charged by major banks. We expressly reserve the right to effect delivery with cash on delivery.

Data Protection

BEROLLKA and the customer are responsible themselves for compliance with the respective data protection rules applicable to them. The customer must not provide BEROLLKA with personal data that is illegal or that the customer is not legally authorised to distribute.

The place of performance and the place of jurisdiction will be the place where the invoice is issued.

Heidelberg, Germany, will be the exclusive place of jurisdiction for any and all disputes arising under the contractual relationship. However, BEROLLKA is entitled to sue the customer at any other legal venue. These General Terms and Conditions of Business are governed by the laws of the Federal Republic of Germany, excluding of the United Nations Convention on Contracts for the International Sale of Goods of Goods (CISG).

Technical modifications and printing errors excepted. A current list of prices can be ordered at any time.