



货物运输条件鉴定书



Certification

for Safe Transport of Chemical Goods

危险品

样品名称 : D11电动轮椅车 (D11锂电池24V 6Ah 144Wh)

Sample name: /

委托单位 : 


生产单位 : 安徽金百合医疗器械有限公司



SICIT 上海化工院检测有限公司

Shanghai Institute of Chemical Industry Testing Co., Ltd



货物运输条件鉴定书

NO. 212200104644

Page

Certification for Safe Transport of Chemical Goods

样品名称 Sample Name	中文 Chinese	D11电动轮椅车 (D11锂电池24V 6Ah 144Wh)				
	英文 English	/				
委托单位 Consignor		[REDACTED]				
生产单位 Manufacturer		/				
检验方法、程序 Inspection method and procedure		国际航空运输协会《危险品规则》63版 IATA Dangerous Goods Regulations (DGR) 63rd Edition				
样品外观 Sample appearance		多种颜色塑胶及金属外壳 multicolor plastics cement and metal shell				
包装件信息 Package information		/				
序号 NO.	电池种类 Battery type	型号 Model	容量Capacity /锂含量Li content	放置方式 Placement	单颗重量kg Unit weight	数量 Quantity
1	可充电锂离子电池组 Rechargeable Li-ion battery	D11	6Ah 144Wh	安装在车辆内 Contained in vehicle	1.44	2
鉴定结论 IDENTIFICATION CONCLUSION	1. 危险性识别 (Hazards identification) 杂项。 Miscellaneous.					
	2. 空运按照国际航空运输协会《危险品规则》办理的类项 (Suggestion according to IATA DGR) Shipping name: Battery-powered vehicle Class or division: 9 UN Number: UN3171					
备注 Comment	3. 包装要求 (Packaging requirements) 按UN3171要求办理。 The goods are packaged according to UN3171.					
	检验日期: Inspection Date: 2021-12-06		签发日期: Issue Date: 2021-12-06		生效日期: Effective Date: 2022-01-01	



批准
Approver: 王军

审核
Checker: 董学胜

主检
Appraiser: 钱玉婷




货物运输条件鉴定书

Certification for Safe Transport of Chemical Goods

NO. 21220010464

Page 2

序号 No.	检验结果及其他事项 Inspection results and other things
1	<p>本报告所述锂电池按照《危险品规则》(63版)[以下简称DGR] 3.9.2.6.1(e)规定的质量管理体系进行制造。 本报告所述锂电池不属于因安全原因召回或损毁的锂电池。 Lithium cells and batteries listed in this report were manufactured under the quality management program described in DGR 63rd 3.9.2.6.1(e). Lithium cells and batteries listed in this report are not the defective cells or batteries returned to the manufacturer for safety reasons.</p>
2	<p>本报告所述锂电池已通过《联合国试验和标准手册》第III部分38.3小节相应测试要求。 Lithium cells and batteries listed in this report are of the types proved to meet the requirements of each applicable test in the UN Manual of Tests and Criteria, Part III, sub-section 38.3. UN38.3试验概要编号 The UN38.3 Test Summary No. (s) 812100700310628 详细信息请扫描右侧二维码。 Please scan the QR code on the right for more information.</p> 
3	<p>电池应安装并固定在车辆的电池仓内，并有适当的防止损坏及短路的保护措施。 The batteries should be installed and securely fastened in the battery holder of the vehicle and be protected in a suitable manner to prevent damage and short circuits.</p>
4	<p>如果车辆可能非直立放置，车辆应固定在坚固的刚性外包装内。车辆应牢固固定在外包装内，以防止任何运输过程中可能改变方向或造成车辆损坏的移动。 Where vehicles could possibly be handled in other than an upright position, the vehicle must be secured in a strong, rigid outer packaging. The vehicle must be secured by means capable of restraining the vehicle in the outer packaging to prevent any movement during transport which would change the orientation or cause the vehicle to be damaged.</p>
5	<p>车辆必须固定(或包装)以防意外启动。 车辆操作或安全所需的危险货物(如灭火器、轮胎充气罐或安全装置)应牢固地安装在车内。 装有防盗装置、无线电通信设备或导航系统的车辆应禁用此类装置、设备或系统。 Vehicles must be packed to prevent accidental activation. Dangerous goods required for the operation or safety of the vehicle, such as fire extinguishers, tyre inflation canisters or safety devices, must be securely mounted in the vehicle. Vehicles equipped with theft-protection devices, installed radio communications equipment or navigational system must have such devices, equipment or system disabled.</p>
6	<p>每个包装件必须贴有9类危险性标签。 Each package must be labelled with the Class 9 hazard label.</p>
7	/



-验证码:860782-

货物运输条件鉴定书

Certification for Safe Transport of Chemical Goods

NO. 21220010464

Page



电池型号: D11

FC 3055 102 55AH
D11 11 4 1000 1100 1000 1000 1000 1000 1000 1000
+3.7V 25.1x200x100

Standard Voltage: 24V
Standard Capacity: 5Ah
Power: 144Wh
Limited Charging Voltage: 29.1V
Limited Charging Current: 2A
Standard: GB31241-2014



报告结束