



Product Code	45688
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2194611
CND	V0899
EAN/UPC	8023279456882

Description

EMERGENCY TROLLEY - small
Supplied with rotating defibrillator/monitor tray, height adjustable I.V. pole with 2 hooks.

- antimicrobial epoxy painted body with seal lock
- top surface with 4-sided rim
- top surface with built-in compartment
- ball-bearing full-extension sliding drawer with self-closing system
- built-in handle
- 4 swivelling castors Ø 125 mm - 2 antistatic and with brake
- size: small 760x714xh 1,022 mm

Made in Italy.