



Product Code	45718
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2476078
CND	V0880
EAN/UPC	8023279457186

Description

NEO PLUS ANESTHESIA TROLLEY

Anesthesia trolley made of ABS and antibacterial materials externally, while the internal structure is composed of steel and aluminium. It is equipped with five big and deep drawers with dividers and four anti-static castors Ø 10 cm (two with brake).
Size: 75x47.5x95 cm
Size with accessories: 90x52x138 cm
Weight with accessories: 44.5 Kg

Drawers inner size:
2 drawers 42.4x37.5x6.8 cm
2 drawers 42.4x37.5x11 cm
1 drawer 42.4x37.5x22 cm

Table board covered by transparent glass
Hand sanitizer holder
Handle on both sides
Pull-out worktop
Foot pedal controls lids opening
Removable container

Standard accessories

Central lock
2 ABS waste bins with single lid
Pull-out worktop
10 anesthesia boxes
Dividers
Hand sanitizer holder
Utility container
Sharp container 2l