

NEO PLUS EMERGENCY TROLLEY



Product Code	45720
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2476077
CND	V0880
EAN/UPC	8023279457209

Description

NEO PLUS EMERGENCY TROLLEY

Emergency trolley made of ABS and antibacterial materials externally, while the internal structure is composed of steel and aluminium. It is equipped with five big and deep drawers with dividers and 4 castors Ø 10 cm (two with brake).

Size: 75x47.5x92 cm

Weight with accessories: 41 Kg

Drawers inner size:

2 drawers 42.4x37.5x6.8 cm

2 drawers 42.4x37.5x11 cm

1 drawer 42.4x37.5x22 cm

Movable IV pole and defibrillator plate

Table board covered by transparent glass

Handle on both sides

Standard accessories

IV pole
Oxygen bottle holder
Defibrillator plate
Utility container
Pull-out worktop
CPR board
Dividers
Hand sanitizer holder
Power outlet
2 sharp containers 2l
Central lock
2 ABS waste bins with lid