



Product Code	45722
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2476060
CND	V0880
EAN/UPC	8023279457223

Description

NEO EMERGENCY TROLLEY

The double locking system (key + seal) is available starting from lot 6 onward

Heavy duty construction emergency trolley size with five drawers.
Four castors (Ø 10 cm) assembled on the vertical column line for better stability.
One piece counter top with a fully guarded handle for comfort and protection.

Size: 75x47.5x92 cm
Weight with accessories: 39 Kg

Multilingual instructions: GB, FR, IT, ES, PT, DE, NL, PL, RO, CZ, HU, SE, FI, LT, HR, GR, BG.

Standard accessories

- Defibrillator plate
- IV pole
- Sliding side shelf
- CPR board
- Sharp container
- 2 waste bins
- Oxygen bottle holder
- Dividers
- Power outlet
- Central lock
- Utility container