



DELPHIN


POOL LIFT

PRACTICAL HELP

For independent bathing



Poollift



Elegant and robust **design** is suitable for public and private pools



Optional: transport wheels and/or folding armrests




DELPHIN PRODUCT HIGHLIGHTS

- Independent operation without an assistant
- Operated with water pressure
- No electrical components
- Choose from two seat colors (blue and white)
- Attractive design in 316 stainless steel
- For indoor and outdoor use
- Convenient operation with control lever




FAST INSTALLATION


No underwater installation is necessary, an attachment point on the pool edge is sufficient. The lift is operated with household water pressure.



Variant A: Recessed console



Variant B: Screwed bracket



Specifications

Load capacity	130 kg
Lifting height standard	1220 mm
Lifting height optional	1067 to 1600 mm
Water pressure min.	4 bar
Water depth min.	1150 mm
Material	Stainless steel AISI 316

Scan here to
view the product



Lehner Lifttechnik GmbH
Mitterweg 4
4722 Peuerbach
Austria
sales@lehnerlifttechnik.at
www.lehner-lifttechnik.com

09/2025