

Adress

Shipping address

☐ see address☐ Order☐ Estimate

Purchaser: _____

Telephone: _____

Com.: _____

Basic Price 2900,- €Price valid from January 1st 2026
(until revoked)

ROOKIE

Foldable mobility base for seating systems

Capacity up to 120kg

user weight + seating system: _____ kg



Standard equipment:

- foldable aluminium frame with abduction
- back with high adjustable handles
- seat angle adjustable
- back angle adjustable, tiltable
- castor and rear wheels

basic measure

ROOKIE foldable frame  abduction frame (10 cm)  Code U301	Seat width (SW) ☐ 28 ☐ 32 ☐ 36 ☐ 40 ☐ 44	Seat depth (SD) Length of the seat frame Gr. 1 frame short ☐ ST 32 Standard Gr. 2 frame long ☐ ST 40 add charge (see U334 page 3) (seat extension 6cm backwards see page 3)	Back height (BH)  45	Lower Leg Length (LLL): LLL: cm
---	---	---	--	---

technical data

T seat depth: 32cm with short frame T seat depth: 40cm (optional) with long frame 		Seat tilt -6 to 35° at straight / flat adjusted frame
A1 = total width front A1 = SB + 18cm A2 = frame width front A2 = SB + 10cm B = seat frame width B = SB C = seat width (= SW) E = max. seat shell width E = C (SW) + 5cm by wheel cover D = total width D = SW + 24cm (without drum brake) D = SW + 28cm (with drum brake)		Incline back 80 to 120°
folding measure D = total width in folding position D = 34cm (with / without drum brake)		Weight: from 10,5 kg without footrest without wheel cover without drum brake All dimensions in cm Dimensions tol. +/- 0,5cm

Technical changes and print mistakes are reserved.

VD240E-26

measurements

seat height H:

Assumption for dimensions: the seat plate is parallel to the frame !

seat height back

(rear wheels) to the seat plate

20"

22"

24"

wheel up	41,0 cm	43,5 cm	46,0 cm
wheel centred	43,5 cm	46,0 cm	48,5 cm
wheel down	46,0 cm	48,5 cm	51,0 cm

seat heights with trapezium adapter +2,0cm

seat height front

(castor / fork)
to the seat platefork
extra short
2564213fork
short
4064201fork
long
4064202fork
45 wide
4064811fork
60 wide
4064811

4"

wheel up	39,2 cm	42,2 cm	44,8 cm
wheel centred	40,5 cm	43,5 cm	46,0 cm
wheel down	41,8 cm	44,8 cm	47,2 cm

5"

wheel up	40,5 cm	43,5 cm	46,0 cm
wheel centred	41,8 cm	44,8 cm	47,2 cm
wheel down	43,0 cm	46,0 cm	49,8 cm

140mm

wheel up	45,5 cm
wheel centred	46,8 cm
wheel down	48,0 cm



6"

wheel up	44,8 cm	47,2 cm
wheel centred	46,0 cm	48,5 cm
wheel down	47,2 cm	51,0 cm

7"

wheel up	47,2 cm
wheel centred	48,5 cm
wheel down	49,8 cm

seat heights with trapezium adapter +2,0cm

seat height
backseat height
front

Note: up to the position of the wheels, if the frame is not straight, you can get other angle of the seat plate !!

Technical changes and print mistakes are reserved.

VD2460E-26

frame colour

Because of the material are different colours possible

shiny

☐  25 RAL9010 white	☐  16 RAL5002 ultramarine blue	☐  19 silvermetallic
☐  24 RAL1028 melonyellow	☐  17 dormantblue	☐  46 RAL 9023 perl darkgrey
☐  27 mandarin met.	☐  64 petrolgreen met.	☐  63 platin met.
☐  06 dormantred	☐  51 applegreen metallic	☐  23 RAL9005 black
☐  40 redmetallic	☐  45 RAL 6018 yellow green	

met. = metalized

frame colour

Because of the material are different colours possible

mat / satined

☐  79 edelweiss matt	☐  70 ultramarine blue satined	☐  77 silvermetallic satined
☐  66 mandarin met. satined	☐  75 petrolgreen met. satined	☐  78 platin met. satined
☐  67 traffic red satined	☐  73 yellowgreen met. satined	☐  58 satinblack

met. = metalized

seat height

seat height in cm
the seat plate is parallel to the frame!

front: _____ cm





back: _____ cm

note



note

frame accessories

armrest

Tip shackle ☐ right ☐ links  Code U311 Piece 94,- €	Anti-tip wheel ☐ right ☐ left  Code U312 Piece 152,- €	Armrest adjustable in height with raster, for mounting on the seating shell ☐ armrest long ☐ armrest short  Code U319 (self installation on the seating shell) Pair 234,- €	Armrest short adjustable in height and angle adjustable with raster, for mounting on the seating shell ☐ armrest long ☐ armrest short  Code U318 (self installation on the seating shell) Pair 272,- €
--	---	---	--

footrest

Single footrest alu-profile, black with knee angle and gear segment adjustable  Code U322 ☐ Piece 356,- €	Single footrest alu. flat-profile with gear segment knee angle  Code U325 ☐ Piece 356,- €	Double footrest alu. flat-profile with gear segment knee angle  Code U326 ☐ Pair 432,- €	Heel edge for single footrest Code U328 ☐ Pair 148,- €	Heel edge for double footrest Code U327 ☐ Pair 159,- €
--	--	---	---	---



back

Back with single push handles, high adjustable, back angle adjustable with stabilization bar back angle displaceable -10/0/10/20/30° ☐ ☒ Code U340	Back guide for seating shells Code U346 ☐ Piece 292,- €	Single push handles angle adjustable (with gear segment) Code U341 ☐ Piece OE 242,- €	Head rest with attachment at the stabilization bar (with black cover) Code U342 ☐ Piece 325,- €	Attachment for head rest (without the head rest) adapter tube 20x20 Code U343 ☐ Piece OE 102,- €	Additional support for the head rest attachment (only in combination with CODE 342 possible) Code U344 ☐ Pc. OE 198,- €
Standard					

seat

Seat standard (short) seat tubes with seat plate (angle adjustable) size 1 SD = 32cm seat frame short Code U330	Seat frame long, seat tube with seat plate (angle adjustable) size 2 SD = 40cm seat frame long (and basic frame long) Code U334	Seat extension 6cm backwards tube of seat is 6cm longer backwards (instead of normal seat tube) size 1 SD = 32 + 6cm = 38cm size 2 SD = 40 + 6cm = 46cm Code U331	Standard ☐ Set OE 219,- €	☐ Set OE 192,- €
--	---	---	--	---

adapter

Trapezium adapter 2-pieces (for seating shell)  Code U335 ☐ Set 324,- €	Trapezium adapter upper part (for seating shell)  Code U336 ☐ Piece 137,- €	Trapezium adapter lower part (for mobility base)  Code U337 ☐ Piece 191,- €	Edition "DUO-Color" ROOKIE colour 1: _____ - frame part right - frame part left - back tube colour 2: _____ - footrest - design-caster fork - hand rim* *only in combination with CODE U359 ☐ Code U305 OE 106,- €
--	--	--	---

Rear wheel with quick release axle and pneumatic tyre ☐ 20" x1 3/8 Code U360 ☐ 22" x1 3/8 Code U361 ☐ 22" x1 Code U362 ☐ 24" x1 3/8 Code U363 ☐ 24" x1 Code U364	Rear wheel with quick release axle and solid tyre ☐ 20" x1 3/8 Code U365 ☐ 22" x1 3/8 Code U366 ☐ 22" x1 black Code U366S ☐ 24" x1 3/8 Code U367 ☐ 24" x1 Code U368 ☐ 24"x1 black Code U368S	Rear wheel with drum brake, hand rims aluminium and quick release axle (for companion) ☐ 20" Code U351 ☐ 22" Code U352 ☐ 24" Code U353	"Black Edition" drum brake rear wheels with black hub, black wheels rim and quick release axle with black cover and with black tyres Marathon Plus Evolution 22"x1 ☐ U352S 24"x1 ☐ Code U353S (only in combination with CODE U355 / U354)	Wheel cover plastic (with holder)  Code U390	Camber 0° ☒
Standard		Pair OE 326,- €	Pair OE 432,- €	☐ Pair 264,- €	

Handrim aluminium black ☐ 20" Code U356 ☐ 22" Code U355 ☐ 24" Code U354 extension narrow ☐ wide ☐	Hand rim stainless steel, polished ☐ 22" Code U358 ☐ 24" Code U357 extension narrow ☐ wide ☐	handrim for rear wheel powder coated  Code 8881 only with hand rim alu.	Spoke cover transparent  ☐ 20" Code U370 ☐ 22" Code U371 ☐ 24" Code U372	Handrim cover black  ☐ 20" Code U394 ☐ 22" Code U395 ☐ 24" Code U396	Wheel extension Code U376	
Standard		☐ Pair OE 194,- €	Pair 180,- €	Pair 138,- €	☐ Pair 190,- €	

Castor solid with screw axle and design fork in colour of the wheelchair ☐ 4" Code U381 ☐ 5" Code U382 ☐ 6" Code U383 ☐ 6" PU soft black Code U383P ☐ 140mm soft Code U384 ☐ 140mm soft black Code U384S ☐ 7" Code U385 *	Castor pneumatic with screw axle and design fork in colour of the wheelchair ☐ 6" Code U386 ☐ 7" Code U387 * * Castor fork only in black possible	Soft alu castor black soft tyre with black aluminium rim ☐ 4" Code U3894S ☐ 5" Code U3895S	Castor fork with lock Lock for driving forwards (only for 4"/5"/6") Code U380	
Standard		Pair OE 102,- €	☐ Pair OE 220,- €	

Knee agnle brake with short lever stroke  Code U374	Brake lever extension straight, foldaway  Code U375		Anchor plate system anchor plates for transport in cars (ISO 7176-19) Assembly kit for mobility bases including mounting  Code U3091
Standard	☐ Pair 84,- €		☐ Set 895,- €

[illegible]

note

VD2460E-26
Technical changes and print mistakes are reserved.

General Terms and Conditions of Business (GTC)

Validity

All sales are subject to our delivery and payment conditions, without any exception. Orders and other agreements will not be valid unless confirmed by us in writing. We expressly object to the customer's purchase conditions.

Quotations and contract conclusion

Quotations are always subject to change, unless they are expressly limited in terms of time. Orders require a written confirmation. After confirmation, prices quoted in orders are fixed but we reserve the right to make adequate price adjustments in the event of price increases for which we are not responsible (e.g. tax increases, wage increases, price changes on the part of our suppliers). Cost estimates for customized products and repairs are always non-binding. We reserve the right to make technical modifications. The contract will not be concluded before the order is confirmed.

Shipment

To the extent the customer has not made any specification, shipment will be made in usual packaging using an adequate method. Insurance will only be taken out upon request and at the expense of the customer.

Delivery will be made carriage forward. We are free to select a method of shipment. If the customer has special requests, e.g. express delivery, express freight, etc., the additional cost will also be borne by the recipient. The goods will always be shipped at the buyer's risk. In case of damage during transport, the recipient must apply for an ascertainment of the facts with the respective forwarding agent.

Prices

All prices quoted by BEROLLKA will be quoted ex works, excluding any applicable statutory value added tax, public charges and levies such as potential customs tariffs, and insurance and packaging costs that will be charged separately.

Minimum quantity surcharge

If a small order is placed with a net value below 20 euros, an additional lump-sum minimum quantity surcharge of 15.00 euros will additionally be charged to cover administrative costs.

Delivery periods

The delivery dates and periods stated are subject to change. Any delay in delivery will not entitle the customer to withdraw from the contract. The assertion of damage claims is ruled out. BEROLLKA can make partial deliveries for justified reasons.

Wrong orders / return of merchandise

Any goods erroneously ordered will only be taken back provided we are responsible for such an error and the merchandise is in impeccable condition. In the event that an item erroneously ordered or merchandise that was returned is taken back in an individual case as a gesture of goodwill, attention must be paid to the following:

1. It is necessary to state the article number, the designation of the article, the date and number of the delivery note and the invoice number.
2. The return shipment must be carriage paid.
3. The merchandise must correspond to the article number(s) listed in the delivery note. In no event will any individual parts removed from an article be taken back.
4. The returned merchandise must be in impeccable condition.
5. Processing fees amounting to 10% of the net merchandise value not, however, less than €35.- will be charged for any instance where goods are taken back.

Deliveries effected more than 3 months earlier as well as sewing pieces and customized products cannot be exchanged. A copy of the delivery note or of the invoice will serve as proof in this regard.

Warranty and notices of defects

Excluding any further claims, the seller is liable for defects of the merchandise in that the seller will mend or replace all those parts free of charge that have become defective or unusable within 6 months due to provable defects of material or workmanship, and of which the seller has been notified in writing immediately after knowledge was gained of the defect. Should longer periods of time be required for statutory reasons, the warranty period is agreed to be 12 months.

The customer must grant the seller the required time and opportunity for effecting delivery of the replacement parts. Replaced parts will become the property of the seller. The buyer must check the merchandise without undue delay after it has been delivered by the seller and, if a defect is found, notify the seller of the defect without undue delay. Should the buyer fail to notify the seller, the merchandise will be deemed accepted unless a defect occurs that could not have been found when the merchandise was checked.

In any event in which a complaint is brought to the notice of BEROLLKA, BEROLLKA invariably has the right to inspect and check the item delivered that is the subject matter of said complaint, for which the customer will grant BEROLLKA the required time and opportunity. BEROLLKA can also demand that the customer sends the item delivered that is the subject matter of a complaint back to BEROLLKA at the expense of BEROLLKA. Should a notice of defects given by the customer turn out to be wilfully or gross negligently unfounded, the customer will be obliged to compensate BEROLLKA for any and all expenses incurred in this connection. BEROLLKA will be entitled to remedy warranty defects at its own discretion either by rectification of the defect or by replacement of the defective part or of the entire item delivered at no cost for the customer. Should the customer itself repair an item that gave rise to a complaint, or should the customer endeavour to rectify the defects itself without prior written approval of the seller, the seller will not grant any claims for reimbursement of expenses (such as transport, road, labour or material costs) or warranty claims.

The following are excluded from the liability for damages:

damage arising from natural wear and tear, damage due to negligent and/or improper treatment, or violent damage. A revocation of the purchase contract, a reduction of the purchase price or damage claims of any kind are ruled out. The warranty obligation will cease if repairs and/or modifications or replacements of individual parts are performed by other parties than the seller or the delivering company.

Reservation of title

With regard to resellers: The merchandise delivered will remain our property until any and all invoices – also for different deliveries – have been fully paid. Our merchandise that is hence still our property may only be resold within the ordinary course of business. Amounts receivable from any resale are assigned to us. We must be notified without undue delay and on the reseller's own initiative of any disposal of said merchandise outside of the ordinary course of business, e.g. through pledging our property or amounts receivable assigned to us.

Terms of payment

Any invoices are payable within 30 days from the date of the invoice strictly net cash. A 2% discount will be granted for payments received within 8 days provided any and all payment claims under earlier deliveries have been settled. We are entitled to charge default interest if the 30-day payment deadline is exceeded. The interest rate will be the customary market rate charged by major banks. We expressly reserve the right to effect delivery with cash on delivery.

Data Protection

BEROLLKA and the customer are responsible themselves for compliance with the respective data protection rules applicable to them. The customer must not provide BEROLLKA with personal data that is illegal or that the customer is not legally authorised to distribute.

The place of performance and the place of jurisdiction will be the place where the invoice is issued.

Heidelberg, Germany, will be the exclusive place of jurisdiction for any and all disputes arising under the contractual relationship. However, BEROLLKA is entitled to sue the customer at any other legal venue. These General Terms and Conditions of Business are governed by the laws of the Federal Republic of Germany, excluding of the United Nations Convention on Contracts for the International Sale of Goods (CISG).

Technical modifications and printing errors excepted. A current list of prices can be ordered at any time.