

Adress

Shipping address

☐ see address

- ☐ Order
☐ Estimate

Purchaser: _____

Telephone: _____

 **MADE IN GERMANY**

 **CRASH TESTED** ISO 7178-10

Com.: _____

Basic Price 2790,- €

Price valid from January 1st 2026
(until revoked)

Tourer

mobility base for seating systems

Standard equipment:

- rigid aluminium with **abduction**
- handbow angle adjustable, height adjustable and tiltable
- castor
- rear wheel with drum brake

capacity up to 120kg

user weight + seating system: _____ kg

max. width
seat hell

basic measure						
frame	frame tube	seat width B	frame width front A	frame width back C	width between wheel D	total width E
frame size 0 ☐ Code U4400	25mm	26 cm	45 cm	32 cm	37,5 cm 40,5cm*	55,5 cm 47cm*
frame size 1 ☐ Code U4401	25mm	32 cm	51 cm	38 cm	43,5 cm 46,5cm*	61,5 cm 53cm*
* frame size 2 ☐ Code U4402	25mm	38 cm	57 cm	44 cm	49,5 cm 52,5cm*	67,5 cm 59cm*
* frame size 3 ☐ Code U4403	25mm	44 cm	63 cm	50 cm	55,5 cm 58,5cm*	73,5 cm 65cm*

* Note: for size 2 and 3 suggest the frame with 2 gas springs for seat angle adjustment (see CODE U4408) !!!

* without wheel cover and without camber

* width with 12" wheel

technical data

Seat tilt

- ☐ -8 to 30°
☐ 0 to 35°
☐ 10 to 45°

min. seat height

45 cm

Back incline
with gas spring

90 to 120°


concerning seat frame
(trapezium adapter +2cm)

Seat about 8cm
height adjustable


Weight:

from 11 kg (without rear wheels)
from 15 kg (without rear wheels 24")

H total height: 95 - 120cm*
L total length: 80 - 90cm
T seat depth: 33 - 40cm



* according to the angle of the back/ hand bow and the size of the rear wheels



frame colour

Because of the material are different colours possible

shiny

☐ 	25 RAL9010 white	☐ 	53 cassis shiny	☐ 	45 RAL 6018 yellow green	☐ 	63 platin met.
☐ 	24 RAL1028 melonyellow	☐ 	59 atlanticblue met.	☐ 	60 mojitogreen met.	☐ 	23 RAL9005 black
☐ 	27 mandarin met.	☐ 	16 RAL5002 ultramarine blue	☐ 	19 silvermetallic		
☐ 	06 dormantred	☐ 	17 dormantblue	☐ 	46 RAL 9023 perl darkgrey		
☐ 	40 redmetallic	☐ 	64 petrolgreen met.				

met. = metalized

frame colour

Because of the material are different colours possible

mat / satined

☐ 	79 edelweiss mat	☐ 	75 petrolgreen met. satined	☐ 	77 silvermetallic satined
☐ 	66 mandarin met. satined	☐ 	73 yellowgreen met. satined	☐ 	78 platin met. satined
☐ 	67 traffic red satined	☐ 	74 mojitogreen met. satined	☐ 	58 satinblack
☐ 	71 atlanticblue met. satined				
☐ 	70 ultramarine blue satined				

met. = metalized

frame accessories

Edition "Duo-Color"

Color 1: _____

Color 2: _____

- Frame- Back frame
- Top backrest bow / top hand bow
- Seat plate- Seat crossbar
- Footrest flat profile (U2422 / U2423 / U2425 / U2426) *
- Footrest holder (for U2425 / U2426)

- Aluminium wheel cover / clothes guard (U4476)
- Steering wheel fork + cover **
- Inlays / frame covers ***
- Aluminium / stainless steel handrim ****

Code U2405

- * only in combination with CODE U2424
- ** only in combination with CODE U2498
- *** only in combination with CODE U2495
- **** only in combination with CODE 0000189 or 0000189E



MP 118,- €

frame accessories

armrest

Tip shackle


- ☐ right
☐ left




Code U2411

Anti-tip wheel
with adaptermax.
2 part s
possible!

- ☐ right
☐ left



Code U2412

Armrest holder with
arm pads fixed at the
back tube


Code U2446

Armrest adjustable in height
with raster, for mounting
on the seating shell

- ☐ armrest long
☐ armrest short

Code U2448

(self installation on the seating shell)

Armrest short adjustable in
height and angle adjustable
with raster, for mounting
on the seating shell

Code U2447

(self installation on the seating shell)



Piece	94,- €	Piece	152,- €	Pair	296,- €	Pair	234,- €	Pair	286,- €
-------	--------	-------	---------	------	---------	------	---------	------	---------

note


Technical changes and print mistakes are reserved.

VD2421E-18

footrest



Single footrest with knee angle adjustable  Code U2420	Single footrest with knee angle and gear segment adjustable  Code U2421	Footrest angle adjustable additional joint for foot angle adjustable (for single footrest with knee angle CODE U2420 or U2421)  Code U2428	Footrest folding upwards (for single footrest with knee angle CODE U2420 and U2421)  Code U2429 not in combination with CODE U2428	Spasticity lock for footrest U2429 height retractable (only in combination with CODE U2429 possible) Code U2429S	
☐ Piece 296,- €	☐ Piece 348,- €	☐ Piece OE 116,- €	☐ Piece OE 182,- €	☐ Set OE 75,- €	

footrest




Single footrest alu. flat-profile, black with gear segment knee angle  Code U2422	Heel edge for single footrest (black)  Code U2422F	Double footrest alu. flat-profile black with gear segment knee angle  Code U2423	Heel edge for double footrest (black) for CODE U2423 Code U2423F	Calf plate (pad) double footrest for CODE U2423 Code U7023 ☐ Pair 172,- € single footrest for CODE U2422 Code U7022 ☐ Piece 156,- €	Coloured footrest: for flat profile footrest for CODE U2422, U2423, U2425 or CODE U2426 Farbe: _____ Code U2424 (See colours from page 2) ☐ OE 85,- €
☐ Piece 356,- €	☐ Piece 148,- €	☐ Pair 432,- €	☐ Pair 159,- €		

footrest

back





Double footrest with footrest holder detachable, alu. flat-profile, black from size 0  Code U2425	Double footrest with footrest holder detachable and height adjustable, alu. flat-profile, black from size 1  Code U2426	admission height adjustable for footrest holder  Code U2427	Back with handbow height adjustable, angle adjustable and tiltable (with 2 single brake levers) ☐ hand bow U-type at size 2 = standard ☐ hand bow trapezium type at size 1 = standard ☒ Code U2440 back angle: Back incline stepless fixed adjustable 90° - 120° (without gas spring)	Back with individual sliding handles, height adjustable  Code U2440-S
☐ Pair 412,- €	☐ Pair 612,- €	☐ Pair OE 302,- €	Standard	☐ Stück MP 298,- €

back


Back angle with gas spring stepless adjustable from 90° - 120° Code U2441 ☐ Pc. OE 132,- €	Back guide for seating shells with slide bush on back tube  Code U2442 ☐ Piece 282,- €	Back guide for seating shells with slide bush on back tube and angle adjustable bracket  Code U2443 ☐ Stück 299,- €	Back guide PD central with guide for seating shells with a physiological centre  Code U2444 ☐ Piece 342,- €	Headrest admission for back guide PD center only in combination with Code U2444 Code U2444KH ☐ Pc. 96,- €
---	---	--	--	--

seat


adapter

Seating plate seat angle adjustable with gas spring Seat depth adjustable from 33 to 40cm  Code U2430	Seat plate long (+10cm) seat angle adjustable with gas spring seat depth adjustable from 43 to 50cm Code U2435 ☐ Pc. OE 219,- €	Trapezium adapter 2-pieces (for seating shell)  Code U2431 ☐ Set 324,- €	Trapezium adapter upper part (for seating shell)  Code U2432 ☐ Piece 137,- €	Trapezium adapter lower part alu. (for mobility bases)  Code U2433 ☐ Piece 191,- €	HC-Trapezium adapter for seating shell, 2-pieces alu. upper and lower part made of aluminium  Code U2439 ☐ Set 424,- €	Seat frame with 2 gas springs for seat adjustment ## Special option## (from 90kg recommended) for size 2 and 3 Code U4408 ☐ Set OE 552,- €
Standard						

rear wheel

SENTINEL	Rear wheel with drum brake, black hub, black rim and hand rims aluminium black, quick release axle, black pneumatic tyre (for companion)	Rear wheel with drum brake, black hub, black rim and hand rims aluminium black, quick release axle, black solid tyre (for companion)	Camber	Wheel cover plastic (with holder and reinforced plate)	Wheel cover aluminium (clothing guard round in frame colour - with holder)
	☐ 12" Code U4461S without push rim / only 0° camber possible	☐ 20"x1 Code U4467S	☐ 0° standard		
	☐ 16"x1 3/4 Code U4465S without push rim	☐ 22"x1 Code U4468S	☐ 4°	☐ 20" Code U4475	for rear wheel 20", 22" or 24"
	☐ 20"x1 Code U4462S	☐ 24"x1 Code U4469S		☐ 22" / 24" Code U4479	Code U4476
	☐ 22"x1 Code U4463S				
	☐ 24"x1 Code U4464S				
Standard				Pair 264,- €	☐ Pair 369,- €

rear wheel





Spoke cover transparent	Handrim cover black	Wheelbase extension (about 10cm) for 16" / 20" / 22" / 24" only with 0° camber possible	Inlay made of metal, coloured powder-coated (plugged in frame) colour: _____ (see colour from page 2)	
				
☐ 20" Code U4470	☐ 20" Code U4494	Code U4477	Code U2495	
☐ 22" Code U4471	☐ 22" Code U4495	combination with tip shackle not possible II	combination mit brake Code U4492 is not possible by 12"	
☐ 24" Code U4472	☐ 24" Code U4496			
Pair 180,- €	Pair 138,- €	☐ Pair 219,- €	☐ Set 319,- €	

castor

Castor with Design aluminium fork and screw axle	Soft castor with Design aluminium fork and screw axle	Alu fork diagonal and screw axle	Design-Alu-fork coloured powder coated
☐ 140mm black soft (solid) Code U2481S	☐ 6" soft castor black soft tyre black plastic rim Code U2488	☐ 8" solid castor black black tyre black plastic tyre Code U4488	colour: _____ (see colour on page 2)
☐ 140mm grey soft (solid) Code U2481			
☐ 6" PU soft black tyre Code U2484P			
☐ 6" Vollgummi grey tyre Code U2484			
☐ 6" luftbereift grey tyre Code U2483			
☐ 7" Vollgummi black tyre Code U2485S			
☐ 7" Vollgummi grey tyre Code U2485			
		- only in combination with 22" or 24" - seat height around 2,5cm higher	Code U2498
Standard	☐ Pair OE 76,- €	☐ Pair OE 99,- €	☐ Pair OE 102,- €

brake

accessoires

Fixing brake for rear wheels for 20", 22" or 24"  Code U4491	for 12" and 16"  Code U4492	Brake lever extension straight, foldaway  Code U4478	Accessory box (box made of aluminium, black powder coated)  Code U2490	Fixation bow for accessory box  Code U2491
Pair OE 174,- €	Pair OE 174,- €	☐ Pair 84,- €	☐ Piece 598,- €	☐ Piece 110,- €

Restraint system for transport in car

Restraint system for hook for transport with a mobility base for seating system in a car (according to ISO 7176-19) (without pelvic restraint and without buckle tongue) 4 piece transport brackets (2 in front and 2 back)  Code U4499 additional you need an extended pelvic strap and matching retractors from specialist shop.	Pelvic restraint for the restraint system (for CODE U4499) Pelvic restraint for connecting the person with the restraint system Code U4498	Anchor plate system anchor plates for transport in cars (ISO 7176-19) Assembly kit for mobility base including mounting and integrated pelvic restraint  Code U4497
☐ Set Net price 272,- €	☐ Piece net 98,- €	☐ Set 895,- €

Tourer

Trike

trakking extension
as steering and pushing aid



Tourer

trakking extension **TRIKE**
as steering and pushing aid

hook closure



capacity up to **60kg** size 0/1
(wheelchair incl. payload)

CODE U2496

special holder lateral at frame
with hook closure
castor 10x2" black alu fork
arrangement with gear segment
with ramp as mounting aid

capacity up to **90kg** size 2/3
(wheelchair incl. payload)

CODE U2497

Note:

- combination TRIKE with anchor plate system
is possible !!

- combination TRIKE with Inlay (Code 2495)
is possible !!

☐ Piece 718,- € ☐ Piece 804,- €

TRIKE colour frame

Edition "...-Color"

Note: the TRIKE frame parts are made of aluminium	TRIKE- frame colour powder coated Colour : frame colour (see page 1) - frame type right / left - cross tube - castor fork Code 9908	Edition "Duo-Color" colour 1: _____ (castor fork) colour 2: _____ (frame type right / left, cross tube) Code 9909		Castor fork with brake for directional stability Code 9909F
without extra charge	☐ OE 136,- €	☐ OE 146,- €		☐ Piece OE 156,- €

Tourer

accessories

trakking extension **TRIKE**
as steering and pushing aid

quick release
(clamp)



capacity up to **60kg**
(base incl. payload)
for size 0 and size 1

CODE U24991

special holder for quick release (clamp)
at the frame
castor 10x2" black alu fork
arrangement of levers with gear segment
with ramp as mounting aid

capacity up to **90kg**
(base incl. payload)
for size 2 and size 3

CODE U24992

☐ Piece 728,- € ☐ Piece 814,- €

note

note

VD2421E-18
Technical changes and print mistakes are reserved.

General Terms and Conditions of Business (GTC)

Validity

All sales are subject to our delivery and payment conditions, without any exception. Orders and other agreements will not be valid unless confirmed by us in writing. We expressly object to the customer's purchase conditions.

Quotations and contract conclusion

Quotations are always subject to change, unless they are expressly limited in terms of time. Orders require a written confirmation. After confirmation, prices quoted in orders are fixed but we reserve the right to make adequate price adjustments in the event of price increases for which we are not responsible (e.g. tax increases, wage increases, price changes on the part of our suppliers). Cost estimates for customized products and repairs are always non-binding. We reserve the right to make technical modifications. The contract will not be concluded before the order is confirmed.

Shipment

To the extent the customer has not made any specification, shipment will be made in usual packaging using an adequate method. Insurance will only be taken out upon request and at the expense of the customer.

Delivery will be made carriage forward. We are free to select a method of shipment. If the customer has special requests, e.g. express delivery, express freight, etc., the additional cost will also be borne by the recipient. The goods will always be shipped at the buyer's risk. In case of damage during transport, the recipient must apply for an ascertainment of the facts with the respective forwarding agent.

Prices

All prices quoted by BEROLLKA will be quoted ex works, excluding any applicable statutory value added tax, public charges and levies such as potential customs tariffs, and insurance and packaging costs that will be charged separately.

Minimum quantity surcharge

If a small order is placed with a net value below 20 euros, an additional lump-sum minimum quantity surcharge of 15.00 euros will additionally be charged to cover administrative costs.

Delivery periods

The delivery dates and periods stated are subject to change. Any delay in delivery will not entitle the customer to withdraw from the contract. The assertion of damage claims is ruled out. BEROLLKA can make partial deliveries for justified reasons.

Wrong orders / return of merchandise

Any goods erroneously ordered will only be taken back provided we are responsible for such an error and the merchandise is in impeccable condition. In the event that an item erroneously ordered or merchandise that was returned is taken back in an individual case as a gesture of goodwill, attention must be paid to the following:

1. It is necessary to state the article number, the designation of the article, the date and number of the delivery note and the invoice number.
2. The return shipment must be carriage paid.
3. The merchandise must correspond to the article number(s) listed in the delivery note. In no event will any individual parts removed from an article be taken back.
4. The returned merchandise must be in impeccable condition.
5. Processing fees amounting to 10% of the net merchandise value not, however, less than €35.- will be charged for any instance where goods are taken back.

Deliveries effected more than 3 months earlier as well as sewing pieces and customized products cannot be exchanged. A copy of the delivery note or of the invoice will serve as proof in this regard.

Warranty and notices of defects

Excluding any further claims, the seller is liable for defects of the merchandise in that the seller will mend or replace all those parts free of charge that have become defective or unusable within 6 months due to provable defects of material or workmanship, and of which the seller has been notified in writing immediately after knowledge was gained of the defect. Should longer periods of time be required for statutory reasons, the warranty period is agreed to be 12 months.

The customer must grant the seller the required time and opportunity for effecting delivery of the replacement parts. Replaced parts will become the property of the seller. The buyer must check the merchandise without undue delay after it has been delivered by the seller and, if a defect is found, notify the seller of the defect without undue delay. Should the buyer fail to notify the seller, the merchandise will be deemed accepted unless a defect occurs that could not have been found when the merchandise was checked.

In any event in which a complaint is brought to the notice of BEROLLKA, BEROLLKA invariably has the right to inspect and check the item delivered that is the subject matter of said complaint, for which the customer will grant BEROLLKA the required time and opportunity. BEROLLKA can also demand that the customer sends the item delivered that is the subject matter of a complaint back to BEROLLKA at the expense of BEROLLKA. Should a notice of defects given by the customer turn out to be wilfully or gross negligently unfounded, the customer will be obliged to compensate BEROLLKA for any and all expenses incurred in this connection. BEROLLKA will be entitled to remedy warranty defects at its own discretion either by rectification of the defect or by replacement of the defective part or of the entire item delivered at no cost for the customer. Should the customer itself repair an item that gave rise to a complaint, or should the customer endeavour to rectify the defects itself without prior written approval of the seller, the seller will not grant any claims for reimbursement of expenses (such as transport, road, labour or material costs) or warranty claims.

The following are excluded from the liability for damages:

damage arising from natural wear and tear, damage due to negligent and/or improper treatment, or violent damage. A revocation of the purchase contract, a reduction of the purchase price or damage claims of any kind are ruled out. The warranty obligation will cease if repairs and/or modifications or replacements of individual parts are performed by other parties than the seller or the delivering company.

Reservation of title

With regard to resellers: The merchandise delivered will remain our property until any and all invoices – also for different deliveries – have been fully paid. Our merchandise that is hence still our property may only be resold within the ordinary course of business. Amounts receivable from any resale are assigned to us. We must be notified without undue delay and on the reseller's own initiative of any disposal of said merchandise outside of the ordinary course of business, e.g. through pledging our property or amounts receivable assigned to us.

Terms of payment

Any invoices are payable within 30 days from the date of the invoice strictly net cash. A 2% discount will be granted for payments received within 8 days provided any and all payment claims under earlier deliveries have been settled. We are entitled to charge default interest if the 30-day payment deadline is exceeded. The interest rate will be the customary market rate charged by major banks. We expressly reserve the right to effect delivery with cash on delivery.

Data Protection

BEROLLKA and the customer are responsible themselves for compliance with the respective data protection rules applicable to them. The customer must not provide BEROLLKA with personal data that is illegal or that the customer is not legally authorised to distribute.

The place of performance and the place of jurisdiction will be the place where the invoice is issued.

Heidelberg, Germany, will be the exclusive place of jurisdiction for any and all disputes arising under the contractual relationship. However, BEROLLKA is entitled to sue the customer at any other legal venue. These General Terms and Conditions of Business are governed by the laws of the Federal Republic of Germany, excluding of the United Nations Convention on Contracts for the International Sale of Goods (CISG).

Technical modifications and printing errors excepted. A current list of prices can be ordered at any time.