

EASYLINER®

ELDT | ELFR

Gel Liner, Reinforced –
With or Without Distal Attachment



DESCRIPTION

ALPS® EasyLiner® is formulated with EasyGel that contains powerful antioxidants proven to be highly beneficial to residual limb skin. Through the EasyGel and high-performance fabric, the EasyLiner offers a perfect combination of softness, strength, elasticity, comfort, and durability.

Refer to the case studies in the product catalog for evidence of the benefits.

PART #	SUSPENSION	L-CODE
ELDT	LOCKING	L5673
ELFR	CUSHION	L5679

Recommended for all activity levels.



TOP FEATURES & BENEFITS



BEIGE FABRIC
HORIZONTAL &
VERTICAL STRETCH



EASYGEL



LOCKING




CUSHION

- ✓ Reinforced seams improves durability and appearance
- ✓ Flat encapsulation provides additional comfort

- ✓ High performance fabric to improve durability
- ✓ Available in 3mm and 6mm uniform, and tapered from 6mm at distal end to 3mm


CROSS SECTION DIAGRAM

3mm Uniform, 6mm Uniform, 6/3mm Tapered



HOW TO SIZE

Measurements taken 6cm from the distal end



SIZING CHART

MEASURED SIZE (CM)	LOCKING 3.0MM, 6MM, 6/3MM, OR 9MM	CUSHION 3.0MM, 6MM, 6/3MM, OR 9MM
16-19 cm	ELDT 16 - (3, 6, or 6/3)	ELFR 16 - (3, 6, or 6/3)
20-23 cm	ELDT 20 - (3, 6, or 6/3)	ELFR 20 - (3, 6, or 6/3)
24-25 cm	ELDT 24 - (3, 6, 6/3, or 9)	ELFR 24 - (3, 6, 6/3, or 9)
26-27 cm	ELDT 26 - (3, 6, or 6/3)	ELRF 26 - (3, 6, or 6/3)
28-31 cm	ELDT 28 - (3, 6, 6/3, or 9)	ELRF 28 - (3, 6, 6/3, or 9)
32-34 cm	ELDT 32 - (3, 6, or 6/3)	ELFR 32 - (3, 6, or 6/3)
35-37 cm	ELDT 35 - (3 or 6)	ELFR 35 - (3 or 6)
38-43 cm	ELDT 38 - (3, 6, or 6/3)	ELFR 38 - (3, 6, or 6/3)
44-53 cm	ELDT 44 - (3, 6, or 6/3)	ELFR 44 - (3, 6, or 6/3)

*Expired patent owned by ALPS and/or one of its subsidiaries



T: 727.528.8566 | 800.574.5426 | F: 727.528.8862
www.easyliner.com | info@easyliner.com
 2895 42nd Avenue North St. Petersburg, Florida 33714

