

ottobock.

Kimba Family



Products for combinations



Kimba seating unit



Kimba outdoor mobility base



Kimba Cross



Kimba Inline



Kimba Home





Kimba indoor mobility base



Robby bathing support

What do you want to order?

- | | |
|--------------------------|------------------------------------|
| <input type="checkbox"/> | Kimba seating unit |
| <input type="checkbox"/> | Kimba outdoor mobility base |
| <input type="checkbox"/> | Kimba Cross |
| <input type="checkbox"/> | Kimba Inline |
| <input type="checkbox"/> | Kimba Home |
| <input type="checkbox"/> | Kimba indoor mobility base |
| <input type="checkbox"/> | Robby bathing support |

	Dimensions kid	Dimensions Kimba size 1	Dimensions Kimba size 2
			
① Back length [mm]		410 - 610	560 - 700
② Seat depth [mm]		190 - 310 (140 shortened)	240 - 400
Seat width [mm]		200 - 300	260 - 350
③ Lower leg length [mm]		190 - 310 (122 shortened)	200 - 370
Load capacity		Seating unit: 40 kg* Outdoor mobility base: 55 kg	
a) User weight			
b) Medical Equipment			
* With option „footrest angle adjustable for short lower leg length“ the maximum load decreases to 35 kg.			

Kimba Seating unit

Maximum load capacity: 40 kg

Article number: 470G71=5_AA02_C

valid from 01.01.2022 to 31.12.2022

Order no.:



Figure includes accessories

☐ **Provide estimate only**

☐ **Order**

☐ **Attachment for custom fabrication enclosed**

Billing address/customer number

Delivery address

Sales person

Name

.....
Date Signature

By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals.










Grundausrüstung












- Seat and back support
- Seat width, seat depth, back support height and lower leg length are continuously adjustable
- Seat inclination and knee angle continuously adjustable
- Back support angle adjustable
- Optimized quick-release mechanism to remove the seat
- Position the seat in/against the direction of travel
- Leg support assembly with flip-up foot plate
- Lateral side supports
- Adjustable pelvic laterals
- Frame colour: black
- Padding colour grey











Notice:

This product can be combined with outdoor mobility bases Kimba 470G71=5_AA03_C, Kimba Cross and Kimba Inline and also with indoor mobility bases Kimba Home and Kimba. Please fill out the corresponding order sheets separately.






Size		
KC01	<input type="checkbox"/>	Size 1 (from approx. 1 to 6 years) -
KC02	<input type="checkbox"/>	Size 2 (from approx. 4 to 10 years) Add. charge

Head support				
KD01	<input type="checkbox"/>		Without head support	-
KD02	<input type="checkbox"/>		With integrated head support with slight lateral support	Add. charge
KD03B/KD03			Head support (laterally adjustable)	<input type="checkbox"/> enclosed <input type="checkbox"/> mounted AC/AC
KL01B	<input type="checkbox"/>		Head support concave, S, enclosed (incl. holder; is delivered in combination with a short back)	Add. charge
KL02B	<input type="checkbox"/>		Head support concave, M, enclosed (incl. holder; is delivered in combination with a short back)	Add. charge
KL03B	<input type="checkbox"/>		Head support lateral, S, enclosed (incl. holder; is delivered in combination with a short back)	Add. charge
KL06B	<input type="checkbox"/>		Head support lateral adjustable, M, enclosed (incl. holder; is delivered in combination with a short back)	Add. charge
KL05B	<input type="checkbox"/>		Head support holder, enclosed (is delivered in combination with a short back)	Add. charge
Padding				
KD00	<input type="checkbox"/>		Without (for individual seat and upholstery adjustment)	-
KD41B	<input type="checkbox"/>		Padding, enclosed	-
KD41	<input type="checkbox"/>		Padding, mounted	AC
Inlays				
KD33B	<input type="checkbox"/>		Inlay grey, enclosed	AC
KD34B	<input type="checkbox"/>		Inlay turquoise, enclosed	AC
KD35B	<input type="checkbox"/>		Inlay pink, enclosed	AC
KD36B	<input type="checkbox"/>		Inlay green, enclosed	AC

Pelvic belts					
KP02B/KP02		Lap belt, padded	<input type="checkbox"/> enclosed	<input type="checkbox"/> mounted	AC/AC
KP18B/KP18		4-Point paediatric lap belt, with plastic buckles (width: 2,5 cm; length: 100 cm)	<input type="checkbox"/> enclosed	<input type="checkbox"/> mounted	AC/AC
KP06B/KP06		Groin strap	<input type="checkbox"/> enclosed	<input type="checkbox"/> mounted	AC/AC
KP15B/KP15		Pediatric lap belt, with buckles (width: 2,5 cm; length: 100 cm)	<input type="checkbox"/> enclosed	<input type="checkbox"/> mounted	AC/AC
Chest and shoulder supports					
KP16B/KP16		Shoulder harness Size S	<input type="checkbox"/> enclosed	<input type="checkbox"/> mounted	AC/AC
KP17B/KP17		Shoulder harness Size M	<input type="checkbox"/> enclosed	<input type="checkbox"/> mounted	AC/AC
KP19B/KP19		Chest and shoulder pad Size S	<input type="checkbox"/> enclosed	<input type="checkbox"/> mounted	AC/AC
KP20B/KP20		Chest and shoulder pad Size M	<input type="checkbox"/> enclosed	<input type="checkbox"/> mounted	AC/AC
KP21B/KP21		Slimline sternum harness	<input type="checkbox"/> enclosed	<input type="checkbox"/> mounted	AC/AC
KP22B/KP22		Dynamic harness Size XS	<input type="checkbox"/> enclosed	<input type="checkbox"/> mounted	AC/AC
KP23B/KP23		Dynamic harness Size S	<input type="checkbox"/> enclosed	<input type="checkbox"/> mounted	AC/AC
KP24B/KP24		Dynamic harness Size M	<input type="checkbox"/> enclosed	<input type="checkbox"/> mounted	AC/AC
KP03B/KP03		4-point belt, padded	<input type="checkbox"/> enclosed	<input type="checkbox"/> mounted	AC/AC
KP05B/KP05		5-Point belt, small buckle (recommended for seating unit size 1)	<input type="checkbox"/> enclosed	<input type="checkbox"/> mounted	AC/AC
KP04B/KP04		5-Point belt, large buckle	<input type="checkbox"/> enclosed	<input type="checkbox"/> mounted	AC/AC
KP10B/KP10		Positioning vest	<input type="checkbox"/> enclosed	<input type="checkbox"/> mounted	AC/AC

Preassembly					
KP14	<input type="checkbox"/>		Mounting kit strap/vest, mounted		Add. charge
Trunk support					
KD04B/KD04			Thoracic support	<input type="checkbox"/> enclosed <input type="checkbox"/> mounted	AC/AC
KD31B	<input type="checkbox"/>		Swing away thoracic supports (laterals), enclosed		Add. charge
KD10B	<input type="checkbox"/>		Lumbar support padding, enclosed		Add. charge
KD11B	<input type="checkbox"/>		Neck support, enclosed		Add. charge
KC10	<input type="checkbox"/>		Seat depth reducer (seat depth possible up to 14 cm)		Add. charge
KC13B	<input type="checkbox"/>		Inlay small, enclosed (spacer fabric grey, only for size 1)		Add. charge
Arm supports					
KE01	<input type="checkbox"/>		Arm support, Standard, mounted		Add. charge
KE02	<input type="checkbox"/>		Arm support high/low, adjustable, mounted (not possible with country of destination Germany)		Add. charge
Tray					
KS01B	<input type="checkbox"/>		Tray, grey, enclosed (only possible in direction of travel and in combination with KE 02 „arm support high/ low“ in lowest setting only)		Add. charge
KS02B	<input type="checkbox"/>		Tray, transparent, enclosed (only possible in direction of travel and in combination with KE 02 „arm support high/ low“ in lowest setting only))		Add. charge

Leg and foot support					
KB03	<input type="checkbox"/>		Leg support, Standard, mounted (from LLL 19 cm) –		
KB05	<input type="checkbox"/>		Leg support, angle-adjustable, mounted (from LLL 19 cm) Add. charge		
KB10	<input type="checkbox"/>		Leg support, angle-adjustable for short lower leg length, mounted (only for size 1, from LLL 14 cm) Add. charge		
Accessories leg and foot support					
KB08B/KB08			Slip resistant footrest cover (not mounted in combination with KB10)	<input type="checkbox"/> enclosed	<input type="checkbox"/> mounted AC/AC
KB09B	<input type="checkbox"/>		Padding, enclosed Add. charge		
KB12B	<input type="checkbox"/>		Foot support – sandals, enclosed Add. charge		
KB07B	<input type="checkbox"/>		Velcro footstraps, instep padded Add. charge		
KC06B/KC06			Abductor size 1, only for seat size 1, (L 9 cm × W 6 cm × H 9 cm)	<input type="checkbox"/> enclosed	<input type="checkbox"/> mounted AC/AC
KC07B/KC07			Abductor size 2 (L 10 cm × W 7,5 cm × H 10 cm)	<input type="checkbox"/> enclosed	<input type="checkbox"/> mounted AC/AC
KC08B/KC08			Abductor size 3, only for seat size 2 (L 12 cm × W 9 cm × H 10 cm)	<input type="checkbox"/> enclosed	<input type="checkbox"/> mounted AC/AC
Canopy					
KW13	<input type="checkbox"/>		Canopy in grey-melange, inside grey, incl. rain cover and mounted adapter Add. charge		

Slip sack				
KW01B	<input type="checkbox"/>		Summer/winter slipsack Size S, enclosed	Add. charge
KW02B	<input type="checkbox"/>		Summer/winter slipsack Size M, enclosed	Add. charge
KW03B	<input type="checkbox"/>		Summer/winter slipsack Size L, enclosed	Add. charge
Accessories seat				
KQ04B	<input type="checkbox"/>		Grip bar with padding, enclosed	Add. charge
KD20	<input type="checkbox"/>		Back angle limitation, mounted	Add. charge
KW04B	<input type="checkbox"/>		Rain cape, grey, enclosed	Add. charge
KD12B	<input type="checkbox"/>		Incontinence cover, enclosed	Add. charge
KC16B	<input type="checkbox"/>		Retrofit kit for Kimba Neo street frame (without tilt; not possible in combination with costum fabrication platform for respirators up to 90 kg/ 198,4 lb)	Add. charge

Kimba Outdoor mobility base

Maximum load capacity: 55 kg
(including seating unit)

Article number: 470G71=5_AA03_C
valid from 01.01.2022 to 31.12.2022

Order no.:



Standard equipment

- Folding outdoor mobility base made of high-grade aluminium
- Optimized rear axle suspension
- Push bar with ratchet joint, height-adjustable
- Seat tilt with tilt activation at the frame incl. locking pin
- Integrated storage space (also usable as a platform for respirators)
- Automatic folding lock
- Optimized seat adapter with automatic tilt limitation
- EVA tyres with frictionbrake
- Front wheels 7", swivelling with integrated swivel lock
- Integrated anchor points for mounting to wheelchair accessible vehicle

☐ **Provide estimate only**

☐ **Order**

☐ **Attachment for custom fabrication enclosed**

Billing address/customer number

Delivery address

Sales person

Name

.....
Date Signature

By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals.

Size

KA01



Size 1 (from approx. 1 to 6 years)

KA05



Size 2 (from approx. 4 to 10 years)

Add. charge

Adapter Sitzschale

KC03



Seating shell interface

Add. charge

KB01



Leg support assistance for seating shell interface

Add. charge

KC14



Universal adapter (seat size 1 to base size 2)
(only for country of destination Brazil)

Add. charge

KC15B



Retrofit-kit: Kimba Neo seat

Add. charge

Frame colour			
013	<input type="checkbox"/>	Jet black RAL 9005 matte	
023	<input type="checkbox"/>	White RAL 9016	




Reflector sticker			
KA13	<input type="checkbox"/>	Classic	Add. charge
KA14	<input type="checkbox"/>	Sporty	Add. charge





Classic









Sporty

Colour reflector sticker			
	<input type="checkbox"/>		Grey
	<input type="checkbox"/>		Green
	<input type="checkbox"/>		Turquoise

Tyres			
KF03	<input type="checkbox"/>		EVA -
KF04	<input type="checkbox"/>		Pneumatic Add. charge

Rear wheel type			
KG01	<input type="checkbox"/>	Wheel lock (actuated by the foot via wheel lock bar)	-
KG02	<input type="checkbox"/>	With frictionbrake (manually actuated via brake levers at the push bar)	Add. charge

Accessories			
KA07B	<input type="checkbox"/>		Frame padding, pair, enclosed Add. charge
KS03B	<input type="checkbox"/>		Buggy board, enclosed Add. charge
KA08	<input type="checkbox"/>		Tip-assist, left Add. charge
KA12	<input type="checkbox"/>		Tip-assist, right Add. charge
KS14B	<input type="checkbox"/>		Tether band, enclosed Add. charge
Accessories general			
KS07B	<input type="checkbox"/>		Cupholder, enclosed Add. charge
KS13B	<input type="checkbox"/>		Diaper bag, grey-melange, enclosed (Contains extensive accessories such as: changing mat, dummy bag, isolation pocket for bottles and a plastic utensil bag, as well as numerous compartments and zipped pockets) Add. charge

Kimba Cross

Maximum load capacity of the outdoor mobility base: 50 kg/110 lbs (including seating unit)

Article number: HR3401=0000_K

valid from 01.01.2022 to 31.12.2022

Order no.:



Standard equipment

- Three-wheeled outdoor mobility base with telescoping frame, for mounting a Kimba neo seating unit, size 1 or size 2, or a custom seating shell
- Can be folded/disassembled
- Elastomer spring system
- Drum brake
- 18" rear wheels with pneumatic tyres
- Rigidly mounted 16" front wheel with pneumatic tyre, includes splash guard
- Height-adjustable push bar
- Frame colour: black

For children from approx. 1 to 10 years

☐ **Provide estimate only**

☐ **Order**

☐ **Attachment for custom fabrication enclosed**

Billing address/customer number

Delivery address

Sales person

Name

.....
Date

Signature

By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals.

Outdoor mobility base

☐

Kimba Cross outdoor mobility base

–

For use in conjunction with

☐

Kimba seating unit, size 1

☐

Kimba seating unit, size 2

Seat adapter for Kimba seating unit

HR34011300

☐

Size 1

–

HR34021300

☐

Size 2

–

HR32x48300

☐

Seating shell interface Kimba

Add. charge

Accessories for outdoor mobility base			
HR34012100	<input type="checkbox"/>	12" swivelling wheel with quick-lock attachment	Add. charge
HR34017000	<input type="checkbox"/>	Anti-tipper, complete	Add. charge
HR34018000	<input type="checkbox"/>	Clothing guard (plug-on wheel covers)	Add. charge
HR34018100	<input type="checkbox"/>	Clothing guard extension around the wheel (not available in combination with spoke protector)	Add. charge
HR34019000	<input type="checkbox"/>	Spoke protector, transparent	Add. charge
HR34015000	<input type="checkbox"/>	Storage bag (e.g. for the swivelling front wheel)	Add. charge
HR34016100	<input type="checkbox"/>	Bicycle draw-bar, complete (includes coupling from Weber company, type E, safety strap)*	Add. charge
HR34016200	<input type="checkbox"/>	Equipment package required according to German Road Traffic Licensing Regulations (StVZO) (lighting/pennant/roll bar)	Add. charge
HR16001102	<input type="checkbox"/>	Buggy board	Add. charge
HR32xx2700-	<input type="checkbox"/>	Foot support assembly for seating shell interface Kimba	Add. charge
* According to a leaflet for bicycle trailers issued by the German Federal Ministry of Transport: (www.pdeleuw.de/fahrrad/fakten) – Recommended overall weight of the bicycle trailer incl. load: 40 kg. Recommended in combination with four-point belt, groin strap, spoke protectors.			

HR34016200 – Equipment package required according to German Road Traffic Licensing Regulations (StVZO)
According to a leaflet for bicycle trailers issued by the German Ministry of Transport: (www.pdeleuw.de/fahrrad/fakten)
Recommended overall weight of the bicycle trailer, including load: 40 kg. Recommended in combination with four-point belt, groin strap, spoke protectors. Seating shell interface and footrest assembly in combination with Squiggles and Mygo not possible.
HR34018000 – Clothing guard is a plug-on wheel cover
HR34018100 – Clothing guard extension around the wheel cannot be used in combination with spoke protector
Accessories only for corresponding size
Footrest assembly only available in combination with seating shell interface

Kimba Inline

Maximum load capacity of the outdoor mobility base: 50 kg/110 lbs (including seating units)

Article number: HR3255=0000_K

valid from 01.01.2022 to 31.12.2022

Order no.:



Standard equipment

Kimba Inline twin outdoor mobility base for two Kimba seating units, size 1

- Pneumatic tyres provide suspension comfort
- Rear wheels with drum brake
- Frame colour: black

Twin buggy for children from approx. 1–6 years

* Notice: Seating unit, size 2, is available in custom fabrication upon request.

☐ Provide estimate only

☐ Order

☐ Attachment for custom fabrication enclosed

Billing address/customer number

Delivery address

Sales person

Name

.....
Date

Signature

By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals.

Inline outdoor mobility base

HR32551000-



Inline twin outdoor mobility base for up to two seating systems positioned one behind the other, with pneumatic front wheel 200 and pneumatic rear wheel 300, complete with drum brake

–

Accessories

HR16001102



Buggy board

Add. charge

HR32548300



Front seating shell interface

Add. charge

HR32548301



Rear seating shell interface

Add. charge

HR32582700



Footrest assembly for front seating shell interface

Add. charge

HR32582701



Footrest assembly for rear seating shell interface

Add. charge

Kimba Home

Max. load capacity (including seating unit):
Size 1: 50 kg
Size 2: 75 kg

valid from 01.01.2022 to 31.12.2022

Order no.:



- Indoor chassis for Kimba seating units size 1 and size 2
- Tilt adjustment of the seat unit via gas pressure spring 500 N (tilt adjustment range -5° till +35°)
- Height adjustment via gas pressure spring 1200 N (height of seat unit adapter min. 19 cm, max. 56 cm)
- 4 castors (75 mm), of these 2 with locking device
- Frame colour: Kimba white
- Push bar incl. release lever for seat angle adjustment

☐ Provide estimate only

☐ Order

Billing address/customer number

Delivery address

Sales person

Name

.....
Date

Signature

By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals.

Basic model

470G71=90001

☐

Kimba Home indoor chassis size 1

470G71=90002

☐

Kimba Home indoor chassis size 2

Accessories

471A71=SK910

☐

Anti-tip mechanism

Add. charge

471C71=SK730

☐

Kimba seating shell interface, size 1

Add. charge

471C71=SK830

☐

Kimba seating shell interface, size 2

Add. charge

HR32692700-000

☐

Foot support for seating shell interface, complete, size 2

Add. charge

471A71=SK921

☐

Gas spring (2 x 300 N) as support for seat lift function, recommended recommended from 30 kg load

Add. charge

471A71=SK922

☐

Gas spring (2 x 500 N) as support for seat lift function, recommended recommended from 40 kg load

Add. charge

Information



Castor with locking device



Foot pedal for height adjustment



Push bar with release lever for seat angle adjustment



Seat angle adjustment



Seat height adjustment



Adapter for Kimba seat unit



Anti-tip mechanism

Kimba Indoor mobility base

Article number: HR2233=0000_K

valid from 01.01.2022 to 31.12.2022

Order no.:



- Height-adjustable steel construction
- Caster with swivel lock, diameter 75 mm
- Angle-adjustable seat interface
- Adjustable wheel base
- Frame colour: anthracite metallic

☐ Provide estimate only ☐ Order

Billing address/customer number

Delivery address

Sales person

Name

.....
Date Signature

By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals.

Standard model

HR22330000-011	<input type="checkbox"/>	Size 1	-
HR22430000-011	<input type="checkbox"/>	Size 2	-

For use in conjunction with

	<input type="checkbox"/>	Kimba Neo seating unit, size 1
	<input type="checkbox"/>	Kimba Neo seating unit, size 2

Caster with swivel lock

HR22321150	<input type="checkbox"/>	Diameter 75 mm	-
HR22321160	<input type="checkbox"/>	Diameter 125 mm	Add. charge

Seat interface

HR32046x00	<input type="checkbox"/>	Parallel	Add. charge
HR32x48300	<input type="checkbox"/>	Kimba	Add. charge

Accessories

HR22327100	<input type="checkbox"/>	Push handle	Add. charge
HR32044x00-	<input type="checkbox"/>	Footrest assembly for seating shell interface	Add. charge

Robby

Bathing support

valid from 01.01.2022 to 31.12.2022

Order no.:



Standard equipment

- Folding bathing and lying support that can be adapted to grow with the user
- Net upholstery (blue)
- Lap belt and leg strap
- Frame colour: silver

☐ Provide estimate only

☐ Order

Billing address/customer number

Delivery address

Sales person

Name

Date

Signature

By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals.

Robby bathing support

HR51040000-027	<input type="checkbox"/>	Robby, size 1 (seat depth 20 – 24 cm, back support height incl. head section 63 – 72.5 cm) Bathing and lying support, includes washable net upholstery, lap belt, leg straps; max. load capacity 30 kg	Add. charge
HR51050000-027	<input type="checkbox"/>	Robby, size 2 (seat depth 31 – 42.5 cm, back support height incl. head section 89 – 94 cm) Bathing and lying support, includes washable net upholstery, lap belt, leg straps; max. load capacity 60 kg	Add. charge

Accessories

HR51045300	<input type="checkbox"/>	Trunk belt	Add. charge
HR51045200-103	<input type="checkbox"/>	Adjustable head support pads (can also be used as trunk support pads)	Add. charge

Information



Adjustable head support pads

Technical data

	Size 1	Size 2
Seat depth (depending on strap tension and expansion of telescoping tubes)	200–240 mm	310–425 mm
Back support height	630–725 mm	890–940 mm
Overall width	360 mm	380 mm
Overall length (support surface)	1150–1330 mm	1470–1730 mm
Max. load capacity	30 kg	60 kg
Weight	3.8 kg	4.3 kg

Notes & Remarks