

# ottobock.

## Motus CS

### Adaptive wheelchair

Article-number: 480F61=2\_AA03\_C

valid from 01.10.2021 to 31.12.2021

Order no.:

-----



**Provide estimate only**

**Order**

**Attachment for custom fabrication enclosed**

Billing address/costumer number

Delivery address

Sales person

Name

.....  
Date

Signature

By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals.

### Standard equipment

- Aluminium tube frame with crossbrace
- Folding
- Integrated leg support holders
- Single-panel foot plates, angle-adjustable
- Back support tubes with different height setting options
- Back support upholstery, adaptable
- Clothing guard side panels
- Drive wheels with hollow rims and quick-release axles, 1" tyres
- Knee lever wheel lock

Max. user weight: 125 kg/275 lbs

Crash-tested according to DIN ISO 7176-19

Wheelchair type	
MB50	<input type="checkbox"/> Self-propelled chair * -
	<input type="checkbox"/> Standard -

\* Risk due to getting caught on obstacles. Please observe safety information in the instructions for use!

Seat measurements	
Seat width SW [cm]	<input type="checkbox"/> 35,5 <input type="checkbox"/> 38 <input type="checkbox"/> 40,5 <input type="checkbox"/> 43 <input type="checkbox"/> 45,5 <input type="checkbox"/> 48
Seat depth SD [cm]	<input type="checkbox"/> 36 <input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 44 <input type="checkbox"/> 46 <input type="checkbox"/> 48 <input type="checkbox"/> 50 <input type="checkbox"/> 52 <input type="checkbox"/> 54

Anterior seat height [cm] (see table page 6)	
<input type="checkbox"/> 38 <input type="checkbox"/> 39 <input type="checkbox"/> 40 <input type="checkbox"/> 41 <input type="checkbox"/> 42 <input type="checkbox"/> 43 <input type="checkbox"/> 44 <input type="checkbox"/> 45 <input type="checkbox"/> 46 <input type="checkbox"/> 47 <input type="checkbox"/> 48	
<input type="checkbox"/> 49 <input type="checkbox"/> 50 <input type="checkbox"/> 51 <input type="checkbox"/> 52 <input type="checkbox"/> 53 <input type="checkbox"/> 54 <input type="checkbox"/> 55	

Lower leg length [cm] (without seat cushion)	
Lower leg length (LL) [cm] ..... (available adjustment range 34 cm – 55 cm)	

Posterior seat height [cm]	
<input type="checkbox"/> 36 <input type="checkbox"/> 37 <input type="checkbox"/> 38 <input type="checkbox"/> 39 <input type="checkbox"/> 40 <input type="checkbox"/> 41 <input type="checkbox"/> 42 <input type="checkbox"/> 43 <input type="checkbox"/> 44 <input type="checkbox"/> 45 <input type="checkbox"/> 46	
<input type="checkbox"/> 47 <input type="checkbox"/> 48 <input type="checkbox"/> 49 <input type="checkbox"/> 50 <input type="checkbox"/> 51 <input type="checkbox"/> 52	

Back support height [cm]	
<input type="checkbox"/> 30 <input type="checkbox"/> 32,5 <input type="checkbox"/> 35 <input type="checkbox"/> 37,5 <input type="checkbox"/> 40 <input type="checkbox"/> 42,5 <input type="checkbox"/> 43 <input type="checkbox"/> 45 <input type="checkbox"/> 47,5 <input type="checkbox"/> 48 <input type="checkbox"/> 50	

Horizontal drive wheel position	
<input type="checkbox"/> A <input type="checkbox"/> B <input type="checkbox"/> C <input type="checkbox"/> D	

Drive wheel camber	
<input type="checkbox"/> 0° <input type="checkbox"/> 2° <input type="checkbox"/> 3° <input type="checkbox"/> 4°	

Drive wheel rim	
MG06	<input type="checkbox"/> Hollow rim -
MG15	<input type="checkbox"/> Infinity Ultralight Add. charge
MG80	<input type="checkbox"/> Drive wheel with drum brake/attendant Add. charge

Tyres	
MG31	<input type="checkbox"/> Schwalbe RightRun 1" -
MG30	<input type="checkbox"/> Pneumatic tyre, 1 3/8" -
MG34	<input type="checkbox"/> PU 1 3/8" -
MG341	<input type="checkbox"/> PU 1", black Add. charge
MG37	<input type="checkbox"/> Schwalbe Marathon Plus Evolution 1" Add. charge
MG371	<input type="checkbox"/> Schwalbe Marathon Plus Evolution 1.1" Add. charge
MG200	<input type="checkbox"/> Schwalbe One 1" Add. charge

Handrims	
MG52	<input type="checkbox"/> Aluminium silver anodised, high -
MG53	<input type="checkbox"/> Aluminium silver anodised, low Add. charge
MG54	<input type="checkbox"/> Rubber-coated, high Add. charge
MG56	<input type="checkbox"/> Stainless steel, polished, high Add. charge
MG65	<input type="checkbox"/> GEKKO, ergonomic handrim Add. charge

Wheelbase extension	
MG86	<input type="checkbox"/> Yes Add. charge

Drive wheel accessories			
MG87	<input type="checkbox"/>	High-pressure air pump, 14 bar	Add. charge
MG88	<input type="checkbox"/>	Airman pump, battery-operated	Add. charge
Caster wheels			
MF30	<input type="checkbox"/>	4" solid rubber, black	-
MF31	<input type="checkbox"/>	5" solid rubber, black	-
MF04	<input type="checkbox"/>	5,5" soft rubber, Aluminium	-
MF32	<input type="checkbox"/>	6" solid rubber, black	-
MF08	<input type="checkbox"/>	7" solid rubber	-
Frame accessories			
MA51	<input type="checkbox"/>	Tip-assist, right	Add. charge
MA52	<input type="checkbox"/>	Tip-assist, left	Add. charge
MA55	<input type="checkbox"/>	Anti-tipper, right	Add. charge
MA56	<input type="checkbox"/>	Anti-tipper, left	Add. charge
MA57	<input type="checkbox"/>	Crutch holder, right	Add. charge
MA58	<input type="checkbox"/>	Crutch holder, left	Add. charge
Wheel lock			
MH01	<input type="checkbox"/>	Push-to-lock knee lever wheel lock	-
MH05	<input type="checkbox"/>	OutFront scissor wheel lock	Add. charge
MH06	<input type="checkbox"/>	Knee lever wheel lock for user and attendant	Add. charge
Wheel lock accessories			
MH10	<input type="checkbox"/>	Wheel lock lever extension	Add. charge
Back support system			
MD01	<input type="checkbox"/>	Back support, straight	-
MD09	<input type="checkbox"/>	Back support with lumbar curve 8°	-
MD10	<input type="checkbox"/>	Back support, angle-adjustable 30° (back support height 43 or 48 cm)	Add. charge
MD11	<input type="checkbox"/>	Foldable	Add. charge
Push handles			
MD50	<input type="checkbox"/>	Short	-
MD51	<input type="checkbox"/>	Long	-
MD53	<input type="checkbox"/>	Telescoping, short	Add. charge
MD52	<input type="checkbox"/>	Telescoping, long	Add. charge
MD55	<input type="checkbox"/>	Height-adjustable/removable	Add. charge
MD56	<input type="checkbox"/>	Foldable	Add. charge

<b>Back support system range of adjustment</b>			
	<input type="checkbox"/>	<b>30 - 32,5</b> (back support straight with push handles height-adjustable/telescoping)	-
	<input type="checkbox"/>	<b>30 - 37,5</b> (back support straight with push handles short/long/folding)	-
	<input type="checkbox"/>	<b>32,5 - 40</b> (back support straight with push handles height-adjustable/telescoping)	-
	<input type="checkbox"/>	<b>35 - 40</b> (back support 8° with push handles height-adjustable/telescoping)	-
	<input type="checkbox"/>	<b>37,5 - 45</b> (back support straight with push handles height-adjustable/telescoping)	-
	<input type="checkbox"/>	<b>40 - 45</b> (back support straight with push handles short/long/folding)	-
	<input type="checkbox"/>	<b>42,5 - 50</b> (back support 8° or back support foldable with push handles long)	-
	<input type="checkbox"/>	<b>45 - 50</b> (back support 8°)	-
	<input type="checkbox"/>	<b>47,5 - 50</b> (back support straight with push handles short/long/folding or back support foldable with push handles long)	-
<b>Back support upholstery</b>			
<b>MD04</b>	<input type="checkbox"/>	<b>Standard</b> (only for back support, angle-adjustable 30°)	-
<b>MD03</b>	<input type="checkbox"/>	<b>Adaptable, padded</b>	-
<b>MD05</b>	<input type="checkbox"/>	<b>Adaptable, breathable</b>	<b>Add. charge</b>
<b>MD06</b>	<input type="checkbox"/>	<b>Without</b>	-
<b>Back support pad widening</b>			
<b>MD80</b>	<input type="checkbox"/>	<b>+ 2 cm</b>	<b>Add. charge</b>
<b>MD81</b>	<input type="checkbox"/>	<b>+ 4 cm</b>	<b>Add. charge</b>
<b>MD82</b>	<input type="checkbox"/>	<b>Without</b>	-
<b>Use of BAXX back support system</b> (include order form, model-specific restrictions are possible)			
	<input type="checkbox"/>	<b>Baxx Aluminium back support</b> (enclosed)	<b>Add. charge</b>
<b>Back support accessories</b>			
<b>MD58</b>	<input type="checkbox"/>	<b>Stabiliser bar</b>	<b>Add. charge</b>
<b>Side panels</b>			
<b>ME01</b>	<input type="checkbox"/>	<b>Aluminium with protection against cold</b>	-
<b>ME11</b>	<input type="checkbox"/>	<b>Plastic plug-on side panel</b>	-
<b>ME02</b>	<input type="checkbox"/>	<b>Aluminium with clothing guard and against cold</b>	-
<b>ME04</b>	<input type="checkbox"/>	<b>Carbon with clothing guard</b>	<b>Add. charge</b>
<b>ME17</b>	<input type="checkbox"/>	<b>Aluminium with arm support long, depth-adjustable</b>	<b>Add. charge</b>
<b>ME05</b>	<input type="checkbox"/>	<b>Aluminium with arm support short</b>	<b>Add. charge</b>
<b>Side panels accessories</b>			
<b>ME50</b>	<input type="checkbox"/>	<b>Arm support padded</b> <input type="checkbox"/> <b>short</b> <input type="checkbox"/> <b>long</b> <input type="checkbox"/> <b>long, extra high</b>	<b>Add. charge</b>
<b>ME60</b>	<input type="checkbox"/>	<b>Therapy tray</b>	<b>Add. charge</b>
<b>Seat upholstery</b>			
<b>MC13</b>	<input type="checkbox"/>	<b>Standard</b>	-
<b>MC14</b>	<input type="checkbox"/>	<b>Readjustable</b>	<b>Add. charge</b>
<b>MC06</b>	<input type="checkbox"/>	<b>Without</b>	-

Seat cushion			
MC01	<input type="checkbox"/>	<b>Standard</b> (seat cushion thickness 5 cm)	<b>Add. charge</b>
MC02	<input type="checkbox"/>	<b>Cubic Foam</b> (seat cushion thickness 5 cm)	<b>Add. charge</b>
MC05	<input type="checkbox"/>	<b>Evolight</b> (seat cushion thickness 6 cm)	<b>Add. charge</b>
MC29	<input type="checkbox"/>	<b>Livo</b> (seat cushion thickness 5 cm)	<b>Add. charge</b>
MC30	<input type="checkbox"/>	<b>Evoflair</b> (seat cushion thickness 8 cm)	<b>Add. charge</b>

Leg supprts			
MB03	<input type="checkbox"/>	<b>Single-panel, angle-adjustable</b>	-

Design of foot support bar			
MB200	<input type="checkbox"/>	<b>Standard</b>	-
MB201	<input type="checkbox"/>	<b>Swing away, left</b>	<b>Add. charge</b>
MB202	<input type="checkbox"/>	<b>Swing away, right</b>	<b>Add. charge</b>

Lenght of foot support bar (without seat cushion)			
MB40	<input type="checkbox"/>	<b>Short</b> (max. 34 cm – 41 cm)	-
MB41	<input type="checkbox"/>	<b>Medium</b> (max. 41 cm – 51 cm)	-
MB42	<input type="checkbox"/>	<b>Long</b> (max. 47 cm – 55 cm)	-

Other accessories			
MZ51	<input type="checkbox"/>	<b>Lap belt with metal buckle</b>	<b>Add. charge</b>
MZ52	<input type="checkbox"/>	<b>Tool kit</b>	<b>Add. charge</b>

Please also use the order form for wheelchair accessories 647F661=\*

Spoke protector			
MG90	<input type="checkbox"/>	<b>Transparent</b>	<b>Add. charge</b>

Frame colours			
002	<input type="checkbox"/>	<b>Anthracite metallic</b>	-
010	<input type="checkbox"/>	<b>Marine blue</b>	<b>Add. charge</b>
012	<input type="checkbox"/>	<b>Sky blue</b>	<b>Add. charge</b>
014	<input type="checkbox"/>	<b>Candy red</b>	<b>Add. charge</b>
016	<input type="checkbox"/>	<b>Silver metallic</b>	<b>Add. charge</b>
021	<input type="checkbox"/>	<b>Sparkle granny-smith</b>	<b>Add. charge</b>
023	<input type="checkbox"/>	<b>White RAL 9016</b>	<b>Add. charge</b>
048	<input type="checkbox"/>	<b>Jet black matte</b> (structure)	<b>Add. charge</b>

Side panes in frame colour*			
ME20	<input type="checkbox"/>	<b>Yes</b>	<b>Add. charge</b>


\*: Standard colour jet black matte (texture); not possible in combination with side panels ME04 & ME11

Special accessories			
MZ58	<input type="checkbox"/>	<b>Bag for drive wheels</b> (for 2 drive wheels)	<b>Add. charge</b>
MZ59	<input type="checkbox"/>	<b>Pocket for mobile phone</b>	<b>Add. charge</b>
MZ89	<input type="checkbox"/>	<b>Protection for quick-release-axle</b>	<b>Add. charge</b>
MZ92	<input type="checkbox"/>	<b>Wheel Blades</b> (mini-skis for caster wheels)	<b>Add. charge</b>
MZ93	<input type="checkbox"/>	<b>Backpack</b>	<b>Add. charge</b>

**Information**

<b>Anterior seat height [mm]</b>					
caster wheel fork	wheel diameter				
	4"	5"	5,5"	6"	7"
short	380 - 460	380 - 470	400 - 480	400 - 480	420 - 490
long	440 - 530	460 - 540	440 - 550	430 - 550	440 - 550

**Positions drive wheel adapter**



Information



Drive wheel adapter



MA51 Tip-assist



MA55 Anti-tipper



MA57 Crutch holder



MH01 Knee lever wheel lock



MH05 Outfront-scissor wheel lock



MH10 Knee lever wheel lock extension




MD10 Back support, angle-adjustable



MD11 Foldable back support



MD55 Push handles height-adjustable, removable



MD56 Push handles foldable

**Information**



MD04 Back support upholstery  
Standard



MD03 Back support upholstery,  
padded



MD05 Back support upholstery,  
breathable



MD58 Stabiliser bar



ME11 Plug-on side panel plastic



ME02 Side panel with clothing  
guard and against cold



ME04 Carbon side panel



ME17 Side panel with arm support  
long, height-adjustable



ME50 Side panel, padded



ME60 Therapy tray



MB29 Livo



MB30 Evoflair



MB03 Single-panel leg supports,  
angle-adjustable



MB201 Swing away foot support  
bar, left



## Information



MZ58 Bag for drive wheels



MZ59 Poket for mobile phone



MZ89 Protection for quick-release-axle



MZ92 Wheel Blades



MZ93 Backpack

**Notes & Remarks**