

ottobock.**Motus CV****Adaptive wheelchair****Article-number: 480F61=2_AA02_C****valid from 01.10.2021 to 31.12.2021**

Order no.:

 Provide estimate only **Order** **Attachment for custom fabrication enclosed**

Billing address/costumer number

Delivery address

Sales person

Name

Date

Signature

By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals.

**Standard equipment**

- Aluminium tube frame with crossbrace
- Folding
- Swing-away leg support holders
- Individual foot plate, angle-adjustable
- Back support tubes with different height setting options
- Back support upholstery, adaptable
- Clothing guard side panels
- Drive wheels with hollow rims and quick-release axles, 1" tyres
- Knee lever wheel lock

Max. user weight: 125 kg/275 lbs

Crash-tested according to DIN ISO 7176-19

Wheelchair type	
MB50	<input type="checkbox"/> Self-propelled chair * -
	<input type="checkbox"/> Standard -

* Risk due to getting caught on obstacles. Please observe safety information in the instructions for use!

Seat measurements	
Seat width SW [cm]	<input type="checkbox"/> 35,5 <input type="checkbox"/> 38 <input type="checkbox"/> 40,5 <input type="checkbox"/> 43 <input type="checkbox"/> 45,5 <input type="checkbox"/> 48
Seat depth SD [cm]	<input type="checkbox"/> 36 <input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 44 <input type="checkbox"/> 46 <input type="checkbox"/> 48 <input type="checkbox"/> 50 <input type="checkbox"/> 52 <input type="checkbox"/> 54

Anterior seat height [cm] (see table page 6)	
<input type="checkbox"/> 38 <input type="checkbox"/> 39 <input type="checkbox"/> 40 <input type="checkbox"/> 41 <input type="checkbox"/> 42 <input type="checkbox"/> 43 <input type="checkbox"/> 44 <input type="checkbox"/> 45 <input type="checkbox"/> 46 <input type="checkbox"/> 47 <input type="checkbox"/> 48	
<input type="checkbox"/> 49 <input type="checkbox"/> 50 <input type="checkbox"/> 51 <input type="checkbox"/> 52 <input type="checkbox"/> 53 <input type="checkbox"/> 54 <input type="checkbox"/> 55	

Lower leg length [cm] (without seat cushion)	
Lower leg length (LL) [cm] (available adjustment range 34 cm – 55 cm)	

Posterior seat height [cm]	
<input type="checkbox"/> 36 <input type="checkbox"/> 37 <input type="checkbox"/> 38 <input type="checkbox"/> 39 <input type="checkbox"/> 40 <input type="checkbox"/> 41 <input type="checkbox"/> 42 <input type="checkbox"/> 43 <input type="checkbox"/> 44 <input type="checkbox"/> 45 <input type="checkbox"/> 46	
<input type="checkbox"/> 47 <input type="checkbox"/> 48 <input type="checkbox"/> 49 <input type="checkbox"/> 50 <input type="checkbox"/> 51 <input type="checkbox"/> 52	

Back support height [cm]	
<input type="checkbox"/> 30 <input type="checkbox"/> 32,5 <input type="checkbox"/> 35 <input type="checkbox"/> 37,5 <input type="checkbox"/> 40 <input type="checkbox"/> 42,5 <input type="checkbox"/> 43 <input type="checkbox"/> 45 <input type="checkbox"/> 47,5 <input type="checkbox"/> 48 <input type="checkbox"/> 50	

Horizontal drive wheel position	
<input type="checkbox"/> A <input type="checkbox"/> B <input type="checkbox"/> C <input type="checkbox"/> D	

Drive wheel camber	
<input type="checkbox"/> 0° <input type="checkbox"/> 2° <input type="checkbox"/> 3° <input type="checkbox"/> 4°	

Drive wheel rim	
MG06	<input type="checkbox"/> Hollow rim -
MG15	<input type="checkbox"/> Infinity Ultralight Add. charge
MG80	<input type="checkbox"/> Drive wheel with drum brake/attendant Add. charge

Tyres	
MG31	<input type="checkbox"/> Schwalbe RightRun 1" -
MG30	<input type="checkbox"/> Pneumatic tyre, 1 3/8" -
MG34	<input type="checkbox"/> PU 1 3/8" -
MG341	<input type="checkbox"/> PU 1", black Add. charge
MG37	<input type="checkbox"/> Schwalbe Marathon Plus Evolution 1" Add. charge
MG371	<input type="checkbox"/> Schwalbe Marathon Plus Evolution 1.1" Add. charge
MG200	<input type="checkbox"/> Schwalbe One 1" Add. charge

Handrims	
MG52	<input type="checkbox"/> Aluminium silver anodised, high -
MG53	<input type="checkbox"/> Aluminium silver anodised, low Add. charge
MG54	<input type="checkbox"/> Rubber-coated, high Add. charge
MG56	<input type="checkbox"/> Stainless steel, polished, high Add. charge
MG65	<input type="checkbox"/> GEKKO, ergonomic handrim Add. charge

Wheelbase extension	
MG86	<input type="checkbox"/> Yes Add. charge

Drive wheel accessories			
MG87	<input type="checkbox"/>	High-pressure air pump, 14 bar	Add. charge
MG88	<input type="checkbox"/>	Airman pump, battery-operated	Add. charge
Caster wheels			
MF30	<input type="checkbox"/>	4" solid rubber, black	-
MF31	<input type="checkbox"/>	5" solid rubber, black	-
MF04	<input type="checkbox"/>	5,5" soft rubber, Aluminium	-
MF32	<input type="checkbox"/>	6" solid rubber, black	-
MF08	<input type="checkbox"/>	7" solid rubber	-
Frame accessories			
MA51	<input type="checkbox"/>	Tip-assist, right	Add. charge
MA52	<input type="checkbox"/>	Tip-assist, left	Add. charge
MA55	<input type="checkbox"/>	Anti-tipper, right	Add. charge
MA56	<input type="checkbox"/>	Anti-tipper, left	Add. charge
MA57	<input type="checkbox"/>	Crutch holder, right	Add. charge
MA58	<input type="checkbox"/>	Crutch holder, left	Add. charge
Wheel lock			
MH01	<input type="checkbox"/>	Push-to-lock knee lever wheel lock	-
MH05	<input type="checkbox"/>	OutFront scissor wheel lock	Add. charge
MH06	<input type="checkbox"/>	Knee lever wheel lock for user and attendant	Add. charge
Wheel lock accessories			
MH10	<input type="checkbox"/>	Wheel lock lever extension	Add. charge
Back support system			
MD01	<input type="checkbox"/>	Back support, straight	-
MD09	<input type="checkbox"/>	Back support with lumbar curve 8°	-
MD10	<input type="checkbox"/>	Back support, angle-adjustable 30° (back support height 43 or 48 cm)	Add. charge
MD11	<input type="checkbox"/>	Foldable	Add. charge
Push handles			
MD50	<input type="checkbox"/>	short	-
MD51	<input type="checkbox"/>	long	-
MD53	<input type="checkbox"/>	telescoping, short	Add. charge
MD52	<input type="checkbox"/>	telescoping, long	Add. charge
MD55	<input type="checkbox"/>	Height-adjustable/removable	Add. charge
MD56	<input type="checkbox"/>	Foldable	Add. charge

Back support system range of adjustment		
<input type="checkbox"/>	30 - 32,5 (back support straight with push handles height-adjustable/telescoping)	-
<input type="checkbox"/>	30 - 37,5 (back support straight with push handles short/long/folding)	-
<input type="checkbox"/>	32,5 - 40 (back support straight with push handles height-adjustable/telescoping)	-
<input type="checkbox"/>	35 - 40 (back support 8° with push handles height-adjustable/telescoping)	-
<input type="checkbox"/>	37,5 - 45 (back support straight with push handles height-adjustable/telescoping)	-
<input type="checkbox"/>	40 - 45 (back support straight with push handles short/long/folding)	-
<input type="checkbox"/>	42,5 - 50 (back support 8° or back support foldable with push handles long)	-
<input type="checkbox"/>	45 - 50 (back support 8°)	-
<input type="checkbox"/>	47,5 - 50 (back support straight with push handles short/long/folding or back support foldable with push handles long)	-

Back support upholstery		
MD04	<input type="checkbox"/> Standard (only for back support, angle-adjustable 30°)	-
MD03	<input type="checkbox"/> adaptable, padded	-
MD05	<input type="checkbox"/> adaptable, breathable	Add. charge
MD06	<input type="checkbox"/> without	-

Back support pad widening		
MD80	<input type="checkbox"/> + 2 cm	Add. charge
MD81	<input type="checkbox"/> + 4 cm	Add. charge
MD82	<input type="checkbox"/> without	-

Use of BAXX back support system (include order form, model-specific restrictions are possible)		
<input type="checkbox"/>	Baxx Aluminium back support (enclosed)	Add. charge

Back support accessories		
MD58	<input type="checkbox"/> Stabiliser bar	Add. charge

Side panels		
ME01	<input type="checkbox"/> Aluminium with protection against cold	-
ME11	<input type="checkbox"/> Plastic plug-on side panel	-
ME02	<input type="checkbox"/> Aluminium with clothing guard and against cold	-
ME04	<input type="checkbox"/> Carbon with clothing guard	Add. charge
ME17	<input type="checkbox"/> Aluminium with arm support long, depth-adjustable	Add. charge
ME05	<input type="checkbox"/> Aluminium with arm support short	Add. charge

Side panels accessories		
ME50	<input type="checkbox"/> Arm support padded <input type="checkbox"/> short <input type="checkbox"/> long <input type="checkbox"/> long, extra high	Add. charge
ME60	<input type="checkbox"/> Therapy tray	Add. charge


Seat upholstery		
MC13	<input type="checkbox"/> Standard	-
MC14	<input type="checkbox"/> readjustable	Add. charge
MC06	<input type="checkbox"/> without	-

Seat cushion			
MC01	<input type="checkbox"/>	Standard (seat cushion thickness 5 cm)	Add. charge
MC02	<input type="checkbox"/>	Cubic Foam (seat cushion thickness 5 cm)	Add. charge
MC05	<input type="checkbox"/>	Evolight (seat cushion thickness 6 cm)	Add. charge
MC29	<input type="checkbox"/>	Livo (seat cushion thickness 5 cm)	Add. charge
MC30	<input type="checkbox"/>	Evoflair (seat cushion thickness 8 cm)	Add. charge
Leg supports			
MB10	<input type="checkbox"/>	Individual, angle-adjustable, flip up	-
MB21	<input type="checkbox"/>	Amputation leg support, right	Add. charge
MB22	<input type="checkbox"/>	Amputation leg support, left	Add. charge
MB23	<input type="checkbox"/>	Amputation leg support, pair	Add. charge
MB26	<input type="checkbox"/>	Elevating leg support, right (plus MB10 left)	Add. charge
MB27	<input type="checkbox"/>	Elevating leg support, left (plus MB10 right)	Add. charge
MB25	<input type="checkbox"/>	Elevating leg support, pair	Add. charge
MB99	<input type="checkbox"/>	without	Add. charge
Length of foot support bar (without seat cushion)			
MB40	<input type="checkbox"/>	short (max. 34 cm – 41 cm)	-
MB41	<input type="checkbox"/>	medium (max. 41 cm – 51 cm)	-
MB42	<input type="checkbox"/>	long (max. 47 cm – 55 cm)	-
Other accessories			
MZ51	<input type="checkbox"/>	Lap belt with metal buckle	Add. charge
MZ52	<input type="checkbox"/>	Tool kit	Add. charge
Please also use the order form for wheelchair accessories 647F661=*			
Spoke protector			
MG90	<input type="checkbox"/>	Transparent	Add. charge
Frame colours			
002	<input type="checkbox"/>	Anthracite metallic	-
010	<input type="checkbox"/>	Marine blue	Add. charge
012	<input type="checkbox"/>	Sky blue	Add. charge
014	<input type="checkbox"/>	Candy red	Add. charge
016	<input type="checkbox"/>	Silver metallic	Add. charge
021	<input type="checkbox"/>	Sparkle granny-smith	Add. charge
023	<input type="checkbox"/>	White RAL 9016	Add. charge
048	<input type="checkbox"/>	Jet black matte (structure)	Add. charge
Side panel in frame colour*			
ME20	<input type="checkbox"/>	Yes	Add. charge
*: Standard colour jet black matte (texture); not possible in combination with side panels ME04 & ME11			
Special accessories			
MZ58	<input type="checkbox"/>	Bag for drive wheels (for 2 drive wheels)	Add. charge
MZ59	<input type="checkbox"/>	Pocket for mobile phone	Add. charge
MZ89	<input type="checkbox"/>	Protection for quick-release-axle	Add. charge
MZ92	<input type="checkbox"/>	Wheel Blades (mini-skis for caster wheels)	Add. charge
MZ93	<input type="checkbox"/>	Backpack	Add. charge

Information

Anterior seat height [mm]					
caster wheel fork	wheel diameter				
	4"	5"	5,5"	6"	7"
short	380 - 460	380 - 470	400 - 480	400 - 480	420 - 490
long	440 - 530	460 - 540	440 - 550	430 - 550	440 - 550

Positions drive wheel adapter



Information



Drive wheel adapter



MA51 Tip-assist



MA55 Anti-tipper



MA57 Crutch holder



MH01 Knee lever wheel lock



MH05 Outfront-scissor wheel lock



MH10 Knee lever wheel lock extension




MD10 Back support, angle-adjustable



MD11 Foldable back support



MD55 Push handles height-adjustable, removable




MD56 Push handles foldable

Information



MD04 Back support upholstery
Standard



MD03 Back support upholstery,
padded



MD05 Back support upholstery,
breathable



MD58 Stabiliser bar



ME11 Plug-on side panel plastic



ME02 Side panel with clothing
guard and against cold



ME04 Carbon side panel



ME17 Side panel with arm support
long, height-adjustable



ME50 Side panel, padded



ME60 Therapy tray



MB29 Livo



MB30 Evoflair



MB10 Leg support, individual, angle-
adjustable, swivel up



MB22 Amputation leg support, left



MB25 Leg supports, elevating

Information



MZ58 Bag for drive wheels



MZ59 Poket for mobile phone



MZ89 Protection for quick-release-axle



MZ92 Wheel Blades



MZ93 Backpack

Notes & Remarks