

Please contact Invacare EU Export Customer Services on :-
Tel : +33 2 47 62 69 80 ; Email : serviceclient_export@invacare.com
Or visit our website : www.invacare.eu.com



Küschall Champion

Active manual scissor folding chair



Description :

- Aluminum frame with hydroformed tubes
- Personalization kit with premium decals
- Angle-adjustable backrest and folding seat with passive seat lock as standard
- Useful seat width up to 500 mm.
- Maximum user weight: 120 kg
- Full weight: from 9.2 kg.
- Transport weight: from 6.7 kg.
- Warranty: 2 years excluding wearing parts

New Champion SL :
Lighter thanks to its welded sides (UHA1260)

Discover the new possibilities of the special fabrication **UNIQUE**. Push back the limits of the personalization with a wide range of possibilities. Create **UNIQUE** wheelchairs.

Discover our new tool to help calculate overall dimensions and find the visual configurator to make your configuration!

Date :

Purchase Order :

Customer Code :

Distributor Stamp:

We reserve the right to deliver a + or - tolerance of +/- 10 mm on all measurements We reserve the right to make technical changes without notice

☐ Please make your selection with a cross in the box

◆ Standard feature

○ Free of charge option

CODE		DESIGNATION	NC : Non Compatible NC : 9999 = not compatible with the reference 9999	Euros	
DHA0031	☐	küschall® Champion®		3982,30	
SPECIAL CONSTRUCTION					
DHA0000	☐	MyKüschall special construction		On Request	
UHA0001	☐	UNIQUE - Personalization		On Request	
POSITIONING PACKAGE					
DHA0220	☐	Posture And Pressure Relief Package Libra/MX2 Include : Matrx Libra Cushion (DHA2720) Matrx MX2 Backrest (DHA2300 To DHA2315) Bodypoint Evoflex Belt (Dha7810 To DHA7830)	NC : SW 360 mm only possible with SD 400/420 mm NC 1350, 1405, 1415, 1425, 7230	1222,23	Save 15% by selecting this PACK, i.e. €277,55
DHA0230	☐	Posture And Pressure Relief Package Libra/Elite E2 Include : Matrx Libra Cushion (DHA2720) Matrx Elite/Elite Deep E2 Backrest (DHA2320-Dha2350) Bodypoint Evoflex Belt (DHA7810 To DHA7830)	NC : SW 360 mm only possible with SD 400/420 mm NC 1350, 1405, 1415, 1425, 7230	1105,59	Save 10% by selecting this PACK, i.e. €106,70
SEATING CHOICE					
FRAME					
DHA1150	☐	Frame Angle 75° Aluminium	NC : 1350	◆	
DHA1160	☐	Frame Angle 90° Aluminium	NC : 1350	○	
DHA1170	☐	Frame Angle 75° Aluminium - Tapered Front Frame	NC : 1250, 1320, 1325, 7010	○	
DHA1180	☐	Frame Angle 90° Aluminium - Tapered Front Frame	NC : 1250, 1320, 1325, 7010	○	
FRAME OPTION					
DHA1250	☐	Champion SK Option : Folding Front Frame - Only With Non Tapered Frame	NC : 1350, 1924 to 1934, 4190, 8700	618,01	
UNIQUE - CHAMPION SL FRAME					
UHA1260	☐	CHAMPION SL : COMPLETE FRAME WITH FRONT & REAR WHEEL SUPPORT WELDED	SELECT THE OPTION UHA0001 TO ACCESS THIS OPTION NC: 3490, 6180	597,16	NEW

SEAT WIDTH (SW)						
DHA1320	☐	360 mm	NC : 2364, 2366, 2368, 2370, 2372, 2710, 2742, 2744, 2746, 2748, 7000			○
DHA1325	☐	380 mm	NC : 2366, 2368, 2370, 2372, 2742, 2744, 7000			○
DHA1330	☐	400 mm	NC : 2360, 2366, 2368, 2370, 2372, 2744			○
DHA1335	☐	420 mm	NC : 2360, 2362, 2368, 2370, 2372			○
DHA1340	☐	440 mm	NC : 2360, 2362, 2364, 2370, 2372 2740, 2742			○
DHA1345	☐	460 mm	NC : 2360, 2362, 2364, 2370, 2372 2740, 2742, 2744			○
DHA1350	☐	480 mm (only with tapered frame)	NC : 2300, 2305, 2310, 2315, 2345, 2360, 2362, 2364, 2366, 2740, 2742, 2744			○
SEAT DEPTH (SD)						
DHA1405	☐	340 mm	NC : 2710, 2720, 2740, 2742, 2744, 2746, 2748, 2860, 7810, 7820, 7830			○
DHA1415	☐	360 mm	NC : 2710, 2720, 2742, 2744, 2746, 2748, 2860, 7810, 7820, 7830			○
DHA1425	☐	380 mm	NC : 2720, 2744			○
DHA1435	☐	400 mm	NC : 2744			◆
DHA1445	☐	420 mm				○
DHA1460	☐	440 mm	NC : 2740, 2742			○
DHA1470	☐	460 mm	NC : 2740, 2742, 2744			○
DHA1480	☐	480 mm	NC : 2710, 2740, 2742, 2744, 2746, 2748			○
BACK HEIGHT (BH) - FOR MATRX BACKRESTS: ONLY DEFINES THE PUSHHANDLES HEIGHT						
DHA1615	☐	300 mm	NC : 2200, 2322, 2325, 2330, 2331, 2345, 2350, 2355, 2357, 2400			○
DHA1620	☐	315 mm	NC : 2200, 2325, 2330, 2331, 2350, 2355, 2357, 2400			○
DHA1625	☐	330 mm	NC : 2325, 2330, 2331, 2350, 2355, 2357			○
DHA1630	☐	345 mm	NC : 2330, 2331, 2355, 2357			○
DHA1635	☐	360 mm				○
DHA1640	☐	375 mm				◆
DHA1645	☐	390 mm				○
DHA1650	☐	405 mm				○
DHA1655	☐	420 mm	DHA1660	☐	435 mm	○
DHA1665	☐	450 mm	DHA1670	☐	465 mm	○
FRONT SEAT-TO-FLOOR-HEIGHT (FSTF) - see table 1 and 2						
DHA1740	☐	450 mm	NC : 1946 to 1960			○
DHA1745	☐	460 mm	NC : 1948 to 1960			○
DHA1750	☐	470 mm	NC : 1950 to 1960			○
DHA1755	☐	480 mm	NC : 1952 to 1960			○
DHA1760	☐	490 mm	NC : 1954 to 1960			○
DHA1765	☐	500 mm	NC : 1956 to 1960			◆
DHA1770	☐	510 mm	NC : 1958, 1960			○
DHA1775	☐	520 mm	NC : 1960			○
DHA1780	☐	530 mm	NC : roues de 3"			○
DHA1785	☐	540 mm	NC : roues de 3" et 4"			○
REAR SEAT-TO-FLOOR HEIGHT (RSTF)						
- SELECT 75° FRAME IF 0 <FSTF - RSTF < 70MM; SELECT 90° FRAME IF 30 <FSTF - RSTF < 100MM						
DHA1815	☐	390 mm	DHA1820	☐	400 mm	NC : 6420, 6430
DHA1825	☐	410 mm	DHA1830	☐	420 mm	NC : 6430
DHA1835	☐	430 mm				◆
DHA1840	☐	440 mm				○
DHA1845	☐	450 mm	DHA1850	☐	460 mm	○
DHA1855	☐	470 mm	DHA1860	☐	480 mm	NC : 6410

For any chair renewal, we advise you to completely redo the measurement. A change in chair generation involves changes in dimensions as well as changes in the positioning of the user.

RSTF / Rear Wheel Size

24" : 390 - 460 mm

25" : 410 - 480 mm

26" : 430 - 480 mm




For any chair renewal, we advise you to completely redo the measurement. A change in chair generation involves changes in dimensions as well as changes in the positioning of the user.



RSTF / Rear Wheel Size
24" : 390 - 460 mm
25" : 410 - 480 mm
26" : 430 - 480 mm

LOWER LEG LENGTH (LLL)							
- SELECT 75° FRAME IF LLL ≤ FSTF- 30 MM ; SELECT 90° FRAME IF LLL ≤ FSTF- 40 MM							
DHA1924	☐	320 mm	DHA1926	☐	330 mm	NC : 1250	☐
DHA1928	☐	340 mm	DHA1930	☐	350 mm	NC : 1250	☐
DHA1932	☐	360 mm	DHA1934	☐	370 mm	NC : 1250	☐
DHA1936	☐	380 mm				WITH SK OPTION (DHA1250) ONLY IN COMBINATION WITH SD ≥ 400MM	☐
DHA1938	☐	390 mm				WITH SK OPTION (DHA1250) ONLY IN COMBINATION WITH SD ≥ 380MM	☐
DHA1940	☐	400 mm	DHA1942	☐	410 mm		☐
DHA1944	☐	420 mm					☐
DHA1946	☐	430 mm					☐
DHA1948	☐	440 mm	DHA1950	☐	450 mm		☐
DHA1952	☐	460 mm	DHA1954	☐	470 mm		☐
DHA1956	☐	480 mm	DHA1958	☐	490 mm		☐
DHA1960	☐	500 mm					☐
BACKREST							
DHA2140	☐	Angle Adjustable, Forward Folding Backrest With Soft Locking - Backrest Remains In Position During Transfer					☐
BACK UPHOLSTERY							
DHA2200	☐	Standard Light Upholstered				NC : BH < 330 mm, 2300 to 2357, 2430	☐
DHA2220	☐	Tension Adjustable Premium Black Back Upholstery With Loop Bands				NC : 2300 to 2357	77,73
DHA2230	☐	Aero Back Tension Adjustable				NC : 2300 to 2357	258,77
DHA2240	☐	Without Backrest Upholstery - Mandatory If Ordering Positioning Backrests Like Matrx				NC : 2430	☐
UNIQUE - BACK UPHOLSTERY STITCHING COLOUR - Only for Aero Back upholstery							
UHA2241	☐	NEON YELLOW					101,42
UHA2242	☐	NEON PINK					101,42
UHA2243	☐	BLUE					101,42
UHA2244	☐	WHITE					101,42
UHA2245	☐	RED					101,42
UHA0002	☐	"Your personnal color": Please contact our customers service					On Request
UNIQUE - BACK UPHOLSTERY - EMBROIDERY - Only for Aero Back upholstery							
UHA2246	☐	SMILEY					218,01
UHA2248	☐	TRAVELER					218,01
UHA2249	☐	SAKURA BEAUTY					218,01
UHA2250	☐	FANNY CACTUS					218,01
UHA2251	☐	KUSCHALL LOGO: Please contact our customers service					On Request
UHA0003	☐	"Your personnal design": Please contact our customers service - Provide the wished design in .ai / .pdf or .png format					On Request
UNIQUE - BACK UPHOLSTERY - FLOCKING - Only for Aero Back upholstery							
UHA2253	☐	SMILEY					118,48
UHA2255	☐	TRAVELER					118,48
UHA2256	☐	SAKURA BEAUTY					118,48
UHA2257	☐	FANNY CACTUS					118,48
UHA2259	☐	PAINT STAIN					118,48
UHA0004	☐	"Your personnal design": Please contact our customers service - Provide the wished design in .ai / .pdf or .png format					On Request
MATRX BACKRESTS							
DHA2300	☐	Matrx MX2, Height 150 mm With Quick Release Hardware - User max. weight 114 kg				NC : 2370, 2372, 3480	927,01
DHA2305	☐	Matrx MX2, Height 230 mm With Quick Release Hardware - User max. weight 114 kg				NC : 2370, 2372, 3480	927,01
DHA2310	☐	Matrx MX2, Height 305 mm With Quick Release Hardware - User max. weight 114 kg				NC : 2370, 2372, 3480	927,01
DHA2315	☐	Matrx MX2, Height 405 mm With Quick Release Hardware - User max. weight 114 kg				NC : 2370, 2372, 3480	927,01
DHA2316	☐	Elite Deep E2 , Height 160 mm With Mini-Set Mounting Hardware				NC : 3480	717,54
DHA2317	☐	Elite Deep E2 , Height 210 mm With Mini-Set Mounting Hardware				NC : 3480	717,54
DHA2318	☐	Elite Deep E2 , Height 260 mm With Mini-Set Mounting Hardware				NC : 3480	717,54
DHA2320	☐	Elite Deep E2 , Height 310 mm With Mini-Set Mounting Hardware				NC : 3480	717,54

User's LLL must be measured taking away the cushion thickness





LLL = LLL patient - Seat cushion



2300, 2305, 2310, 2315

MX2: only 140g more than a tension-adjustable backrest!

DHA2322	☐	Elite Deep E2, Height 360 mm With Mini-Set Mounting Hardware	NC : 3480	717,54	 2316 to 2331  2332 to 2357	
DHA2325	☐	Elite Deep E2, Height 410 mm With Mini-Set Mounting Hardware	NC : 3480	717,54		
DHA2330	☐	Elite Deep E2, Height 460 mm With Easy-Set Mounting Hardware	NC : 3480	717,54		
DHA2331	☐	Elite Deep E2, Height 510 mm With Easy-Set Mounting Hardware	NC : 3480	717,54		
DHA2332	☐	Elite E2, Height 160 mm With Mini-Set Mounting Hardware	NC : 3480	717,54		
DHA2334	☐	Elite E2, Height 210 mm With Mini-Set Mounting Hardware	NC : 3480	717,54		
DHA2335	☐	Elite E2, Height 260 mm With Mini-Set Mounting Hardware	NC : 3480	717,54		
DHA2340	☐	Elite E2, Height 310 mm With Mini-Set Mounting Hardware	NC : 3480	717,54		
DHA2345	☐	Elite E2, Height 360 mm With Mini-Set Mounting Hardware	NC : 3480	717,54		
DHA2350	☐	Elite E2, Height 410 mm With Mini-Set Mounting Hardware	NC : 3480	717,54		
DHA2355	☐	Elite E2, Height 460 mm With Easy-Set Mounting Hardware	NC : 3480	717,54		
DHA2357	☐	Elite E2, Height 510 mm With Easy-Set Mounting Hardware	NC : 3480	717,54		
MATRX BACKRESTS WIDTH CHOICE					 2430  2450  Backrest Angle BA	
DHA2360	☐	360 mm - Indicative Bust Width 265 mm		☐		
DHA2362	☐	380 mm - Indicative Bust Width 290 mm		☐		
DHA2364	☐	410 mm - Indicative Bust Width 315 mm		☐		
DHA2366	☐	430 mm - Indicative Bust Width 340 mm		☐		
DHA2368	☐	460 mm - Indicative Bust Width 370 mm		☐		
DHA2370	☐	480 mm - Indicative Bust Width 390 mm		☐		
DHA2372	☐	510 mm - Indicative Bust Width 420 mm		☐		
PUSH HANDLES					 2600	
DHA2400	☐	Without Push Handles	NC : BH < 330 mm	☐		
DHA2410	☐	Push Handles Fixed, Mini		☐		
DHA2420	☐	Push Handles Fixed, Long		☒		
DHA2430	☐	Height Adjustable, Rear Set	NC : 3480	212,32		
DHA2440	☐	Height Adjustable, Integrated		254,03		
DHA2450	☐	Foldable Push Handles		213,27		
BACKREST ANGLE (BA)					 2600	
DHA2510	☐	76°	DHA2570 ☐	89,5°		☐
DHA2530	☐	80,5°	DHA2590 ☐	94°		☐
DHA2550	☐	85°				☒
SEAT UPHOLSTERY					 2600	
DHA2600	☐	STANDARD - CENTRAL LOOP BAND		☒		
UNIQUE - SEAT UPHOLSTERY STITCHING COLOUR						
UHA2602	☐	NEON YELLOW		118,48		
UHA2603	☐	NEON PINK		118,48		
UHA2604	☐	BLUE		118,48		
UHA2605	☐	WHITE		118,48		
UHA2606	☐	RED		118,48		
UHA0005	☐	"Your personnal color": Please contact our customers service		On Request		
CUSHION & CUSHION OPTIONS					 2700  2710	
Cushions						
DHA2700	☐	STANDARD SEAT CUSHION, 60 MM DOUBLE FOAM WITH VISCO TOP	NC : 2810, 2820, 2830, 2840, 2850	101,42		
DHA2710	☐	FLO-TECH LITE	NC : 2810, 2820, 2830, 2840, 2850	118,48		
DHA2720	☐	MATRX LIBRA CUSHION, BLACK SMOOTH COVER - NOT POSSIBLE WITH SD < 400 MM ; SW 360 MM ONLY POSSIBLE WITH SD 400/420 MM	NC: Depth < to 400 mm, width 360 mm possible only with depth 400 and 420 mm	293,84		

DHA2740	☐	CLASSIC VISCO CUSHION 375 X 400 MM - ONLY FOR USERS WEIGHTING UP TO 50 KG	NC : 2810, 2820, 2830, 2840, 2850	76,78	 2720  2740 / 2742 2744 / 2746  2748  2800
DHA2742	☐	CLASSIC VISCO CUSHION 405 X 400 MM - ONLY FOR USERS WEIGHTING UP TO 100 KG	NC : 2810, 2820, 2830, 2840, 2850	65,40	
DHA2744	☐	CLASSIC VISCO CUSHION 430 X 430 MM - ONLY FOR USERS WEIGHTING UP TO 100 KG	NC : 2810, 2820, 2830, 2840, 2850	65,40	
DHA2746	☐	CLASSIC VISCO CUSHION 455 X 400 MM - ONLY FOR USERS WEIGHTING UP TO 100 KG	NC : 2810, 2820, 2830, 2840, 2850	76,78	
DHA2748	☐	CONTOUR VISCO NG CUSHION	NC : 2810, 2820, 2830, 2840, 2850	76,78	
DHA2750	☐	VICAIR ACTIVE O2 CUSHION		683,73	
SEAT OPTION					
DHA2800	☐	Seat Bag, Click In And Removable, Includes The Fixation Clips (7010)		81,52	 2810  2820  2830  2840  2850  2860
MATRX LIBRA CUSHION ACCESSORIES					
DHA2810	☐	Fluid Supplement Sac (Pair)		91,00	
DHA2820	☐	Tapered Wedge Pair		41,71	
DHA2830	☐	Abductor		47,39	
DHA2840	☐	Lateral Side Wedge		41,71	
DHA2850	☐	Anterior Front Wedge		54,98	
SEAT RIGIDIZER					
DHA2860	☐	Seat Rigidizer	(Mandatory: choose a seat cushion)	56,87	
FOOTREST					
Footplate / SW : Straight frame: High mounted footplate not possible for SW 360 mm Tapered frame: High mounted footplate not possible for SW 360 mm - 400 mm					 3240
DHA3240	☐	Flip-Up Footplate, Angle Adjustable, Set To The Front		◆	
CALF STRAP					
DHA3320	☐	MOUNTED ON SUPPORTER		○	 3440
DHA3330	☐	MOUNTED ON FRAME		◆	
SIDE GUARDS					
SIDE GUARDS / CLOTHES GUARDS					 3440
DHA3405	☐	Economy Clothes Guards, Flip Back	NC: 6430	○	
DHA3420	☐	Carbon Clothes Guards, Flip Back		◆	
DHA3430	☐	Carbon Mudguards, Fixed	NC : 3490	223,70	
DHA3440	☐	Carbon Mudguards, Removable	NC : 3490	336,49	
UNIQUE - CUSTOM TRANSPARENT COLORS - Only for economy clothes guards					
UHA3422	☐	TRANSPARENT		520,38	 3480 3490
UHA3423	☐	TRANSPARENT BLUE		745,02	
UHA3424	☐	NEON YELLOW		745,02	
UHA3425	☐	NEON ORANGE		745,02	
UHA3426	☐	NEON RED		745,02	
UHA0006	☐	"Your personnal color": Please contact our customers service		On Request	
UHA0007	☐	"Your personnal design": Please contact our customers service - Provide the wished design in .ai / .pdf or .eps format		On Request	
ARMREST					
DHA3480	☐	Tubular Armrest, Non-Locking, Height Adjustable, Removable, Swing-Away, Comfort Padded		264,45	
DHA3490	☐	Siderest, Height Adjustable, Short Armpad (T-Arm), Non-Locking		253,08	
COLOUR OFFERING & CUSTOMIZATION					
COLOUR OFFERING & CUSTOMIZATION					EXTRA LEAD TIME +15 DAYS
DHA4000	☐	SPACE BLACK (RAL 9005)		○	
DHA4010	☐	ANTHRACITE GREY (RAL 7016)		○	
DHA4020	☐	PURE WHITE (RAL 9016)	NC : 4220	○	
DHA4030	☐	CARMINE RED (RAL 3002)		○	
DHA4040	☐	AZURE BLUE (RAL 5017)		○	

DHA4060	☐	FLASHY YELLOW (RAL 1016)		216,11	EXTRA LEAD TIME +15 DAYS
DHA4100	☐	BLACK GLOSSY (RAL 9005)		216,11	
DHA4130	☐	SILVER GREY (RAL 9006)		216,11	
DHA4145	☐	PEARL NIGHT BLUE	NEW	216,11	
DHA4150	☐	RUBY RED (RAL 3003)		216,11	
DHA4170	☐	DYNAMITE ORANGE - specific color mix, no RAL definition		216,11	
DHA4180	☐	FOREST GREEN (RAL 6029)		216,11	
UHA4101	☐	DIAMOND SILVER		696,68	
UHA4102	☐	PASTEL GREEN		750,71	
UHA4103	☐	SKY BLUE		750,71	
UHA4104	☐	PLUM PURPLE		750,71	
UHA4105	☐	SALMON		750,71	
DHA4070	☐	ATOMIC VIOLET - Specific Color Mix, No Ral Definition		On quote	
UHA0008	☐	SPECIAL FRAME COLOUR, ON DEMAND ONLY -Manufacturing time + 10 days: for a color already offered on the other ranges of Invacare manual wheelchairs. WHERE - Manufacturing time + 15 days: for any other request (subject to feasibility) - specify RAL number		750,71	


TAG

CERAMIC GLAZE FINISH				
UHA4066	☐	CERAMIC GLAZE FINISH (ONLY FOR UHA**** COLOURS)		643,60

STYLING				
DHA4190	☐	Bold	NC : 4300, 4310	64,45
DHA4195	☐	Tag	NC : 4200 to 4250	◆


BOLD

BOLD COLOUR				
DHA4200	☐	BLACK		○
DHA4210	☐	DARK GREY		○
DHA4220	☐	LIGHT GREY		○
DHA4230	☐	RED		○
DHA4240	☐	BLUE		○
DHA4250	☐	GOLD		○
UHA4252	☐	NEON YELLOW		161,14
UHA4253	☐	NEON GREEN		161,14
UHA4254	☐	NEON ORANGE		161,14
UHA4255	☐	NEON PINK		161,14
UHA4256	☐	NEON RED		161,14
UHA0009	☐	"Your personnal color": Please contact customers service		On Request

TAG COLOUR				
DHA4300	☐	BLACK		○
DHA4310	☐	RED		○

FRONT WHEEL				
DHA5000	☐	3" SUPERSPORT		○
DHA5050	☐	4" EVERYDAY WHEEL WITH BLACK TYRE		○
DHA5100	☐	4" COMFORT WHEEL WITH BLACK TYRE		○
DHA5150	☐	4" STAREC BLACK HUB, BLACK TYRE		75,83
DHA5150	☐	4" FROGLEGS FRONT WHEEL WITH BLACK HUB, BLACK TYRE	Select the Froglegs front wheels with black hub (DHA5150) to access UNIQUE anodizations	121,33
DHA5250	☐	4" FROGLEGS FRONT WHEEL WITH SILVER HUB, BLACK TYRE		121,33
DHA5300	☐	5" EVERYDAY WITH BLACK TYRE		○
DHA5350	☐	5" COMFORT WHEEL WITH BLACK TYRE		◆
DHA5400	☐	5" STAREC BLACK HUB, BLACK TYRE		103,32
DHA5450	☐	5" FROGLEGS FRONT WHEEL WITH BLACK HUB, BLACK TYRE	Select the Froglegs front wheels with black hub (DHA5450) to access UNIQUE anodizations	147,87
DHA5500	☐	5" FROGLEGS FRONT WHEEL WITH SILVER HUB, BLACK TYRE		147,87

 5000

5200 / 5250 / 5450 / 5500
 5350

UNIQUE - FROGLEGS HUB COLOURS				
UHA5502	☐	PINK		98,58
UHA5504	☐	ORANGE		98,58
UHA5506	☐	BLUE		98,58
UHA5507	☐	VIOLET		98,58
UHA5508	☐	RED		98,58
UHA5509	☐	GOLD		98,58


5050 / 5300

UNIQUE - CERAMIC BEARING				
UHA5510	☐	CERAMIC BEARING		98,58



TAG

TAG



BOLD



5000



5200 / 5250 / 5450 / 5500




5350



5050 / 5300

REAR WHEEL POSITION			
DHA6100	☐	1ST FROM FRONT (THE MOST ACTIVE)	○
DHA6110	☐	2ND FROM FRONT (POSITION 1 ON THE CHAMPION 2008)	○
DHA6120	☐	3RD FROM FRONT	○
DHA6130	☐	4TH FROM FRONT (POSITION 2 ON THE CHAMPION 2008)	○
DHA6140	☐	5TH FROM FRONT (POSITION 3 ON THE CHAMPION 2008)	○
DHA6150	☐	6TH FROM FRONT	○
DHA6160	☐	7TH FROM FRONT (POSITION 2 ON THE CHAMPION 2008)	○
DHA6170	☐	8TH FROM FRONT	○
DHA6180	☐	9TH FROM FRONT (THE MOST STABLE)	○
QUAD QUICK RELEASE AXLE			
DHA6300	☐	TETRA-CLIP	NC : 6550, 6950 143,13
REAR WHEEL SIZE			
DHA6410	☐	24"	◆
DHA6420	☐	25"	○
DHA6430	☐	26"	517,54
REAR WHEEL TYPE			
DHA6500	☐	24" ACTIVE WHEELS (BLACK RIM, HUB AND RADIAL SPOKING INCL. ALU HANDRIMS)	24" NC : All Handrims, 6900, 6910 ◆
DHA6510	☐	KUSCHALL STAREC HUB (BLACK RIM AND HUB, BLACK CROSS SPOKING)	24", 26" ○
DHA6520	☐	HIGH PERFORMANCE (BLACK RIM, BLACK STAREC HUB, BLACK CROSS SPOKING)	24", 25" NC : 6660 532,70
DHA6550	☐	SPINERGY LXX (BLACK RIM, HUB AND RADIAL SPOKING)	24", 25" NC : 6660, 6910 771,56
UHA6551	☐	SPINERGY LITE EXTREME LX - 680g per wheel	24", 25" 1060,66
UHA6552	☐	SPINERGY LITE EXTREME X-LACED XLX - 660g per wheel	24", 25" 1140,28
UHA6553	☐	SPINERGY X-LACED SPORT LITE EXTREME XSLX - 757g per wheel	24", 25" 1206,64
UHA6554	☐	SPINERGY CARBON FIBER CLX - 570g per wheel	24", 25" 1358,29
UNIQUE - RIM ANODIZED COLOURS			
UHA6557	☐	BLUE	162,09
UHA6558	☐	GOLD	162,09
UHA6559	☐	GREEN	162,09
UHA6560	☐	ORANGE	162,09
UHA6561	☐	RED	162,09
UHA6562	☐	VIOLET	162,09
UNIQUE - HUB ANODIZED COLOURS			
UHA6563	☐	BLUE	162,09
UHA6564	☐	GOLD	162,09
UHA6565	☐	SILVER GREY (RAL 9006)	162,09
UHA6566	☐	GREEN	162,09
UHA6567	☐	ORANGE	162,09
UHA6568	☐	RED	162,09
UHA6569	☐	VIOLET	162,09
UNIQUE - SPOKE COLOURS			
UHA6570	☐	BLUE	162,09
UHA6571	☐	GREEN	162,09
UHA6572	☐	ORANGE	162,09
UHA6573	☐	PINK	162,09
UHA6574	☐	RED	162,09
UHA6575	☐	VIOLET	162,09
UHA6576	☐	WHITE	162,09
UHA6577	☐	YELLOW	162,09
UNIQUE - CERAMIC BEARING			
UHA6579	☐	CERAMIC BEARING	332,70



6300

REAR WHEEL TYRE				
DHA6600	☐	MAXXIS LITE GUARD BLACK TYRE (puncture protection)	24"	◆
DHA6625	☐	MAXXIS ALPHA PROTECT BLACK TYRE (puncture protection)	24"	72,04
DHA6630	☐	MBL SPEEDLITE (puncture protection)	24"	72,04
DHA6660	☐	SOLID TYRE BLACK	24" NC : 7230	72,04
HANDRIM				
DHA6700	☐	ANTI-SLIP DIST. 3 CM ADJUSTABLE 4CM	24", 26"	142,18
DHA6710	☐	STEEL, DIST. 3 CM ADJUSTABLE 4CM	24", 26"	113,74
DHA6720	☐	TITANIUM, DIST. 3 CM ADJUSTABLE 4CM	24"	277,73
DHA6740	☐	CARBOLIFE GEKKO, WITH ANTI-SLIP DIST. 3 CM ADJUSTABLE 4CM	24", 25"	373,46
DHA6750	☐	CARBOLIFE GEKKO S, WITH ANTI-SLIP DIST. 3 CM ADJUSTABLE 4CM	24", 25"	373,46
DHA6760	☐	CARBOLIFE CURVE L WITH ERGONOMIC PROFILE DIST. 3 CM ADJUSTABLE 4CM	24"	282,46
DHA6770	☐	CARBOLIFE QUADRO, WITH ERGONOMIC PROFILE DIST. 3 CM ADJUSTABLE 4CM	24"	373,46
DHA6780	☐	LIGHT ALU, GREY, DIST. 3 CM ADJUSTABLE 4CM	24", 25", 26"	◆
DHA6790	☐	LIGHT ALU, BLACK, DIST. 3 CM ADJUSTABLE 4CM	24", 25", 26"	○
REAR WHEEL CAMBER				
DHA6820	☐	0°		◆
DHA6830	☐	1°		○
DHA6840	☐	3°	NC : 7230	○
SPOKE GUARDS				
DHA6900	☐	Transparent	24", 26"	200,95
DHA6910	☐	Black	24", 26"	200,95
ALBER ADDS ON KIT				
DHA6950	☐	E-MOTION & TWION ADAPTATION PLATE 2 ALBER ANTITIPPERS INCLUDED	Only with 24" wheels NC : 7500, 7510, 7520, 8200	579,15
<i>Tip : for perfect compatibility of the brakes with the manual wheels, we advise you to equip them with the same coating as the motorized wheels.</i>				
FIXATION CLIPS				
DHA7000	☐	FIXATION STRAP		○
DHA7010	☐	FIXATION CLIP - HOLDS FRAME TOGETHER IN FOLDED POSITION		◆
UNFOLDING HELP				
DHA7100	☐	UNFOLDING HELP - TO OPEN FOLDING MECHANISM FROM FRONT		◆
UNIQUE - CUSTOMIZATION				
UHA6960	☐	UNIQUE - CUSTOM COLOURS FORKS AND BEARING COVER		197,16
UHA6961	☐	UNIQUE - CUSTOM COLOURS FORKS, BEARING COVER, REAR WHEEL SUPPORTS, LATERAL SUPPORTS, BRAKE CLAMPS		551,66



6700 6710



6720



6740



6790 6760




7010

UNIQUE - COLOURS PREMIUM ANODIZATION				
UHA6963	☐	RED		○
UHA6964	☐	VIOLET		○
UHA6965	☐	COPPER		○
UHA6966	☐	GREEN		○
UHA6967	☐	ORANGE		○
UHA6968	☐	GOLD		○
UHA6969	☐	CHAMPAGNE		○
UHA6970	☐	PINK		○
UHA6971	☐	BLUE		○
UHA6972	☐	POLY MIRROR FINISH		○
BRAKES				
DHA7200	☐	PULL BRAKES		○
DHA7210	☐	STANDARD BRAKES, LEVER BENT		◆
DHA7220	☐	STANDARD BRAKES, LEVER STRAIGHT		○
DHA7230	☐	PERFORMANCE BRAKES	NC : 7300, 7810, 7820, 7830	76,78
DHA7250	☐	LIGHT ACTIVE BRAKES		145,02
BRAKE EXTENSION LEVER				
DHA7300	☐	FOR STRAIGHT LEVER BRAKES		76,78
SAFETY OPTIONS				
Assembly of a maximum of 2 elements per side from the following list: anti-tippers, tipping tubes, cane holder or transit wheels				
ANTI-TIPPER				
DHA7500	☐	ANTI-TIPPER RIGHT	NC : 8300	131,75
DHA7510	☐	ANTI-TIPPER LEFT	NC : 8300	131,75
DHA7520	☐	ANTI-TIPPER BOTH SIDES		263,51
BELTS				
DHA7800	☐	STANDARD BELT - PLASTIC BUCKLE		42,65
DHA7810	☐	BODYPOINT POSTURE BELT EVOFLEX 2 POINTS, SMALL, FLAT-MOUNT WITH PUSH BUTTON		217,06
DHA7820	☐	BODYPOINT POSTURE BELT EVOFLEX 2 POINTS, MEDIUM, FLAT-MOUNT WITH REHAB LATCH		217,06
DHA7830	☐	BODYPOINT POSTURE BELT EVOFLEX, 2 POINTS, LARGE, FLAT-MOUNT WITH REHAB LATCH		217,06
PASSIVE LIGHTS				
DHA7900	☐	PASSIVE LIGHTS ASSEMBLED (MOUNTED ON SPOKES)		22,75
ATTENDANT OPTIONS				
FRAME PROTECTION				
DHA8010	☐	FRONT FRAME PROTECTORS; 2 PC. - REMOVABLE, WITH LOOP BANDS		25,59
TIPPER AID				
DHA8100	☐	TIPPER-AID, RIGHT		66,35
CANE HOLDER Max. 2 adaptations per side (adaptations = anti tipper, cane holder, transit wheels)				
DHA8200	☐	CANE HOLDER, MOUNTED RIGHT		66,35
CANE HOLDER, MOUNTED RIGHT				
DHA8300	☐	TRANSIT WHEELS 3" ONLY IN COMBINATION WITH ANTITIPPER BOTH SIDES (DHA7520)		301,45
OTHERS OPTIONS				
DHA8500	☐	AIR PUMP - NEW		33,18
DHA8550	☐	Trousse to outils		22,75
DHA8680	☐	24" MOUNTAIN BIKE WHEELS WITH BAG - MOUNTAIN BIKE WHEELS AND TIRES WITH ALU GREY HANDRIMS ; QUICK RELEASE AXLES INCLUDED		646,45
DHA8700	☐	TRAVEL BAG FOR CHAMPION		432,23
DHA8650	☐	WHEELS BAG	(22", 24", 25", 26")	65,40
DHA8750	☐	TRAVEL BAG FOR CHAMPION WITH SK OPTION (FOLDABLE FRONT FRAME)		499,53
DHA8800	☐	BODYPOINT CALF STRAP		104,27
DHA8850	☐	BODYPOINT MOBILITY BAG		104,27


7210 7230


7250






7810 / 7820 / 7830


8550


8680


8700

FRAME TYPE AND ANGLE, DIFFERENCE BETWEEN FSTF AND RSTF, FOOTPLATE MOUNTING, SEAT DEPTH AND FRONT WHEEL SIZE													FRAME TYPE AND ANGLE, FSTF, FRAME LENGTH, SEAT DEPTH AND FRONT WHEEL SIZE													
		FOOTPLATE SET TO THE FRONT						FOOTPLATE SET TO THE REAR						Long frame						Short frame						
		SD340 - SD380			SD ≥ 400			SD340 - SD380			SD ≥ 400			SD340 - ST380			SD ≥ 400			SD340 - ST380			SD ≥ 400			
		3"	4"	5"	3"	4"	5"	3"	4"	5"	3"	4"	5"	3"	4"	5"	3"	4"	5"	3"	4"	5"	3"	4"	5"	
75° Frame	Difference between FSTF and RSTF)	3"	4"	5"	3"	4"	5"	3"	4"	5"	3"	4"	5"	FSTF	3"	4"	5"	3"	4"	5"	3"	4"	5"	3"	4"	5"
	0 mm								x	x			x	x												
	10 mm									x																
	20 mm									x																
	30 mm									x																
	40 mm																									
	50 mm																									
	60 mm																									
	70 mm																									
90° Frame	30 mm			x				x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	
	40 mm			x				x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x		
	50 mm						x		x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x		
	60 mm								x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x		
	70 mm									x	x	x	x	x	x	x	x	x	x	x	x	x	x	x		
	80 mm										x	x	x	x	x	x	x	x	x	x	x	x	x	x		
	90 mm										x	x	x	x	x	x	x	x	x	x	x	x	x	x		
	100 mm																							x		
	There is no conflict if the lower leg length (LLL) is smaller than "FSTF minus 100 mm"																									
MATRIX BACKRESTS WIDTH AVAILABLE RANGE																										
SEAT WIDTH		Matrix Elite E2 / Elite E2 Deep backrest width choice (mm)							SEAT WIDTH		Matrix MX2 backrest width choice (mm)															
		360	380	410	430	460	480	510			360	380	410	430	460	480	510									
360																										
380																										
400																										
420																										
440																										
460																										
480																										
Ramps																										
1440530	☐	Fiberglass Folding Mono Ramp Folded dimensions: 2000 x 230 mm 16.3 kg - maximum capacity: 300 kg											1 787,68													
1654128	☐	Basic ramp 2 telescopic parts in aluminum and non-slip coating Folded dimensions: 1210 x 194 mm 8.8 kg per pair - maximum capacity: 270 kg											402,84													
1654129	☐	Quick2go Compact 3-part aluminum telescopic ramp with non-slip coating Folded dimensions: 870 x 237 mm 13.6 kg per pair - maximum capacity: 260 kg											805,69													
1654130	☐	Aluminum Folding Mono Ramp Folded dimensions: 1824 x 714 mm 13.9 kg - maximum capacity: 300 kg											426,54													
													1440530													
													1654128													
													1654129													
													1654130													

Ramps

1440530	☐	Fiberglass Folding Mono Ramp Folded dimensions: 2000 x 230 mm 16.3 kg - maximum capacity: 300 kg		1 787,68	
1654128	☐	Basic ramp 2 telescopic parts in aluminum and non-slip coating Folded dimensions: 1210 x 194 mm 8.8 kg per pair - maximum capacity: 270 kg		402,84	
1654129	☐	Quick2go Compact 3-part aluminum telescopic ramp with non-slip coating Folded dimensions: 870 x 237 mm 13.6 kg per pair - maximum capacity: 260 kg		805,69	
1654130	☐	Aluminum Folding Mono Ramp Folded dimensions: 1824 x 714 mm 13.9 kg - maximum capacity: 300 kg		426,54	

*Στις παραπάνω τιμές δεν συμπεριλαμβάνεται ο ΦΠΑ 13%