

Küschall® Compact

Please contact Invacare EU Export Customer Services on :-
Tel : +33 247 626980 ; Email : serviceclient_export@invacare.com
Or visit our website : www.invacare.eu.com



Küschall Compact



Active aluminum manual wheelchair folding by brace



Description :

- Aluminum chassis with hydroformed tubes
- Customization kit with premium decals
- Available in RETRACTABLE BRACKETS and FIXED BRACKETS
- Rear frame protection as standard
- Max user weight: 130 kg
- Complete weight: from 10.5 kg (retractable) and 9.6 kg (fixed)
- Transport weight: from 7.8 kg (retractable) and 6.7 kg (fixed)
- Warranty: 2 years excluding wearing parts

Discover the new possibilities of the special fabrication **UNIQUE**. Push back the limits of the personalization with a wide range of possibilities. Create **UNIQUE** wheelchairs.

Discover our new tool to help calculate overall dimensions and find the visual configurator to make your configuration!

Date :

Purchase Order :

Customer Code :

Distributor Stamp :

We reserve the right to deliver a + or - tolerance of +/- 10 mm on all measurements We reserve the right to make technical changes without notice

☐ Please make your selection with a cross in the box

◆ Standard feature

○ Free of charge option

CODE		DESCRIPTION	NC : Non Compatible NC : 9999 = not compatible with the reference 9999	Euros	
DCA0100	☐	küschall® Compact® - Fixed backrest		2850	 0120
DCA0110	☐	küschall® Compact® - Angle adjustable backrest	NC : 2120	2965	
DCA0120	☐	küschall® Compact® - Foldable and angle adjustable backrest	NC : 1615, 1620, 1625, 2120, 2180, 2430	3050	

SPECIAL CONSTRUCTION

DCA0000	☐	Myküschall special construction		On Request	
UCA0001	☐	UNIQUE - Personalization		On Request	


POSITIONING PACKAGE

DCA0220	☐	Posture and pressure relief package Libra Include : Matrix Libra cushion (DCA2720) Matrix MX2 backrest (DCA2300 to DCA2315) Bodypoint Evoflex belt (DCA7810 to DCA7830) Please select size of back and belt below in this order form!	SW 360 mm only possible with SD 400/420 mm (NC : 1300, 1305, 1310, 1315, 1350, 1355, 1400, 1405, 1415, 1425, 7230)	1222,23	Save 15% by selecting this PACK, i.e. €277,53
DCA0230	☐	Posture and pressure relief package Libra/Elite E2 include : Matrix Libra cushion (DCA2720) Matrix Elite/Elite Deep E2 backrest (DCA2320-DCA2350) Bodypoint Evoflex belt (DCA7810 to DCA7830) Please select size of back and belt below in this order form!	SW 360 mm only possible with SD 400/420 mm (NC : 1300, 1305, 1310, 1315, 1400, 1405, 1415, 1425, 7230)	1105,59	Save 10% by selecting this PACK, i.e. €106,70

SEATING CHOICE


Front Frame					
DCA1100	☐	Swing away standard	NC : 1900, 1902, 3040, 3140	◆	 Retractable
DCA1110	☐	Fixed frame active, 70°	NC : 1700 to 1730, 1805	○	
DCA1120	☐	Fixed frame dynamic, 80° Minimum seat angle (fstf to rstf) is 40 mm with 5" front wheels	NC : 1700 to 1730, 1805, 1962, 5550, 5600, 5650	117,54	
DCA1130	☐	Fixed frame tapered front frame 70° Additionnal delivery time + 10 days	NC : 1300, 1305, 1700 to 1730, 1805	132,70	

DCA1140	☐	Fixed frame tapered front frame 80° Minimum seat angle (fstf to rstf) is 30 mm - Additional delivery time + 10 days			NC : 1300, 1305, 1700 to 1730, 1805, 1962, 5550, 5600, 5650	132,70	 Fixed
Seat width (SW) (Width 280 to 420 mm: single cross / width from 440 to 500 mm: double cross)							
DCA1300	☐	280 mm	DCA1305	☐	300 mm	NC : 2210, Matrx backrest, Matrx and Classic Visco cushions, 2490, 2860 for L.280 mm, 3020, 3030, 3120, 3130, 3260, 6000, 6010, 8450	○
DCA1310	☐	320 mm	DCA1315	☐	340 mm	NC : Dossiers Matrx, coussins Matrx and Classic Visco, 2490, 3020, 3030, 3120, 3130, 3260, 6000, 6010, 8450	○
DCA1320	☐	360 mm			NC : 2364, 2366, 2368, 2370, 2372, 2490, 2710, 2742, 2744, 2746, 2748, 3020, 3030, 3120, 3130	○	
DCA1325	☐	380 mm			NC : 2366, 2368, 2370, 2372, 2490 2742, 2744, 3020, 3030, 3120, 3130	○	
DCA1330	☐	400 mm			NC : 2360, 2366, 2368, 2370, 2372, 2744	○	
DCA1335	☐	420 mm			NC : 2360, 2362, 2368, 2370, 2372, 2746	○	
DCA1340	☐	440 mm			NC : 2360, 2362, 2364, 2370, 2372, 2740, 2742	○	
DCA1345	☐	460 mm			NC : 2360, 2362, 2364, 2372, 2740, 2742, 2744	○	
DCA1350	☐	480 mm			NC : 2300, 2305, 2310, 2315, 2345, 2360, 2362, 2364, 2366, 2740, 2742, 2744	○	
DCA1355	☐	500 mm			NC : 2210, 2300, 2305, 2310, 2315, 2345, 2360, 2362, 2364, 2366, 2368, 2740, 2742, 2744, 2746	○	
Seat depth (SD) (Chassis adjusted to seat depth from depth 380 mm (single cross) and from depth 400 mm (double cross))							
DCA1400	☐	320 mm	DCA1405	☐	340 mm	NC : 2155, 2160, Matrx cushions and Classic Visco, 2860, 2910, 5650	○
DCA1415	☐	360 mm			NC : 2710, 2720, 2742, 2744, 2746, 2748, 2860, 5650	○	
DCA1425	☐	380 mm			NC : 2720, 2744, 5650	○	
DCA1435	☐	400 mm			NC : 2744	◆	
DCA1445	☐	420 mm				○	
DCA1460	☐	440 mm			NC : 2740, 2742	○	
DCA1470	☐	460 mm			NC : 2740, 2742, 2744	○	
DCA1480	☐	480 mm			NC : 2740, 2742, 2744, 2746, 2748	○	
DCA1490	☐	500 mm			NC : 2710, 2740, 2742, 2744, 2746, 2748	○	
Backrest Height (BH) (Height 300 to 330 mm NC angle adjustable and folding backrest)							
DCA1615	☐	300 mm	DCA1620	☐	315 mm	NC : 2240, 2595, 7600	○
DCA1625	☐	330 mm	DCA1630	☐	345 mm	NC : 2240, 2595, 7600	○
DCA1635	☐	360 mm			NC : 2240, 2595, 7600	○	
DCA1640	☐	375 mm			NC : 2240, 2595, 7600	◆	
DCA1645	☐	390 mm			NC : 2240, 2595, 7600	○	
DCA1650	☐	405 mm	DCA1655	☐	420 mm	NC : 2240	○
DCA1660	☐	435 mm	DCA1665	☐	450 mm	NC : 2240	○
DCA1670	☐	465 mm	DCA1675	☐	480 mm	NC : 2240	○
DCA1680	☐	495 mm	DCA1685	☐	510 mm	NC : 2240	○



Front Seat to Floor Height (FSTF) - see table 3 at the end of the measurement sheet						
DCA1700	☐	370 mm		NC : 1810 to 1870, 1932 to 1962, 3030, 3130, 6410, 6420, 7600		☐
DCA1705	☐	380 mm		NC : 1815 to 1870, 1934 to 1962, 3030, 3130, 6410, 6420, 7600		☐
DCA1710	☐	390 mm		NC : 1820 to 1870, 1936 to 1962, 3030, 3130, 6410, 6420, 7600		☐
DCA1715	☐	400 mm		NC : 1825 to 1870, 1938 to 1962, 3030, 3130, 6410, 6420, 7600		☐
DCA1720	☐	410 mm		NC : 1830 to 1870, 1940 to 1962, 3030, 3130		☐
DCA1725	☐	420 mm		NC : 1835 to 1870, 1942 to 1962		☐
DCA1730	☐	430 mm		NC : 1840 to 1870, 1944 to 1962		☐
DCA1735	☐	440 mm		NC : 1845 to 1870, 1946 to 1962		☐
DCA1740	☐	450 mm		NC : 1850 to 1870, 1948 to 1962		☐
DCA1745	☐	460 mm		NC : 1855 to 1870, 1950 to 1962		☐
DCA1750	☐	470 mm		NC : 1860 to 1870, 1952 to 1962		☐
DCA1755	☐	480 mm		NC : 1805, 1865, 1870, 1954 to 1962		☐
DCA1760	☐	490 mm		NC : 1805, 1810, 1956 to 1962		☐
DCA1765	☐	500 mm		NC : 1805, 1810, 1815, 1958 to 1962		☒
DCA1770	☐	510 mm		NC : 1805 to 1820, 1960, 1962		☐
DCA1775	☐	520 mm		NC : 1805 to 1825, 1962		☐
DCA1780	☐	530 mm		NC : 1805 to 1830		☐
Rear Seat to Floor Height (RSTF) With fixed brackets at 70°, FSTF - 100 mm < RSTF < FSTF - 0 mm With fixed stems at 80°, FSTF - 100 mm < RSTF < FSTF - 30 mm With 70° and 80° retractable brackets, FSTF - 100 mm < RSTF < FSTF - 0 mm						
DCA1805	☐	370 mm		Only with 22" wheels and swing away frame		☐
DCA1810	☐	380 mm		NC : 6410, 6420, 6430		☐
DCA1815	☐	390 mm	DCA1820	☐	400 mm	☐
DCA1825	☐	410 mm		NC : 6420, 6430		☐
DCA1830	☐	420 mm	DCA1835	☐	430 mm	☐
DCA1840	☐	440 mm	DCA1845	☐	450 mm	☐
DCA1850	☐	460 mm				☒
DCA1855	☐	470 mm				☐
DCA1860	☐	480 mm	DCA1865	☐	490 mm	☐
DCA1870	☐	500 mm		NC : 6400, 6410		☐
Lower Leg Length (LLL) (Depending on the configuration, the pallets can be mounted high) - see table 1. With retractable and fixed brackets at 70°, LLL max = FSTF - 20 mm With retractable and fixed brackets at 80°, LLL max = FSTF - 30 mm With legrest in 75° position, LLL max = FSTF - 40 mm and LLL min = 380 mm With legrest in 60° position, LLL max = FSTF and LLL min = 380 mm "						
DCA1900	☐	200 mm		NC : 3030, 3130		☐
DCA1902	☐	210 mm		NC : 3030, 3130		☐
DCA1904	☐	220 mm	DCA1906	☐	230 mm	☐
DCA1908	☐	240 mm	DCA1910	☐	250 mm	☐
DCA1912	☐	260 mm	DCA1914	☐	270 mm	☐
DCA1916	☐	280 mm	DCA1918	☐	290 mm	☐
DCA1920	☐	300 mm	DCA1922	☐	310 mm	☐
DCA1924	☐	320 mm	DCA1926	☐	330 mm	☐
DCA1928	☐	340 mm	DCA1930	☐	350 mm	☐
DCA1932	☐	360 mm	DCA1934	☐	370 mm	☐
DCA1936	☐	380 mm	DCA1938	☐	390 mm	☐
DCA1940	☐	400 mm	DCA1942	☐	410 mm	☐
DCA1944	☐	420 mm				☐
DCA1946	☐	430 mm				☒
DCA1948	☐	440 mm	DCA1950	☐	450 mm	☐
DCA1952	☐	460 mm	DCA1954	☐	470 mm	☐
DCA1956	☐	480 mm	DCA1958	☐	490 mm	☐
DCA1960	☐	500 mm				☐
DCA1962	☐	510 mm		NC : 3000, 3100, 3010, 3110		☐

IMPORTANT:
Depending on the seat height selected and the diameter of the front wheels: assembly with a hemi or standard chassis and with/without high mounting - refer to table 3 at the end of the measurement sheet




Hemi Frame



High Mount

Rear wheel diameter / RSTF 22":
370 - 470mm 24":
390 - 490mm 25":
420 - 500mm 26":
420 - 500mm



LLL = LLL patient ; Seat cushion

BACKREST TYPE					
Backrest Tube					
DCA2120	☐	Fixed backrest with lumbal bend 7° (degree)		NC : Matrx Backrests	52,13
Backrest angle adjustable by screws in relation to the seat (BA) (only if 0110 or 0120 selected / with folding backrest only 2 adjustment ranges: 82° to 90° or 94° to 98°) (Angle measured in relation to the seat))					
DCA2155	☐	82°	DCA2160	☐	86°
DCA2165	☐	90°			◆
DCA2170	☐	94°	DCA2175	☐	98°
				NC : dossier rabattable + dossiers Matrx	○
DCA2180	☐	102°			NC : Dossier rabattable
Backrest Upholstery					
DCA2200	☐	Standard upholstery			NC : 2300 to 2357, 2400, 2430, 2595
DCA2210	☐	Light back (hand-sewn, airy, anti-perspirant)			NC : 2300 to 2357, 2430, 2595
DCA2220	☐	Tension adjustable premium black back upholstery with loop bands			NC : 2300 to 2357
DCA2230	☐	Aero back, tension adjustable with loop bands			NC : 2300 to 2357
DCA2240	☐	Without back upholstery, mandatory if Matrx backrest ordered			NC : 2430, 2595
UNIQUE - BACK UPHOLSTERY STITCHING COLOUR - Only for Aero Back upholstery					
UCA2241	☐	NEON YELLOW			101,42
UCA2242	☐	NEON PINK			101,42
UCA2243	☐	BLUE			101,42
UCA2244	☐	WHITE			101,42
UCA2245	☐	RED			101,42
UCA0002	☐	"Your personal color": Please contact our customers service			On Request
UNIQUE - BACK UPHOLSTERY - EMBROIDERY - Only for Aero Back upholstery					
UCA2246	☐	SMILEY			218,01
UCA2248	☐	TRAVELER			218,01
UCA2249	☐	SAKURA BEAUTY			218,01
UCA2250	☐	FANNY CACTUS			218,01
UCA2251	☐	KUSCHALL LOGO: Please contact our customers service			On Request
UCA0003	☐	"Your personal design": Please contact our customers service - Provide the wished design in .ai / .pdf or .png format			On Request
UNIQUE - BACK UPHOLSTERY - FLOCKING - Only for Aero Back upholstery					
UCA2253	☐	SMILEY			118,48
UCA2255	☐	TRAVELER			118,48
UCA2256	☐	SAKURA BEAUTY			118,48
UCA2257	☐	FANNY CACTUS			118,48
UCA2259	☐	PAINT STAIN			118,48
UCA0004	☐	"Your personal design": Please contact our customers service - Provide the wished design in .ai / .pdf or .png format			On Request
Matrx Backrests - The height shown is that of the backrest. Find the trunk supports and all the accessories on the Matrx measurement sheet!					
DCA2300	☐	Matrx MX2 , height 150 mm with quick release hardware - User max. weight 114 kg			NC : SW > 460 mm and SW < 360 mm, 2370, 2372, 2490
DCA2305	☐	Matrx MX2 , height 230 mm with quick release hardware - User max. weight 114 kg			NC : SW > 460 mm and SW < 360 mm, 2370, 2372, 2490
DCA2310	☐	Matrx MX2 , height 305 mm with quick release hardware - User max. weight 114 kg			NC : SW > 460 mm and SW < 360 mm, 2370, 2372, 2490
DCA2315	☐	Matrx MX2 , height 405 mm with quick release hardware - User max. weight 114 kg			NC : SW > 460 mm and SW < 360 mm, 2370, 2372, 2490
DCA2316	☐	Elite Deep E2 , height 160 mm with mini-set mounting hardware			NC : SW < 360 mm, 2490
DCA2317	☐	Elite Deep E2 , height 210 mm with mini-set mounting hardware			NC : SW < 360 mm, 2490
DCA2318	☐	Elite Deep E2 , height 260 mm with mini-set mounting hardware			NC : SW < 360 mm, 2490



2120

MX2: only 140g more than a tension-adjustable backrest!



2300, 2305, 2310, 2315





2120

MX2: only 140g more than a tension-adjustable backrest!



2300, 2305, 2310, 2315

DCA2320	☐	Elite Deep E2 , height 310 mm with mini-set mounting hardware	NC : SW < 360 mm, 2490	717,54	 2316 to 2331
DCA2322	☐	Elite Deep E2 , height 360 mm with mini-set mounting hardware	NC : SW < 360 mm, 2490	717,54	
DCA2325	☐	Elite Deep E2 , height 410 mm with mini-set mounting hardware	NC : SW < 360 mm, 2490	717,54	
DCA2330	☐	Elite Deep E2 , height 460 mm with mini-set mounting hardware	NC : SW < 360 mm, 2490	717,54	
DCA2331	☐	Elite Deep E2 , height 510 mm with mini-set mounting hardware	NC : SW < 360 mm, 2490	717,54	
DCA2332	☐	Elite E2 , height 160 mm with mini-set mounting hardware	NC : SW < 360 mm, 2490	717,54	 2332 to 2357
DCA2334	☐	Elite E2 , height 210 mm with mini-set mounting hardware	NC : SW < 360 mm, 2490	717,54	
DCA2335	☐	Elite E2 , height 260 mm with mini-set mounting hardware	NC : SW < 360 mm, 2490	717,54	
DCA2340	☐	Elite E2 , height 310 mm with mini-set mounting hardware	NC : SW < 360 mm, 2490	717,54	
DCA2345	☐	Elite E2 , height 360 mm with mini-set mounting hardware	NC : SW > 460 mm and SW < 360 mm, 2490	717,54	
DCA2350	☐	Elite E2 , height 410 mm with mini-set mounting hardware	NC : SW < 360 mm, 2490	717,54	
DCA2355	☐	Elite E2 , height 460 mm with mini-set mounting hardware	NC : SW < 360 mm, 2490	717,54	 2430 2450
DCA2357	☐	Elite E2 , height 510 mm with mini-set mounting hardware	NC : SW < 360 mm, 2490	717,54	
Matrx Backrest Width - see tables 4 at the end of the measurement sheet					
DCA2360	☐	360 mm - indicative bust width 265 mm		○	
DCA2362	☐	380 mm - indicative bust width 290 mm		○	
DCA2364	☐	410 mm - indicative bust width 315 mm		○	
DCA2366	☐	430 mm - indicative bust width 340 mm		○	
DCA2368	☐	460 mm - indicative bust width 370 mm		○	
DCA2370	☐	480 mm - indicative bust width 390 mm		○	 2490
DCA2372	☐	510 mm - indicative bust width 420 mm		○	
Push Handles					
DCA2400	☐	Without push handles	NC : 2490	○	 2430 2450
DCA2410	☐	Push handles fixed, mini	NC : 2490, 6580	○	
DCA2420	☐	Push handles fixed, long		◆	
DCA2430	☐	Height adjustable rear set	NC : 2490, 3520, 3620	212,32	
DCA2440	☐	Height adjustable, integrated	NC : 2490, 2595	254,03	
DCA2450	☐	Foldable push handles	NC : 2490, 2595, 6580	213,27	
Headrest					 2490
DCA2490	☐	Headrest with lateral adjustment and adjustable ear cushions - supplied with backrest adjuster	NC : 2595, 6580	168,72	
Backrest Stabilization Bar					 2610
DCA2595	☐	Backrest stabilization bar, swing away and removable	NC : BH ≤ 390 mm	151,66	
SEAT UPHOLSTERY					
DCA2610	☐	Standard seat upholstery - central loop band		◆	 2610
DCA2620	☐	Seat upholstery with 1 central (soft) and 2 side straps fasteners bands (rough)		37,91	

UNIQUE - SEAT UPHOLSTERY STITCHING COLOUR				
UCA2602	☐	NEON YELLOW		118,48
UCA2603	☐	NEON PINK		118,48
UCA2604	☐	BLUE		118,48
UCA2605	☐	WHITE		118,48
UCA2606	☐	RED		118,48
UCA0005	☐	"Your personnel color": <i>Please contact our customers service</i>		On Request
SEATING				
Seat Cushion				
DCA2700	☐	STANDARD SEAT CUSHION, 60 MM DOUBLE FOAM WITH VISCO TOP	NC : 2810, 2820, 2830, 2840, 2850	101,42
DCA2710	☐	MATRX FLO-TECH LITE	NC : 2810, 2820, 2830, 2840, 2850	118,48
DCA2720	☐	MATRX LIBRA CUSHION, BLACK SMOOTH COVER	NOT POSSIBLE WITH SW ≤ 340 MM & SD ≤ 380 MM SW 360 MM ONLY POSSIBLE WITH SD 400/420 MM	293,84
DCA2740	☐	CLASSIC VISCO CUSHION 375 X 400 MM ONLY FOR USERS WEIGHTING UP TO 50 KG	NC : 2810, 2820, 2830, 2840, 2850	76,78
DCA2742	☐	CLASSIC VISCO CUSHION 405 X 400 MM ONLY FOR USERS WEIGHTING UP TO 100 KG	NC : 2810, 2820, 2830, 2840, 2850	65,40
DCA2744	☐	CLASSIC VISCO CUSHION 430 X 430 MM ONLY FOR USERS WEIGHTING UP TO 100 KG	NC : 2810, 2820, 2830, 2840, 2850	65,40
DCA2746	☐	CLASSIC VISCO CUSHION 455 X 400 MM ONLY FOR USERS WEIGHTING UP TO 100 KG	NC : 2810, 2820, 2830, 2840, 2850	76,78
DCA2748	☐	CONTOUR VISCO NG CUSHION	NC : 2810, 2820, 2830, 2840, 2850	76,78
DCA2750	☐	VICAIR ACTIVE O2 CUSHION		683,73
SEAT OPTION				
DCA2800	☐	SEAT BAG, CLICK IN AND REMOVABLE		81,52
Accessoires coussin Matrx Libra				
DCA2810	☐	FLUID SUPPLEMENT SAC (PAIR)		91,00
DCA2820	☐	TAPERED WEDGE PAIR		41,71
DCA2830	☐	ABDUCTOR		47,39
DCA2840	☐	LATERAL SIDE WEDGE		41,71
DCA2850	☐	ANTERIOR FRONT WEDGE		54,98
Base rigide				
DCA2860	☐	SEAT RIGIDIZER	(Mandatory: choose a seat cushion)	56,87
SEAT BAGS				
DCA2910	☐	BIG BAG (TO FIX UNDER SEAT)	NC : SD ≤ 340 mm	61,61
LEGRESTS AND FOOTRESTS				
LEGRESTS FOR SWING AWAY				
RIGHT LEGREST				
DCA3000	☐	LEGREST, ANGLE OPEN, 70° - RIGHT FOOTREST MANDATORY	NC : 3110	◆
DCA3010	☐	LEGREST, ANGLE MID-CLOSED, 80° - RIGHT FOOTREST MANDATORY	NC : 3100, 5650	○
DCA3020	☐	AMPUTEE LEGREST - RIGHT	NC : 3280	274,88
DCA3030	☐	ANGLE ADJUSTABLE LEGREST - RIGHT INCLUDING ANGLE, WIDTH AND DEPTH ADJUSTABLE PLASTIC FOOTPLATE	NC : 3280	100,25
DCA3040	☐	WITHOUT LEGRESTS ON THE RIGHT SIDE ONLY FIXED FRAME		○

 2700

 2710

 2720


2740 / 2742
2744 / 2746

 2748

 2800

 2810

 2820

 2830

 2840

 2850

 2860


3000 / 3100

LEFT LEGREST				 3020 / 3120  3030 / 3130
DCA3100	☐	LEGREST, ANGLE OPEN, 70° - LEFT FOOTREST MANDATORY	◆	
DCA3110	☐	LEGREST, ANGLE MID-CLOSED, 80° - LEFT FOOTREST MANDATORY	NC : 5650 ○	
DCA3120	☐	AMPUTEE LEGREST - LEFT	NC : 3280 274,88	
DCA3130	☐	ANGLE ADJUSTABLE LEGREST - LEFT INCLUDING ANGLE, WIDTH AND DEPTH ADJUSTABLE PLASTIC FOOTPLATE	NC : 3280 100,25	
DCA3140	☐	WITHOUT LEGRESTS ON THE LEFT SIDE ONLY FIXED FRAME	○	
FOOTREST				
DCA3260	☐	2 PC. FLIP UP, PLASTIC - NEW PART	NOT POSSIBLE WITH SW < 340 MM, LLL < 260 MM, TAPERED FRAME SW < 380 MM	◆
DCA3270	☐	2 PC. - ANGLE AND DEPTH ADJUSTABLE, CARBON	NC : 3310	311,85
DCA3280	☐	1 PC. FLIP-UP CARBON FOOTPLATE, ANGLE AND DEPTH ADJUSTABLE	NOT POSSIBLE WITH LLL < 300 MM, NC : 3310	329,86
FOOT AND LEG SUPPORT				
DCA3300	☐	CALF STRAP MOUNTED ON THE FRAME		◆
DCA3305	☐	CALF STRAP WITH BRACKET SUPPORT		76,78
DCA3310	☐	HEEL LOOPS		72,99
SIDE GUARDS				
CLOTHES-GUARDS, MUDGUARDS				
DCA3400	☐	ECONOMY CLOTHES GUARDS, FIXED <i>Increase seat width by 0 to +10 mm</i>	NC : 3500, 3510, 3530, 3540, 3600, 3610, 3630, 3640, 7600, 8450	◆
DCA3405	☐	ECONOMY CLOTHES GUARDS, FLIP BACK <i>Decrease seat width by 0 to +10 mm</i>	NC : 3500, 3510, 3530, 3540, 3550, 3600, 3610, 3630, 3640, 7600, 8450	○
DCA3410	☐	CARBON CLOTHES GUARDS, FIXED <i>Increase seat width by 0 to +10 mm</i>	NC : 3500, 3510, 3530, 3540, 3600, 3610, 3630, 3640, 7600, 8450	167,77
DCA3425	☐	ECONOMY MUDGUARDS, FIXED <i>Increase seat width by 0 to +10 mm</i>	NC : 3500, 3510, 3530, 3540, 3600, 3610, 3630, 3640, 6000, 6010, 6140, 6150, 6580, 7600, 8450	○
DCA3430	☐	CARBON MUDGUARDS, FIXED NOT POSSIBLE WITH FSTF < 410 MM <i>Increase seat width by 0 to +10 mm</i>	NC : 3500, 3510, 3530, 3540, 3600, 3610, 3630, 3640, 6000, 6010, 6140, 6150, 6580, 6950, 7600, 8450, NC FSTF < 410 mm	223,70
DCA3440	☐	CARBON MUDGUARDS, REMOVABLE <i>Decrease seat width by 0 to +10 mm</i>	NC : 3500, 3510, 3530, 3540, 3600, 3610, 3630, 3640, 6140, 6150, 8450	336,49
DCA3450	☐	ECONOMY MUDGUARDS, REMOVABLE <i>Decrease seat width by 0 to +10 mm</i>	NC : 3500, 3510, 3530, 3540, 3600, 3610, 3630, 3640, 6140, 8450	○
UNIQUE - CUSTOM TRANSPARENT COLORS - Only for economy clothes guards				
UCA3422	☐	TRANSPARENT		520,38
UCA3423	☐	TRANSPARENT BLUE		745,02
UCA3424	☐	NEON YELLOW		745,02
UCA3425	☐	NEON ORANGE		745,02
UCA3426	☐	NEON RED		745,02
UCA0006	☐	"Your personal color": Please contact our customers service		On Request
UCA0007	☐	"Your personal design": Please contact our customers service - <i>Provide the wished design in .ai / .pdf or .eps format</i>		On Request

	3500 3510 3600 3610
	3530 3630
	3540 3640

UCA4101	☐	DIAMOND SILVER		696,68
UCA4102	☐	PASTEL GREEN		750,71
UCA4103	☐	SKY BLUE		750,71
UCA4104	☐	PLUM PURPLE		750,71
UCA4105	☐	SALMON		750,71
DCA4070	☐	ATOMIC VIOLET- specific color mix, no RAL definition		On quote
DCA4050	☐	Special Frame Colour (on demand only): - Manufacturing time + 10 days: for a color already offered on the other ranges of Invacare manual wheelchairs. - Manufacturing time + 15 days: for any other request (subject to feasibility) - specify RAL number		750,71
UNIQUE - CERAMIC GLAZE FINISH				
UCA4066	☐	CERAMIC GLAZE FINISH (ONLY FOR UKA**** COLOURS)	NC : 4300, 4310	643,60
STYLING				
DCA4190	☐	Bold	NC : 4300, 4310	64,45
DCA4195	☐	Tag	NC : 4200 to 4250	◆
BOLD COLOR				
DCA4200	☐	BLACK		○
DCA4210	☐	DARK GREY		○
DCA4220	☐	LIGHT GREY		○
DCA4230	☐	RED		○
DCA4240	☐	RED		○
DCA4250	☐	GOLD		○
UCA4252	☐	NEON YELLOW		161,14
UCA4253	☐	NEON GREEN		161,14
UCA4254	☐	NEON ORANGE		161,14
UCA4255	☐	NEON PINK		161,14
UCA4256	☐	NEON RED		161,14
UCA0009	☐	"Your personnal color": Please contact customers service		On Request
TAG COLOR				
DCA4300	☐	BLACK		○
DCA4310	☐	RED		○
FRONT WHEEL				
DCA5050	☐	4" EVERYDAY WHEEL WITH BLACK TYRE		○
DCA5100	☐	4" COMFORT WHEEL WITH BLACK TYRE		○
DCA5150	☐	4" COMFORT WHEEL WITH BLACK TYRE		75,83
DCA5200	☐	4" FROGLEGS FRONT WHEEL WITH BLACK HUB, BLACK TYRE	Select the Froglegs front wheels with black hub (DCA5200) to access UNIQUE anodizations	121,33
DCA5250	☐	4" FROGLEGS FRONT WHEEL WITH SILVER HUB, BLACK TYRE		121,33
DCA5300	☐	5" EVERYDAY WITH BLACK TYRE		○
DCA5350	☐	5" COMFORT WHEEL BLACK TYRE		◆
DCA5400	☐	5" STAREC BLACK HUB, BLACK TYRE		103,32
DCA5450	☐	5" FROGLEGS FRONT WHEEL WITH BLACK HUB, BLACK TYRE	Select the Froglegs front wheels with black hub (DCA5450) to access UNIQUE anodizations	147,87
DCA5500	☐	5" FROGLEGS FRONT WHEEL WITH SILVER HUB, BLACK TYRE		147,87
DCA5550	☐	6" COMFORT WITH BLACK TYRE		○
DCA5600	☐	6" PNEUMATIC WHEEL		38,86
DCA5650	☐	7" COMFORT WITH BLACK TYRE		38,86
UNIQUE - FROGLEGS HUB COLOURS				
UCA5502	☐	PINK		98,58
UCA5504	☐	ORANGE		98,58
UCA5506	☐	BLUE		98,58
UCA5507	☐	VIOLET		98,58
UCA5508	☐	RED		98,58
UCA5509	☐	GOLD		98,58
UNIQUE - CERAMIC BEARING				
UCA5510	☐	CERAMIC BEARING		98,58



BOLD



TAG



5200 / 5250 / 5450 / 5500



5350 / 5550



5050 / 5300




BOLD



TAG



5200 / 5250 / 5450 / 5500



5350 / 5550



5050 / 530

REAR WHEEL CHOICE				
ONE-ARM-DRIVE SYSTEM				
DCA6000	☐	ONE ARM DRIVE SYSTEM RIGHT SIDE (WITH LIGHT WHEEL 7 BAR TYRES, ALU HANDRIMS)	NC : 6040, 6130, 6140, 6150, 6300, 6400, 6420, 6430, 6500, 6510, 6520, 6550, 6580, 6680, 6840, 6900, 6910, 6950, 8650, 8680	690,05
DCA6010	☐	ONE ARM DRIVE SYSTEM LEFT SIDE (WITH LIGHT WHEEL 7 BAR TYRES, ALU HANDRIMS)	NC : 6040, 6130, 6140, 6150, 6300, 6400, 6420, 6430, 6500, 6510, 6520, 6550, 6580, 6680, 6840, 6900, 6910, 6950, 8650, 8680	690,05
REAR WHEEL FIXATION				
DCA6020	☐	STANDARD FIXATION	NC : 6140, 6150	◆
DCA6040	☐	REAR WHEEL EXTENSION	NC : 6100, 6580, 6950	134,60
REAR WHEEL POSITION				
DCA6100	☐	1ST POSITION FROM THE FRONT	NC : 6580, 6840	○
DCA6110	☐	2ND POSITION FROM THE FRONT	NC : 6950	○
DCA6120	☐	3RD POSITION FROM THE FRONT		◆
DCA6130	☐	4TH POSITION FROM THE FRONT	NC : 6580, 6950	○
DCA6140	☐	5TH POSITION FROM THE FRONT		○
DCA6150	☐	6TH POSITION FROM THE FRONT		○
QUAD QUICK RELEASE AXLE				
DCA6300	☐	TETRA-CLIP	NC : 6550, 6580, 6950	143,13
REAR WHEEL SIZE				
DCA6400	☐	22"	NC : 6500, 6520, 6550, 6560, 6600, 6630, 6660, 6760, 6710, 6720, 6740, 6750, 6760, 6770, 6950	336,49
DCA6410	☐	24"		◆
DCA6420	☐	25"	NC : 6500, 6510, 6580, 6600, 6630, 6660, 6700, 6710, 6720, 6760, 6770, 6900, 6910, 6950	○
DCA6430	☐	26"	NC : 6500, 6520, 6550, 6600, 6630, 6660, 6720, 6740, 6750, 6760, 6770, 6840, 6950	517,54
REAR WHEEL TYPE				
DCA6500	☐	24" ACTIVE WHEELS (BLACK RIM, HUB AND RADIAL SPOKING INCL. ALU HANDRIMS)	24" (NC : All Handrims, All Spokeguards)	◆
DCA6510	☐	KUSCHALL STAREC HUB (BLACK RIM AND HUB, BLACK CROSS SPOKING)	22", 24", 26"	○
DCA6520	☐	HIGH PERFORMANCE (BLACK RIM, BLACK STAREC HUB, BLACK CROSS SPOKING)	24", 25" NC : 6660	533,65
DCA6550	☐	SPINERGY LXX (BLACK RIM, HUB AND RADIAL SPOKING)	24", 25" NC : 6660, 6910	771,56
UCA6551	☐	SPINERGY LITE EXTREME LX - 680g per wheel	24", 25"	1060,66
UCA6552	☐	SPINERGY LITE EXTREME X-LACED XLX - 660g per wheel	24", 25"	1140,28


6300

Tyre size
MAXXIS LIGHT GUARD /
ALPHA PROTECT :
24" : 25-540

MBL SPEEDLITE :
24" : 25-540

 6600  6625

UCA6553	☐	SPINERGY X-LACED SPORT LITE EXTREME XSLX - 757g per wheel	24", 25"		1206,64
UCA6554	☐	SPINERGY CARBON FIBER CLX 570g per wheel	24", 25"		1358,29
DCA6580	☐	22", 24", 26" WHEELS WITH ATTENDANT DRUM BRAKE	22", 24", 25"	NC : 6950, 8680	518,48
UNIQUE - RIM ANODIZED COLOURS					
UCA6557	☐	BLUE			162,09
UCA6558	☐	GOLD			162,09
UCA6559	☐	GREEN			162,09
UCA6560	☐	ORANGE			162,09
UCA6561	☐	RED			162,09
UCA6562	☐	VIOLET			162,09
UNIQUE - HUB ANODIZED COLOURS					
UCA6563	☐	BLUE			162,09
UCA6564	☐	GOLD			162,09
UCA6565	☐	SILVER GREY (RAL 9006)			162,09
UCA6566	☐	GREEN			162,09
UCA6567	☐	ORANGE			162,09
UCA6568	☐	RED			162,09
UCA6569	☐	VIOLET			162,09
UNIQUE - SPOKE COLOURS					
UCA6570	☐	BLUE			162,09
UCA6571	☐	GREEN			162,09
UCA6572	☐	ORANGE			162,09
UCA6573	☐	PINK			162,09
UCA6574	☐	RED			162,09
UCA6575	☐	VIOLET			162,09
UCA6576	☐	WHITE			162,09
UCA6577	☐	YELLOW			162,09
UNIQUE - CERAMIC BEARING					
UCA6579	☐	CERAMIC BEARING			332,70

REAR WHEEL TYRE					 6630  6660
DCA6600	☐	MAXXIS LITE GUARD BLACK TYRE (puncture protection)	24"	◆	
DCA6625	☐	MAXXIS ALPHA PROTECT BLACK TYRE (puncture protection)	24"	72,04	
DCA6630	☐	MBL SPEEDLITE (puncture protection)	24"	72,04	
DCA6660	☐	24" SOLID TYRE BLACK	24" NC : 7230	72,04	
HANDRIM					 6700 6710  6720  6740  6790 6760
DCA6700	☐	ANTI-SLIP DIST. 3 CM ADJUSTABLE 4CM		142,18	
DCA6710	☐	STEEL, DIST. 3 CM ADJUSTABLE 4CM		113,74	
DCA6720	☐	TITANIUM, DIST. 3 CM ADJUSTABLE 4CM		277,73	
DCA6740	☐	CARBOLIFE GEKKO, WITH ANTI-SLIP DIST. 3 CM ADJUSTABLE 4CM		373,46	
DCA6750	☐	CARBOLIFE GEKKO S, WITH ANTI-SLIP DIST. 3 CM ADJUSTABLE 4CM		373,46	
DCA6760	☐	CARBOLIFE CURVE L WITH ERGONOMIC PROFILE DIST. 3 CM ADJUSTABLE 4CM		282,46	
DCA6770	☐	CARBOLIFE QUADRO, WITH ERGONOMIC PROFILE DIST. 3 CM ADJUSTABLE 4CM		373,46	
DCA6780	☐	LIGHT ALU, GREY, DIST. 3 CM ADJUSTABLE 4CM		◆	
DCA6790	☐	LIGHT ALU, BLACK, DIST. 3 CM ADJUSTABLE 4CM		○	
REAR WHEEL CAMBER					
DCA6820	☐	0°		◆	
DCA6840	☐	3°	NC : 7230	○	
SPOKE GUARDS					
DCA6900	☐	22", 24", 26" TRANSPARENT		200,95	
DCA6910	☐	22", 24", 26" BLACK		200,95	
ALBER ADDS ON KIT					
DCA6950	☐	TWION AND E-MOTION M25 ADAPTER PLATE, 2 ALBER ANTI-TIPPERS INCLUDED	Only 24" wheels. NC : 7500, 7510, 7520, 8100, 8110, 8120, 8200, 8210	579,15	
Tip : for perfect compatibility of the brakes with the manual wheels, we advise you to equip them with the same coating as the motorized wheels.					
UNIQUE - CUSTOMIZATION					
UCA6960	☐	UNIQUE - CUSTOM COLOURS FORKS AND BEARING COVER		197,16	
UCA6961	☐	UNIQUE - CUSTOM COLOURS FORKS, BEARING COVER, FORKS HOLDERS, REAR WHEEL SUPPORTS, BRAKE CLAMPS, HANGER RELEASE LEVER		551,66	

UNIQUE - COLOURS PREMIUM ANODIZATION				
UCA6963	☐	RED		○
UCA6964	☐	VIOLET		○
UCA6965	☐	COPPER		○
UCA6966	☐	GREEN		○
UCA6967	☐	ORANGE		○
UCA6968	☐	GOLD		○
UCA6969	☐	CHAMPAGNE		○
UCA6970	☐	PINK		○
UCA6971	☐	BLUE		○
UCA6972	☐	POLY MIRROR FINISH		○
Brakes - see table 2 at the end				
DCA7200	☐	PULL BRAKES		◆
DCA7210	☐	PUSH BRAKES, BENT		○
DCA7220	☐	PUSH BRAKES LEVER STRAIGHT		○
DCA7230	☐	PERFORMANCE BRAKES	NC : 7300, 7810, 7820, 7830	76,78
DCA7250	☐	LIGHT ACTIVE BRAKE	NC : 7300	145,02
BRAKE EXTENSION LEVER				
DCA7300	☐	BRAKE EXTENSION LEVER		76,78
SAFETY OPTIONS				
Assembly of a maximum of 2 elements per side from the following list: anti-tippers, tipping tubes, cane holder, tube armrests or transit wheels				
ANTITIPPER				
DCA7500	☐	ANTITIPPER RIGHT		130,81
DCA7510	☐	ANTITIPPER LEFT		130,81
DCA7520	☐	ANTITIPPER BOTH SIDES		261,61
TRANSPORTATION KIT				
DCA7600	☐	TRANSPORTATION KIT (LABELS)	NC : 7800, BH ≤ 400 mm	57,82
POSTURE BELT				
DCA7800	☐	STANDARD POSTURE BELT - PLASTIC BUCKLE		42,65
DCA7810	☐	BODYPOINT POSTURE BELT EVOFLEX, 2 POINTS, SMALL, FLAT-MOUNT WITH PUSH BUTTON		217,06
DCA7820	☐	BODYPOINT POSTURE BELT EVOFLEX, 2 POINTS, MEDIUM, FLAT-MOUNT WITH REHAB LATCH		217,06
DCA7830	☐	BODYPOINT POSTURE BELT EVOFLEX, 2 POINTS, LARGE, FLAT-MOUNT WITH REHAB LATCH		217,06
PASSIVE LIGHTS				
DCA7900	☐	PASSIVE LIGHTS ASSEMBLED (MOUNTED ON SPOKES)		22,75
ATTENDANT OPTIONS				
FRAME PROTECTION				
DCA8000	☐	REAR FRAME PROTECTION		◆
DCA8010	☐	FRONT FRAME PROTECTION		44,55
TIPPER AID				
DCA8100	☐	TIPPER-AID, RIGHT	NC : 8200	66,35
DCA8110	☐	TIPPER-AID, LEFT	NC : 8210	66,35
DCA8120	☐	TIPPER-AID, BOTH SIDES	NC : 8200, 8210	132,70
CANE HOLDER				
DCA8200	☐	CANE HOLDER, RIGHT		66,35
DCA8210	☐	CANE HOLDER, LEFT		66,35
TRANSIT WHEELS				
DCA8300	☐	TRANSIT WHEELS 3"		301,42



7250



7500 / 7520



7810 / 7820 / 7830



8010



8100
8200



8300

OTHERS OPTIONS

DCA8450	☐	TABLE TRAY		391,47	 8550  8650  8850
DCA8500	☐	AIR PUMP - NEW		33,18	
DCA8550	☐	TOOL KIT		22,75	
DCA8650	☐	WHEELS BAG		65,40	
DCA8680	☐	24" MOUNTAIN BIKE WHEELS WITH BAG MOUNTAIN BIKE WHEELS AND TIRES WITH ALU GREY HANDRIMS ; QUICK RELEASE AXLES INCLUDED		646,45	
DCA8800	☐	BODYPOINT CALF STRAP		104,27	
DCA8850	☐	BODYPOINT MOBILITY BAG		104,27	

Table 1 :
HIGH OR STANDARD MOUNTED SITUATION (Foot plate, legrest angle and LLU)

	plastic footplate standard		Carbon footplate 2 parts				Carbon footplate 1 piece			
	Swing Away Frame		Fixed Frame		Swing Away Frame		Fixed Frame		Swing Away Frame	
	70°	80°	70°	80°	70°	80°	70°	80°	70°	80°
LLL 200	STD	STD	HIGH	HIGH	STD	STD	HIGH	HIGH	STD	STD
LLL 210	STD	STD	HIGH	HIGH	STD	STD	HIGH	HIGH	STD	STD
LLL 220	STD	STD	HIGH	HIGH	STD	STD	HIGH	HIGH	STD	STD
LLL 230	STD	STD	HIGH	HIGH	STD	STD	HIGH	HIGH	STD	STD
LLL 240	STD	STD	HIGH	HIGH	STD	STD	HIGH	HIGH	STD	STD
LLL 250	STD	STD	HIGH	HIGH	STD	STD	HIGH	HIGH	STD	STD
LLL 260	STD	STD	HIGH	HIGH	STD	STD	HIGH	HIGH	STD	STD
LLL 270	STD	STD	HIGH	HIGH	STD	STD	HIGH	HIGH	STD	STD
LLL 280	STD	STD	HIGH	HIGH	STD	STD	HIGH	HIGH	STD	STD
LLL 290	STD	STD	HIGH	HIGH	STD	STD	HIGH	HIGH	STD	STD
LLL 300	STD	STD	HIGH	HIGH	STD	STD	HIGH	HIGH	STD	STD
LLL 310	STD	STD	HIGH	HIGH	STD	STD	HIGH	HIGH	STD	STD
LLL 320	STD	STD	HIGH	HIGH	STD	STD	HIGH	HIGH	STD	STD
LLL 330	STD	STD	HIGH	HIGH	STD	STD	HIGH	HIGH	STD	STD
LLL 340	STD	STD	HIGH	HIGH	STD	STD	HIGH	HIGH	STD	STD
LLL 350	STD	STD	HIGH	HIGH	STD	STD	HIGH	HIGH	STD	STD
LLL 360	STD	STD	HIGH	HIGH	STD	STD	HIGH	HIGH	STD	STD
LLL 370	STD	STD	HIGH	HIGH	STD	STD	HIGH	HIGH	STD	STD
LLL 380	STD	STD	HIGH	HIGH	STD	STD	HIGH	HIGH	STD	STD
LLL 390	STD	STD	HIGH	HIGH	STD	STD	HIGH	HIGH	STD	STD
LLL 400	STD	STD	HIGH	HIGH	STD	STD	HIGH	HIGH	STD	STD
LLL 410	STD	STD	HIGH	HIGH	STD	STD	HIGH	HIGH	STD	STD
LLL 420	STD	STD	HIGH	HIGH	STD	STD	HIGH	HIGH	STD	STD
LLL 430	STD	STD	HIGH	HIGH	STD	STD	HIGH	HIGH	STD	STD
LLL 440	STD	STD	HIGH	HIGH	STD	STD	HIGH	HIGH	STD	STD
LLL 450	STD	STD	HIGH	HIGH	STD	STD	HIGH	HIGH	STD	STD
LLL 460	STD	STD	HIGH	HIGH	STD	STD	HIGH	HIGH	STD	STD
LLL 470	STD	STD	HIGH	HIGH	STD	STD	HIGH	HIGH	STD	STD
LLL 480	STD	STD	HIGH	HIGH	STD	STD	HIGH	HIGH	STD	STD
LLL 490	STD	STD	HIGH	HIGH	STD	STD	HIGH	HIGH	STD	STD
LLL 500	STD	STD	HIGH	HIGH	STD	STD	HIGH	HIGH	STD	STD
LLL 510	STD	STD	HIGH	HIGH	STD	STD	HIGH	HIGH	STD	STD

Table 2 :
Brakes in relationship with rear wheel position, rear wheel size and SD

	REAR WHEEL POS.	22" REAR WHEEL						24" REAR WHEEL						25" REAR WHEEL						26" REAR WHEEL					
		1	2	3	4	5	6	1	2	3	4	5	6	1	2	3	4	5	6	1	2	3	4	5	6
push brakes	SD 320-380																								
	SD 400																								
	SD 420																								
	SD 440																								
	SD 460																								
pull brakes	SD 480																								
	SD 500																								
	SD 320-380																								
	SD 400																								
	SD 420																								
performance brakes	SD 440																								
	SD 460																								
	SD 480																								
	SD 500																								
active brakes	SD 320-380																								
	SD 400																								
	SD 420																								
	SD 440																								
	SD 460																								
light active brakes	SD 480																								
	SD 500																								
	SD 320-380																								
	SD 400																								
	SD 420																								

NC
NC with Swing away frames

Table 3 :
SWING AWAY HEMI FRAME CONDITION

Front wheel	4"	5"	6"	7"
FSTF 370	Hemi*			
FSTF 380	Hemi	Hemi*		
FSTF 390	Hemi	Hemi*		
FSTF 400	Hemi	Hemi*	Hemi*	
FSTF 410	Std*	Hemi	Hemi*	
FSTF 420	Std*	Hemi	Hemi*	
FSTF 430	Std*	Hemi	Hemi	
FSTF 440	Std	Std*	Hemi	
FSTF 450	Std	Std*	Hemi	
FSTF 460	Std	Std*	Std*	Hemi
FSTF 470	Std	Std	Std*	Hemi
FSTF 480	Std	Std	Std*	Hemi
FSTF 490	Std	Std	Std	Std*
FSTF 500	Std	Std	Std	Std*
FSTF 510	Std	Std	Std	Std*
FSTF 520	Std	Std	Std	Std
FSTF 530	Std	Std	Std	Std

NC
FSTF NOT IN THE EPF
WITH HOUSING CONNECTOR PLATE




Table 4 :
MATRIX BACKREST WIDTH AVAILABLE RANGE RANGE

SEAT WIDTH (SW)	ELITE E2 / ELITE DEEP E2 MATRIX BACKREST WIDTH CHOICE						
	360	380	410	430	460	480	510
280							
300							
320							
340							
360							
380							
400							
420							
440							
460							
480							
500							

SEAT WIDTH (SW)	MX2 MATRIX BACKREST WIDTH CHOICE						
	360	380	410	430	460	480	510
280							
300							
320							
340							
360							
380							
400							
420							
440							
460							
480							
500							

Accessory Ramps

1440530	☐	Fiberglass Folding Mono Ramp Folded dimensions: 2000 x 230 mm 16.3 kg - maximum capacity: 300 kg		1 787,68	 1440530
1654128	☐	Basic ramp 2 telescopic parts in aluminum and non-slip coating Folded dimensions: 1210 x 194 mm 8.8 kg per pair - maximum capacity: 270 kg		402,84	 1654128
1654129	☐	Quick2go Compact 3-part aluminum telescopic ramp with non-slip coating Folded dimensions: 870 x 237 mm 13.6 kg per pair - max. : 260 kg		805,69	 1654129
1654130	☐	Aluminium folding mono ramp Folded dimensions: 1824 x 714 mm 13.9 kg - maximum capacity: 300 kg		426,54	 1654130

*Στις παραπάνω τιμές δεν συμπεριλαμβάνεται ο ΦΠΑ13%