

Levitop Comfort - Mobile patient lift

Product short description:

The **Levitop Comfort** lift makes it easier to lift and transport patients with disabilities and reduced mobility.



Product gallery:



Product description:

Levitop Comfort - Mobile patient lift

The Levitop Comfort Lift is a mobile lift designed to assist in lifting and transporting people with reduced mobility due to illness or dysfunction. It can be used to lift patients from a rehabilitation bed, bathtub, shower tray, or even directly from the floor with a suitable suspension.

The battery has a large capacity, allowing for approximately 20 full lifting and lowering cycles on a single charge with a maximum load of 230 kg.

Features:

The lift is equipped with a standard four-point sling.

The sling rotates 360° with the option to lock in four positions (90°).

Wheels allow for maneuvering even in small and narrow spaces.

Rear wheels are equipped with a brake.

A remote control is included.

The width of the leg span can be adjusted with the remote control.

An emergency lowering mechanism enhances patient safety.

Scale - optional.

VertiUp eFlow drive system

The system consists of two main parts: the drive system and the steering handle. The drive wheel, located between the two castors at the bottom of the lift, is powered by LINAK's lithium-ion battery and allows the lift to move in any direction with minimal effort. Control of the lift is intuitive and independent of the weight of the patient or the lift itself. The sensor integrated into the steering handle ensures smooth control of the lift without increasing the operator's effort, even in tight spaces, while removing pressure from the handle causes the device to stop via the electromagnetic brake built into the drive wheel.

Drive system

Steering handle

Levitop Comfort will be an excellent choice for:

Hospital

Developed for use in hospitals, the innovative Levitop Comfort mobile patient lift is an indispensable tool to support the patient care process. The use of this advanced device brings numerous benefits to both medical personnel and patients.

Care home

The Levitop Comfort mobile patient lift is an invaluable aid for medical staff and patients in care home environments. With its advanced features and ergonomic design, it enables the effective and safe movement of individuals with various health needs during rehabilitation therapy

Private usage

Levitop Comfort - Mobile patient lift is an indispensable tool for home use, providing support and safety during daily activities. Its use enables the effective and comfortable lifting and moving of individuals with limited mobility.

Sanatorium

The Levitop Comfort mobile patient lift is an invaluable aid for medical staff and patients in the sanatorium environment. With its advanced features and ergonomic design, it enables the effective and safe movement of individuals with various health needs during rehabilitation therapy.

Advantages of Levitop Standard

Sling strap.

Rotatable 360°.

Strong construction.

Frame made of durable metal.

Robust suspension attachment.

For lifting arms.

Remote controller.

Ergonomic and easy to use.

Braked wheels.

Brakes in the rear wheels prevent the unit from moving during use.

Suspension lock.

Blocking the rotation of the beam.

Compatible slings:

Universal

Head support: Yes Material: Polyester Sizes: M, L, XL, Bariatric Lifting capacity: 318 kg / 454 kg Attachment: Loops

Bath & Toilet

Head support: YesMaterial: Polyester meshSizes: M, L, XLLifting capacity: 250 kgAttachment: Loops

Bath & Toilet

Head support: NoMaterial: Polyester meshSizes: M, L, XLLifting capacity: 250 kgAttachment: Clips

Technical data:

Base legs spreading:	electrical
Max. safe working load (SWL)[kg]:	230
Length [cm]:	130
Width [cm]:	67
Height [cm]:	135
Lifting range [cm]:	68 - 196
Min. turning diameter [cm]:	144
Min. base legs width (int./ext.) [cm]:	52 / 67
Max. base legs width (int./ext.) [cm]:	103 / 118
Base legs height [cm]:	11,5
Wheels diameter (with / without lock) [cm]:	10 / 12,5
Weight [kg]:	66
Battery capacity [Ah]:	2,9
Charger (integrated with control box) [V/Hz]:	100-240 / 50 / 60
Remote control functions:	• up / down • el. leg spread: outside / inside • overload indicator • battery charge status • service indicator
Lifting speed (raising / lowering) [s]:	36 / 34
Packed dimensions (L x W x H) [cm]:	150 x 80 x 155

Video presentation:

Product codes:

Reference: 01-WPP-0110-001

EAN13: -

UPC: -

Product features:

Product attributes: