


PROSTHETIC LINER SIZING CHARTS

HOW TO SIZE

Measurements taken
6cm from the distal end



T: 727.528.8566 | 800.574.5426 | F: 727.528.8862

www.easyliner.com | infous@easyliner.com

2895 42nd Avenue North St. Petersburg, Florida 33714

CAN #10248 Form #100770 Rev H

PART #	DESCRIPTION	THICKNESS (MM)	SIZE									
AKDT AKDTHD AKFR AKFRHD	Extreme Liner	3 or 6	16	20	24	26	28	32	35	38	44	
APDT-TP APFR-TP	Anterior Posterior Tapered Liner	6/3		20	24	26	28	32				
ECDTHD ECFRHD	ECO Liner	3 or 6		20	24	26	28	32	35	38	44	
ELDT ELFR	EasyLiner	3, 6, 6/3* or 9*	16	20	24*	26	28*	32	35*	38	44	
GPDT GPFR	General Purpose	3 or 6	16	20	24	26	28	32	35	38	44	
SPDT SPDTHD SPFR SPFRHD	Superior Performance	3 or 6	16	20	24	26	28	32	35	38	44	
SADT SAFR SBFR	Smart Seal SADT	3 or 6		20	24	26	28	32		38		
	Smart Seal SAFR	3 or 6				26	28	32	35	38		
	Smart Seal SBFR	3 or 6		20	24	26	28					
EZF	EasyFlex	3 or 6	16	20	24	26	28	32	35	38	44	
VSOT	Liberty Liner	6		20	24	26	28	32		38		
PSI	Proximal Seal Interface	6 or 6/3	16	20	24	26	28	32		38	44	
PFGFR PFGDT	ProFit Liner	3 or 6	16	20	24	26	28	32	35	38	44	
USHGL USHGC USLGL USLGC USUGL USUGC	UltraSeal Gel Liner	3		20	24	26	28	32	35	38	44*	
CIRCUMFERENCE (CM)			16-19	20-23	24-25	26-27	28-31	32-34	35-37	38-44	44-53	

*EasyLiner 9mm only available in sizes 24 and 28; 3mm and 6mm only available for size 35

*UltraSeal size 44 circumference ranges 44-47cm

PART #	DESCRIPTION	THICKNESS (MM)	SIZE										
LSC LSL	Silicone Pro	2	20	22	24	26	28	30	32	35	38	41	44
SSGSC SSGSL	SoftSil®	2	20	22	24	26	28	30	32	35	38	41	44
USHSL USHSC USLSL USLSC USUSL USUSC	UltraSeal Silicone	2	20	22	24	26	28	30	32	35	38	41	44
CIRCUMFERENCE (CM)			20-22	22-24	24-26	26-28	28-30	30-32	32-35	35-38	38-41	41-44	44-47


PART #	DESCRIPTION	THICKNESS (MM)	SIZE									
ELPX	EasyLiner SuperStretch	6	10	16	20	24	26	28	32	38		
CIRCUMFERENCE (CM)			10-15	16-19	20-23	24-25	26-27	28-31	32-37	38-43		

PROSTHETIC SLEEVES

SIZING CHART

PART #	DESCRIPTION	THICKNESS (MM)	SIZE							
SFB	Superior Suspension	2		20	24	26	28	32		
SFP	Superior Performance	3 or 6	16	20	24	26	28	32	38	44
SFK	Extreme	3 or 6	16	20	24	26	28	32	38	44
SFR	EasySleeve	3 or 6	16	20	24	26	28	32	38	44
SLE	Easyleeve SuperStretch	6	16	20	24	26	28	32		
SFS	EasySleeve Fabric Reinforced	3 or 6	16	20	24	26	28	32	38	44
RFS	Reinforced Straight Sleeve	3		20	24	26	28	32	38	44
CIRCUMFERENCE (CM)			20-24	25-30	31-34	35-37	38-40	41-45	46-60	61-70
PART #	DESCRIPTION	THICKNESS (MM)	SIZE							
SFX	Flex Sleeve	3		20	23	26	29	32		
RFX	Reinforced Flex Sleeve	3		20	23	26	29	32		
CIRCUMFERENCE (CM)				34-37	36-39	38-42	41-46	45-50		

Measurements taken 10cm above the Patella



T: 727.528.8566 | 800.574.5426 | F: 727.528.8862

www.easyliner.com | infous@easyliner.com

2895 42nd Avenue North St. Petersburg, Florida 33714

