



53 CD 2 NE 19 25 ED

NEMO ATIVA GROWTHABLE 19/22/25 (550 spokes air/150 solid)

NEMO

Color - structure	pink, red, blue, white
Color - upholstery	black
Front wheels - solid	Ø150/Ø200mm
Solid / pneumatic rear wheels	Ø550mm
Rear wheels	with spokes
Adjustable height and reclining backrest	15°
Backrest position	85.5°/90.5°/96°/101°
Footrest adjustment	-10° / 10°

NEMO, aluminum active paediatric growth wheelchair. Adjustable tension back and seat, with padded seat and backrest for maximum comfort. Seat height adjustable. Backrest height and angle adjustable. Height adjustable handles. Side panels. Unique footplate height, length and angle adjustable. 550mm pneumatic rear wheels and 100mm solid front wheels. Quick release rear wheels and foldable backrest for easy transport. See accessories.



MEASURES

in cm / kg - tolerance +/- 1

MEASURES													
in cm / kg - tolerance +/- 1													
													
NEMO 19 25	19-25	39-45	26	69.5	44-47	22-32	69	65	22-34	85.5°-101°	100	10.6	
NEMO 25 31	25-31	45-51	26	74.5	44-47	22-32	69	65	22-34	85.5°-101°	100	10.6	
NEMO 31 37	31-37	51-57	26	82	46-50	22-32	90.5	75	33-43	85.5°-101°	100	10.6	