





Aqua Whirl WKG | Whirlpool bath tub for upper limbs



Tubs are perfect for hydro massage treatment of limbs and post-injury conditions, in nervous system disorders or muscular and nervous system fatigue.

Standard accessories:

- Swivel chair



Features:

- Basin available in two colours: white or calypso green
- Basin made of high quality acrylic (fibreglass mat reinforced)
- Economical use of water – min. 30 l
- Effective hydro massage therapy via 44 nozzles
- „Dry-run“ pump protection system
- Adjustable air supply system for the water nozzles
- Electronically controlled descaling system
- Timer: 1-9 h:59 min
- Electronic control panel
- Short filling and draining time
- Big diameter drain and filling system
- Semi-automatic draining system
- Temperature sensor
- Water level sensor
- Shower
- Smooth, profiled treatment chamber
- Ergonomic design
- Automatic filling (option)

Technical data:

Treatment capacity [l]:	min. 30
Volume up to overflow [l]:	45
Filling/drainage [min.]:	1,3/0,7
Number of aeration nozzles:	44
Power consumption [A]:	5
Supply voltage [V/Hz]:	230/50
Dimensions (L x W x H) [cm]:	88 x 93 x 90
Weight [kg]:	60

