



Aquai | Device for dry hydromassage



AQUAI has wide range of uses: it is used in the wellness and fitness areas. It perfectly suits rehabilitation centers. It is also a very good option for relaxation rooms in companies.

A patient in a rehabilitation center or a guest in a SPA center can experience the beneficial power of hydrojet massage without getting wet. The massage is done through a layer of plastic coat.

 Packing dimensions: **259 x 139 x 105 cm** Weight: **+100 kg**

 Length: **237 cm**

 Height: **70 cm**

 Width: **114 cm**

 Weight: **190/410 kg**

 Capacity: **220 l**

Advantages of the treatment:

- Improvement in blood circulation and metabolism
- Increase in oxygen transport in muscles
- Improvement in movement range of extremities
- Tension relief
- Stress reduction and relaxation
- Well-being


Features:

- 6 programmable and independent zones
- Jets can work simultaneously or in zones
- Dry water massage nozzles with increased durability
- 6 standard programs
- 1 user program
- Membrane surface temperature setting: temperature adjustable from 30 to 40°C
- Adjustable treatment time (1-30 min)
- Delayed start at programmed number of hours
- Adjustable frequency and sequence of zones
- Integrated heating system
- Internal water cooling system maintains water temperature during preset hours regardless of the number of treatments
- Air vent located on the diaphragm as part of deaeration system
- No need to connect to water tap
- Easily accessible modern colour touch panel
- LED lighting at the base

Technical data:

Dimensions (L x W x H) [cm]:	237 x 114 x 70
Capacity [l]:	220
Weight (empty/filled with water) [kg]:	190/410
Max. load [kg]:	135
Power supply [V/Hz]:	230/50 or 400/50

Massage zones:



Additional options:

AUDIO: sound input from an external device such as Hi-Fi, mp3 etc.

CARD MANAGER: a system supporting the sale of treatment sessions using stored-value cards; together with top-up software.

