

Customer ID: \_\_\_\_\_

# ORDER FORM

 ORDER  QUOTE 

## ICON 60

LIGHTWEIGHT, RIGID FRAME  
WHEELCHAIR

### Base model

- Max user weight - 120kg
- Series 7000 aluminium rigid active wheelchair
- Multiple adjustments, seat depth, front and rear seat heights and centre of gravity
- Angle adjustable front forks
- 24" Daylight rear wheel with back anodised hand rim, QR axle and pneumatic tyre
- 4" Alu core castor
- Lockdown folding, angle, height and tension adjustable backrest with 3D breathable fabric and seat cushion.
- Low profile scissor lock brakes
- Aluminium side guards
- Calf strap and under seat zipped pouch included



SUITABLE FOR USE WITH POWERED  
ASSISTED WHEELCHAIR SYSTEMS




**Help to save paper and use the interactive PDF.**  
Please order by e-mail to [info@rehasense.com](mailto:info@rehasense.com). Thank's a lot.

|                 |                  |
|-----------------|------------------|
| Name:           | Date:            |
| Customer ID:    | Order number:    |
| Contact person: | Delivery address |
| E-Mail:         | Str. & house nr: |
| Tel./Fax:       | City & postcode: |

### PLEASE MARK TO SELECT AN OPTION


Hereby I consent to the processing of my personal data for the purpose of the order fulfilment for a rehabilitation equipment in accordance with Article 4 Paragraph 11 of the Regulation (EU) 2016/679 of the European Parliament and of the Council of 27 April 2016 on the protection of natural persons with regard to the processing of personal data and on the free movement of such data, and repealing Directive 95/46/EC (General Data Protection Regulation)

I am aware of the fact that all dimensions can be provided within the dimensional tolerance of  $\pm 10$  mm. A dimensional deviation within this tolerance range shall not constitute a reason to claim the wheelchair and deem it non-compliant with this contract.

| SEAT WIDTH            |             |                   |                         |       |   |
|-----------------------|-------------|-------------------|-------------------------|-------|---|
| Tick to order         | Part number | Description       | Standard or cost option | Price |   |
| <input type="radio"/> | SW34        | Seat width 340 mm | Standard                |       | Seat width range fixed from 340 mm to 480 mm<br> |
| <input type="radio"/> | SW36        | Seat width 360 mm | Standard                |       |   |
| <input type="radio"/> | SW38        | Seat width 380 mm | Standard                |       |   |
| <input type="radio"/> | SW40        | Seat width 400 mm | Standard                |       |   |
| <input type="radio"/> | SW42        | Seat width 420 mm | Standard                |       |   |
| <input type="radio"/> | SW44        | Seat width 440 mm | Standard                |       |   |
| <input type="radio"/> | SW46        | Seat width 460 mm | Standard                |       |   |
| <input type="radio"/> | SW48        | Seat width 480 mm | Standard                |       |   |

| FRAME SEAT DEPTH      |        |                   |          |  |   |
|-----------------------|--------|-------------------|----------|--|---|
| <input type="radio"/> | SD40   | Seat depth 400 mm | Standard |  | Seat depth range adjustable between 400 mm - 490 mm<br> |
| <input type="radio"/> | SD41.5 | Seat depth 415 mm | Standard |  |   |
| <input type="radio"/> | SD43   | Seat depth 430 mm | Standard |  |   |
| <input type="radio"/> | SD44.5 | Seat depth 445 mm | Standard |  |   |
| <input type="radio"/> | SD46   | Seat depth 460 mm | Standard |  |   |
| <input type="radio"/> | SD47.5 | Seat depth 475 mm | Standard |  |   |
| <input type="radio"/> | SD49   | Seat depth 490 mm | Standard |  |   |

| CUSHION SEAT DEPTH    |         |                   |          |  |  |
|-----------------------|---------|-------------------|----------|--|--|
| <input type="radio"/> | CSD40   | Seat depth 400 mm | Standard |  | Seat depth range adjustable between 400 mm - 490 mm<br> |
| <input type="radio"/> | CSD41.5 | Seat depth 415 mm | Standard |  |  |
| <input type="radio"/> | CSD43   | Seat depth 430 mm | Standard |  |  |
| <input type="radio"/> | CSD44.5 | Seat depth 445 mm | Standard |  |  |
| <input type="radio"/> | CSD46   | Seat depth 460 mm | Standard |  |  |
| <input type="radio"/> | CSD47.5 | Seat depth 475 mm | Standard |  |  |
| <input type="radio"/> | CSD49   | Seat depth 490 mm | Standard |  |  |

| FRONT SEAT HEIGHT     |      |   |          |  |  |
|-----------------------|------|---|----------|--|--|
|                       | FSHA | Seat height front 430 mm (4" castor wheel only / std. fork) | Option   |  | Front seat height range adjustable from 430 mm - 560 mm<br> |
|                       | FSHB | Seat height front 445 mm (4" castor wheel only / std. fork) | Option   |  |  |
|                       | FSHC | Seat height front 455 mm (5" castor wheel only / std. fork) | Option   |  |  |
|                       | FSHD | Seat height front 455 mm (4" castor wheel only / jazz fork) | Option   |  |  |
|                       | FSHE | Seat height front 470 mm (4" castor wheel only / jazz fork) | Option   |  |  |
|                       | FSHF | Seat height front 480 mm (5" castor wheel only / jazz fork) | Option   |  |  |
| <input type="radio"/> | FSH1 | Seat height front 495 mm (4" castor wheel only / std. fork) | Standard |  |  |
| <input type="radio"/> | FSH2 | Seat height front 505 mm (4" & 5" castor wheel / std. fork) | Standard |  |  |
| <input type="radio"/> | FSH3 | Seat height front 515 mm (4" & 5" castor wheel / std. fork) | Standard |  |  |
| <input type="radio"/> | FSH4 | Seat height front 525 mm (5" castor wheel only / std. fork) | Option   |  |  |
| <input type="radio"/> | FSH5 | Seat height front 530 mm (4" castor wheel only / jazz fork) | Option   |  |  |
| <input type="radio"/> | FSH6 | Seat height front 545 mm (4" castor wheel only / jazz fork) | Option   |  |  |
| <input type="radio"/> | FSH7 | Seat height front 560 mm (5" castor wheel only / jazz fork) | Option   |  |  |



### BUSHING TYPE

|                       |     |              |          |  |
|-----------------------|-----|--------------|----------|--|
| <input type="radio"/> | BT1 | Inch 1/2"    | Standard |  |
| <input type="radio"/> | BT2 | Metric 12 mm | Option   |  |

### BACKREST ANGLE

|                       |      |   |          |  |
|-----------------------|------|---|----------|--|
| <input type="radio"/> | BRA1 | Forward (from +1° to +10°)<br><small>enter the value</small>  | Option   |  |
| <input type="radio"/> | BRA0 | Vertical backrest (0°)  | Standard |  |
| <input type="radio"/> | BRA2 | Backward (from -1° to -10°)<br><small>enter the value</small> | Option   |  |

### FRAME COLOUR

|                       |     |             |          |  |
|-----------------------|-----|-------------|----------|--|
| <input type="radio"/> | FC1 | White       | Option   |  |
| <input type="radio"/> | FC2 | Black Matte | Standard |  |
| <input type="radio"/> | FC3 | Silver      | Option   |  |
| <input type="radio"/> | FC4 | Gold        | Option   |  |
| <input type="radio"/> | FC5 | Blue        | Option   |  |
| <input type="radio"/> | FC6 | Red         | Option   |  |
| <input type="radio"/> | FC7 | Green       | Option   |  |

### SEAT CUSHION

|                       |     |                       |          |  |
|-----------------------|-----|-----------------------|----------|--|
| <input type="radio"/> | SC1 | Standard seat cushion | Standard |  |
| <input type="radio"/> | SC0 | Without seat cushion  | Option   |  |

### BACKREST OPTIONS

|                       |             |  |          |  |
|-----------------------|-------------|--|----------|--|
| <input type="radio"/> | BRO1        | Height adjustable backrest without push handles  | Standard |  |
| <input type="radio"/> | BRO2        | Height adjustable backrest fitted with height adjustable push handles                              | Option   |  |
| <input type="radio"/> | BRO3 i BRO5 | Height adjustable backrest fitted with ergonomic backrest frame                                    | Option   |  |
| <input type="radio"/> | BRO4 i BRO6 | Height adjustable backrest fitted with ergonomic backrest frame and height adjustable push handles | Option   |  |

|   |     |  |        |  |
|---|-----|--|--------|--|
| <b>BACKREST AD-ONS (available with BRO1 - BRO4)</b> |     |  |        |  |
| <input type="radio"/>                               | EVI | Big radius, expandable in length and width direction, up to the narrowing of extension elements (from -30 to +30 mm) | Option | <br><b>BRO1</b><br>Height adjustable backrest without push handles |





|   |  |  |   |
|---|--|--|---|
| <br>Option with adjustable handles in the backrest frame. | <br>Option with height adjustable backrest fitted with ergonomic backrest frame. | <br>Option with adjustable handles in the backrest frame. Crossbar curved, two radius. | <br>Big radius, expandable in length and width direction, up to the narrowing of extension elements (from -30 to +30 mm). |
|---|--|--|---|





### BACKREST HEIGHT (mm) (available with BRO1 - BRO4)










|                                    |                           |                           |                           |                           |                           |
|------------------------------------|---------------------------|---------------------------|---------------------------|---------------------------|---------------------------|
| <input type="radio"/> 320 standard |                           |                           |                           |                           |                           |
| <input type="radio"/> 330          | <input type="radio"/> 340 | <input type="radio"/> 350 | <input type="radio"/> 360 | <input type="radio"/> 370 | <input type="radio"/> 380 |


### BACKREST HEIGHT (mm) (available with BRO5 - BRO6)

|                                    |                           |                           |                           |                           |                           |
|------------------------------------|---------------------------|---------------------------|---------------------------|---------------------------|---------------------------|
| <input type="radio"/> 390 standard |                           |                           |                           |                           |                           |
| <input type="radio"/> 400          | <input type="radio"/> 410 | <input type="radio"/> 420 | <input type="radio"/> 430 | <input type="radio"/> 440 | <input type="radio"/> 450 |


| BACKREST OPTIONS (not available with EV1) |      |  |        |  |   |
|---|------|--|--------|--|---|
| <input type="radio"/>                     | BRH1 | Backrest with foldable push handles, single bended crossbar, frame height <b>250mm + 60 mm with handle (310 mm)</b>  | Option |  |   <p>Backrest height non-extendable, for folding handles in the backrest frame.</p>                               |
| <input type="radio"/>                     | BRH2 | Backrest with foldable push handles, single bended crossbar, frame height <b>300 mm + 60 mm with handle (360 mm)</b> | Option |  |   |
| <input type="radio"/>                     | BRH3 | Backrest with foldable push handles, single bended crossbar, frame height <b>350 mm + 60 mm with handle (410 mm)</b> | Option |  |   |
| <input type="radio"/>                     | BRH4 | Backrest with foldable push handles, double bended crossbar, frame height <b>250 mm + 60 mm with handle (310 mm)</b> | Option |  |   <p>Backrest height non-extendable, for folding handles in the backrest frame. Crossbar curved, two radiuses</p> |
| <input type="radio"/>                     | BRH5 | Backrest with foldable push handles, double bended crossbar, frame height <b>300 mm + 60 mm with handle (360 mm)</b> | Option |  |   |
| <input type="radio"/>                     | BRH6 | Backrest with foldable push handles, double bended crossbar, frame height <b>350 mm + 60 mm with handle (410 mm)</b> | Option |  |   |
| <input type="radio"/>                     | BRH7 | Backrest with foldable push handles, double bended crossbar, frame height <b>400 mm + 60 mm with handle (460 mm)</b> | Option |  |   |


| FOOTREST              |     |  |          |  |  |
|-----------------------|-----|--|----------|--|--|
| <input type="radio"/> | FR1 | Standard footrest                                  | Standard |  |    <p>FR1      FR2      FR3</p> |
| <input type="radio"/> | FR2 | Reinforced footrest for track wheel                | Option   |  |  |
| <input type="radio"/> | FR3 | Elevated footrest                                  | Option   |  |  |
| FOOTREST OPTIONS      |     |  |          |  |  |
| <input type="radio"/> | FC1 | Footrest cover aluminium (available for FR1 & FR3) | Option   |  |   |


| SIDE PANELS           |     |   |          |  |   |
|-----------------------|-----|---|----------|--|---|
| <input type="radio"/> | SP1 | Aluminium side panels, flat, swing away                                 | Standard |  |  |
| <input type="radio"/> | SP2 | Carbon side panels, flat, swing away                                    | Option   |  |  |
| <input type="radio"/> | SP3 | Aluminium fenders with radius/removable black structured gloss          | Option   |  |  |
| <input type="radio"/> | SP4 | Aluminium fender with radius/removable black structured matte           | Option   |  |  |
| <input type="radio"/> | SP5 | Carbon fenders with radius/removable                                    | Option   |  |  |
| <input type="radio"/> | SP6 | Aluminium fenders with radius drop-in/removable, black structured gloss | Option   |  |  |
| <input type="radio"/> | SP7 | Aluminium fenders with radius drop-in/removable, black structured matte | Option   |  |  |
| <input type="radio"/> | SP8 | Carbon fenders with radius drop-in/removable                            | Option   |  |  |
| <input type="radio"/> | SP9 | Aluminium side panels with armrest cushion                              | Option   |  |  |

| BRAKES                |    |   |          |  |   |
|-----------------------|----|---|----------|--|---|
| <input type="radio"/> | B1 | Horizontal (left, right) (only with long clamp BC2) | Option   |  |  |
| <input type="radio"/> | B2 | Horizontal Premium Lock-P (L+R)                     | Standard |  |   |
| <input type="radio"/> | B3 | Horizontal Lock-A in aluminium                      | Option   |  |   |
| <input type="radio"/> | B4 | Vertical - push to lock (only with long clamp BC2)  | Option   |  |   |




| BRAKE MOUNTING CLAMPS |     |                                   |          |  |  |
|-----------------------|-----|-----------------------------------|----------|--|--|
| <input type="radio"/> | BC1 | Standard short clamp (only B2/B3) | Standard |  |  <p><b>BC1</b> <b>BC2</b></p> |
| <input type="radio"/> | BC2 | Long clamp (all brakes)           | Option   |  |  |

| FORKS                 |    |                |          |  |   |
|-----------------------|----|----------------|----------|--|---|
| <input type="radio"/> | F1 | Standard Light | Standard |  |  <p><b>F1</b> <b>F2</b></p> |
| <input type="radio"/> | F2 | Jazz-Shock     | Option   |  |   |

| FRONT WHEELS          |      |  |          |  |  |
|-----------------------|------|--|----------|--|--|
| <input type="radio"/> | FW1  | Alucore - 4" 100x34 mm rim anodised, Silver coloured | Standard |  |  <p><b>FW1</b><br/>Alucore - 4" 100x34 mm rim anodised, Silver coloured</p> |
| <input type="radio"/> | FW2  | Alucore - 4" 100x34 mm rim anodised, black coloured  | Standard |  |  |
| <input type="radio"/> | FW3  | Alucore - 5" 125x34 mm rim anodised, silver coloured | Option   |  |  |
| <input type="radio"/> | FW4  | Alucore - 5" 125x34 mm rim anodised, black coloured  | Option   |  |  |
| <input type="radio"/> | FW5  | Tristar - 4" 100x34 mm rim anodised, silver coloured | Option   |  |  |
| <input type="radio"/> | FW6  | Tristar - 4" 100x34 mm rim anodised, black coloured  | Option   |  |  |
| <input type="radio"/> | FW7  | Tristar - 5" 125x34 mm rim anodised, silver coloured | Option   |  |  |
| <input type="radio"/> | FW8  | Tristar - 5" 125x34 mm rim anodised, black coloured  | Option   |  |  |
| <input type="radio"/> | FW9  | Lotus fibercore - 4" 100x34 mm rim plastic, black    | Option   |  |  |
| <input type="radio"/> | FW10 | Lotus fibercore - 5" 125x34 mm rim plastic, black    | Option   |  |  |





## REAR WHEELS

|                       |      |   |          |   |
|-----------------------|------|---|----------|---|
| <input type="radio"/> | RW1  | DAYLITE - 24" x 1" wheel<br>(Pushrim: PR1, PR2, PR3, PR4, PR5, PR6, PR7, PR8 Tyre: T1, T2, T3)            | Standard |  <p><b>RW1</b><br/>DAYLITE 930S 24" wheel/2 pcs.</p> |
| <input type="radio"/> | RW2  | SKYLITE - 24" x 1" wheel<br>(Pushrim: PR1, PR2, Tyre: T1, T3)   | Option   |   |
| <input type="radio"/> | RW3  | TENSIO RAY 12S - 24" wheel<br>(Pushrim: PR1, PR2, PR3, PR4, PR5, PR6, PR7, PR8 Tyre: T1, T2, T3)          | Option   |   |
| <input type="radio"/> | RW4  | TENSIO RAY 12S - BLACK 24" wheel<br>(Pushrim: PR1, PR2, PR3, PR4, PR5, PR6, PR7, PR8 Tyre: T1, T2, T3)    | Option   |   |
| <input type="radio"/> | RW5  | HURRICANE 24" x 1" wheel<br>(Pushrim: PR1, PR2, PR3, PR4, PR5, PR6, PR7, PR8 Tyre: T1, T2, T3)            | Option   |   |
| <input type="radio"/> | RW6  | CYCLONE 24" wheel<br>(Pushrim: PR1, PR2, PR3, PR4, PR5, PR6, PR7, PR8, bushing type: BT Tyre: T1, T2, T3) | Option   |   |
| <input type="radio"/> | RW7  | SMARTLITE 24" wheel<br>(Pushrim: PR1, PR2 Tyre: T1)   | Option   |   |
| <input type="radio"/> | RW8  | OFF ROAD 24" x 2" buddy wheel<br>(Pushrim: PR1, PR7 Tyre: T4)   | Option   |   |
| <input type="radio"/> | RW9  | CHALLENGER 24" wheel<br>(Pushrim: PR1 Tyre: T5)   | Option   |   |
| <input type="radio"/> | RW10 | DAYLITE - 25" wheel<br>(Pushrim: PR1, PR5, PR8 Tyre: T3, T2, T3)  | Option   |   |
| <input type="radio"/> | RW11 | SKYLITE - 25" wheel<br>(Pushrim: PR1, PR5, PR8 Tyre: T3, T2, T3)  | Option   |   |



## PUSH RIMS


|                       |     |                                       |          |  |
|-----------------------|-----|---------------------------------------|----------|--|
| <input type="radio"/> | PR1 | Anodised aluminium, Black - coloured  | Standard |  <p><b>PR1</b><br/>Anodised aluminium,<br/>Black - coloured</p> |
| <input type="radio"/> | PR2 | Anodised aluminium, Silver - coloured | Option   |  |
| <input type="radio"/> | PR3 | Steel                                 | Option   |  |
| <input type="radio"/> | PR4 | Titanium                              | Option   |  |
| <input type="radio"/> | PR5 | Tetra Optimum F                       | Option   |  |
| <input type="radio"/> | PR6 | Tetra silicone, Black                 | Option   |  |
| <input type="radio"/> | PR7 | Nova H                                | Option   |  |



| TYRES                 |    |  |          |  |
|-----------------------|----|--|----------|--|
| <input type="radio"/> | T1 | Speedlite 24" x 1"                                   | Standard |  |
| <input type="radio"/> | T2 | Gazelle 24" x 1"                                     | Option   |  |
| <input type="radio"/> | T3 | Trailblazer 24" x 1"                                 | Option   |  |
| <input type="radio"/> | T4 | Schwalbe Big Apple (only to be used with 24" Buddy)  | Option   |  |
| <input type="radio"/> | T5 | Innova IA-2003 (only to be used with 24" Challenger) | Option   |  |



| QR AXLES              |     |                                       |          |  |
|-----------------------|-----|---------------------------------------|----------|--|
| <input type="radio"/> | QR1 | Quick release axle / 2pcs             | Standard |  |
| <input type="radio"/> | QR2 | Quick release tetra handle axle / 2pc | Option   |  |



| OPTIONS               |    |   |        |  |
|-----------------------|----|---|--------|--|
| <input type="radio"/> | O1 | Anti-Tip Left<br>(Available alone: O1 or together with O2 or O4)  | Option |  |
| <input type="radio"/> | O2 | Anti-Tip Right<br>(Available alone: O2 or together with O1 or O3) | Option |  |
| <input type="radio"/> | O3 | Tip-on Left<br>(Available alone: O3 or together with O2 or O4)    | Option |  |
| <input type="radio"/> | O4 | Tip-on Right<br>(Available alone: O4 or together with O1 or O3)   | Option |  |



| DOUBLE REAR AXLE EXTENSION |     |   |        |  |
|----------------------------|-----|---|--------|--|
| <input type="radio"/>      | RAE | Double rear axle extension 6° bushing angle<br>(Use only with additional set of Challenger rear wheels, not available with Anti-tip and Tip-on) | Option |  |



| TRANSPORT CLAMPS (Suitable for passenger transport only in combination with transportation clamps) |     |                              |        |  |
|--|-----|------------------------------|--------|--|
| <input type="radio"/>  | TC1 | Transport clamps for vehicle | Option |  |



| TRACK WHEEL           |      |                                  |        |
|-----------------------|------|----------------------------------|--------|
| <input type="radio"/> | STWA | Single arm track wheel Aluminium | Option |
| <input type="radio"/> | DTWA | Double arm track wheel Aluminium | Option |
| <input type="radio"/> | STWC | Single arm track wheel Carbon    | Option |
| <input type="radio"/> | DTWC | Double arm track wheel Carbon    | Option |



Your notes

Technical data will vary according to the frame dimensions chosen and the way the wheelchair is set up (front and rear wheels). The data in the table refers to a seat width of 400 mm and a seat depth of 420 mm. All measurements are in millimetres (mm) for distance and kilograms (kg) for weight, unless otherwise stated.