



SUNRISE
MEDICAL

zippie

Edition 11.24
Valid till 11.25

Voyage

Deliver to: Name: Street: Town: Postal code: E-Mail: Tel: Fax:	CE: Invoice to: Name: Street: Town: Postal code: E-Mail: Tel: Fax:
---	---

Order date: Purchase order: Number required:	Marked for: Sunrise order no:
--	----------------------------------

D Order D Calculation D Built 4 me

Reimbursement code:

Art. no.: 99000601 (IC) / 99000600 (GER)
EIZ18 (USA)

• = Standard D = Cost option
0 = No cost option

Max. user weight: 34 kg

USER WEIGHT kg

Mandatory information to ensure correct configuration

Base Price 4300 €

BASE UNIT

199BS9	- Frame, silver - Foot Operated Wheel Lock 12" Rear Tires & 7" Front Casters - Front Caster Suspension & Swivel Locks - Quick Release Wheels - Infinitely variable tilt; range 35° (-5° to 30°) - Angle Adjustable Stroller Handle - Canopy with breathable window/mesh pockets - Fwd & Rwd facing seating capability - Under carriage storage (may be limited when some options chosen)	-	-
199KJ11	D Transit option	285 €	

TYRES - with quick release rear and front wheels, and front wheel suspension with caster swivel locks

199WQ07	12" Pneumatic wheels with airless inserts and 7" casters	-	
---------	--	---	--

SEATING SYSTEM

199KM01	- Comfort Seating - Seat width, seat depth, backheight and lower leg length adjustable - Angle adjustable footrest	-	-
199KM02	- Recline backrest from 80° - 180° - Canopy with mesh window	-	-

UPHOLSTERY

199KN01	- Seat upholstery with incontinence protection, removable, breathable, spacer fabric, black with dots 0 Blue trim on upholstery	-	-
199KN02	0 Red trim on upholstery	-	-
199KN03	0 Beige trim on upholstery	-	-
199KN04	0 Grey trim on upholstery	-	-

FOOTREST

199F329	Footboard, platform, composite, angle-adjustable, flip-up, heel loop, black	-	
---------	---	---	--

For more information
www.SunriseMedical.eu
 © Sunrise Medical 2025

ZI PPIE Voyage

TECHNICAL DATA

Max. user weight:	34 kg
Seat width:	240 - 300 mm, with insert up to 155 mm (excl. upholstery)
Seat depth:	200 - 330 mm (excl. upholstery)
Seat angle:	-5° to 30°
Seat height:	575 - 640 mm (3 steps possible 575/ 600 / 664 mm)
Footplate angle range:	50°
Lower leg length:	155 - 310 mm
Backrest height:	410 - 610 mm (excl. upholstery)
Backrest angle:	80° - 180°
Total width:	615 mm
Total length:	1100 mm
Total height:	810 - 1220 mm
Folded length:	800 mm
Folded height:	310 mm
Max. safe slope for use of parking brake:	8°
Rear wheel:	12"
Front wheel:	7"
Seat weight:	9,8 kg
Base weight:	12,7 kg
Removable parts:	Seat and wheels - for transport only
Colours:	4 colours
Crash tested:	Yes - tested in accordance with Annex A of ISO 7176-19 (2008)
Transport without user:	Fixe the product safe in the vehicle (e.g. car or plane)
Intended use and environment:	For children with limited mobility: indoor and outdoor use, operated by caregiver

ZIPPIE Voyage			
● = Standard ○ = No cost option D = Cost option			Price in €
BELLY BAR			
199KR01	D	Belly bar, padded, removable	€ 120
LATERALS			
199K001	D	Laterals, padded, spacer fabric, height and width-adjustable, small (105 x 150 mm, pair)	€ 330
199K002	D	Laterals, padded, spacer fabric, height and width-adjustable, large (110 x 200 mm, pair)	€ 330
HIP GUIDES			
199KP01	D	Hip guides, padded, spacer fabric, width-adjustable (pair)	€ 330
199KP03	D	Seat insert (Hip guides required, reduces SW up to 155 mm)	€ 205
HEAD SUPPORT			
199VJ01	D	Head support, padded, spacer fabric, height and width adjustable, small (pair)	€ 330
199VJ02	D	Head support, padded, spacer fabric, height and width adjustable, large (pair)	€ 330
POMMEL			
199KQ01	D	Pommel, padded, spacer fabric, square (piece)	€ 175
199KQ02	D	Pommel, padded, spacer fabric, triangular (piece)	€ 175
FOOTREST ACCESSORIES			
199A68	D	Shoe holders, small (150 mm, pair)	€ 240
199A69	D	Shoe holders, medium (200 mm, pair)	€ 240
199A79	D	Foot straps with covers, small (400 mm; pair; only in combination with shoe holders)	€ 56
199AB0	D	Foot straps with covers, medium (450 mm; pair; only in combination with shoe holders)	€ 56
JAY ANTERIOR TRUNK SUPPORTS			
Anterior Trunk Support - Style			
199JY01	D	Chest strap, lateral opening, dynamic	€ 240
199JY02	D	Chest strap, lateral opening, classic (not available with size: early intervention)	€ 220
199JY06	D	Chest strap, central opening, dynamic	€ 240
199JY05	D	Chest strap, central opening, classic (not available with size: early intervention)	€ 220
Anterior Trunk Support - Size			
199JZ05	0	Early Intervention (measurements see backpage)	-
199JZ06	0	Xtra-Small (measurements see backpage)	-
199JZ01	0	Small (measurements see backpage)	-
JAY ANTERIOR PELVIC BELTS			
			Must use when ordering Anterior Trunk Supports
Standard 2 Point Belt - Push Button, Small (25 mm)			
199ZZ10	D	no pad	€ 140
199ZZ11	D	with pad, small (110 mm)	€ 150
199ZZ12	D	with pad, medium (140 mm)	€ 150
Standard 2 Point Belt - Push Button Large (38 mm)			
199ZZ14	D	no pad	€ 140
199ZZ15	D	with pad, large (190 mm)	€ 150
V-Contour Pelvic 4 Point Belt - Push Button Varies			
199ZZ23	D	with pad, small, push button small (25 mm)	€ 205
199ZZ24	D	with pad, medium, push button large (38 mm)	€ 205

JAY ANTERIOR TRUNK SUPPORTS


Chest strap, lateral opening

	A	B
El	170 mm	250 mm
XS	190 mm	320 mm
S	220 mm	380 mm



Chest strap, central opening

	A	B
El	125 mm	250 mm
XS	150 mm	280 mm
S	170 mm	310 mm



JAY ANTERIOR PELVIC BELTS

Standard 2 Point Belt - Push Button 25 mm or 38 mm

199KJ11



V-Contour Pelvic 4 Point Belt
Push Button Varies



ZI PPIE Voyage				Price in €	
• = Standard		O = No cost option		D = Cost option	
ACCESSORIES				Below accessories not available with JCM Hi-Lo Base	
199LF1	D	Hardware for LTV ventilators, dovetail, vertical frame mounted	(limit 7 Kg)	485 €	
199LF11	D	Hardware for Trilogy ventilators, dovetail, vertical frame mounted	(limit 7 Kg)	485 €	
199LF2	D	Battery and vent/utility underframe tray, central mounted, incl. 1 utility hook	(limit 23 Kg)	520 €	
199LF3	D	Utility platform, front mounted	(limit 23 Kg)	505 €	
199LF4	D	Rear add-on accessory platform, incl. 1 hook	(remove before folding base; limit 20 kg)	250 €	
199A74	D	IV support pole	(limit 1 kg)	240 €	
199A713	D	Bottle holder, oxygen tank	(max. 4,5 kg)	240 €	
199A714	D	IV/O2 combination		345 €	
199LF5	D	Cup holder	(max. load 0,5 kg)	62€	
199LF7	D	Hook utility	(Attention: risk of tipping)	56 €	
199A462	D	Rain cover		95 €	
	D	Tray, push-fit, height-adjustable	(angle adjustable hardware included)		
199A91	D	Light Grey		550 €	
199A92	D	Clear Poly		690 €	
COMMENTS FOR SPECIAL BUILD:					
					
Οι τιμές περιλαμβάνουν 13% ΦΠΑ.					
Zippie Voyage - Base Price					
Option price					
Packaging and transport					
Exclusive VAT					
VAT.					
Total Price					
Sunrise Medical Customer Service Address: Kahlbachring 2-4, D-69254 Malsch/Heidelberg Tel. +497253 980-0 Fax. +497253 980-222				For more information visit our webpage www.SunriseMedical.eu © Sunrise Medical 2025	
				E-mail: export@sunrisemedical.de Intl-Mobility@sunmed.com Website: www.sunrisemedical.eu	

ACCESSORIES

HUMOM1J5{1127



Dovetail Vertical Vent Mount
(LTV/Trilogy)



Battery and VenU Utility
Underframe Tray



Front Mount Utility Platform



Rear Add On Accessory
Platform



IV Support Pole



O2 Tank Holder



Cup holder



Utility Hook



Upper Extremity Support
Tray



JCM Hi-Lo Base

199KJ11