

# Order form Sigvaris Roundknit Custom Made Stockings

## Leg measurements



Country / Distributor \_\_\_\_\_ Contact person \_\_\_\_\_

Patient name \_\_\_\_\_ Sex  female  male

Order No. \_\_\_\_\_ Your Reference No. \_\_\_\_\_

Re-order No. \_\_\_\_\_ Order date \_\_\_\_\_

Remarks \_\_\_\_\_

**SIGVARIS AG**  
 9014 St. Gallen  
 Switzerland  
 Tel. +41 71 272 40 00  
 Fax +41 71 272 40 01  
**custom-made@sigvaris.com**  
 www.sigvaris.com



Product code	New product name	Old product name
SSE	Style SEMITRSPARENT	MAGIC
ECF	Essential COMFORTABLE	COMFORT
ETH	Essential THERMOREGULATING	COTTON
ECL	Essential CLASSICAL	TOP FINE
TRA	Specialities TRADITIONAL	TRADITIONAL

**Quality:**

	SSE	ECF	ETH	ECL	TRA
<b>CCL:</b>	1 2	1 2	1 2 3	1 2	2 3 4
<b>CCL:</b>	<input type="checkbox"/> <input type="checkbox"/>	<input type="checkbox"/> <input type="checkbox"/>	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<input type="checkbox"/> <input type="checkbox"/>	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>

For **Specialities TRADITIONAL** compression up to the cG level (without soft-part)

**Models:**

A-D Calf	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
A-D with grip top	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
A-F Half-thigh	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
A-F Half-thigh with grip top*	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
A-G Thigh	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
A-G with grip top	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
A-G Thigh with waist attachment:					
<input type="checkbox"/> left <input type="checkbox"/> right	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
A-T Pantyhose	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
A-T Bodyform	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
A-T Maternity	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
A-T for men	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
A-T One leg Pantyhose	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

\* for SSE only available in CCL 2

**Colors:**

beige	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
caramel	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
crispa	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
nature	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
savanna	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
black	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
skin	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
grey	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
marine	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Pieces  
 Pairs

**TREND COLORS:**

salted caramel	<input type="checkbox"/>
blueberry dream	<input type="checkbox"/>
fresh basil	<input type="checkbox"/>
hot chili	<input type="checkbox"/>

**Grip tops:**

- Sensinova grip top** neutral design 5.0 cm
- Sensinova grip top** feminine design 5.5 cm (only available for SSE)
- Knobbed grip top** 5.0 cm
- Narrow knobbed grip top** 2.5 cm (for A-D calf stockings available only)
- Micro knobbed grip top** 5.0 cm (only available for ETH A-G thigh stockings in CCL1 and 2)

# Order form Sigvaris Roundknit Custom Made Stockings

## Leg measurements

### Taking measurements

Please use a separate order form for each model.

### When to measure

The measurements should be taken in the morning, on the edema-free leg, in a standing position.

The measurements should be made directly on the skin. Both legs must be measured, if garments for both legs are required.

### How should be measured

**Circumference:** The tape measure should be tight but not cut in when measuring.


**Length:** Measure the length (with the exception of  $\ell T$ ) on the inside of the extended leg in a perpendicular position.

### NOTE

We cannot be held responsible for incorrect measurements given. Custom-made compression stockings cannot be returned.

Product liability cannot be accepted if the label has been removed.

## Leg measurements:



### Toe:

Open toe

Closed toe