

Στη φόρμα παραγγελίας δεν περιλαμβάνονται έξοδα αποστολής και ΦΠΑ 13%

invoice / delivery address

Reference:

☐ **Order**

☐ **Offer**

Rehab Car Seat **HERO** **NXT**

Exportprices

Release: 2025

BASIC

Details:

	HERO-NXT WITH STANDARD CUSHION	HERO-NXT WITHOUT STANDARD CUSHION
SEAT WIDTH (pelvis)	21	25 cm
SEAT WIDTH (knee)	28	28 cm
SEAT DEPTH	25	28 cm
BACK HEIGHT	21 (24) - 34	27 - 35 cm
CROWN HEIGHT	41 - 53	44 - 56 cm
DISTANCE LATERALS	12 - 21 cm	
total height	66 - 73 cm (without feet 60 cm)	
total width	43 cm	
total depth	63 cm	
capacity	18 kg	

common information (standard):

E4 08025600
UN Regulation
No. 129/03

- approved according to ECE-R 129-00 (i.Size)
- 5-point-harness with upholstery, height adjustable
- belt mechanics installed into the front of the seat
- height adjustable back and head rest
- additional upholstery for the seat, back and head rest
- large seat inclination
- swivel 360° and adjustable stilt
- ISOFIX
- covers washable (30 °C)

Note: ISOFIX adaptation in the car is absolutely necessary!!!



☐ **HERO-NXT
red-black**
(470000-E1)
€ 2.250,00



☐ **HERO-NXT
grey-black**
(470000-E2)
€ 2.250,00

accessory for stabilization

abduction block

☐ (470400-E)
[€ 170,--]



security tray

☐ (470600-E)
[€ 375,--]



side supports made of
aluminum fully adjustable

☐ (470800-E)
[€ 415,--]



longer 5-point-harness
(+10 cm)

☐ (470320-E)
[€ 105,--]



further accessories

protection mat
(to protect original seats of the car)

☐ (470780-E)
[€ 85,--]



head fixation
(stabilizes the head if there's no control)

☐ (470630-E)
[€ 375,--]



security pad for the buckle
with snap fastener

☐ (470650-E), [€ 165,--]



belt to connect the chest straps

☐ (470640-E)
[€ 125,--]



notes