

Στη φόρμα παραγγελίας δεν περιλαμβάνονται έξοδα αποστολής και ΦΠΑ 13%

invoice / delivery address

Reference:

☐ **Order**

☐ **Offer**

Rehab Car Seat **HERO** **XL**

Exportprices

Release: 2025

BASIC

Details:

	HERO-XL
SEAT WIDTH (pelvis)	24 cm
SEAT WIDTH (knee)	27 cm
SEAT DEPTH	25 - 34 cm
BACK HEIGHT	28 - 45 cm
CROWN HEIGHT	48 - 68 cm
HEIGHT LATERALS	13 - 20 cm
total height	65 - 82 cm
total width	43 cm
total depth	55 cm
capacity	36 kg

common information (standard):

- homologated according to ECE-R 129/03
- rear-facing car seat
- 5-point-harness with upholstery
- height adjustable head rest
- additional upholstery for the seat, back and head rest
- large seat inclination
- tie-down straps
- side impact protection
- covers washable (30 °C)



Note: ISOFIX latches in the car are NOT necessary!!!



☐ **HERO-XL red-black**
(490000-E1)



☐ **HERO-XL grey-black**
(490000-E2)

accessories for stabilization

- ☐ **MOUNTING FRAME**
(necessary for the options on the right!)
(490100-E)

- ☐ abduction block
(490400-E)



- ☐
- seatdepth adapter
-
- (490850-E)




further accessories

- ☐ security tray
(490600-E)



- ☐ side supports made of aluminum, fully adjustable (490800-E)



- ☐ security pad for the buckle
with snap fastener
(490650-E)



- ☐ belt to connect the chest straps (490640-E)



- ☐ protection mat
(to protect the seats of the car)
(490780-E)



- ☐ head fixation
(in case of poor head control)
(490630-E)



- ☐ Tether straps
(490710-E)



- ☐ Attachment loops, 70 mm
(490720-E)



- ☐ Lower Tether Anchorages
120 mm
(490730-E)



- ☐ Lower Tether Anchorages
160 mm
(490740-E)



notes