

Στη φόρμα παραγγελίας δεν περιλαμβάνονται έξοδα αποστολής και ΦΠΑ 13%

invoice / delivery address

Kom.:

☐ **Order**

☐ **Offer**

Rehab Car Seat **KIDSFLEX^{PRO}**

Exportprices

Release: 2026


BASIC

	size 1	size 2	size XL
crown height	56 - 71 cm	58 - 73 cm	59 - 80 cm
back height	33 - 48 cm	33 - 51 cm	34 - 60 cm
seath width	28 cm	33 cm	40 cm
seat depth (adjustable)	31 - 37 cm	36 - 42 cm	39 - 46 cm
height trunk supports, small	19 - 38 cm	24 - 43 cm	24 - 44 cm
height trunk supports, large	23 - 42 cm	28 - 47 cm	30 - 50 cm
distance trunk supports	14 - 25 cm	17 - 28 cm	17 - 35 cm
body height	100 - 125 cm (max. 50 kg)	100 - 150 cm (max. 50 kg)	ca. 125 - 175 cm (max. 80 kg)
total height	71 - 86 cm	72 - 87 cm	72 - 92 cm

common information (standard):

- Sizes 1 and 2 approved according to ECE R129/04
- Size XL crash-tested based on ECE R129*
- 5-point-harness with upholstery, height adjustable
- head support „standard“ inclusive
- base board with ISOFIX adapter inclusive
- stepless back angle adjustment
- adjustable seat depth
- shoulder supports inclusive
- trunk supports, small inclusive
- adjustable back plate

* ECE R129 applies only to children up to 150 cm in height. Size XL is designed for taller children, for whom no approval standard exists. The seat has therefore been voluntarily crash-tested using the ECE R129 test procedures.



☐ size 1
grey-black
(360000-E1)

☐ size 1
red-black
(360000-E2)

3.650,00 €



☐ size 2
grey-black
(370000-E1)

☐ size 2
red-black
(370000-E2)

3.750,00 €



☐ size XL
grey-black
(380000-E1)

3.850,00 €

belt guide (please ALWAYS choose one of them)

☐ LEFT
(300000-L)

☐ RIGHT
(300000-R)



☐ seat cushion
standard

(350355-E)
[€ 0,-]

☐ seat cushion PRO (ergonomically shaped
seat cushion including abduction pads and
width-adjustable leg supports)

size 1 (360350-E), [€ 415,-]
size 2 (370350-E), [€ 435,-]
size XL (380350-E), [€ 455,-]

fittings for different needs

head support, standard
(free of charge)



☐ head support PRO
(350025-E), [€ 295,-]



☐ head fixation
(350630-E), [€ 375,-]



☐ wider hip cushions
(reduces the seat width by 4 cm)

(350850-E), [€ 90,-]

abduction block, removable

☐ size 1
(360410-E)
[€ 275,-]

☐ size 2
(370410-E)
[€ 285,-]

☐ size XL
(380410-E)
[€ 295,-]



foot rest, height and angle adjustable


☐ size 1 and 2
(350450-E)
[€ 335,-]

☐ swing up
(350452-E)
[€ 190,-]



foot roll, height and angle adjustable

☐ size 1 and 2
(350455-E)
[€ 335,-]



trunk supports, small
(height: 11 cm, depth: 15 cm)

☒ fix
(350500-E)
[€ 0,-]

☐ swing away
(350505-E)
[€ 240,-]



trunk supports, large
(height: 21 cm, depth: 16 cm)

☐ fix
(350510-E)
[€ 80,-]

☐ swing away
(350515-E)
[€ 290,-]



shoulder supports
(free of charge)


 fittings frame

swivel base

☐ (350200-E)
[€ 435,-]

☐ RIGHT

☐ LEFT



☐ angular adjustment (5° - 20°)

size 1 (360100-E), [€ 335,-]

size 2 (370100-E), [€ 355,-]

size XL (380100-E), [€ 395,-]



☐ conversion set

☐ from ISOFIX adapter
to clamp bow

(350020-E)

[€ 265,00]

☐ from clamp bow
to ISOFIX adapter

(350030-E)

[€ 265,00]

fittings for car installation

Fixation (one option included in delivery)

ISOFIX adapter

☐ size 1
(360575-E)

☐ size 2 / XL
(370575-E)



clamp bow

☐ size 1
(360550-E)

☐ size 2 / XL
(370550-E)


outrigger*

☐ (350453-E), [€ 190,-]

*350450-E
or
350450-E
necessary!

additional belt guide
(on the other side)

☐ (350075-E), [€ 135,-]

 fittings security

security tray

(in height and depth adjustable)

☐ size 1
(360600-E)
[€ 455,-]

☐ size 2
(370600-E)
[€ 475,-]

☐ size XL
(380600-E)
[€ 495,-]


security pad for the buckle

Cover with snap fastener

☐ (350650-E)
[€ 165,-]

security pad for the buckle
with HERNIK Security Key*

☐ (350655-E)
[€ 210,-]

NEW!

* incl. belt cutter

protection mat

(to protect the original seats of the car)

☐ (350780-E)
[€ 85,-]


individual colors for the cushions

(300800-E), [€ 330,-]

(one X means single-colored, two X means bicolored)

	inside	outside
☐ red	☐	☐
☐ black	☐	☐
☐ grey	☐	☐

belt to connect the chest straps

☐ (350640-E)
[€ 125,-]

 notes
