

KIDSFLEX R129 Edition

Hernik

Crash-tested and secure



Size 1: for ages 3 to 9 years

Size 2: for ages 6 to 14 years

Size 2XL: suitable for ages 14 and up



Comfortable & secure

We rely on the latest **technology** to make our products safe.

✓ **Fast switch** between active and passive **seating position**



Trust in experience

20 years of experience developing and manufacturing the best products.

✓ **European certification** in accordance with ECE R129/04.



Made in Germany

All our products are developed and **manufactured in Germany**.

✓ **Certified according** to **ISO 9001**

Further refined to meet individual requirements

The KIDSFLEX R129 Edition merges the established adaptability of the KIDSFLEX with improved features. The seat has undergone crash testing in compliance with the R129/03 standard and provides critical components for the secure and stable positioning of the child, even in its fundamental configuration.

Carefully considered essential equipment included:

- 5-point harness featuring a quick-release mechanism
- ISOFIX-Adapter
- Headrest featuring kick-back protection
- Padded shoulder support
- Thoracic supports



Customizable – for your individual needs

Thanks to the customizable **backrest (1)** and by **tilting the seat (2)**, it can be precisely adjusted to align with the car seat or the child's requirements.

The modification of the **seat depth (3)**, along with the adjustable belt guide in the backrest, provides ample options to customize the car seat as the child grows up.

Even if one of three different sizes of the KIDSFLEX R129 Edition do not correspond to the child's measurements, we provide a tailored solution for this, by customizing the dimensions to meet your individual needs!



We are very proud of the many positive customer ratings. On our website you can find out, why Hernik regularly achieves top ratings on comparison portals.

Markus Jacksch,
Managing Director

Accessory Options

A range of accessory options accommodates the unique needs of each child, ensuring that no accessory is mutually exclusive.

A highly sought-after feature is the **swivel base (1)**, which significantly facilitates placing the child in and removing them from the seat.

The **kick-back headrest (2)** provided in the standard equipment can be adjusted to accommodate the child's needs and absorbs energy during a side impact through its lateral support.



Additional Accesories

For optimal positioning, supplementary accessories such as **angular adjustment (3)**, **leg guides (4)**, an **abduction block (5)**, and **swing-away thoracic supports (6)** are available. Kindly consult our order form for a comprehensive list of other equipment options.



We are certain that we can only develop the best seats together with our customers. That is why we are very happy to receive your suggestions for improvement and will incorporate them into our development.

Please send us a message:
office@hernik.de



Further products from HERNIK




Child car seats for every size

HERNIK accompanies your child from infancy to adulthood (up to 75 kg). For more complex requirements and individual requests we also offer customised solutions.

Head protection helmets and seat shells

With our wide range of head protection helmets we offer every user the necessary protection against injuries. Our CAPOVARIO in particular is a real highlight that makes everyday life much easier thanks to its special design.

We also supply you with individual seat shells, such as the ANATOMIC UNICO and ANATOMIC FLEX aluminium shells.



Hernik

**Discover all our products
on hernik.de**

HERNIK GmbH

Gewerbering 8,
83543 Rott am Inn

+49 (08039) 9093930
office@hernik.de
www.hernik.de

