

PowerPack to Wheelchair Fitting Chart - May 2018

Wheelchair	Crossbar*	Drop Plate *	End Brackets*	Fitting Kit (126xx = with new P'pack 125xx = Spares orders)	Extras*
Ability 2	As Requested	As Requested	As Requested	Ability 2 wheelbase (All types)	None
Action NG	Short/Standard	100mm	Black Hinged	None	None
Aktiv X1	Standard	Standard	black hinged	n/a	
Aktiv X2	Standard	50mm	Black hinged	n/a	
Aktiv X3	Standard	50mm	Black hinged	n/a	
Aktiv X4	Long	50mm	Black hinged	n/a	
Aktiv X5	Long	50mm	Black hinged	n/a	
Aktiv X6	Short	75mm	Black hinged	n/a	
Aktiv X7	n/a	n/a	n/a	1266014 or 1252006 (Handivipp 25mm)	Alternate mounting bracket
Aktiv X8	Short	100mm	Black hinged	n/a	
CHUNC					Refer to Specmat (CHUNC)
Cirrus 4	Short	50/75mm on 50mm setting	Cirrus Jubilee	1266001 or 1252009 (Cirrus Jubilee)	Alternate mounting bracket
Cirrus Classic	Short	50/75mm on 50mm setting	Cirrus Jubilee	1266001 or 1252009 (Cirrus Jubilee)	
Cirrus Jubilee	Short	50/75mm on 50mm setting	Cirrus Jubilee	1266001 or 1252009 (Cirrus Jubilee)	Alternate mounting bracket
Days Escape	Standard	50/75mm on 50mm setting (rise position)	Black Hinged	None	None
Days Escape Lite	Short	50mm (rise position)	Black Hinged	None	None
Days Swift	Short/Standard	50/75mm on 50mm setting	63mm Forward/25mm offset	None	None
Drive Enigma	Standard	50/75mm on 50mm setting	Black Hinged	None	None
Drive Medical Expedition	Standard	Standard	Blk hinged	N/A	N/A
Drive Medical I.D. Soft	Standard	100mm	Dahlia	N/A	Alternate mounting bracket
Drive Medical K-Chair	Standard	100mm	Blk hinged	N/A	N/A
Drive Medical Phantom	Standard	100mm	Blk hinged	N/A	N/A
Drive Medical S1	Standard	Standard	Blk hinged	N/A	N/A
Drive Medical S4	Standard	Standard	Blk hinged	N/A	N/A
Drive Medical SD2	Standard	100mm	Blk hinged	N/A	N/A
Drive Medical Sentra	Heavy Duty	Standard	black hinged	n/a	Alternate mounting bracket
Drive Medical Silver Sport	Standard	Standard	Blk hinged	N/A	N/A
Drive Medical Ultra Light	Standard	75mm	Blk hinged	N/A	N/A
Drive Medical XS	Standard	75mm	Square Tubing	N/A	N/A
Drive Medical XS2	Standard	100mm	Blk hinged	N/A	N/A
Enigma Spirit	Long	50/75mm on 50mm setting	Black Hinged	None	Alternate mounting bracket
Enigma XS	Long	50/75mm on 50mm setting	Black Hinged	None	Alternate mounting bracket
ETAC Cross	Standard	75mm	ETAC brkts	n/a	Alternate mounting bracket
ETAC cross XL	Heavy Duty	75mm	ETAC	N/A	Alternate mounting bracket
Excel Vanos Xeryusc	Long	50/75mm on 75mm setting	Black Hinged	None	None
Freedom CGX	As Requested	50/75mm on 50mm setting	Neo Brackets	None	Alternate mounting bracket
G Lite Pro	Standard	100mm	Black Hinged	None	None
Gillingham Tilt	As Requested	50/75mm on 50mm setting	Gillingham Tilt	1266002 or 1252004 (Gillingham Tilt)	Pedal
Greencare	Heavy Duty	50/75mm drop plate on 50mm setting	63mm/25mm offset	None	Alternate mounting bracket
Greencare DB1	As requested	50/75mm on 75mm setting	63mm/25mm offset	None	Alternate mounting bracket
Handicare Emineo	Standard	50/75mm on 50mm setting	63mm Forward/25mm offset	None	2 x Alternate mounting bracket
Handicare Ibis	Short	Standard	Ibis XA	None	None
Handivipp 20mm	n/a	n/a	n/a	1266013 or 1252005 (Handivipp 20mm)	Alternate mounting bracket
Handivipp 25mm	n/a	n/a	n/a	1266014 or 1252006 (Handivipp 25mm)	Alternate mounting bracket
Handivipp 27.5mm	n/a	n/a	n/a	1266015 or 1252007 (Handivipp 27.5mm)	Alternate mounting bracket
Hoggi Wheelbase	Short	75mm	Neo Brackets	None	Alternate mounting bracket
IBIS Pro	Standard	50/75mm rise	IBIS XA Brackets	None	Alternate mounting bracket
Invacare Action 2000	As Requested	50/75mm on 50mm setting	Black Hinged	None	None
Invacare Action 3	Standard	100mm	Black Hinged	None	Alternate mounting bracket
Invacare Action 3 Junior	Short	75mm	Blk hinged	N/A	Alternate mounting bracket
Invacare Action 4 (NG)	Long	100mm	Standard and 63mm Forward/25mm offset	None	None
Invacare Action 5	Standard	100 mm	63mm fwd	n/a	Alternate mounting bracket
Invacare Action Express	Standard	100mm	63mm Forward/25mm offset	None	None

*Where a fitting kit exists these parts should be included in the fitting kit. Please confirm at the time of order.

PowerPack to Wheelchair Fitting Chart - May 2018

Wheelchair	Crossbar*	Drop Plate *	End Brackets*	Fitting Kit (126xx = with new P'pack 125xx = Spares orders)	Extras*
Invacare Alu Lite	Standard	50mm (rise position)	Black hinged	n/a	
Invacare Assist	n/a	n/a	n/a	1266014 or 1252006 (Handivipp 25mm)	
Invacare Azalea	n/a	n/a	n/a	1266014 or 1252006 (Handivipp 25mm)	None
Invacare Ben 9	Long	50/75mm on 75mm setting	63mm Forward	None	None
Invacare BEN NG	Long	75mm	Blk hinged	N/A	N/A
Invacare Blade	short/std	75mm	black hinged	n/a	n/a
Invacare Clematis	N/A	N/A	N/A	Handivipp 25mm	N/A
Invacare Echo 2000	As Requested	100mm	Black Hinged	None	None
Invacare Imperial	Standard	Standard	Black Hinged	None	None
Invacare NXT	Standard	50mm	Blk hinged	N/A	Alt Mtg Brkt
Invacare Rea Assist (S/P)	n/a	n/a	n/a	1266014 or 1252006 (Handivipp 25mm)	None
Invacare Rea Assist (T)	n/a	n/a	n/a	None	None
Invacare Rea Dhalia	Standard	50/75mm on 50mm setting	Dhalia Brackets	None	None
Invacare Topaz	Heavy duty	50mm	63x25mm offset	n/a	Alternate mounting bracket
Karma Advantage	Heavy duty	50mm	Blk hinged	N/A	Alternate mounting bracket
Karma Condor 24" seat	Heavy Duty	50mm (rise position)	Square tubing	None	Alternate mounting bracket
Karma Dove	Standard	50/75mm on 50mm setting	63mm Forward/25mm offset	None	None
Karma Ergo 125	Standard	50/75mm on 50mm setting	Black Hinged	None	None
Karma Ergo 305	Standard	50/75mm on 50mm setting	63mm Forward/25mm offset	None	Alternate mounting bracket
Karma Ergo Lite	Standard	75mm	Blk hinged	N/A	Alternate mounting bracket
Karma Ergo/ Karma Ergo 115	Standard	50/75mm on 75mm setting	Black Hinged	None	None
Karma Flex	Standard	100mm	Black Hinged	None	None
Karma i- Explore	Standard	75mm	Blk hinged	N/A	Alternate mounting bracket
Karma i-lite	Standard	Standard	Blk hinged	N/A	N/A
Karma Martin	Short	50/75mm on 50mm setting	Black Hinged	None	None
Karma MVP	Standard	100mm	Blk hinged	N/A	Alternate mounting bracket
Karma Robin	Long	50/75mm on 50mm setting	Screw Z	None	None
Karma Sparrow	Standard	50/75mm on 75mm setting	Standard and 63mm Forward/25mm offset	None	None
Karma VIP	Standard	50mm	Blk hinged	N/A	Alternate mounting bracket
Karma Wren	Short/Standard	100mm	Black Hinged	None	None
Karma Wren 2	Short	Standard	Standard and 63mm Forward/25mm offset	None	None
Lomax Standard	Short	Standard	Black Hinged	None	None
Lomax Standard Tapper	As Requested	As Requested	63mm Forward/25mm offset	None	None
Lomax Uni8 /9	Standard	50/75mm on 75mm setting	63mm Forward	None	None
Meyra Eurochair	Standard	75mm	Blk hinged	N/A	Alternate mounting bracket
Meyra Service chair	Standard	Standard	Blk hinged	N/A	N/A
Neba Pram	Special	Standard	Black Hinged	None	Pram handle alternate mounting bracket
Neo wheelbase	n/a	n/a	n/a	1266003 or 1252016 (Neo fitting kit)	
OttoBock Avantgarde	Standard	100 mm	Blk hinged	N/A	Alternate mounting bracket
Ottobock Dino	n/a	n/a	n/a	1266014 or 1252006 (Handivipp 25mm)	
Ottobock Discovery	n/a	n/a	n/a	1266014 or 1252006 (Handivipp 25mm)	
Ottobock Motus	Std/H'Duty	100mm	black hinged	n/a	Alternate mounting bracket
OttoBock Motus	Standard	100mm	Blk hinged	N/A	Alternate mounting bracket
Ottobock Start	Standard	Standard	Black Hinged	None	None
Patterson Medical Whirl	Standard	Standard	black hinged	n/a	n/a
PDG Fuze	short/std	75mm	Fuze brackets	n/a	Alternate mounting bracket
PDG Stellar	N/A	N/A	N/A	NEO	Alternate mounting bracket
R82 Cougar	Short	75mm	Gillingham Tilt	None	None
Radcliffe Netti 25mm	n/a	n/a	n/a	1266014 or 1252006 (Handivipp 25mm)	Long lifting strap
Radcliffe Shadow 15" Single Gas Strut (MK4 On)	Extra Short	Radcliffe Shadow 15"	Radcliffe Shadow	1266005 or 1252002 (Radcliffe Shadow 15")	Alternate mounting bracket
Radcliffe Shadow 18" Single Gas Strut (MK4 On)	Short	Radcliffe Shadow 18"	Radcliffe Shadow	1266006 or 1252001 (Radcliffe Shadow 18")	Alternate mounting bracket
Remploy (R Healthcare) 8TRL	Standard	75mm	63mm fwd	n/a	

*Where a fitting kit exists these parts should be included in the fitting kit. Please confirm at the time of order.

PowerPack to Wheelchair Fitting Chart - May 2018

Wheelchair	Crossbar*	Drop Plate *	End Brackets*	Fitting Kit (126xx = with new P'pack 125xx = Spares orders)	Extras*
Remploy (R Healthcare) Accent	Standard	75mm	63mm fwd	n/a	
Remploy (R Healthcare) Access	Standard	75mm	63mm fwd	n/a	
Remploy (R Healthcare) Dash Life	Long	75mm	63mm x 25mm offset	n/a	Alternate mounting bracket
Remploy (R Healthcare) Dash Lite	Standard	75mm	63mm fwd	n/a	
Remploy 9TRL	Long	75mm	63mm Forward	None	None
Remploy Aquila	Long	50/75mm on 50mm setting	Black Hinged	None	None
Remploy GL	Long	50/75mm on 50mm setting	Black Hinged	None	None
RMS Juditta	N/A	N/A	N/A	Handivipp 25mm	N/A
Roma 1500	Standard	75mm	Blk hinged	N/A	N/A
Roma 1530	Standard	75mm	Blk hinged	N/A	N/A
Roma 1630	Standard	Standard	Blk hinged	N/A	N/A
Roma Avant	Standard	Standard	Blk hinged	N/A	N/A
Roma Elite 1485	Long	50/75mm on 50mm setting	63mm Forward/25mm offset	None	None
Roma Medical 1235	Standard	Standard	Blk hinged	N/A	N/A
ROMA Medical steel range	Standard	Standard	black hinged	n/a	n/a
Roma Orbit	Standard	50mm	Blk hinged	N/A	N/A
SEABASS Wheelbase	Standard	30mm drop 75mm forward	Neo Brackets	n/a	Alternate mounting bracket
S.O.S Mojo	Standard	n/a	n/a	1266009 or 1252023 (Mojo fitting kit)	
Strongback att prop 40cm	Short	50mm	n/a	1266012 or Please call	Alt mtg brkt (DO NOT SUPPLY/FIT Handlebar Extensions)
Strongback att prop 45cm	Standard	50mm	n/a	1266012 or Please call	Alt mtg brkt (DO NOT SUPPLY/FIT Handlebar Extensions)
Strongback self prop 40cm	Short	100 mm	n/a	1266012 or Please call	Alt mtg brkt (DO NOT SUPPLY/FIT Handlebar Extensions)
Strongback self prop 45cm	Standard	100 mm	n/a	1266012 or Please call	Alt mtg brkt (DO NOT SUPPLY/FIT Handlebar Extensions)
SUMED Goliath	Heavy duty	75mm	Goliath brkts	n/a	Alternate mounting bracket
SUMED Netti4U	Standard	Standard	black hinged	n/a	Alternate mounting bracket
SUMED Phonix	Heavy duty	100mm	Invacare Dhalia	n/a	Alternate mounting bracket
Sunrise Breezy 300	As Requested	50/75mm on 50mm setting	63mm Forward/25mm offset	None	None
Sunrise Breezy Basix	Standard	100mm	Blk hinged	N/A	N/A
Sunrise Breezy Elegance	Standard	50/75mm on 50mm setting	63mm Forward/25mm offset	None	None
Sunrise Breezy RelaX	N/A	N/A	N/A	Handivipp 25mm	N/A
Sunrise Breezy Rubix	Standard	100mm	Blk hinged	N/A	N/A
Sunrise Breezy SL	Standard	100mm	Blk hinged	N/A	N/A
Sunrise Breezy TL	Standard	100mm	Blk hinged	N/A	N/A
Sunrise Moonlite	Standard	100mm	Black hinged	n/a	
Sunrise Quickie Easy 200	Standard	75mm	Blk hinged	N/A	Alternate mounting bracket
Sunrise Quickie Helix	Standard	100mm	Blk hinged	N/A	Alternate mounting bracket
Sunrise Quickie Iris	N/A	N/A	N/A	Handivipp 25mm	Alternate mounting bracket
Sunrise Quickie Life	Standard	100mm	Blk hinged	N/A	Alternate mounting bracket
Sunrise QuicKie Life T	Short	100mm	Blk hinged	N/A	Alternate mounting bracket
Sunrise Quickie RXS	Standard	100mm	Blk hinged	N/A	Alternate mounting bracket
Sunrise Quickie RXS	Standard	100mm	Blk hinged	N/A	Alternate mounting bracket
Sunrise Zippie 2	Standard	50/75mm on 75mm setting	63mm Forward/25mm offset	None	None
Sunrise Zippie Special	Extra Short	50mm	Black Hinged	None	None
Sunrise Zippie TS	Standard	Standard	Black Hinged	None	None
TGA Elite	Standard	30mm drop 75mm forward	None	None	None
Ultimate Healthcare Drift	Standard	75mm	Black hinged	n/a	Alternate mounting bracket
Ultimate Healthcare Gravity	n/a	n/a	n/a	1266014 or 1252006 (Handivipp 25mm)	Alternate mounting bracket
Ultimate Healthcare Impulse	Standard	75mm	63mm fwd	n/a	Alternate mounting bracket
Van Os Excel G5	Standard	50mm	Blk hinged	N/A	Alternate mounting bracket
Van Os Excel G7	N/A	N/A	N/A	Handivipp 25mm	N/A
Van Os Excel G-Lite	Standard	100mm	Blk hinged	N/A	N/A
Van Os Excel Globe Traveller	Standard	Standard	Blk hinged	N/A	N/A
Wheeltech Enigma	Short	Standard (rise position)	63mm Forward/25mm offset	None	None
Wolturnus A-Run	N/A	N/A	N/A	Handivipp 25mm	N/A
Wren 2	Standard	Standard	Black Hinged	None	None
Z-Tec 610A Transit	Standard	50mm (rise position)	63mm/25	None	None
Z-Tec Di-Lite	Standard	75mm	Blk hinged	N/A	Alternate mounting bracket
Z-Tec H/D Steel	Heavy duty	Standard	Blk hinged	N/A	Alternate mounting bracket
Z-Tec Heavy Duty	Heavy duty	Standard	Blk hinged	N/A	Alternate mounting bracket
Z-Tec Hi Line	Standard	75mm	Blk hinged	N/A	Alternate mounting bracket

*Where a fitting kit exists these parts should be included in the fitting kit. Please confirm at the time of order.

PowerPack to Wheelchair Fitting Chart - May 2018

Wheelchair	Crossbar*	Drop Plate *	End Brackets*	Fitting Kit (126xx = with new P'pack 125xx = Spares orders)	Extras*
Z-Tec T-Line	Standard	75mm	Blk hinged	N/A	Alternate mounting bracket
Z-Tec Trusty	Standard	Standard	Blk hinged	N/A	N/A