

Address

Fax +49 7261-7351-10

Shipping address

Berollka-aktiv Rollstuhltechnik GmbH
 Jahnstraße 16
 74889 Sinsheim - Germany
Phone: +49 7261-7351-0
 info@berollka.de
 www.berollka.de



Com.:

- Order
 Estimate

Basic Price 2298,- €

Price valid from February 1st 2022
 (until revoked)

Racer

Όλες οι αναγραφόμενες τιμές
 δεν περιλαμβάνουν μεταφορικά έξοδα και Φ.Π.Α

Mobility base with short frame for seating systems

Capacity up to 90kg

Standard equipment:

- Rigid short aluminium frame with abduction
- Castor
- Back with hand bow, angle adjustable, high adjustable and tiltable
- Rear wheels with drum brake (and camber)

User weight + Seating system: _____ kg

max. width
 seat shell

basic measure

frame	frame pipe	seat frame width B	frame width front A	frame width back C	width between wheels D	total width E
Frame size 0 <input type="checkbox"/> Code U5000	30mm	26 cm	44 cm	34 cm	37,5 cm 40,5cm*	55,5 cm 52cm*
Frame size 1 <input type="checkbox"/> Code U5001	30mm	32 cm	50 cm	40 cm	43,5 cm 46,5cm*	61,5 cm 58cm*
Frame size 2 <input type="checkbox"/> Code U5002	30mm	38 cm	56 cm	46 cm	49,5 cm 52,5cm*	67,5 cm 64cm*

* Attention: for size 2 we suggest the frame with two gas-springs for seat adjustment (see CODE U5008) !!!

* without wheel cover and camber

* width with 12" wheel

frame colour

- 25 RAL 9010 white
- 24 RAL 1028 melon yellow
- 27 mandarin met.
- 06 dormant red
- 40 red metallic
- 53 dormant cassis
- 59 atlantic blue
- 16 RAL 5002 light blue
- 17 dormant blue
- 45 RAL 6018 yellow green
- 60 mojito green
- 19 silver metallic
- 46 RAL 9023 perl dark grey
- 20 grey metallic
- 22 RAL 9011 black mat
- 58 satin black
- 23 RAL 9005 black

Because of the material are different colours possible

Edition "Duo-Color"

Colour 1: _____
 - Seat frame
 - Footrest holder
 - Footrest

Colour 2: _____
 - Frame
 - Gas spring bar
 - Back down
 - Hand bow

Code U5005

OE 68,- €

met. = metallized

technical data

- Seat tilt
 -10 to 30°
 -2 to 35°



min. seat height with rear wheel
 12" a. 20" **42 cm**
 22" **43 cm**
 24" **44 cm**

Incline back with gas spring
90 to 120°

concerning seat frame (trapezium adapter +2cm)
 seat can be changed 9cm in height

Weight:
 from 11 kg (without rear wheel)
 from 14 kg (with rear wheel 24")

H total height: 95 - 117cm*
 L total length: 70 - 80cm
 T seat depth: size 0 = 32cm
 size 1 / 2 = 35cm



* according to the angle of the back / hand and the size of the rear wheels

frame accessories		armrest		
Tip shackle with adapter <input type="checkbox"/> right <input type="checkbox"/> left Code U5011	Anti-tip wheel with adapter <input type="checkbox"/> right <input type="checkbox"/> left Code U5012	Armrest holder with arm pads, fixed at the back tube (arm pads swing up) Code U5046	Armrest adjustable in height with raster, for mounting on the seating shell <input type="checkbox"/> armrest long <input type="checkbox"/> armrest short Code U5048 (self installation on the seating shell)	Armrest short adjustable in height and angle adjustable with raster, for mounting on the seating shell Code U5047 (self installation on the seating shell)
<input type="checkbox"/> Piece 78,- €	<input type="checkbox"/> Piece 130,- €	<input type="checkbox"/> Pair 268,- €	<input type="checkbox"/> Pair 198,- €	<input type="checkbox"/> Pair 248,- €

footrest					
Single footrest with knee angle adjustable Code U5022	Single footrest with knee angle and gear segment adjustable Code U5023	Footrest angle adjustable additional joint for angle adjustable (for single footrest with knee angle CODE U5022 or U5023) Code U5028	Footrest folding upwards (for single footrest with knee angle CODE U5022 or U5023) Footrest without heel edge from size 1 Code U725 (not in combination with CODE U5028)	Single footrest alu-profile, with holder at the seat tube angle adjustable, folding upwards and LLL adjustable LLL from upper edge seat frame to footrest <input type="checkbox"/> LLL 15 - 32cm (mounted in frame) <input type="checkbox"/> LLL 30 - 40cm* (Footrest holder long with footrest bow) Code U5026 * Note: only possible at a very high seat height (because of the castor)	Footrest holder high adjustable only in combination with CODE U5026 or U5027 additional gear segment Code U5029
<input type="checkbox"/> Piece 258,- €	<input type="checkbox"/> Piece 318,- €	<input type="checkbox"/> Pc. OE 95,- €	<input type="checkbox"/> Pc. OE 158,- €	<input type="checkbox"/> Piece 332,- €	<input type="checkbox"/> Pair OE 170,- €

footrest					
Spasticity lock for footrest folding upwards for CODE U725 <input type="checkbox"/> Code U725S single footrest for CODE U5026 <input type="checkbox"/> Code U5024 double footrest for CODE U5027 <input type="checkbox"/> Code U723	Double footrest alu-profile, with holder on the seat tube angle adjustable, folding upwards and LLL adjustable LLL from upper edge seat frame to footrest <input type="checkbox"/> LLL 15 - 32cm (mounted in frame) <input type="checkbox"/> LLL 30 - 40cm* (Footrest holder long with footrest bow) Code U5027 * Note: only possible at a very high seat height (because of the castor)	Single footrest alu. flat-profile gear segment knee angle Code U7025	Double footrest alu. flat-profile gear segment knee angle Code U7026	Heel edge for double footrest for CODE U7026 <input type="checkbox"/> Pair Code U7028 for single footrest for CODE U7025 <input type="checkbox"/> Piece Code U7027	Calves pads for double footrest for CODE U7026 <input type="checkbox"/> Pair Code U7023 for single footrest for CODE U7025, U5021, U5023 <input type="checkbox"/> Piece Code U7022
OE 68,- €	OE 78,- €	OE 142,- €	<input type="checkbox"/> Pair 351,- €	<input type="checkbox"/> Piece 326,- €	<input type="checkbox"/> Pair 406,- €
				each 128,- €	each 130,- €

back						
Back with hand bow high adjustable, angle adjustable and tiltable without one hand brake (2 single brake lever) <input type="checkbox"/> Hand bow U form by size 2 = standard <input type="checkbox"/> Hand bow trapezium form by size 0 a. 1 = standard Code U5044	Back with single push handles, high adjustable (2 single brake lever) <input type="checkbox"/> Code U5040	Back tubes 10cm higher <input type="checkbox"/> Code U5045	Back guide in the centre with fork for seating shell Code U5042	Back guide for seating shells with slide bush on back tube and angle adjustable holder <input type="checkbox"/> Code U742R with Code U5040 not possible!! <input type="checkbox"/> Code U742RW with Code U5040 not possible!!	Back angular adjustment fixed without gas spring <input type="checkbox"/> Code U5043	Back with 2 gas spring for back angle adjustable # Special option # for size 2 <input type="checkbox"/> Code U5041
Standard		OE 98,- €	<input type="checkbox"/> Pc. 240,- €	Pc. 246,- €	Pc. 268,- €	without OE
						Set OE 440,- €

adapter				seat		
Trapezium adapter 2-pieces (for seating shell) Code U5031	Trapezium adapter upper part (for seating shell) Code U5032	Trapezium adapter alu (for mobility base) Code U5033	HC-Trapezium adapter for seating shell 2-pieces alu Code U5039	Seat 10cm longer seat frame long (forwards) Code U5034	Seat extension 7cm (backwards) Not with Scalamobil Code U5035	Seat frame with 2 gas spring for seat angle adjustment ## Special option ## for size 2 Code U5008
<input type="checkbox"/> Set 286,- €	<input type="checkbox"/> Piece 118,- €	<input type="checkbox"/> Piece 169,- €	<input type="checkbox"/> Set 374,- €	<input type="checkbox"/> Set OE 170,- €	<input type="checkbox"/> Set OE 154,- €	<input type="checkbox"/> Set OE 495,- €

Sum

Footrest swing awa from the wheel SPRINT is in soecial construction possible !!

Technical changes and print mistakes are reserved.

VD2425E-30

rear wheel					
Rear wheel with drum brake, hand rims alu and quick release axle, pneumatic tyre (for companion) <input type="checkbox"/> 12 inch Code U5061 <small>without handrim / without crew axle</small> <input type="checkbox"/> 20 inch Code U5062 <small>only with 140mm castor possible</small> <input type="checkbox"/> 22 inch Code U5063 <input type="checkbox"/> 24 inch Code U5064	Rear wheel with drum brake, hand rims alu. and quick release axle, solid tyre (for companion) <input type="checkbox"/> 20 inch Code U5067 <input type="checkbox"/> 22 inch Code U5068 <input type="checkbox"/> 24 inch Code U5069	Camber for rear wheel 20", 22" a. 24" <input type="checkbox"/> 0° standard <input type="checkbox"/> 4° * <input type="checkbox"/> 7° * * not possible: - attachment for scalamobil - wheel base long	Wheel cover plastic (with holder) (for rear wheel 20", 22" or 24") Code U5079	Wheel cover aluminium round, coloured coated (with holder) Code U5076	Wheel cover PC-Carbon (Clothing cover round with holder) (for rear wheel OFF-ROAD 24") Code U776PC
Standard			<input type="checkbox"/> Pair 237,- €	<input type="checkbox"/> Pair 333,- €	<input type="checkbox"/> Pair 274,- €

rear wheel					
SAFETY-GRIP wheel rim with integrated hand rim <input type="checkbox"/> 22 inch Code U7602 <input type="checkbox"/> 24 inch Code U7604 <small>only in combination with wheel CODE U5063, U5064, U5068 or U5069</small>	"Black Edition" Rear wheel with black hub, black rim and quick release axle with black cap and black tyre Marathon Plus Evolution <input type="checkbox"/> 22"x1 luft Code U763S <input type="checkbox"/> 24"x1 luft Code U764S	Spoke cover transparent <input type="checkbox"/> 20 inch Code U5070 <input type="checkbox"/> 22 inch Code U5071 <input type="checkbox"/> 24 inch Code U5072	handrim cover black <input type="checkbox"/> 20" Code U5094 <input type="checkbox"/> 22" Code U5095 <input type="checkbox"/> 24" Code U5096	Wheel base long (extension by 9cm) (In combination with Scalamobil not possible) Code U5077	
<input type="checkbox"/> Pair OE 76,- €	<input type="checkbox"/> Pair OE 99,- €	<input type="checkbox"/> Pair 166,- €	<input type="checkbox"/> Pair 118,- €	<input type="checkbox"/> Pair 154,- €	

rear wheel		accessories	
Rear wheel 24" with OFF-ROAD tyre (with special size of wheel 507) with aluminium hand rim black <small>instead of normal rear wheels -> only in combination with U5064 IIIII</small> <input type="checkbox"/> OFF ROAD 24" x 1,75 (47-507) Mountainbike Profil Code U759 <small>Note: only with 0° camber possible II</small>	Rear wheel OFF-ROAD 24" additional (special size of wheel 507) with aluminium hand rim black and quick release axle additional to the normal rear wheel 24" (to change the wheel) <input type="checkbox"/> Code U759BB OFF ROAD 24" x 1,70 (47-507) Mountainbike profile <small>Note: only with 0° camber possible II</small>		Anchor plate system anchor plates for transport in cars (ISO 7176-19) Assembly kit for mobility base systems including mounting Code U5097 <small>Note: not in combination with E-Fix and not in combination with wheel base long</small>
<input type="checkbox"/> Pair OE 298,- €	<input type="checkbox"/> Pair 720,- €		<input type="checkbox"/> Set net price 440,- €

castor	accessories	brake	
Castor with screw axle <input type="checkbox"/> 140mm black soft (solid) Code U5081S <input type="checkbox"/> 140mm soft (solid) Code U5081 <input type="checkbox"/> 6 Zoll solid * Code U5084 <input type="checkbox"/> 6 Zoll pneumatic * Code U5083 <input type="checkbox"/> 7 Zoll solid * Code U5082 <input type="checkbox"/> 7 Zoll pneumatic * Code U5088 <small>*6 inch or 7 inch not by 20" and 22" rear wheel</small>	Attachment for Scalamobil onyl by rear wheel 24" # from size 1 Code U5009 <small>Combination with tip shackle, anti-tip wheel, wheel base long and seat 7cm longer to the back is not possible!!</small>	Parking brake for rear wheel <input type="checkbox"/> 20 inch Code U5091 <input type="checkbox"/> 22 inch Code U5092 <input type="checkbox"/> 24 inch Code U5093	Brake lever extension straight, foldaway Code U5078
Standard	<input type="checkbox"/> Set 284,- €	<input type="checkbox"/> Pair OE 146,- €	<input type="checkbox"/> Pair 68,- €

Racer	Mobility base for seating systems	basic price	€ 2298,-
For all orders apply the general terms and conditions.		Sum optinal extras	€
		Shipping and handling	€
		Sum	€
		Tax	€
Purchaser _____	Telephone _____	Total price	€

Attention: Scalamobil-extension: Wheel base long, anti-tip wheel, tip shackle and seat 7cm longer to the back not possible !!!

Racer

accessories

Trike


traking extension
as steering and pushing aid



Racer

trekking extension **TRIKE**
as steering and pushing aid

hook closure



CODE U5019

special holder lateral at frame
with hook closure
castor 10x2" black alu fork
arrangement with gear segment
with ramp as mounting aid

capacity up to **90kg**
(wheelchair incl. payload)

Piece

670,- €

TRIKE colour frame Edition "...-Color"

<p>Note: the TRIKE frame parts are made of aluminium</p>	<p>TRIKE- frame colour powder coated</p> <p>Colour : frame colour (see page 1)</p> <ul style="list-style-type: none"> - frame type right / left - cross tube - castor fork <p>Code 9908</p>	<p>Edition "Duo-Color"</p> <p>colour 1: _____ (castor fork)</p> <p>colour 2: _____ (frame type right / left, cross tube)</p> <p>Code 9909</p>	<p>Castor fork with brake for directional stability</p> <p>Code 9909F</p>
	<p>without extra charge <input type="checkbox"/></p> <p>OE 115,- €</p>	<p><input type="checkbox"/></p> <p>OE 125,- €</p>	<p><input type="checkbox"/> Piece</p> <p>OE 144,- €</p>

accessories

Note:

- Racer: - combination TRIKE with anchor plate system is **not possible !!**
 - combination TRIKE mit parking brake is **not possible !!**
 - combination TRIKE only with 12", 20" or 22" possible !!

General Terms and Conditions of Business (GTC)

Validity

All sales are subject to our delivery and payment conditions, without any exception. Orders and other agreements will not be valid unless confirmed by us in writing. We expressly object to the customer's purchase conditions.

Quotations and contract conclusion

Quotations are always subject to change, unless they are expressly limited in terms of time. Orders require a written confirmation. After confirmation, prices quoted in orders are fixed but we reserve the right to make adequate price adjustments in the event of price increases for which we are not responsible (e.g. tax increases, wage increases, price changes on the part of our suppliers). Cost estimates for customized products and repairs are always non-binding. We reserve the right to make technical modifications. The contract will not be concluded before the order is confirmed.

Shipment

To the extent the customer has not made any specification, shipment will be made in usual packaging using an adequate method. Insurance will only be taken out upon request and at the expense of the customer.

Delivery will be made carriage forward. We are free to select a method of shipment. If the customer has special requests, e.g. express delivery, express freight, etc., the additional cost will also be borne by the recipient. The goods will always be shipped at the buyer's risk. In case of damage during transport, the recipient must apply for an ascertainment of the facts with the respective forwarding agent.

Prices

All prices quoted by BEROLLKA will be quoted ex works, excluding any applicable statutory value added tax, public charges and levies such as potential customs tariffs, and insurance and packaging costs that will be charged separately.

Minimum quantity surcharge

If a small order is placed with a net value below 20 euros, an additional lump-sum minimum quantity surcharge of 15.00 euros will additionally be charged to cover administrative costs.

Delivery periods

The delivery dates and periods stated are subject to change. Any delay in delivery will not entitle the customer to withdraw from the contract. The assertion of damage claims is ruled out. BEROLLKA can make partial deliveries for justified reasons.

Wrong orders / return of merchandise

Any goods erroneously ordered will only be taken back provided we are responsible for such an error and the merchandise is in impeccable condition. In the event that an item erroneously ordered or merchandise that was returned is taken back in an individual case as a gesture of goodwill, attention must be paid to the following:

1. It is necessary to state the article number, the designation of the article, the date and number of the delivery note and the invoice number.
2. The return shipment must be carriage paid.
3. The merchandise must correspond to the article number(s) listed in the delivery note. In no event will any individual parts removed from an article be taken back.
4. The returned merchandise must be in impeccable condition.
5. Processing fees amounting to 10% of the net merchandise value not, however, less than €35.- will be charged for any instance where goods are taken back.

Deliveries effected more than 3 months earlier as well as sewing pieces and customized products cannot be exchanged. A copy of the delivery note or of the invoice will serve as proof in this regard.

Warranty and notices of defects

Excluding any further claims, the seller is liable for defects of the merchandise in that the seller will mend or replace all those parts free of charge that have become defective or unusable within 6 months due to provable defects of material or workmanship, and of which the seller has been notified in writing immediately after knowledge was gained of the defect. Should longer periods of time be required for statutory reasons, the warranty period is agreed to be 12 months.

The customer must grant the seller the required time and opportunity for effecting delivery of the replacement parts. Replaced parts will become the property of the seller. The buyer must check the merchandise without undue delay after it has been delivered by the seller and, if a defect is found, notify the seller of the defect without undue delay. Should the buyer fail to notify the seller, the merchandise will be deemed accepted unless a defect occurs that could not have been found when the merchandise was checked.

In any event in which a complaint is brought to the notice of BEROLLKA, BEROLLKA invariably has the right to inspect and check the item delivered that is the subject matter of said complaint, for which the customer will grant BEROLLKA the required time and opportunity. BEROLLKA can also demand that the customer sends the item delivered that is the subject matter of a complaint back to BEROLLKA at the expense of BEROLLKA. Should a notice of defects given by the customer turn out to be wilfully or gross negligently unfounded, the customer will be obliged to compensate BEROLLKA for any and all expenses incurred in this connection. BEROLLKA will be entitled to remedy warranty defects at its own discretion either by rectification of the defect or by replacement of the defective part or of the entire item delivered at no cost for the customer. Should the customer itself repair an item that gave rise to a complaint, or should the customer endeavour to rectify the defects itself without prior written approval of the seller, the seller will not grant any claims for reimbursement of expenses (such as transport, road, labour or material costs) or warranty claims.

The following are excluded from the liability for damages:

damage arising from natural wear and tear, damage due to negligent and/or improper treatment, or violent damage. A revocation of the purchase contract, a reduction of the purchase price or damage claims of any kind are ruled out. The warranty obligation will cease if repairs and/or modifications or replacements of individual parts are performed by other parties than the seller or the delivering company.

Reservation of title

With regard to resellers: The merchandise delivered will remain our property until any and all invoices – also for different deliveries – have been fully paid. Our merchandise that is hence still our property may only be resold within the ordinary course of business. Amounts receivable from any resale are assigned to us. We must be notified without undue delay and on the reseller's own initiative of any disposal of said merchandise outside of the ordinary course of business, e.g. through pledging our property or amounts receivable assigned to us.

Terms of payment

Any invoices are payable within 30 days from the date of the invoice strictly net cash. A 2% discount will be granted for payments received within 8 days provided any and all payment claims under earlier deliveries have been settled. We are entitled to charge default interest if the 30-day payment deadline is exceeded. The interest rate will be the customary market rate charged by major banks. We expressly reserve the right to effect delivery with cash on delivery.

Data Protection

BEROLLKA and the customer are responsible themselves for compliance with the respective data protection rules applicable to them. The customer must not provide BEROLLKA with personal data that is illegal or that the customer is not legally authorised to distribute.

The place of performance and the place of jurisdiction will be the place where the invoice is issued.

Heidelberg, Germany, will be the exclusive place of jurisdiction for any and all disputes arising under the contractual relationship. However, BEROLLKA is entitled to sue the customer at any other legal venue. These General Terms and Conditions of Business are governed by the laws of the Federal Republic of Germany, excluding of the United Nations Convention on Contracts for the International Sale of Goods of Goods (CISG).

Technical modifications and printing errors excepted. A current list of prices can be ordered at any time.