

Please contact DHCare EU Export Customer Services on :-
Tel : +33 247 626980 ; Email : serviceclient_export@dhcare.com
 Or visit our website : www.dhcare.eu.com



Rea® Azalea® Assist Minor

Manual folding steel chair



Description:

- Drum brakes.
- Seat tilt (+1°/24°) and backrest (0°/30°) by pneumatic cylinders.
- Seated tilt with translation.
- Width-adjustable legrests/stems.
- Heel straps.
- Frame colour: Pearl gray (RAL 0515).
- Weight: from 32 kg.
- Maximum user weight: 75 kg.
- 2-year warranty period excluding wear parts



Date :

Purchase Order :

Customer Code :

Distributor Stamp::

DHCare Unique Customization Service

Please specify the quote or special production number :

☐ Please make your selection with a cross in the box

☒ Standard feature

☐ Free of charge option

CODE		DESCRIPTION	NC : Non Compatible NC : 9999 = not compatible with the reference 9999	Euros
BNI0013	☐	Rea® Azalea® Assist Minor (< 16 ans)		

SEAT

SEAT WIDTH

BNI0038	☐	340 mm - adjustable from 340 to 440 mm		◆
---------	--------------------------	--	--	---

SEAT WIDTH EXTENSION (PRE ADJUSTMENT)

BNI0065	☐	Comes with adjustment + 50 mm (at 390 mm)		○
BNI0075	☐	Comes with adjustment + 100 mm (at 440 mm)	NC : 0110	○

Seat depth

BNI0095	☐	Adjustable from 380 to 450 mm		◆
---------	--------------------------	-------------------------------	--	---

SEAT HEIGHT

BNI0110	☐	400 mm		○
BNI0120	☐	450 mm		◆

SEAT ANGLE

BNI0127	☐	NEGATIVE SEAT ANGLE -1 TO +23° (STD +1 TO +24°)		○
---------	--------------------------	---	--	---

SEAT TYPE

BNI0130	☐	DEPTH ADJUSTABLE		◆
---------	--------------------------	------------------	--	---

SEAT CUSHION

BNI0165	☐	CUSHION FLOSHAPE		◆
---------	--------------------------	------------------	--	---

SEAT UPHOLSTERY

BNI0191	☐	BLACK PLUSH TR35		
BNI0201	☐	BLACK PU TR26		◆
BNI0202	☐	BLACK SPACER TR36		○

BACKREST

BACKREST TYPE

BNI0217	☐	FLEX 2 BACKREST (TENSION ADJ. AZALEA MINOR)		◆
BNI0218	☐	FLEX 2 BACKREST (FIXED, AZALEA MINOR)	NC : 0240	○

BACKREST COVER/ CUSHION

BNI0250	☐	CUSHION, LAGUNA (LATERAL SUPPORT)		◆
BNI0260	☐	CUSHION, MISTRAL 2 (WAIST SUPPORT)	NC : 0322	

BACKREST UPHOLSTERY

BNI0311	☐	BLACK PLUSH TR35		
BNI0321	☐	BLACK PU TR26		◆
BNI0322	☐	BLACK SPACER TR36		○



0165



0217



0218



0250



0260

PUSHHANDLES					 0340
BNI0340	☐	PUSH HANDLES BRACED		◆	
BNI0345	☐	PUSH BAR ANGLE ADJUSTABLE			
SEAT TILT AND BACKREST REGULATION TYPE					 0345
BNI0350	☐	ASSISTANT OPERATED GAS PISTON SEAT AND BACKREST TILT		◆	
BNI1263	☐	TILT LOCK (LOCKING OF PRE-SET ANGLE)		◆	
LEGRESTS / FOOTRESTS					
LEGRESTS / HANGERS LH TYPE					 0380 / 0430
BNI0380	☐	HANGER ANGLE ADJUSTABLE LH. WITH KNEE PAD AND CALF SUPPORT	NC : 0490, 0942		
BNI0390	☐	HANGER FIXED 80 DEG, SWING IN SWING OUT - ALU 2003, LH		◆	
BNI0410	☐	LEGREST AMPUTEE - ALU 2008 LH	NC : 0490, 0510, 0520, 0941, 0942		 0390 / 0440
LEGRESTS / HANGERS RH TYPE					
BNI0430	☐	HANGER ANGLE ADJUSTABLE RH. WITH KNEE PAD AND CALF SUPPORT	NC : 0550, 0942		
BNI0440	☐	HANGER FIXED 80 DEG, SWING IN SWING OUT - ALU 2003, RH		◆	 0410 / 0460
BNI0460	☐	LEGREST AMPUTEE - ALU 2008 RH	NC : 0550, 0570, 0580, 0941, 0942		
FOOTRESTS					
FOOTRESTS ONE PIECE					 0475 / 0476
BNI0475	☐	FOOTREST ONE PIECE DEPTH AND ANGLE ADJUSTABLE SHORT TUBE - NO WIDTH ADJUSTMENT	NC : 0490, 0510, 0520, 0550, 0570, 0580, 0954, 0956		
BNI0476	☐	FOOTREST ONE PIECE DEPTH AND ANGLE ADJUSTABLE LONG TUBE - NO WIDTH ADJUSTMENT	NC : 0490, 0510, 0520, 0550, 0570, 0580, 0954, 0956		
FOOTRESTS LH TYPE					 0510 / 0570
BNI0490	☐	FOOTREST LOWER PART, DEPTH AND ANGLE ADJUSTABLE WITH SHORT TUBE AND HEEL STRAP, LH		○	
BNI0510	☐	FOOTREST LOWER PART, FIXED WITH LONG TUBE AND HEEL STRAP, LH	NC : 0954, 0956	○	
BNI0520	☐	FOOTREST LOWER PART, DEPTH AND ANGLE ADJUSTABLE WITH LONG TUBE AND HEEL STRAP, LH		◆	 0520 / 0580
FOOTRESTS RH TYPE					
BNI0550	☐	FOOTREST LOWER PART, DEPTH AND ANGLE ADJUSTABLE WITH SHORT TUBE AND HEEL STRAP, RH		○	
BNI0570	☐	FOOTREST LOWER PART, FIXED WITH LONG TUBE AND HEEL STRAP, RH.	NC : 0954, 0956	○	 0520 / 0580
BNI0580	☐	FOOTREST LOWER PART, DEPTH AND ANGLE ADJUSTABLE WITH LONG TUBE AND HEEL STRAP, RH		◆	
ARMRESTS (Height adjustable, Autolock attachment included)					
ARMRESTS LH TYPE					 0600 / 0630
BNI0600	☐	STANDARD PAD, LH	NC : 1233, 1240	○	
BNI0605	☐	WIDE PAD, LH	NC : 1233	○	
BNI0606	☐	CONCAVE GUTTER ADJUSTABLE IN HEIGHT AND DEPTH OVE 70 MM WITHOUT TOOLS (COMPATIBLE WITH SHELF)	NC : 1232		 0610 / 0640
BNI0610	☐	PADDED, LH.	NC : 0680, 0690, 1230, 1231, 1233	◆	
BNI0615	☐	HEMIPLEGIC ARMRESTS, LH.	NC : 1230, 1231, 1233, 1240, 1245		
BNI0620	☐	WITHOUT ARMREST, LH	NC : 0690, 1230, 1231, 1240, 1245	○	 0615 / 0645
ARMRESTS RH TYPE					
BNI0630	☐	STANDARD PAD, RH	NC : 1232,1240	○	
BNI0635	☐	WIDE PAD, RH	NC : 1232	○	
BNI6036	☐	CONCAVE GUTTER ADJUSTABLE IN HEIGHT AND DEPTH OVE 70 MM WITHOUT TOOLS (COMPATIBLE WITH SHELF)	NC : 1232		 0615 / 0645
BNI0640	☐	PADDED, RH.	NC : 0680, 0690, 1230, 1231,1232, 1240, 1245	◆	
BNI0645	☐	HEMIPLEGIC ARMRESTS, RH.	NC : 1230, 1231, 1232, 1240, 1245		
BNI0650	☐	WITHOUT ARMREST. RH	NC : 0690, 1230, 1231, 1240, 1245	○	

SEATING OPTIONS				
ARMREST OPTIONS				
BNI0680	☐	SIDE SUPPORT COVER		
BNI0690	☐	CUSHION TO SIDE REST COVER		
SEAT OPTIONS				
BNI0710	☐	INCONTINENCE COVER		
BNI0720	☐	ABDUCTION CUSHION		
BACKREST OPTIONS				
BNI0730	☐	HEADREST, SIDEWISE ADJUSTABLE	NC : 0750, 0770	○
BNI0750	☐	HEADREST WITH CHEEK SUPPORT, SIDEWISE ADJUSTABLE	NC : 0740, 0770, 0780	
BNI0740	☐	COVER FOR HEADREST		
BNI0770	☐	NECKREST		◆
BNI0782	☐	UNIVERSAL HEAD/NECKREST HOLDER		
BNI0790	☐	HANDLES CLINIC KIT FOR HEAD-/NECKREST		
TRUNK SUPPORT LH				
BNI0800	☐	TRUNK SUPPORT MULTI ADJUSTABLE WITH STANDARD ARM (LH)		
BNI0880	☐	TRUNK SUPPORT SWINGAWAY (LH)		
TRUNK SUPPORT RH				
BNI0810	☐	TRUNK SUPPORT MULTI ADJUSTABLE WITH STANDARD ARM (RH)		
BNI0890	☐	TRUNK SUPPORT SWINGAWAY (RH)		
BACKREST OPTIONS (Following)				
BNI0915	☐	TILT SCALE, BACKREST RECLINE		
BNI0920	☐	WEDGES (ONE PAIR)		
BNI0954	☐	Bodypoint Open Buckle Ankle Supports, Side Adjustable, Medium (pair) - Comes loose with bindings		
BNI0956	☐	Bodypoint ankle braces open with buckle, adjustable on the side, Small (pair) - delivered loose with bindings		
LEGREST OPTIONS				
BNI0941	☐	COVER FOOTPLATE		
BNI0942	☐	CALF STRAP PADDED		
CHASSIS				
CHASSIS COLOUR				
BNI0944	☐	PEARL GREY		◆
BNI0946	☐	AZURITE		○
BNI0947	☐	SAND		○
BNI0948	☐	ELECTRIC GREEN		○
BNI0949	☐	HAPPY RED		○
CASTOR ATTACHMENT TYPE				
BNI0960	☐	Fixed		◆
FRONT CASTORS				
BNI0982	☐	5,5"x1.3/4" - 140X45 BLACK MEDIUM HARD		
BNI0992	☐	6" x 1.1/4" - 150X30 BLACK MEDIUM HARD		◆
REAR WHEELS				
REAR WHEELS SIZE				
BNI1047	☐	16" - 405 mm - with fixed rear wheel attachment		◆
REAR WHEELS TYPE				
BNI1076	☐	DRUM BRAKE (22/24" : COLORLESS HUB, CLEAR ANODIZED RIM & STAINLESS STEEL SPOKES / 16" : BLACK PLASTIC HUB & RIM)		◆
TYRE TYPE				
BNI1092	☐	TYRES SOLID BLACK 1" 3/8		◆
SAFETY				
BNI1170	☐	USER MANOUVERED TYRE BRAKES		
BNI1200	☐	PELVIC BELT - PLASTIC BUCKLE		◆
POSITIONING OPTIONS				
Think of Bodypoint positioning accessories: belts, harnesses, ankle supports.				




0720



07300750



07700782



0800 / 0810



0880 / 0890



0915



0954/0956



09440946



09470948



0949



1200



OPTIONS					
TABLE TRAY					
BNI1230	☐	TABLE TRAY INCL MOUNTING KIT	NC : 1237		  1230 1231
BNI1231	☐	TABLE TRAY INCL MOUNTING KIT, WITH ELBOW PADS	NC : 1237		
BNI1232	☐	HALF TABLE TRAY, CLEAR, FLIP UP, LH (ARMREST WITH STRAIGHT NARROW PAD INCL.)	NC : 1235		 1232 / 1233
BNI1233	☐	HALF TABLE TRAY, CLEAR, FLIP UP, RH (ARMREST WITH STRAIGHT NARROW PAD INCL.)	NC : 1235		 1237
BNI1240	☐	SWING AWAY TABLE	NC : 1235, 1237		
TABLE TRAY CUSHION					
BNI1235	☐	TABLE TRAY CUSHION			 1240
BNI1237	☐	PADDING			
OTHER OPTIONS					
BNI1300	☐	HARNESS ATTACHMENT			 1300
BNI1305	☐	LOCKING KIT CASTOR INCL. CAR FIXATION			
BNI1310	☐	ATTACHMENT WHEN USING THE PELVIC BELT FOR			
RAMPS					
1654128	☐	Basic ramp 2 telescopic parts in aluminum and non-slip coating - Folded dimensions: 1210 x 194 mm. 8.8 kg per pair - maximum capacity: 270 kg			 1654128
1654129	☐	Quick2go Compact 3-part aluminum telescopic ramp with non-slip coating Folded dimensions: 870 x 237 mm 13.6 kg per pair - maximum capacity: 260 kg			 1654129
1654130	☐	Aluminium folding mono ramp Folded dimensions: 1824 x 714 mm 13.9 kg - maximum capacity: 300 kg			 1654130