

Please contact DHCare EU Export Customer Services on :-
Tel : +33 247 626980 ; Email : serviceclient_export@dhcare.com
 Or visit our website : www.dhcare.eu.com



Rea® Azalea® Assist Max

Manual steel wheelchair with folding frame



Description:

- Seat tilt with maintenance of the center of gravity and the correct position of the leg rests: 0°/20° and of the backrest: 0°/30° by pneumatic jacks.
- Adjustable legrests from 0° to 80° and width-adjustable armrests.
- Coating: Black Dartex.
- Frame colour: Pearl gray (RAL 0515).
- Weight: from 54 kg.
- Maximum user weight: 180 kg.
- 2-year warranty period excluding wear parts



Date :

Purchase Order :

Customer Code :

Distributor Stamp:

DHCare Unique Customization Service

Please specify the quote or special production number :



☐ Please make your selection with a cross in the box

◆ Standard feature

○ Free of charge option

CODE		DESCRIPTION	NC : Non Compatible NC : 9999 = not compatible with the reference 9999	Euros		
BSB0015	☐	Rea [®] Azalea [®] Assist Max				
PROPULSION TYPE						
BSB0020	☐	TRANSIT: 16" WHEELS WITH HUB BRAKES		◆		
SEAT DIMENSIONS						
SEAT WIDTH						
BSB0050	☐	550 MM		○		
BSB0060	☐	610 MM		○		
SEAT WIDTH EXTENSION (PRE ADJUSTMENT)						
BSB0080	☐	+ 100 MM (ONLY 16" REAR WHEELS)		○		
SEAT DEPTH						
BSB0090	☐	530-600 MM		◆		
SEAT HEIGHT						
BSB0110	☐	450 MM		◆		
SEAT TYPE						
BSB0120	☐	DEPTH ADJUSTABLE		◆		
SEAT CUSHION						
BSB0130	☐	MAX SEAT CUSHION - SEAT UPHOLSTERY MANDATORY		◆		
SEAT UPHOLSTERY						
BSB0155	☐	BLACK PU TR26 - MUTLI STRETCH PU COATED FABRIC		◆		
BACKREST						
BACKREST TYPE						
BSB0160	☐	BACKREST, TENSION ADJUSTABLE - PRE-ADJUSTED SPACE UPWARDS		◆	 0160	
BACKREST COVER/ CUSHION						
BSB0170	☐	COVER, LATERAL - BACKREST UPHOLSTERY MANDATORY		◆	 0170	
BSB0180	☐	CUSHION, LAGUNA (LATERAL SUPPORT) - BACKREST UPHOLSTERY MANDATORY	NC : 0415	○		
BACKREST UPHOLSTERY						
BSB0215	☐	BLACK PU TR26 - MUTLI STRETCH PU COATED FABRIC		◆	 0180	
PUSHHANDLES						
BSB0220	☐	PUSH HANDLES BRACED		◆	 0220	
BSB0230	☐	HIGH PUSH BAR				
SEAT TILT AND BACKREST REGULATION TYPE						
BSB0240	☐	ASSISTANT OPERATED GAS PISTON SEAT AND BACKREST TILT. - TILT LOCK MANDATORY		◆		
BSB0605	☐	TILT LOCK (LOCKING OF PRE-SET ANGLE)		◆		

LEGRESTS / FOOTRESTS				
LEGRESTS / HANGERS TYPE				
BSB0260	☐	HANGER ANGLE ADJUSTABLE WITH KNEE PAD AND REINFORCED, SEPARATE CALF SUPPORT, PAIR		◆
FOOTRESTS TYPE				
BSB0280	☐	FOOTREST LOWER PART, DEPTH AND ANGLE ADJUSTABLE WITH SHORT TUBE AND FOOTBOARD CONVERTER, PAIR		○
BSB0290	☐	FOOTREST LOWER PART, DEPTH AND ANGLE ADJUSTABLE WITH LONG TUBE AND FOOTBOARD CONVERTER, PAIR		◆
ARMRESTS (Height adjustable, Autolock attachmentement included)				
ARMRESTS LH TYPE				
BSB0300	☐	REINFORCED, WIDE PAD, LH		◆
ARMRESTS RH TYPE				
BSB0320	☐	REINFORCED, WIDE PAD, RH		◆
SEATING OPTIONS				
ARMREST OPTIONS				
BSB0340	☐	SIDE SUPPORT COVER		
BSB0350	☐	CUSHION TO SIDE SUPPORT COVER		
BACKREST OPTIONS				
BSB0370	☐	HEADREST, SIDEWISE ADJUSTABLE	NC : 0400	○
BSB0390	☐	HEADREST WITH CHEEK SUPPORT, SIDEWISE ADJUSTABLE	NC : 0380, 0400	
BSB0380	☐	COVER FOR HEADREST		
BSB0400	☐	NECKREST		◆
OTHERS BACKREST OPTIONS				
BSB0415	☐	LATERAL POSITIONING PADS (WITH LOOP BAND BETWEEN COVER AND BACKREST UPHOLSTERY), PAIR		
BSB0420	☐	HANDLES CLINIC KIT FOR HEAD-/NECKREST		
BSB0430	☐	LUMBAR SUPPORT CUSHION		
BSB0440	☐	WEDGES (ONE PAIR)		
CHASSIS				
CHASSIS COLOUR				
BSB0441	☐	PEARL GREY		◆
BSB0442	☐	AZURITE		○
BSB0443	☐	SAND		○
BSB0444	☐	ELECTRIC GREEN		○
BSB0445	☐	HAPPY RED		○
CHASSIS SAFETY				
BSB0450	☐	CHASSIS EXTENDER, CASTORS 50 MM FORWARD.		
CASTOR ATTACHMENT TYPE				
BSB0460	☐	FIXED		◆
FRONT CASTORS				
BSB0495	☐	8"x2" - 200X45 BLACK SOLID TYRE		◆
FRONT CASTORS				
REAR WHEELS SIZE				
BSB0512	☐	16" / 405 MM - FIXED REAR WHEEL ATTACHMENT		◆
REAR WHEELS TYPE				
BSB0524	☐	DRUM BRAKE (22/24" : COLORLESS HUB, CLEAR ANODIZED RIM & STAINLESS STEEL SPOKES / 16" : BLACK PLASTIC HUB & RIM)		◆
TYRE TYPE				
BSB0535	☐	TYRES SOLID BLACK 1" 3/8		◆
SAFETY				
BRAKE TYPE				
BSB0560	☐	USER MANOUEVERED TYRE BRAKES		◆
SAFETY OPTIONS				
BSB0570	☐	PELVIC BELT - PLASTIC BUCKLE		◆
BSB0590	☐	REFLECTORS KIT (WHEELS AND CHASSIS)		
OPTIONS				
Think of Bodypoint positioning accessories: belts, harnesses, ankle supports.				



0370039004300441044204430444044505700590



OPTIONS				
BSB0600	☐	TABLE WITH ELBOW SUPPORTS		
RAMPS				
1654128	☐	Basic ramp 2 telescopic parts in aluminum and non-slip coating. Folded dimensions: 1210 x 194 mm. 8.8 kg per pair - maximum capacity: 270 kg		 0600  1654128
1654129	☐	Quick2go Compact 3-part aluminum telescopic ramp with non-slip coating Folded dimensions: 870 x 237 mm 13.6 kg per pair - maximum capacity: 260 kg		 1654129
1654130	☐	Aluminium folding mono ramp Folded dimensions: 1824 x 714 mm 13.9 kg - maximum capacity: 300 kg		 1654130